




# MaTRx Flo-Tech® Cushion- Fit Kit Instructions



Model #	Accessory	Description
FTSW	Seat Wedge (2pcs)	<ul style="list-style-type: none"> <li>Seat wedges attach under the front of the cushion base and may require trimming to fit the cushion width (see <i>trimming instructions</i>);</li> <li>Increases degree of pre-ischial contour angle to help prevent sliding</li> </ul>
FTHW	Hemi Wedge (1pc)	<ul style="list-style-type: none"> <li>Single wedge attaches under one side of the front of the cushion base and may require trimming to fit the cushion width (see <i>trimming instructions</i>);</li> <li>May be used to provide increased support to one lower extremity - useful for clients who use foot propulsion for mobility or hip flexion discrepancies</li> </ul>
FTTW-LFRR	Tapered Wedge -LFRR	<ul style="list-style-type: none"> <li>Tapered wedge attaches under left front or right rear of cushion base;</li> <li>Accommodates mild pelvic rotation, mild hip flexion discrepancies or increases adductor contour</li> </ul>
FTTW-RFLR	Tapered Wedge -RFLR	<ul style="list-style-type: none"> <li>Tapered wedge attaches under right front or left rear of cushion base;</li> <li>Accommodates mild pelvic rotation, mild hip flexion discrepancies or increases adductor contour</li> </ul>
FTAD-F	Adductor-Front	<ul style="list-style-type: none"> <li>Adductor build-up attaches to front lateral edge of cushion on top of cushion base;</li> <li>Provides increased adductor contour to thighs</li> </ul>
FTAD-R	Adductor-Rear	<ul style="list-style-type: none"> <li>Adductor build-up attaches to rear lateral edge of cushion on top of cushion base;</li> <li>Provides increased adductor contour to hips</li> </ul>
FTAB-L	Abductor-Large (2")	<ul style="list-style-type: none"> <li>Abductor build-up attaches over existing abductor contour on top of cushion base;</li> <li>Provides 2" increased height of abductor contour</li> </ul>
FTAB-S	Abductor-Small (3/4")	<ul style="list-style-type: none"> <li>Abductor build-up attaches over existing abductor contour on top of cushion base;</li> <li>Provides 3/4" increased height of abductor contour</li> </ul>
FTPO	Pelvic Obliquity (1")	<ul style="list-style-type: none"> <li>Pelvic Obliquity build-up attaches on top of cushion base into either side of ischial well;</li> <li>Accommodates obliquities up to 1" - for correctable obliquities use under low side of pelvis and for fixed obliquities use under high side of pelvis;</li> </ul> <p><b>Note:</b> Fluid Supplement may also be used to accommodate mild obliquities.</p>
FTFS	Fluid Supplement	<ul style="list-style-type: none"> <li>Fluid Supplements attach on top of cushion base into desired location in the ischial well;</li> <li>May be placed laterally to accommodate mild pelvic obliquities, sacral placement for additional sacral support/pressure reduction, bilaterally for increased ischial immersion</li> </ul>



# MaTRx Flo-Tech® Cushion- Fit Kit Instructions


The component parts of the modular MaTRx Flo-tech Fit Kit system can generate several variations of the cushion to promote active correction and management of posture where needed. All the add-on Flo-tech components are manufactured from high-density CMHR foam, and have been fitted with Velcro® to ensure that the parts remain in the prescribed position on the cushion even during repeated use.

## **Trimming Instructions:**

Components of the MaTRx Flo-tech Fit Kit can be trimmed to provide the optimal seating position for the end user. The seat wedges (front) and tapered wedges (side) may require trimming to fit the base properly. These parts are indented to indicate the size of the cushion and serve as a guide for trimming the part accurately. Use a utility knife to trim the components- see photos below for examples of trimming.

**⚠** - Please exercise caution when cutting or trimming foam to prevent personal injury.

The base of the Flo-tech cushion may also be trimmed to accommodate leg length discrepancies, or to provide additional hamstring space for knee flexion/contractures, or amputations.



## **Cleaning Instructions:**

Ideally the Matrix Flo-tech Cushion and Fit Kit Components should be stored in a dry, clean environment. The foam cushion/ components and fluid supplement may be wiped down carefully using a dampened cloth from a solution of hot water and mild neutral detergent.

**⚠** - Cleaning agents containing phenols or alcohols should not be used.



700 Ensminger Rd., Suite 112, Tonawanda, NY, USA 14150 TOLL FREE TEL: 1.888.433.6818 FAX: 1.888.433.6834  
84 Citation Dr., Units 1-7, Concord, ON, Canada L4K 3C1 TOLL FREE TEL: 1.866.748.7943 FAX: 905.695.0138