

Pediatric system:



Nicholas, an amazingly active 12 year old, was diagnosed with SMA Type 1 at birth.

Power mobility allows children with physical disabilities to enjoy experiences, activities and interactions that contribute significantly to their development — often for the first time. The modular, growable TRx Pediatric Power Positioning System provides positioning and pressure reduction, plus enhances independence, confidence and motivation. Features include:

- ▶ Small seat sizes growable from 11" to 15", smaller sizes available via RFQ*
- ▶ Pediatric armrests
- ▶ UltraLow/Low-Pro or TRx CG System
- ▶ Low seat-to-floor heights
- ▶ Lifetime structural warranty

Pediatric system

standard features & options

STANDARD FEATURES

- ▶ All systems come with seat assembly, seat pan, rigid back pan, height adjustable armrest and standard arm pads
- ▶ Weight capacities up to 175 lbs, depending on power base
- ▶ Growable seat sizes 11" - 15" (width and depth)
- ▶ Lifetime structural warranty

POSITIONING OPTIONS

- ▶ UltraLow/Low-Pro CG Tilt System (55°)
- ▶ Recline only System (174°)
- ▶ UltraLow/Low-Pro CG Tilt (55°) and Recline (174°) System
- ▶ Power Elevate module available with 6" - 7" elevation
- ▶ Optional Precline module available with Recline
- ▶ Extended Shear Reduction (ESR)

SEATING OPTIONS

- ▶ Kid*ab*ra cushions
- ▶ Headrests

OPTIONAL ACCESSORIES

- ▶ Pediatric Adjustable Height Arms - three sizes

ELECTRONICS

All electronics allow users to control their system through either a separate switch, toggle and/or their wheelchair drive controls


- ▶ Standard electronics included for Single Function system & Tilt/Recline systems
- ▶ Enhanced electronics available
- ▶ Electronics options (ie. Lazarus Jump Starter, M16 Splitter, M06 Attendant Power Switch, Stealth Egg Switch)

OTHER OPTIONS

- ▶ O2 Holders - Articulating and Non-Articulating
- ▶ Vent Trays - Articulating, Gimbal and Compact models - through RFQ

Available on:

3G Storm
TDX3,4,5
M91
M71
Power Tiger
S626/646
Freestyle F11/M11
Zippie 500
P222/220
1113/ATS
Dynamo/ATS
1121 - 1122
1103 - 1107
Q610
Q6000
Q600
Vibe
Other Base Manufacturers
Omega Track



Options



Available Product Feature Reference

LOW STF	TILT	LATERAL TILT	RECLINE	PRECLINE	POWER SLIDING BACK	ELEVATE	SPECIALIZED SYSTEM	SEATING
ARMRESTS	ACCESSORIES	VENTILATOR OPTIONS	O ₂ HOLDER	ELECTRONICS	ANTERIOR ASSIST			