

TerraTrek manual:

TerraTrek



Paul, an active and charismatic founder of Midland in Motions, was injured in an auto accident in 1987.

Take the road less traveled. Experience the freedom and agility of TerraTrek — the rugged, all terrain, manual wheelchair. TerraTrek is built to take whatever you dish out ... from the urban wilderness to the roads off the beaten track. Whether your adventures take you offroad, to the boardroom or down to the farm, the TerraTrek Urban and TerraTrek Outlook help you get there with pizzazz.

Features include:

- ▶ Rugged, durable, versatile and lightweight
- ▶ 8" (Urban) and 12.5" (Outlook) pneumatic front casters
- ▶ Available in seat sizes from 12" - 18" width and 12" - 17" depth

TerraTrek manual

standard features & options

STANDARD FEATURES

- ▶ Chair weight 36 lbs
- ▶ Available in standard red or green
- ▶ Seat widths from 12" to 18"
- ▶ Seat depths from 12" to 17"
- ▶ Sling Seat
- ▶ Adjustable Back Heights 8" to 12" or 13" to 17" in Black Ballistic Nylon
- ▶ 26" All-Terrain Wheels
- ▶ Anodized Aluminum Push Rims
- ▶ Quick Release Axles
- ▶ Wheel Cambers 0°, 3°, 6°
- ▶ Retractable Wheel Locks
- ▶ Bolt On Casters
- ▶ Bolt On Clothing Guards
- ▶ Standard, Wide or Pediatric Footrest

SEATING OPTIONS

- ▶ 6 optional colors
- ▶ Seat widths from 19" to 22"
- ▶ Seat depths from 18" to 21"
- ▶ 1" Horizontal Axle Extenders
- ▶ Lumbar Support Back Upholstery
- ▶ Push Handles
- ▶ Stroller Handles
- ▶ L-Style armrest
- ▶ Spoke Guards
- ▶ Vinyl Coated Push Rims
- ▶ 9° Rear Wheel Camber
- ▶ Quick Release Casters
- ▶ Anti Tipper Bars
- ▶ 2" Foam Cushion
- ▶ Seat Belts, Velcro
- ▶ Fender Style Clothing Guards

OTHER ACCESSORIES

- ▶ Fishing Rod Holder
- ▶ Rifle/Shot Gun Holder
- ▶ Nylon Back Pack
- ▶ Nylon Under Seat Pouch
- ▶ Crutch/Cane Holder
- ▶ Black Neoprene Frame Tube Protectors
- ▶ Two Position Legrest

TerraTrek
Urban



8" front casters

Lifetime structural warranty

TerraTrek
Outlook



12" front casters

Options



Camouflage cover



Spoke Guards



Two-position Footrest



Fishing Rod Holder