

UltraLow system:



Anita, a Research Assistant at a spinal cord rehab hospital, was injured in an auto accident at age 24.

Designed especially for the latest compact power bases, the innovative, modular UltraLow combines stylish good looks with some of the lowest possible seat-to-floor heights on a CG Tilt. Long on both form and function, this entire system can be moved from one power base to another with a simple interface change. Its distinctive slotted design facilitates quick and easy alteration of seat depth as well as fore and aft adjustment of the system on the power base. Features include:

- ▶ Up to 55° tilt
- ▶ Seat-to-floor heights available from 17.25" (depending on base)
- ▶ Increased CG shift of up to 7" for repositioning weight to remain centered on base while tilting
- ▶ Steel seat pan provides thinner, sleek look and 3" of on-chair depth adjustability
- ▶ Steel, angle cut back provides infinite up & down adjustability
- ▶ Legrest receivers adjust independently of system.

UltraLow system *standard features & options*

STANDARD FEATURES

- ▶ All systems come with a depth adjustable seat pan, height adjustable rigid back pan, height adjustable armrest, standard arm pads and Telescoping 70° Center Mount Front Rigging
- ▶ Weight capacities up to 300 lbs
- ▶ Seat Sizes 16" - 21" W x 15" - 21" D; Pediatric sizes from 11" x 15" W x 11" x 14" D
- ▶ Adjustable in 3" range
- ▶ Lifetime structural warranty

POSITIONING OPTIONS

- ▶ UltraLow System (55° tilt)
- ▶ Choice of several low Seat-to-Floor heights from 17.25" - 20.25", depending on base
- ▶ Power Recline available 174° with ESR
- ▶ Optional Precline module available with Recline
- ▶ Power Elevate module available with 6" - 7" elevation, depending on weight
- ▶ Optional Anterior Lift or Assist available

SEATING OPTIONS

- ▶ MaTRx-V cushions
- ▶ MaTRx PB Backs
- ▶ Contour Backs
- ▶ Onyx Package
- ▶ Headrests

NO CHARGE OPTIONS

- ▶ Single Post Back (Not available with Recline.)
- ▶ Armrest Storage Compartment - Glove Box
- ▶ Endomorph or mesomorph back canes can be added to achieve up to 4" difference between seat and back width at no charge.

OPTIONAL ACCESSORIES

- ▶ Fixed and Swing Away Laterals
- ▶ Contour Arm Troughs and Hand Pads
- ▶ Arm Pad Storage - Cup Holder

FRONT RIGGING OPTIONS

All power front rigging may increase STF by approximately 1"

- ▶ Center Mount Power Elevating & Articulating Foot Platform
- ▶ Pivot Plus Power Foot platform
- ▶ Pivot Plus Elevating Legrests - Power & Manual
- ▶ Fixed Center Mount Foot Platform
- ▶ Heavy Duty 70° Lift & Swing Away Footrests
- ▶ Split Power Elevating Articulating Foot Platform

ELECTRONICS

All electronics allow users to control their system through either a separate switch, toggle and/or their wheelchair drive controls

- ▶ Standard electronics included for Tilt and Recline systems
- ▶ Enhanced electronics available
- ▶ Electronics options (ie. Lazarus Jump Starter, M16 Splitter, M06 Attendant Power Switch, Stealth Egg Switch)

OTHER OPTIONS

- ▶ O₂ Holders - Articulating and Non-Articulating
- ▶ Vent Trays available through RFQ

Available on:

Storm (all)
TDX (all)
M71
M11
M91
F11
P220
P222
Q600
Q610
Q6000
1103
1107
1113 ATS
Dynamo ATS
1122
1402



No Charge Options



Available Product Feature Reference

ANTERIOR ASSIST	ANTERIOR LIFT	TILT	LATERAL TILT	RECLINE	PRECLINE	POWER SLIDING BACK	ELEVATE	SPECIALIZED SYSTEM
SEATING	ARMRESTS	FRONT RIGGING	ACCESSORIES	VENTILATOR OPTIONS	O ₂ HOLDER	LOW STF	ELECTRONICS	