

X3 maxX

the COMPLETE SOLUTION

The all new ROVI X3 Power Base is the perfect complement to Motion Concepts' full line of innovative power positioning products! Available exclusively through Motion Concepts, the innovative and exciting ROVI X3 boasts the narrowest wheelbase in the industry and all the function and durability you have come to expect from Motion Concepts.



Follow us on . . .




Toll Free Tel: 1.888.433.6818 ► Toll Free Fax: 1.888.433.6834 ► www.motionconcepts.com

ULTRA LOW *maxX* SYSTEM VITALS

- Industry's most modular system, featuring up to 55° of CG Tilt, 170° of Recline with Extended Shear Reduction (ESR) and up to 12" of nearly silent Power Seat Elevation
- Available with the full range of Motion Concepts' accessories and Invacare® Matrix® Seating products, mounted via our all new Ultra Rail and Ultra Track slotted interfaces
- 300 lb. weight capacity
- Seat-to-Floor Heights as low as 17" with 12" of Power Seat Elevation
- 16" - 24" of seat width and 15" - 22" of seat depth adjustment

WIDTH ADJUSTMENT RANGES	DEPTH ADJUSTMENT RANGES
16" - 20", 19" - 22", 23" - 24"	15" - 19", 17" - 21", 19" - 22"



X3 BASE VITALS

HCPCS CODES	K0856, K0861
SUSPENSION SYSTEM	Active Ride Control (ARC)
BASE WIDTH	23.25"
BASE LENGTH	36.5"
DRIVE WHEEL SIZE ⁽¹⁾	14"
CASTER WHEEL SIZE (front/rear)	6"
TURNING RADIUS	20"
DRIVETRAIN	4-pole EAD Motors
BATTERIES	2 x M34 3 x M40 option for Vents
BATTERY CHARGER	8A, Off Board
BATTERIES WEIGHT (each)	42 lbs
RANGE ⁽²⁾	17.7 miles
WEIGHT CAPACITY	300 lbs
BASE WEIGHT (w/o batteries)	157 lbs
GROUND CLEARANCE	3"
MAXIMUM SPEED	5 mph

⁽¹⁾ All X3 wheels come with solid inserts

⁽²⁾ Actual driving range and speed may vary due to factors such as user weight, type and grade of terrain, battery condition & charge level; type and condition of the drive wheels, and variations in the drive and control systems.

X3



X3 BASE SHROUD COLORS*



BLIZZARD WHITE MATTE BLACK SOLAR YELLOW MOTION GREEN ORIGINAL ORANGE BRIGHT BLUE RED DELICIOUS PARTY PINK

* Wheel rims are always black