



## ELEVATED THINKING . . .

Our **NEW F-Series** takes access to a higher level! Totally redesigned for 2013, the affordable F-Series is proof that you don't need to sacrifice function for economy.

Full featured performance without compromise, each F- Series is now available with up to 12" of seat elevation with no increase in seat to floor height. Our patented, adjustable width, drop seat design ensures you get all the flexibility you are accustomed to with Motion Concepts.

- Seamlessly integrates with virtually any power base in the industry
- Fully modular, featuring up to 50° of CG Tilt, 168° of Recline and 12" of nearly silent Power Seat Elevation
- Now fully width and depth adjustable
- Available with the full range of Motion Concepts' accessories and Invacare® Matrix® Seating
- 250 lb. weight capacity
- 14" - 20" of seat width and 16" - 20" seat depth adjustment
- Seat-to-Floor Heights as low as 16¾" depending on power base.

POWER BASE	MINIMUM STF HEIGHT TILT ONLY	MINIMUM STF HEIGHT TILT/RECLINE	MINIMUM STF HEIGHT TILT/RECLINE /12" ELEVATING LIFT SEAT
Invacare® TDX SP*	16¾"	16¾"	16¾"
Quickie® QM710	17¼"	17¼"	17¼"
Pride® Quantum Q6 Edge	16¾"	16¾"	16¾"

\*Without Transport Ready Option (TRO)



[www.motionconcepts.com](http://www.motionconcepts.com)

*usa*  
 700 Ensminger Rd. Suite 112  
 Tonawanda, NY 14150  
 Toll Free Tel: 1.888.433.6818  
 Toll Free Fax: 1.888.433.6834

*canada*  
 84 Citation Dr.,  
 Concord, ON L4K 3C1  
 Toll Free Tel: 1.866.748.7943  
 Tel: 905.695.0134  
 Fax: 905.695.0138