

Back into a perfect fit

At Motion Concepts, we know that clients come in all shapes and sizes. We understand how frustrating it can be to try to fit uniquely shaped clients into standard seat sizes. We also know that in the real world, not everyone fits into a standard mold — but we believe this shouldn't prevent anyone from enjoying secure, comfortable and proper seating in their wheelchair. And that's why we're thrilled to introduce our new **BodyFit** system.

Our **BodyFit** system is specially designed for clients who need a wide seat pan, combined with a more narrow back. It's the ideal solution for heavier clients with narrower shoulders and larger hips — often referred to as endomorphs. The **BodyFit** accommodates this specific body type by securing the position of clients' arms inside the seat and providing continued lateral support.

How? The **BodyFit** is available in back widths 2" or 4" narrower than the seat width. For example, a 16" back can be combined with an 18" or 20" seat to achieve a more precise client fit. As well, it allows clients to shift the back pan laterally to accommodate orthopedic abnormalities.

But that's not all! This unique system picks up where others leave off. With standard wide or bariatric systems, the armrests are usually left hanging outside of the wheelchair, leaving seats over 24" unable to fit through most normal doors. With the **BodyFit** system, not only are clients' arms securely and naturally tucked inside, they no longer need large, bulky offset brackets on lateral supports and/or other custom modifications. Plus, placement of the armrests closer to clients' bodies and in proper position not only assists in pressure reduction, but simplifies repositioning and transfers. These revised seat side frames also bring the back-canes in towards the center, maintaining a smaller profile.

New recline armrests are available, as are Single Post Flip back arms which can be added to a tilt system and provide optimum upper extremity support — even with a seat pan up to 6" wider than the backrest. Single Post Flip back arms are ideally suited to the **BodyFit** system because they don't interfere with clients' legs when they recline, and they eliminate the need for fancy I-Back cutouts when fixed lateral are required. Armrest spacers are also available on all tilt arms, which allow up to 6" difference between the back and seat.

So, help your clients of all shapes and sizes get a perfect fit today with the new **BodyFit** system from Motion Concepts.

✓ Perfect Fit



Motion Concept

US: 700 Ensminger Rd.,
Suite 112,
Tonawanda, NY 14150
Toll Free Tel: 1.888.433.6818
Toll Free Fax: 1.888.433.6834

Canada: 84 Citation Dr.,
Concord, ON L4K 3C1
Tel: 905.695.0134
Fax: 905.695.0138

www.motionconcepts.com