

August 27, 2004

Motion Concepts

ON A ROLL
A weekly e-bulletin from Motion Concepts

www.motionconcepts.com
U.S.
700 Ensminger Rd.,
Suite 112
Tonawanda, NY 14150
TEL 1.888.433.6818
FAX 1.888.433.6834
Canada
84 Citation Dr.,
Concord, ON
L4K 3C1
TEL 905.695.0134
FAX 905.695.0138

Sizes for the real world



Despite predominant media images, most of us know that in the real world, a significant number of people are larger than average size. For those that also need a power positioning system, weight capacity is always an issue. At Motion Concepts, our bariatric tilt system has always accommodated the needs of bariatric clients. But now, our tried and true bariatric system is even better!

Our new, Low-Pro Bariatric Power Positioning System provides a lower seat-to-floor height, along with better mechanical advantage. This means it can tilt easier — with less stress on the actuator and seat frame. Best of all, it accommodates up to 500 pounds standard, and even more via our Request For Quotation (RFQ) services.

With standard seat widths up to 26" on selected bases, the bariatric Low-Pro seating system from Motion Concepts provides flip-back, removable arms for added accessibility. In fact, we recently installed a customized 30" wide system and even at this size, the removable arms allowed it to fit easily through standard size doorways!

This new bariatric system also boasts a full spectrum of reinforced parts and front rigging designed especially to manage the additional weight requirements of bariatric clients, including an exclusive bariatric Center Mount Foot Plate. It provides a lifetime limited warranty so you can be sure it will stand up to tough situations.

A Single Pivot Low-Pro system is also available for even lower seat-to-floor heights and provides 30° tilt.

So, give your clients of all sizes the best power positioning system on the market. Check out the new bariatric Low-Pro from Motion Concepts today.



Front rigging designed especially to manage the additional size and weight requirements of bariatric clients.