

May 21, 2004

**Motion Concepts**

**ON A ROLL**

A weekly e-bulletin from Motion Concepts

## Now ... North America's best tilt is easier than ever to get funded!

With recent changes to funding codes, the reality is that these days, therapists must sometimes navigate through precarious and challenging paperwork in order to ensure their clients' positioning and mobility needs are properly met. For clients who need tilt, securing appropriate funding can sometimes feel like a job in itself. And that's where Motion Concepts can help.

Motion Concepts creates world class power positioning systems. From our classic TRx and SuperLow systems, to our bariatric and pediatric solutions, every power positioning system we create delivers unparalleled features, functionality and visual appeal. Sometimes however, in certain challenging funding climates, our entire system may not be feasible. And it's precisely for these situations that we've come up with a new way for your clients to get all the benefits of our amazing CG tilt. It's our new **Low-Pro CG Tilt Interface Kit**.

The Low-Pro CG Tilt Interface Kit utilizes most of the major base manufacturers' K-11 base/rehab seat kits, **but** adds our superior, low STF, CG tilt mechanism. So, if your base manufacturer's seat works for your client's application, but they also need a quality tilt that offers fore and aft adjustability in order to balance their weight on the base, you now have a uniquely affordable alternative designed especially to be funding code applicable! Our Low-Pro CG Tilt Interface Kit uses your base manufacturer's rehab seat to provides a sturdy, CG tilt for most of the industry's K-11 power bases including:

- ▶ **Quickie** Freestyle M11 and F11
- ▶ **Quantum** 1121 and 1122
- ▶ **Invacare** TDX-3 and M91

Each kit includes an actuator, standard electronics and all the accompanying hardware needed to add tilt to your application. Best of all, it can be either stocked at your facility or easily installed at our office. So, it's never been easier to enhance the quality of your clients' lives.

For many of today's cost conscious situations, it's easy to see how the Low-Pro CG Tilt Interface Kit can truly mean the difference between pain and pressure relief for your clients.

**U.S.**  
700 Ensminger Rd.,  
Suite 112  
Tonawanda, NY 14150  
TEL 1.888.433.6818  
FAX 1.888.433.6834

**Canada**  
84 Citation Dr.,  
Concord, ON  
L4K 3C1  
TEL 905.695.0134  
FAX 905.695.0138

[www.motionconcepts.com](http://www.motionconcepts.com)

### Each Low-Pro CG Tilt Interface Kit provides:

- ▶ Better, sturdier CG tilt that fits most manufacturer's K-11 base/rehab seat kits
- ▶ 55° tilt
- ▶ Low seat-to-floor heights
- ▶ Funding code applicable
- ▶ 250 lb weight capacity
- ▶ Accommodate seat sizes from 15"-22"\*
- ▶ fore and aft adjustability
- ▶ Better mechanical advantage
- ▶ Lifetime structural warranty
- ▶ Hassle-free and installed by our trained technicians at our facility
- ▶ North American built
- ▶ Can be stocked

\*dependant upon base

