

Clean design eliminates the cluttered look of dual back canes.



Back angle adjustment (closed).



Back angle adjustment (open).



Single and loving it

Although it may seem like it happened suddenly, Motion Concepts' foray into the single scene was actually a long time coming. And now that we've finally taken the plunge, we've never looked better — thanks to our new Single Post Back.

Make no mistake about it — our Single Post Back is hot! In fact, no matter which side you see your clients from, the dreary rehab look is all but gone and forgotten. Available on most tilt and/or elevate systems, our new automotive-style Single Post Back is available in 16" - 20" seat widths, with a Motion Concepts Contour Back.

It can accommodate up to a 250 lb weight capacity, dependent upon the particular seating system it's mounted to. Plus, our Single Post Back can even work with really low seat-to-floor height systems. Contact your local rep or dealer for more information.

Singled out

But we all know that beauty is only skin deep. As far as substance goes, the Single Post Back is much more than just a pretty face. It manually releases to adjust to various different angles — in 6° increments, from vertical to 24° forward and 18° back — so your clients are always comfortable and secure. The Single Post Back also presents a better turning radius with the absence of protruding back canes and features a novel integrated adjustment handle that makes back angle adjustments a snap. As well, it can actually be adjusted while clients are in their chair.

The height of the single post back is even adjustable too, from 21" - 23", in 1" increments. It folds forward to less than 11" off the seat pan for easier than ever transport and fits both regular and Onyx style headrests.

As you can see, our Single Post Back is truly one hot number. So, go ahead and see for yourself. Venture out into the single scene today with a snazzy, new Single Post Back from Motion Concepts.

U.S.
700 Ensminger Rd.,
Suite 112
Tonawanda, NY 14150
TEL 1.888.433.6818
FAX 1.888.433.6834

Canada
84 Citation Dr.,
Concord, ON
L4K 3C1
TEL 905.695.0134
FAX 905.695.0138

www.motionconcepts.com

Back folds down for easy transportation.

