


FOOT PLATE & FOOT PLATFORM OPTIONS:

INDIVIDUAL FOOT PLATE OPTIONS FOR INDEPENDENT ELEVATING & ARTICULATING LEGRESTS (CENTER MOUNTED)



- > Foot shaped with heel loop
- Size: 6" W x 11" D (each)

MULTI-AXIS FOOTPLATE MOUNT



- > Solid clamping mechanism securely locks footplate in virtually any position
- > Low-profile clamp is adjusted with one tool
- > Able to capture complex positions of foot and ankle
- > Interfaces with all Motion Concepts' footrests

FLIP-UP FOOT PLATFORM OPTIONS FOR POWER & MANUAL CENTER MOUNT FRONT RIGGING



Small: 12" W x 8.5" D



Medium: 12" W x 10.75" D



Large: 17" W x 14.75" D



Extra Large: 21" W x 14.75" D

INDIVIDUAL, CAST ALUMINUM FOOT PLATE OPTIONS FOR CMB, CMX, LNX AND FIXED & ADJUSTABLE 70-90° CENTERMOUNT PLATFORMS



- > Built in heel cups & rubber mat
- Small: 5" W x 7.5" D (each)




- > Built in heel cups & rubber mat
- Medium: 5.5" W x 9.5" D (each)




- > Built in heel cups & rubber mat
- Large: 6" W x 11.5" D (each)

FOOT PLATE OPTIONS FOR ELEVATING AND TILT FRONT RIGGING



Small: 4.25" W x 8" D (each)




Medium: 5.25" W x 8" D (each)




Large: 6.25" W x 8" D (each)




Extra Large: 9" W x 11" D (each)



Small: 7.25" W x 9" D (each)



Medium: 7.5" W x 11" D (each)



Large: 10" W x 14" D (each)