

**REQUIRED INFORMATION**


Quote  Order

Date of Order: \_\_\_\_\_ Dealer Account #: \_\_\_\_\_  
 Dealer Name: \_\_\_\_\_  
 PO #: \_\_\_\_\_ Tag \_\_\_\_\_  
 Purchasing Contact: \_\_\_\_\_  
 Phone: \_\_\_\_\_ Fax: \_\_\_\_\_  
 E-mail: \_\_\_\_\_  
 RTS/Therapist: \_\_\_\_\_  
 Ship to Address: \_\_\_\_\_  
 City: \_\_\_\_\_ Province: \_\_\_\_\_  
 Postal Code: \_\_\_\_\_ Client Gender:  M  F  
 Special Client Conditions: \_\_\_\_\_

**REQUIRED MEASUREMENTS**

To ensure system is accurately configured please fill in all required

Client Height:	
Client Weight (lbs):	
A. Seat to Shoulder:	
B. Trunk Depth:	
C. Chest Width:	
D. Knee to Back:	
E. Seat to Top of Head:	
F. Elbow to Hand:	
G. Seat to Elbow:	
H. Hip Width:	
I. Knee to Heel:	
Cushion Thickness:	



Customer Service: 1.866.748.7943 | Fax: 1.905.695.0583 | [orderscdn@motionconcepts.com](mailto:orderscdn@motionconcepts.com)

**ROVI**  
**BASE STATISTICS**

ADP CODES	WP0001287
SUSPENSION SYSTEM	Active Ride Control (ARC)
BASE WIDTH X <sub>3</sub>	23.25"
BASE WIDTH A <sub>3</sub>	25.25"
BASE LENGTH	36.5"
DRIVE WHEEL SIZE <sup>(1)</sup>	14"
CASTER WHEEL SIZE (front/rear)	6"
TURNING RADIUS X <sub>3</sub>	20"
TURNING RADIUS A <sub>3</sub>	20.5"
DRIVETRAIN	4-pole EAD Motors

BATTERIES	2 x M34
BATTERY CHARGER	8A, Off Board
BATTERIES WEIGHT (each)	38.5 lbs (M34) / 17.5 kg
RANGE <sup>(2)</sup>	17.7 miles / 28.5 km/h
WEIGHT CAPACITY	300 lbs / 136 kg
BASE WEIGHT X <sub>3</sub> (w/o batteries)	157 lbs / 71 kg
BASE WEIGHT A <sub>3</sub> (w/o batteries)	159 lbs / 72 kg
GROUND CLEARANCE	3"
MAXIMUM SPEED	6.2 mph

**X<sub>3</sub>**



**A<sub>3</sub>**



\*Shown with Motion Concepts Single Post Back for Contour Backrest. See page 7

<sup>(1)</sup> All X3 and A3 wheels come with solid inserts

<sup>(2)</sup> Actual driving range and speed may vary due to factors such as user weight, type and grade of terrain, battery condition & charge level; type and condition of the drive wheels, and variations in the drive and control systems.

**BASE OPTIONS**

- 888WNLEX3 **ROVI X3<sup>(1)(2)</sup>** . . . . . \$6,800
- 888WNLEA3 **ROVI A3<sup>(1)(2)</sup>** . . . . . \$6,800
- HSMTR **6.2 MPH High Speed Motors<sup>(1)(2)</sup>** . . . . . \$995  
*Must order OnTraxx Enhanced Driving Module*
- OEDM **OnTraxx Enhanced Driving Module** . . . . . \$1,175

1. Price of the powerbase includes: Active Ride Control suspension, tie-down hardware for 4-point strap system, 8 amp off-board charger, flat free inserts and positioning belt. Weight capacity is 300 lbs.  
2. When selecting an increased speed motor package there is a reduction in torque output by roughly 15%.

**FENDER COLOUR**

- 800914-00200-03  **Blizzard White**
- 800914-00200-69  **Solar Yellow**
- 800914-00200-23  **Original Orange**
- 800914-00201-21  **Red Delicious**
- 800914-00201-2C  **Party Pink**
- BCF  **Black Carbon (Matte)**
- BPF  **Black Panther (Gloss)**
- PRF  **Purple Reign**
- IBF  **Indigo Blue**
- MGF  **Motion Green**

**FENDER COLOUR W/ LIGHTS<sup>(1)</sup>**

**SPECIAL OFFER**

- 800914-00201-03  **Blizzard White** . . . . . **Now \$700**
- 800914-00201-69  **Solar Yellow** . . . . . **Now \$700**
- 800914-00201-23  **Original Orange** . . . . . **Now \$700**
- 800914-00201-21  **Red Delicious** . . . . . **Now \$700**
- 800914-00201-2C  **Party Pink** . . . . . **Now \$700**
- MGFL  **Motion Green** . . . . . \$995
- BCFL  **Black Carbon (Matte)** . . . . . \$995
- BPFL  **Black Panther (Gloss)** . . . . . \$995
- IBFL  **Indigo Blue** . . . . . \$995
- PRFL  **Purple Reign** . . . . . \$995

1. Must select compatible Joystick (D51108, D50982, D51711)

**WHEELCHAIR TRANSPORT BRACKETS**

- XBKTS **Wheelchair Transport Brackets(1)** . . . . . STD
1. For unoccupied use only

**TIRE OPTIONS**

- STIRES **Standard Solid Tires** . . . . . STD
- PTIRES **Pneumatic Tires** . . . . . N/C

**BATTERIES**

- M34 SLD G **MK Group Battery x 2** . . . . . \$650

**JOYSTICK SELECTION**

- D51122 **R-Net LED Joystick<sup>(1)</sup>** . . . . . N/C
- D51108 **R-Net LED Joystick for Lights** . . . . . \$500
- D51165 **R-Net Color Remote Joystick with Monojacks<sup>(1)</sup>** . . . . . \$875
- D50982 **R-Net Color Remote Joystick with Monojacks for Lights** . . . . . \$975
- D51711 **R-Net Advanced Joystick<sup>(2)</sup>** . . . . . \$1,595

1. Light package not available with this joystick  
2. 3.5" high-resolution VGA colour display, IR Control as standard, includes integrated Bluetooth Module.

**JOYSTICK MOUNT**

- 800914-001-16-ML **Fixed Joystick Mounting Left** . . . . . N/C
- 800914-001-16-MR **Fixed Joystick Mounting Right** . . . . . N/C
- SAQL **Height Adjustable Swing-Away Quad Link Left** . . . . . \$344
- SAQR **Height Adjustable Swing-Away Quad Link Right** . . . . . \$344

**ELECTRONIC OPTIONS**

- SPH **Smart Phone Holder** . . . . . \$50
- MUSB **Motion Dual USB Charger** . . . . . N/C
- D50612 **R-Net PC Programmer with Dongle** . . . . . \$603
- D50996 **Diagnostic Test Tool (DTT) with VR2 & RNet Cables** . . . . . \$1,876

**ADDITIONAL ACCESSORIES**

- UCH **Universal Cup Holder<sup>(1)</sup>** . . . . . \$50

1. Not available on the joystick side.

**SEAT-TO-FLOOR HEIGHT<sup>(1)</sup>**

- STF17 **17" Seat-To-Floor** . . . . . N/C
- STF18 **18" Seat-To-Floor** . . . . . N/C
- STF19 **19" Seat-To-Floor** . . . . . N/C
- STF20 **20" Seat-To-Floor** . . . . . N/C

1. Seat-to-floor height is measured at the front of the seat pan based on a standard seat size. Seat-to-floor height may vary +/- .25" dependent on condition of the power base and seat depths greater than 18"

**REHAB SEAT**

RSMAXX Ultra Low Maxx Rehab Seat . . . . . \$2,495

**ELEVATING LIFT SEAT MODULE**

MD12L 12 " Elevating Lift for Rehab Seat<sup>(1)</sup> . . . . . \$1,995  
Weight capacity up to 250lbs

SFCB Single Function Electronics<sup>(2)</sup> . . . . . STD

Select Switch:

- SPBO Single Push Button
- STO Single Toggle
- DPBO Dual Push Button
- SOS Omit Switch
- JOY Through Joystick

Select Switch Mounting:

- ESML Left Armpad
- ESMR Right Armpad
- ESMJL Left Joystick (Mounts to Joystick Bracket)
- ESMJR Right Joystick (Mounts to Joystick Bracket)
- ESMD Do not mount

Optional:

RDE Reduced Drive for Elevating Seat<sup>(3)</sup> . . . . . \$400

1. Not available on seat sizes greater than 20" wide by 20" deep.
2. Required when selecting 12" Elevating Lift. These electronics include drive lock-out.
3. If reduced drive is not selected the system will be in drive lock-out when elevated.

**ADD: MODULES**

PRETILT 5° Fixed Pre-Tilt<sup>(1)</sup> . . . . . N/C

1. This Fixed Pre-tilt option adds approximately 2" to the front seat to floor height

**ADULT SEAT WIDTH SETTINGS**

- SW116 16" Wide (Adjustable from 16" to 20") . . . . . N/C
- SW117 17" Wide (Adjustable from 16" to 20") . . . . . N/C
- SW118 18" Wide (Adjustable from 16" to 20") . . . . . N/C
- SW119 19" Wide (Adjustable from 16" to 20") . . . . . N/C
- SW120 20" Wide (Adjustable from 16" to 20") . . . . . N/C
- W219 19" Wide (Adjustable from 19" to 22") . . . . . \$460
- W220 20" Wide (Adjustable from 19" to 22") . . . . . \$460
- W221 21" Wide (Adjustable from 19" to 22") . . . . . \$460
- W222 22" Wide (Adjustable from 19" to 22") . . . . . \$460
- BW123 23" Wide (Adjustable from 20" to 24")<sup>(1)</sup> . . . . . \$460
- BW124 24" Wide (Adjustable from 20" to 24")<sup>(1)</sup> . . . . . \$460

1. This option is only available on the Rovi A3 Base.

**ADULT ULTRA RAIL SEAT DEPTH SETTINGS**

- SD116 16" Deep (Adjustable from 16" to 21") . . . . . N/C
- SD117 17" Deep (Adjustable from 16" to 21") . . . . . N/C
- SD118 18" Deep (Adjustable from 16" to 21") . . . . . N/C
- SD119 19" Deep (Adjustable from 16" to 21") . . . . . N/C
- SD120 20" Deep (Adjustable from 16" to 21") . . . . . N/C
- D219 19" Deep (Adjustable from 19" to 23") . . . . . \$460
- D220 20" Deep (Adjustable from 19" to 23") . . . . . \$460
- D221 21" Deep (Adjustable from 19" to 23") . . . . . \$460
- D222 22" Deep (Adjustable from 19" to 23") . . . . . \$460
- D223 23" Deep (Adjustable from 19" to 23") . . . . . \$460

**LAP BELTS**

- BELT48 Push Button Style Lap Belt, Length 48" . . . . . STD
- BELT60 Push Button Style Lap Belt, Length 60" . . . . . N/C
- BELT71 Push Button Style Lap Belt, Length 71" . . . . . N/C
- PLB2 Bodypoint Padded Lap Belt 2 Point<sup>(1)</sup> . . . . . \$250
- PLB4 Bodypoint Padded Lap Belt 4 Point<sup>(1)</sup> . . . . . \$315

1. 67" length only option

**FINISHED BACK HEIGHT - NON RECLINE SYSTEMS<sup>(1)(2)</sup>**

- 18" FBT18
- 19" FBT19
- 20" FBT20
- 21" FBT21
- 22" FBT22
- 23" FBT23
- 24" FBT24
- 25" FBT25

1. 4" Gap from the seat pan to the bottom of back pan is recommended.
2. Finished back height is measured from seat pan to top of backrest.

**BACK CANES**

- SBC18 18" Straight Back Canes . . . . . STD
- SBC20 20" Straight Back Canes . . . . . STD
- SBC22 22" Straight Back Canes . . . . . STD
- SBC24 24" Straight Back Canes . . . . . STD

Select Cane Angle

- 75° BA75
- 79° BA79
- 85° BA85
- 90° BA90
- 95° BA95
- 101° BA101
- 106° BA106
- 112° BA112
- 116° BA116

- ABC20 20" Angled 5° Back Canes . . . . . N/C
- ABC22 22" Angled 5° Back Canes . . . . . N/C
- ABC24 24" Angled 5° Back Canes . . . . . N/C

Select Cane Angle

- 81° BA81
- 84° BA84
- 91° BA91
- 96° BA96
- 101° BA101
- 106° BA106
- 111° BA111
- 117° BA117
- 121° BA121

**BACK OPTIONS**

- SRB Standard Rehab Back & Cushion . . . . . \$125
- ORCB Omit Rigid Back & Cushion . . . . . N/C

**ARMRESTS**

- DP1 Adult Dual Post 9.75" - 12.75" . . . . . STD
- DP2 Adult Dual Post 12.75" - 15.75" . . . . . STD
- DP3 Adult Dual Post 15.75" - 18.75" . . . . . STD

Select:

- FP Flat Pouch . . . . . STD
- GB Glove Boxes . . . . . N/C

Optional:

- OAL Outback Arm Option Left<sup>(1)</sup> each . . . . . \$195
- OAR Outback Arm Option Right<sup>(1)</sup> each . . . . . \$195

- RMCAN Ultra Rail Seat Mnt Flip Back Cantilever  
MAXX Tilt Arm 9" - 12.5"<sup>(2)</sup> . . . . . \$488

Optional:

- DPM Dual Post Module to RMCAN (pair) . . . . . \$250

- OCTA1 Outback Cantilever Tilt Armrest 9"-12"<sup>(3)</sup> . . . \$488
- OCTA2 Outback Cantilever Tilt Armrest 12"-15"<sup>(3)</sup> . . \$488

1. Only available with Dual Post Armrests.  
2. Mounts to Seat Rail. Slotted upper tube for infinite fore/aft adjustment of arm pad.

**ARMREST TUBE LENGTH - PLEASE SPECIFY**

- ATFL Left Full . . . . . STD
- ATDL Left Desk . . . . . STD
- ATFR Right Full . . . . . STD
- ATDR Right Desk . . . . . STD

**ARMPADS**

**Standard**

- SAPFL Left Full (2.25" x 14") . . . . . STD
- SAPDL Left Desk (2.25" x 10") . . . . . STD
- SAPFR Right Full (2.25" x 14") . . . . . STD
- SAPDR Right Desk (2.25" x 10") . . . . . STD

**Waterfall Arm Pads<sup>(1)</sup>**

- WAPFL Left Full (13") . . . . . N/C
- WAPDL Left Desk (9") . . . . . N/C
- WAPFR Right Full (13") . . . . . N/C
- WAPDR Right Desk (9") . . . . . N/C

**Gel Pads**

- GAPNL Narrow Left (2" x 12") . . . . . \$88
- GAPNR Narrow Right (2" x 12") . . . . . \$88
- GAPWL Wide Left (3.5" x 12") . . . . . \$96
- GAPWR Wide Right (3.5" x 12") . . . . . \$96
- GAPLL Long Left (3.5" x 14") . . . . . \$103
- GAPLR Long Right (3.5" x 14") . . . . . \$103

1. Full length pad requires full length arm tube if using Ultra Rail mounted Tilt Arms.

**MANUAL CENTER MOUNT FRONT RIGGING**

MUST BE COMPLETED	Measurement in Inches
Knee to Heel:	
Cushion Thickness: (Subtract)	
Cushion Compression:(Add)	
Seat Pan to Footplate:(Total)	

- MSFCMFP Maxx Style Fixed Center Mount  
Foot Platform<sup>(1)(2)</sup> . . . . . \$395

- IFP Medium Rubber Coated Footplate  
(11.5" W x 10") . . . . . N/C

- FBIFP Large Rubber Coated Footplate  
(17" w x 13"d) . . . . . N/C


- ICP Individual calf pads . . . . . \$200

Select seat pan to footplate:

- SPTF1 9" - 13" . . . . . N/C
- SPTF14 13" - 19" . . . . . N/C


Set at:

- SA70 Set at 70° . . . . . N/C
- SA90 Set at 90° . . . . . N/C
- SA97 Set at 97° . . . . . N/C



1. Depth adjustable, and easy slide adjustment for seat pan to foot plate extension.  
2. 3" ground clearance required.

**FOOT PLATFORM OPTIONS FOR MSFCMFP<sup>(1)</sup>**



- SFFFP Small Platform 12"w x 8.5"d . . . . . \$180
- MFFFP Medium Platform 12"w x 10.75"d . . . . . \$180
- LFFFP Large Platform 14.5"w x 13"d . . . . . \$250
- XLFFFP Extra Large Platform 17"w x 13"d . . . . . \$250
- XXLFFFP Extra Extra Large Platform 21"w x 13"d . . . \$250

1. May need to adjust leg position to avoid interference with front castors

**MANUAL ELEVATING CENTER MOUNT LEGREST**

MECML Manual Elevating Center Mount Legrest  
With Plantar Flex ..... \$595

Select legrest set up position:

- 90SU Set at 90° (elevating range 90-40°) ..... N/C
- 110TU Tucked at 110° (elevating range 110-60°)<sup>(1)(2)</sup> N/C

Rubber Coated Footplate:

MRCP Rubber Coated Footplate (11.5" W x 10") ..... \$165

Select seat pan to footplate for Rubber Coated Footplate

- CSPTF1 11.5" - 15" ..... N/C
- CSPTF2 14" - 19" ..... N/C

Individual Footplates:


- IFSL Small - (5" x 7.5" deep) - Left ..... \$165
- IFSR Small - (5" x 7.5" deep) - Right ..... \$165
- IFML Medium - (5.5" x 9.5" deep) - Left ..... \$165
- IFMR Medium - (5.5" x 9.5" deep) - Right ..... \$165
- IFLL Large - (6" x 11.5" deep) - Left ..... \$165
- IFLR Large - (6" x 11.5" deep) - Right ..... \$165

Select seat pan to footplate for Individual Footplates

- CSPTF3 9.5" - 14.5" ..... N/C
- CSPTF4 12" - 19" ..... N/C

Optional:

CICP Individual calf pads ..... \$200



1. Not available on seat depths 16", 17".

**MAXX STYLE MANUAL ELEVATING LEGRESTS**

MMPPSA MAXX Swing Away Manual Elevating<sup>(1)(2)</sup> . . . \$680

Select Receivers:

- RS Straight Receivers ..... N/C
- RF Flared 1" Outward Receivers ..... N/C

Select Seat Pan to Footplate:

- SPTF5 Length: 10"-14" ..... N/C
- SPTF2 Length: 13"-17" ..... N/C
- SPTF7 Length: 16"-20" ..... N/C

Options:

- FCPMS Flat Calf Pad Option ..... \$275
- CCPMS Curved Calf Pad Option<sup>(3)</sup> ..... \$275

1. Must choose footplate on page 8
2. Comes standard with moulded i-skin calf pads
3. Not available on 14 and 15" seat widths. For 16" and 17" seat widths you must select flared receivers.

**SWING AWAY LEGREST**

HD70 Heavy Duty 70° Swing Away Footrest<sup>(1)(2)</sup> . . . \$495

Select Receivers:

- RS Straight Receivers ..... N/C
- RF Flared 1" Outward Receivers ..... N/C

Select Seat Pan to Footplate:

- SPTF11 Length: 6"-10" ..... N/C
- SPTF8 Length: 11"-15" ..... N/C
- SPTF2 Length: 13"-17" ..... N/C
- SPTF7 Length: 16"-20" ..... N/C

Options:

- HDFCP Flat Calf Pad Option ..... \$275
- HD70CCP Curved Calf Pad Option<sup>(3)</sup> ..... \$275
- DCS Dynamic Coil Springs ..... \$190

I70SAF Invacare 70° Swing-away Footrests ..... \$325  
 I70-CF Comes with composite foot plates ..... N/C

70HDAIH Invacare 70° Non-Tapered HD Swing-away  
Footrests with AT5543 Adjustable Angle  
Footplates, Heel Loops and Impact Guards . \$457

70TAPAIH Invacare Tapered 70° Swing-away  
Footrests with AT5543 Adjustable Angle  
Footplates, Heel Loops and Impact Guards . \$497

ALPT Longer Pivot and Slide Tube for  
70° Swingaways (Pair) ..... \$132

RLSL Residual Limb Support - Left<sup>(4)</sup>  
6.75"W x 9.75"L x 2"H (each) ..... \$495

RLSR Residual Limb Support - Right<sup>(4)</sup>  
6.75"W x 9.75"L x 2"H (each) ..... \$495

1. Must choose footplate on page 6
2. Heavy duty durable construction (4.5lbs ea.).
3. For 16" and 17" seat widths you must select flared receivers.
4. Swing away and removable, 50° of angle adjustment, 3.5" height adjustment and 2" of for/aft adjustment.





Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**MATRIX SEAT CUSHION OPTIONS**

- Matrix Libra Seat Cushion** ..... \$745  
 LC1616    LC1618    LC1620    LC1716    LC1718  
 LC1720    LC1816    LC1818    LC1820    LC1916  
 LC1918    LC1920    LC2016    LC2018    LC2020
- Matrix Libra Seat Cushion** ..... \$1,025  
 LC1622    LC1722    LC1822    LC1922    LC2022
- Matrix Libra HD Seat Cushion** ..... \$1,025  
 LC2118    LC2120    LC2122    LC2218    LC2220  
 LC2320    LC2322    LC2418    LC2420    LC2422
- Matrix PS Seat Cushion** ..... \$569  
 PS1616    PS1618    PS1620    PS1716    PS1718  
 PS1720    PS1816    PS1818    PS1820    PS2016  
 PS2018    PS2020
- Matrix PS HD Seat Cushion** ..... \$975  
 PS2218    PS2220    PS2222    PS2418    PS2420  
 PS2422    PS2618    PS2620    PS2622
- Matrix PSP Seat Cushion** ..... \$645  
 PSP1616    PSP1618    PSP1620    PSP1716  
 PSP1718    PSP1720    PSP1816    PSP1818    PSP1820  
 PSP2016    PSP2018    PSP2020
- Matrix PSP HD Seat Cushion** ..... \$995  
 PSP2218    PSP2220    PSP2222    PSP2418    PSP2420  
 PSP2422    PSP2618    PSP2620    PSP2622
- Matrix Multi Cushion with reversable cover** ..... \$745  
 MM1616    MM1618    MM1620    MM1716    MM1718  
 MM1720    MM1816    MM1818    MM1820    MM2016  
 MM2018    MM2020
- Matrix Multi Cushion with reversable cover** ..... \$995  
 MM2216    MM2218    MM2220    MM2222

**UPGRADABLE BACK OPTIONS**

- Matrix Elite Back (3" Contour)** ..... \$766  
 PBE1614    PBE1616    PBE1618    PBE1620  
 PBE1714    PBE1716    PBE1718    PBE1720  
 PBE1814    PBE1816    PBE1818    PBE1820  
 PBE1916    PBE1918    PBE1920    PBE2014  
 PBE2016    PBE2018    PBE2020
- Matrix Elite HD Back (3" Contour)** ..... \$1,095  
 PBE2116-HD    PBE2120-HD    PBE2216-HD    PBE2220-HD  
 PBE2316-HD    PBE2320-HD    PBE2416-HD    PBE2420-HD
- Matrix Elite HD Back (3" Contour)** ..... \$1,666  
 PBE2616-HD    PBE2620-HD
- Matrix Elite Deep Back (6" Contour)** ..... \$918  
 EDB1614    EDB1616    EDB1618    EDB1620  
 EDB1714    EDB1716    EDB1718    EDB1720  
 EDB1814    EDB1816    EDB1818    EDB1820  
 EDB1914    EDB1916    EDB1918    EDB1920  
 EDB2014    EDB2016    EDB2018    EDB2020

**UPGRADABLE BACK OPTIONS CONT.**

- Matrix Elite Deep HD Back (6" Contour)** ..... \$1,095  
 EDB2116-HD    EDB2118-HD    EDB2120-HD    EDB2216-HD  
 EDB2218-HD    EDB2220-HD    EDB2316-HD    EDB2318-HD  
 EDB2320-HD    EDB2416-HD    EDB2418-HD    EDB2420-HD
- Matrix Elite E2 Back (3" Contour)** ..... \$965  
 E2S1614    E2S1616    E2S1714    E2S1716  
 E2S1814    E2S1816    E2S1914    E2S1916  
 E2S2014    E2S2016
- Matrix Elite E2 Back (3" Contour)** ..... \$995  
 E2S1618    E2S1620    E2S1718    E2S1720  
 E2S1818    E2S1820    E2S1918    E2S1920  
 E2S2018    E2S2020
- Matrix Elite E2 Deep Back (6" Contour)** ..... \$1,125  
 E2D1614    E2D1616    E2D1714    E2D1716  
 E2D1814    E2D1816    E2D1914    E2D1916  
 E2D2014    E2D2016
- Matrix Elite E2 Deep Back (6" Contour)** ..... \$1,155  
 E2D1618    E2D1620    E2D1718    E2D1720  
 E2D1818    E2D1820    E2D1918    E2D1920  
 E2D2018    E2D2020
- Matrix Elite TR Back (3" Contour)** ..... \$866  
 PBE1616-TR    PBE1618-TR    PBE1620-TR    PBE1716-TR  
 PBE1718-TR    PBE1720-TR    PBE1816-TR    PBE1818-TR  
 PBE1820-TR    PBE1916-TR    PBE1918-TR    PBE1920-TR  
 PBE2016-TR    PBE2018-TR    PBE2020-TR
- Matrix Elite TR HD Back 21" - 22" W (3" Contour)** ..... \$1,151  
 PBE2120-TRHD    PBE2220-TRHD    PBE2320-TRHD  
 PBE2420-TRHD
- Matrix MAC Back** ..... \$1,495  
 MAC1616    MAC1618    MAC1620    MAC1622  
 MAC1716    MAC1718    MAC1720    MAC1722  
 MAC1816    MAC1818    MAC1820    MAC1822  
 MAC1916    MAC1918    MAC1920    MAC1922  
 MAC2016    MAC2018    MAC2020    MAC2022
- Must Select Pad Contour:*  
 MLCSL   **Left 3.5"**    MLCSR   **Right 3.5"**  
 MLCDL   **Left 6"**    MLCDR   **Right 6"**  
 MLCXDL   **Left 9"**    MLCXDR   **Right 9"**
- Matrix MAC Back w/ HUG** ..... \$1,645  
 MAC1616-HUG    MAC1618-HUG    MAC1620-HUG  
 MAC1622-HUG    MAC1716-HUG    MAC1718-HUG  
 MAC1720-HUG    MAC1722-HUG    MAC1816-HUG  
 MAC1818-HUG    MAC1820-HUG    MAC1822-HUG  
 MAC1916-HUG    MAC1918-HUG    MAC1920-HUG  
 MAC1922-HUG    MAC2016-HUG    MAC2018-HUG  
 MAC2020-HUG    MAC2022-HUG
- Must Select Pad Contour:*  
 MLCSL   **Left 3.5"**    MLCSR   **Right 3.5"**  
 MLCDL   **Left 6"**    MLCDR   **Right 6"**  
 MLCXDL   **Left 9"**    MLCXDR   **Right 9"**



Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**MATRIX ELAN HEADREST PAD<sup>(1)</sup>**

<input type="checkbox"/> ESP6	Matrx Elan Standard Pad 6" . . . . .	\$310
<input type="checkbox"/> ESP10	Matrx Elan Standard Pad 10" . . . . .	\$310
<input type="checkbox"/> ESP14	Matrx Elan Standard Pad 14" . . . . .	\$310
<input type="checkbox"/> ESP6-HUG	Matrx Elan Hug Standard Pad 6" . . . . .	\$310
<input type="checkbox"/> ESP10-HUG	Matrx Elan Hug Standard Pad 10" . . . . .	\$310
<input type="checkbox"/> ESP14-HUG	Matrx Elan Hug Standard Pad 14" . . . . .	\$310

1. If mounting on PB Elite must order removable hardware

**HEADREST MOUNTING HARDWARE OPTIONS**

<input type="checkbox"/> LXSHW	Loxx Standard Hardware 15" . . . . .	\$275
<input type="checkbox"/> MEHWRE	Elan Hardware 15" . . . . .	\$285
<input type="checkbox"/> FPLXSHW	FLIPP Loxx Standard Hardware 15" . . . . .	\$705
<input type="checkbox"/> FPLXSHW-EX	FLIPP Loxx Extended Hardware . . . . .	\$705
<input type="checkbox"/> FPMEHW	FLIPP Elan Hardware 15" . . . . .	\$715

**MATRIX FIXED LATERALS - FOR MATRIX BACKS<sup>(1)</sup>**

<input type="checkbox"/> PBSSL	Matrx Small Lateral Fixed-Left (3 x 5.5") . . .	\$113
<input type="checkbox"/> PBSSR	Matrx Small Lateral Fixed-Right (3 x 5.5") . .	\$113
<input type="checkbox"/> PBLSL	Matrx Standard Lateral Fixed-Left (7.5 x 5") .	\$116
<input type="checkbox"/> PBLSR	Matrx Standard Lateral Fixed-Right (7.5 x 5")	\$116
<input type="checkbox"/> PBOLL	Matrx Offset Lateral Fixed Left . . . . .	\$125
<input type="checkbox"/> PBOLR	Matrx Offset Lateral Fixed Right . . . . .	\$125

1. Not available on the MAC backrests.

**MATRIX SWING AWAY LATERALS - FOR MATRIX BACKS<sup>(1)</sup>**

<input type="checkbox"/> SAEHBL	Matrx Elite Swing Away Hardware Left . . . . .	\$250
<input type="checkbox"/> OSAHL	Matrx Offset Elite Swing Away Hardware Left . . . . .	\$250

Select Left Pad:

<input type="checkbox"/> EPADXL	Size XS (3.25 x 4.25") . . . . .	\$115
<input type="checkbox"/> EPADSL	Size S (3.5 x 5.25") . . . . .	\$115
<input type="checkbox"/> EPADML	Size M (3.75 x 5.75") . . . . .	\$115
<input type="checkbox"/> EPADLL	Size L (4.25 x 6.5") . . . . .	\$115

<input type="checkbox"/> SAHEBR	Matrx Elite Swing Away Hardware Right . . . . .	\$250
---------------------------------	--	-------

<input type="checkbox"/> OSAHR	Matrx Offset Elite Swing Away Hardware Right . . . . .	\$250
--------------------------------	---	-------

Select Right Pad:

<input type="checkbox"/> EPADX	Size XS (3.25 x 4.25") . . . . .	\$115
<input type="checkbox"/> EPADSR	Size S (3.5 x 5.25") . . . . .	\$115
<input type="checkbox"/> EPADMR	Size M (3.75 x 5.75") . . . . .	\$115
<input type="checkbox"/> EPADLR	Size L (4.25 x 6.5") . . . . .	\$115

1. Not available on the MAC backrests.

**HIP SUPPORT PADS WITH FIXED HARDWARE**

**Lateral Hip Support Left**

<input type="checkbox"/> HSPXLH	Size XS (5.5" W x 4" H) . . . . .	\$120
<input type="checkbox"/> HSPSMLH	Size S/M (7" W x 4" H) . . . . .	\$120
<input type="checkbox"/> HSPSLH	Size S (7"W x 6"H) . . . . .	\$120
<input type="checkbox"/> HSPMLH	Size M (9" W x 4"H) . . . . .	\$120
<input type="checkbox"/> HSPLLH	Size L (13" W x 4"H) . . . . .	\$120

**Lateral Hip Support Right**

<input type="checkbox"/> HSPXRH	Size XS (5.5"W x 4"H) . . . . .	\$120
<input type="checkbox"/> HSPSMRH	Size S/M (7"W x 4"H) . . . . .	\$120
<input type="checkbox"/> HSPSRH	Size S (7"W x 6"H) . . . . .	\$120
<input type="checkbox"/> HSPMRH	Size M (9"W x 4"H) . . . . .	\$120
<input type="checkbox"/> HSPLRH	Size L (13"W x 4"H) . . . . .	\$120

**KNEE SUPPORT PADS WITH FIXED HARDWARE**

**Lateral Knee Support Left**

<input type="checkbox"/> HSPXLK	Size XS (5.5"W x 4" H) . . . . .	\$120
<input type="checkbox"/> HSPSMLK	Size S/M (7"W x 4" H) . . . . .	\$120
<input type="checkbox"/> HSPSLK	Size S (7"W x 6" H) . . . . .	\$120
<input type="checkbox"/> HSPMLK	Size M (9"W x 4"H) . . . . .	\$120
<input type="checkbox"/> HSPLLK	Size L (13"W x 4"H) . . . . .	\$120

**Lateral Knee Support Right**

<input type="checkbox"/> HSPXRK	Size XS (5.5"W x 4"H) . . . . .	\$120
<input type="checkbox"/> HSPSMRK	Size S/M (7"W x 4"H) . . . . .	\$120
<input type="checkbox"/> HSPSRK	Size S (7"W x 6"H) . . . . .	\$120
<input type="checkbox"/> HSPMRK	Size M (9"W x 4"H) . . . . .	\$120
<input type="checkbox"/> HSPLRK	Size L (13"W x 4"H) . . . . .	\$120

**UPGRADED HIP/ KNEE SUPPORT HARDWARE<sup>(1)</sup>**

**MAXX Style Quick Release, Removable, Multi-Axis Hardware**

Select mounting position:

<input type="checkbox"/> MHSHLH	Left mounted at Hip . . . . .	\$215
<input type="checkbox"/> MHSHRH	Right mounted at Hip . . . . .	\$215
<input type="checkbox"/> MHSHLK	Left mounted at Knee . . . . .	\$215
<input type="checkbox"/> MHSHRK	Right mounted at Knee . . . . .	\$215

**Extended Maxx Style Quick Release, Removable, Multi-Axis Mounting Hardware<sup>(2)</sup>**

Select mounting position:

<input type="checkbox"/> EMQHLH	Left mounted at Hip . . . . .	\$215
<input type="checkbox"/> EMQHRH	Right mounted at Hip . . . . .	\$215
<input type="checkbox"/> EMQHLK	Left mounted at Knee . . . . .	\$215
<input type="checkbox"/> EMQHRK	Right mounted at Knee . . . . .	\$215

1. Must select Hip and/or Knee Support Pad size

2. Recommended when finished seat cushion thickness will be greater than 3.5" high.



