

TRx ULTRA-Low Bariatric Power Positioning System

US PRICE LIST AND ORDER FORM


Effective July 2008



Date of Order: _____
 Dealer Name: _____
 Dealer Account #: _____
 PO #: _____
 Base MFR. PO #: _____
 Tag: _____
 Purchasing Contact: _____
 Phone: _____ Fax: _____
 RTS/Therapist: _____
 Ship to Address: _____
 City: _____ State: _____ Zip: _____
 Client Info: Gender: _____
 Special Client Conditions: _____
300-500 lbs weight limit. Dependent on weight capacity of the power base.

To ensure system is accurately configured please fill in all required measurements.

Client Height:	
Client Weight (lbs):	
A. Seat to Shoulder:	
B. Trunk Depth:	
C. Chest Width:	
D. Knee to Back:	
E. Seat to Top of Head:	
F. Elbow to Hand:	
G. Seat to Elbow:	
H. Hip Width:	
I. Knee to Heel:	
Seat Cushion Thickness:	



INVACARE

BASE TYPE SELECTION

TDX SR HD Select STF⁽¹⁾: 18.25"
400 lbs weight capacity

TDX SP* *must order base with 400 lb weight capacity
Select STF⁽¹⁾: 18.25"
400 lbs weight capacity

3G Arrow Select STF⁽¹⁾: 18.5"

3G Torque Select STF⁽¹⁾: 18.5"

3G RangerX Select STF⁽¹⁾: 18.5"
400 lbs weight capacity

M94 Select STF⁽¹⁾: 18.5"
450 lbs weight capacity

ELECTRONICS

Mark VI PSS with ACC SPJ MPJ ...N/C

Mark V NX™80 with SPJ-80 ...N/C

Select: EX Controller SPJ DPJ MPJ

Mark IVR2 Mark IV Mark IVA ...N/C

OTHER INFORMATION

Short Base Long Base

4 Pole Motors GB Motors

BASES MUST BE ORDERED "TRx" READY FROM MANUFACTURER

QUANTUM

BASE TYPE SELECTION

1170 XL Select STF⁽¹⁾: 19.25"
400 lbs weight capacity

6000 XL Select STF⁽¹⁾: 17.5"
400 lbs weight capacity

1420* Select STF⁽¹⁾: 21"
450 lbs weight capacity

1650* Select STF⁽¹⁾: 19.5"
(Power Pivot Plus legs will increase STF to 21")
500 lbs weight capacity

614 HD Select STF⁽¹⁾: 18.25" 19.25"
400 lbs weight capacity

*Batteries are required at the time of installation.

ELECTRONICS

VSI 6-Key* VSI 7-Key P&G ...N/C

VSI Integral 4-key R-Net ...N/C

Q-Logic ...N/C

(* Interface is required, see electronic options)

*** ON BOARD CHARGERS ARE NOT RECOMMENDED FOR REHAB APPLICATIONS AS THE MAIN CHARGING SOURCE***

QUICKIE

BASE TYPE SELECTION

Rhythm Select STF⁽¹⁾: 18.75"
400 lbs weight capacity

S646 Select STF⁽¹⁾: 19.25"
400 lbs weight capacity

ELECTRONICS

Qtronix VSI 6-key ...N/C

Delphi Select Joystick: QR* QC ...N/C

* You must select QR to run the functions through the joystick.

OTHER INFORMATION

For S646 Bases select:

Standard Width(18") Wide Width (20")

SYSTEMS⁽²⁾ MUST SELECT ONE SYSTEM

(Includes seat pan, back pan, armrests, armpads, electronics and standard switch controls)

TRx Ultra-Low Bariatric CG Tilt Only System (45°Tilt)⁽³⁾⁽⁴⁾\$7,395

Note: Tilt range is 45° up to 400 lbs, Tilt range is 30° for 400 lbs and greater

TRx Ultra-Low Bariatric CG Tilt (45°Tilt)⁽³⁾⁽⁴⁾ and Power Recline (135°) System\$11,390

Note: Tilt range is 45° up to 400 lbs, Tilt range is 30° for 400 lbs and greater

TRx Ultra-Low Bariatric CG Tilt (45°Tilt)⁽³⁾⁽⁴⁾ and Power Recline (135°) System with Extended Shear Reduction (ESR)\$11,955

Note: Tilt range is 45° up to 400 lbs, Tilt range is 30° for 400 lbs and greater

TRx Bariatric Recline Only System (135°)\$6,995

TRx Bariatric Recline Only System (135°) with Extended Shear Reduction\$7,560

MODULES - ONLY AVAILABLE ON SYSTEMS WITH RECLINE

TRx 30° Precline Module⁽⁵⁾\$500

TRx 8° Precline Module⁽⁵⁾N/C

STANDARD SEAT SPECIFICATIONS

Seat Width

19" 20" 21" 22" 23" 24" 25"N/C

Seat Depth

16" 17" 18" 19" 20"N/C

21" 22" 23" 24"N/C

CUSTOM SEAT SPECIFICATIONS

Seat Width

26" and overCall for RFQ

SEAT CUSHION OPTIONS

MaTRx-V Seat Cushion 20" w-30" w\$625

MaTRx-Flo-Tech Cushion 20" w-24" w\$595

"ON-CHAIR" Cushion Credit(\$90)

FABRIC OPTIONS FOR MATRX-VI CUSHION

Onyx Fabric Upgrade for MaTRx-Vi Cushion .. \$100

Vinyl Fabric Upgrade for MaTRx-Vi Cushion ... \$50

RECLINE BACK OPTIONS

TRx Recline Body Fit System (Minimum back widths of 16") ... \$995

(This option allows the seat back to be narrower than the seat pan. For example an 18" back pan width with a 20" seat width is accomplished by selecting "2" Narrower" option.)

2" Narrower 4" Narrower

TILT ONLY BACK CANE SPECIFICATIONS

Select Back Cane Height - Measured from Seat pan to top of back cane

22" 24"N/C

Select Back Cane Angle 76° 79° 86° 96° 101° 112° 116°

Custom Cane Height Specify _____\$250

Angled (5°) Tilt Back Canes⁽⁶⁾N/C

Select Back Cane Angle 81° 84° 91° 101° 106° 117° 121°

Endomorph Back Canes Select: 1" 2" 4" Narrower N/C

FINISHED BACK HEIGHT*

Measured from seat pan to top of back shell

18"⁽⁷⁾ 19"⁽⁷⁾ 20" 21" 22" 23" 24" 25"

26" 27" 28" 29" 30" 31" 32" 33"

34" 35" * Comes with standard cane height

BACK OPTIONS

Standard Rehab Back and Back cushion
Finished height from 18" to 25"STD
Finished height from 26" to 35"\$450

MaTRx Fabric Upcharge -for Standard rehab back\$100

Recessed planar interface plate⁽⁸⁾ N/C

OMIT: Rigid Back Pan & Cushion⁽⁹⁾(\$75)

OMIT: Back Cushion Only(\$25)

UPGRADEABLE BACK OPTIONS

MaTRx Contour Back (HCPCS Code E2613)\$495
(For Systems 15"-21" wide under 300lbs)

Specify actual back height dimension: 12" 15" 17"

MaTRx Contour Back (HCPCS Code E2614)\$595
(For Systems 22" wide under 300lbs)

Specify actual back height dimension: 15" 17"

MaTRx High Back (HCPCS Code E2620)\$595
(For Systems 16"-21" wide under 300lbs)

Specify actual back height dimension: 20" 25" 30"

MaTRx PB Heavy Duty Back (5" Contour)\$835
(For Systems 20"-22" wide and over 500lbs)

Specify actual back height dimension: 16" 20"

MaTRx PB Heavy Duty Back (5" Contour)\$1,580
(For Systems 26"-30" wide and over 500lbs)

Specify actual back height dimension: 16" 20"

MaTRx PB Heavy Duty Deep Back (7" Contour)\$935
(For Systems 20"-22" wide and over 500lbs)

Specify actual back pan dimensions: 16" 20"

MaTRx PB Heavy Duty Deep Back (7" Contour) ...\$1,435
(For Systems 26"-30" wide and over 500lbs.)

Specify actual back pan dimensions: 16" 20"

"ON-CHAIR" Matrx Back Credit(\$125)

Onyx Fabric Upgrade\$100

Select Contour Back High Back

NOTE WHEN CHOOSING MaTRx Backs
4" Gap recommended for Tilt systems to allow for the seat cushion.
6" Gap required for systems with Recline/ESR or PSB
* PB Backs are not available with the Precline Module.

HEADRESTS & MOUNTING HARDWARE

- Headrests with Fixed Mounting Hardware** (HCPCS code E0955)
- Concepts Standard Headrest \$295
 - Combination Head/Neck - Small \$295
 - Combination Head/Neck - Large \$295
 - Stealth Comfort Plus 10" (CP250) \$242
 - Stealth Comfort Plus 14" (CP150) \$242
 - Onyx Auto Style Headrest - Straight Mounting available only . . \$350
- Headrest Swing Aways and Mounts - Does not include the Headrest** (HCPCS code E1028)
- Concepts Detachable Mounting Hardware \$200
 - Stealth's World's Best Headrest Hardware (TWB485) . . \$220
 - OMIT: Headrest** N/C

ARMRESTS

- Adult Dual Post Adjustable Height (Single Button Release) . . STD
 - Add Quick Height Adjustment Lever \$250
 - Select Size: (9.75"-12.75") (12.75"-15.75") (15.75"-18.75")
 - Select: Flat pouch Glove Boxes*
 - (* Not Available with Mesomorph Back Canes)
- Reclining Adjustable Height STD
 - Select Size: (9.5"-13") (13.5"-18.5")
- 1" Offset Armrest Spacer (pair) \$250
 - (Allows seat width to increase by 2")
- 1/2" Offset Armrest Spacer (pair) \$250
 - (Allows seat width to increase by 1")

ARMPADS

- Standard Pads-** STD
 - Select: Full (2.25"x14") Left Right
 - Desk (2.25"x10") Left Right
- Waterfall Arm Pads-(Pair)** \$35
 - Select: Full (2.25"x14") Left Right
 - Desk (2.25"x10") Left Right
- Flat Pads- (Pair)** \$55
 - Select: Full (3"x14") Left Right
 - Desk (3"x10") Left Right
- Gel Pads- (Pair)** \$150
 - Select: Narrow (2"x12") Left Right
 - Wide (3.5"x12") Left Right
- Ergonomic Arm Trough - LEFT** (4.5"x13") \$85
 - (HCPCS Code K0106)
- Ergonomic Arm Trough - RIGHT** (4.5"x13") \$85
 - (HCPCS Code K0106)
- Flat Hand Pad - LEFT** (4.5"x5") \$55
 - (HCPCS Code K0108)
- Flat Hand Pad - RIGHT** (4.5"x5") \$55
 - (HCPCS Code K0108)
- OMIT: Armpads** (7/8" diameter tube) (\$25)
 - Select: Full size (14") Left Right
 - Desk size(10") Left Right
- Horizontal Armpad Extensions (each)** Quantity () . . \$75
 - Select: 3" 4.375" 7"

ARMPADS ACCESSORIES

- TRx Arm Pad Storage/Cup Holder LEFT** * \$175
- TRx Arm Pad Storage/Cup Holder RIGHT** * \$175
 - *Not available on the same side as the joystick is mounted
- TRx Elbow Block LEFT** (prevents arm from slipping off) \$195
 - Select Pad Size:
 - XS(3.25"x4.25") S(3.5"x5.25") M(3.75"x5.75") L(4.25"x6.5")
 - Select Vertical Extension: S(3") L(5")
 - Horizontal Extension: S(5") L(7")
- TRx Elbow Block RIGHT** (prevents arm from slipping off) \$195
 - Select Pad Size:
 - XS(3.25"x4.25") S(3.5"x5.25") M(3.75"x5.75") L(4.25"x6.5")
 - Select Vertical Extension: S(3") L(5")
 - Horizontal Extension: S(5") L(7")
- Multi Purpose Arm Pad Adapter Block (each)** Quantity () \$35
 - (Allows hard mounting & positioning of standard pads or other manufacturer's arm pads. Adapter allows square mounting on round tubes)
- Flip & Fold Away 1/2 Lap Tray** (14"W x 10"D) \$795
 - (Only available with a flat arm pad)
 - Select Mounting*: Left Right
 - *Not available on the same side as the joystick is mounted

LATERALS - FOR ALL OTHER BACKS (EXCLUDING MATRIX PB ELITE AND PB BACKS)

- TRx Lateral Trunk Support with Fixed Mounting LEFT⁽¹⁰⁾**
 - Pad Size: S(3.5"x5.25") M(3.75"x5.75") L(4.25"x6.5")
 - Select Offset 0"-3" 3"-6"
 - (HCPCS code E0956) \$100
- TRx Lateral Trunk Support with Fixed Mounting RIGHT⁽¹⁰⁾**
 - Pad Size: S(3.5"x5.25") M(3.75"x5.75") L(4.25"x6.5")
 - Select Offset 0"-3" 3"-6"
 - (HCPCS code E0956) \$100
- Onyx Fabric Upgrade** \$100

**SELECT SWING AWAY HARDWARE TO CONVERT ABOVE:
NOT AVAILABLE ON MATRIX PB ELITE, PB ELITE FULL, AND MATRIX PB BACKS**

- TRx Swingaway Mounting Hardware⁽¹⁰⁾**
 - LEFT (HCPCS code E1028) \$195
 - RIGHT (HCPCS code E1028) \$195
 - Select Lateral Latch plate length:
 - Right: Short (3.5") Long(5") Left: Short (3.5") Long(5")

STEALTH LATERAL OPTIONS

- TWB-L Stealth Laterals⁽¹⁰⁾ (pair)** \$590
 - Pad Size: S(3"X5") L(4"X6")

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

HIP SUPPORTS

- Lateral Hip Support Pad with Fixed Mounting LEFT . . \$120
 Select Size:
 (5.5"W x 4.5"H) (7"W x 5.75"H) (9"W x 4.5"H) (13"W x 4.5"H)
 (HCPCS code E0956)
- Lateral Hip Support Pad with Fixed Mounting RIGHT . \$120
 Select Size:
 (5.5"W x 4.5"H) (7"W x 5.75"H) (9"W x 4.5"H) (13"W x 4.5"H)
 (HCPCS code E0956)
- Onyx Fabric Upgrade \$100

HARDWARE TO CONVERT ABOVE HIP SUPPORT TO REMOVABLE/SWINGAWAY HIP SUPPORTS

- TRx Lift Off, Removable Hardware for Hip Supports
 - LEFT (HCPCS code E1028) \$195
 - RIGHT (HCPCS code E1028) \$195
- TRx Swingaway Removable Hardware for Hip Supports
 - LEFT (HCPCS code E1028) \$215
 - RIGHT (HCPCS code E1028) \$215
 Select Mounting Position:
 Hip Knee or Both... Add \$215 per side

STANDARD FRONT RIGGING

CALL TO CONFIRM ON SEAT SIZES BELOW 16"W X 16"D

- Legrest Receivers Only STD

****Please Specify Front rigging to be used on TRx System ⁽¹¹⁾**

POWER CENTER MOUNT FOOT PLATFORM OPTIONS

MUST BE COMPLETED	Measurement in Inches
Knee to Heel:	
Cushion Thickness: (Subtract)	
Cushion Immersion: (Add)	
Seat Pan to Footplate: (Total)	

*****CENTER MOUNTS ARE RESTRICTED TO SEAT DEPTHS 16" AND GREATER******

- CMB Belt Drive Center Mount Articulating Foot Platform⁽¹²⁾⁽¹³⁾⁽¹⁴⁾
 (HCPCS code K0108) \$2,595
 Seat pan to footplate: - select one:
 Short 7"-13"- 3"-5" adjustable articulation
 Long 14"-21"- 5"-7.5" adjustable articulation
 - Standard with one piece calf pad

POWER FRONT RIGGING OPTIONS

- TRx Power Elevating Pivot Plus Legrests ⁽¹⁵⁾ \$2,490
 (Includes Pivot Plus Legrests & Actuator)(HCPCS code E1010)
 Select Function: Individual⁽¹⁶⁾ Combined Synchronized W/Recline
 Seat pan to footplate: S (12"-16") M (15"-19") L (18"-21")
 Select Footplate: Regular (7.5"W x 9"D) Large (9"W x 11"D)
 ADD: Power Actuator Quick Release \$350

ELEVATING LEGRESTS

- TRx Pivot Plus Manual Elevating Legrests (Pair) \$720
 Seat pan to footplate: S (12"-16") M (15"-19") L (18"-21")
 Select Footplate: Regular (7.5"W x 9"D) Large (9"W x 11"D)
- AT5544 Standard Elevating Legrests (Pair) \$495
- Jazzy Elevating Legrests (Pair) \$495
- Quickie Standard Elevating Legrests (Pair) \$495

ADDITIONAL FRONT RIGGING OPTIONS

90° FRONT RIGGING/CENTER MOUNTED OPTIONS NOT AVAILABLE ON STORM SERIES BASES

- TRx-Heavy Duty 70° Swingaway Footrests (Pair)⁽¹⁷⁾ . . . \$420
 Select Receivers: Straight Flared inward 1" Flared outward 1"
 Seat pan to footplate: S (12"-16") M (15"-19") L (18"-22")
 Select Footplate: S(4.25"W x 8"D) M(5.25"Wx8"D) L(6.25"Wx8"D)
- Dynamic Coil Springs (Pair) for TRx-Heavy Duty 70° Swingaway Footrests \$165
- Action 70° Swingaway Footrests (Pair) \$395
- Jazzy 70° Swingaway Footrests (Pair) \$395
- Quickie 70° Swingaway Footrests (Pair) \$395
- Adjustable 70°-90° Fixed Center Mount Foot platform⁽¹⁸⁾ \$595
 (Depth and Height Adjustable -- Comes with a 12" x 10" footplate made of durable aluminum and calf pads are included)
 Seat pan to footplate: S (12"-16") M (15"-19") L (18"-22")
- Invacare Split Foot Platform Interface Bracket \$595
 (Invacare CTMT-T Foot Platform is included)

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

FOOT PLATFORM & FOOTPLATE OPTIONS

Upgraded Flip-up Foot Platform* (12"W x 10"D) . . . \$180

Bariatric Flip-up Foot Platform*(15"W x 10"D) \$250

Split Footplates* (6" x 10" each) pair \$275

(*For use with Adjustable 70°-90° Fixed Center Mount Foot platform -Made of durable aluminum 1/4" thick with flexible heel loops.)

Adjustable Angle Footplates(Each)(HCPCS code K0040)

Left:

S(4.25"W x 8"D) M(5.25"Wx8"D) L(6.25"Wx8"D) . . . \$125

Right:

S(4.25"W x 8"D) M(5.25"Wx8"D) L(6.25"Wx8"D) . . . \$125

Adjustable Angle Footplates with Cutouts for Caster Clearance (Each)(HCPCS code K0040)

Left Regular (7.5"W x 9"D) Large (9"W x 11"D) . . . \$125

Right Regular (7.5"W x 9"D) Large (9"W x 11"D) . . . \$125

Heel Loops (Each) QUANTITY () \$20
(HCPCS code E0951)

CALF PAD AND STRAP OPTIONS

Flat Calf pads \$275

Select: TRx HD Swing Away Footrests

Fixed Center Mount Foot platform

Curved Calf pads \$275

Select: TRx HD Swing Away Footrests

TRx Pivot Plus Legrests

Fixed Center Mount Foot platform

One piece calf pad for Fixed Center Mount Footplate . \$275

Custom Bariatric Calf Pads Call for RFQ

Airmesh Calf Strap for chairs 18"-22" wide-3.5" pad height
(HCPCS code E0978) BP-CS-L \$50

SWING AWAY FOOTREST MECHANISM

- Assists user to swing away 70° footrests
- Power version requires enhanced electronics
- Compatible with Motion Concepts TRx HD 70° Swing away footrests or Invacare Action Swing away footrests
- Footrests NOT included

TRx Power Swing Away Mechanism . .Suggested Retail \$1,995

.RFQ Required

Specify : Switch or Joystick (See page 7 for electronics)

TRx Manual Swing Away Mechanism . . .Suggested Retail \$995

.RFQ Required

Specify Lever Mounting Right or Left

ADDITIONAL OPTIONS AND ACCESSORIES

TRx Concepts Back Pack - Single \$125

TRx Concepts Back Pack - Double \$195

Transfer Handles (Height Adjustable) \$350

Select: 4"-6.5"

9"-11.5"

D-Loops (set of 2)* \$250

-D-Loops increase the seat to floor height by 3/8"

BODY POINT LAP BELTS

Padded Lap Belt 2 point - BP-PLB2 (HCPCS code E0978) . . \$120

Padded Lap Belt 4 point - BP-PLB4 (HCPCS code E0978) . . \$150

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

SPECIFY DRIVE

Select: Left hand Right hand

ELECTRONICS SWITCH MOUNTING

Select: Left Right Do Not Mount

Specialty Controls Please call for specification sheet 1-888-433-6818

STANDARD ELECTRONICS - N/A WITH VENT TRAYS OR HEAD ARRAYS

(These electronics are equipped with drive lock-out. They do not include attendant control, tilt limits, and variable speed control(s).)

TRx Single-Function Control BoxSTD

Please Select One: Single Push Button Dual Push Button

Single Toggle Omit Switch*

(*Single-Function Control Box uses a phono jack connection. Compatible with other phono jack switches.)

HIGH END ELECTRONICS

(These electronics are equipped with drive lock-out, attendant control, tilt limits and variable speed control(s).)

SINGLE-FUNCTION SYSTEM - (IE: TILT ONLY)

TRx Enhanced Single-Function Control Box \$500

Select: Single Push Button Dual Push Button⁽¹⁹⁾

Single Toggle

Enhanced Single Function through Joystick* \$1,150

(HCPCS code E2310)

Additional Switch: Single Push Button Dual Push Button⁽¹⁹⁾ \$97

Single Toggle

DUAL-FUNCTION SYSTEM - (IE: TILT/RECLINE)

TRx Enhanced Dual-Function Control Box \$700

Select: Dual Push Button Quad Push Button⁽¹⁹⁾

Dual Toggle Four-Way Toggle

Enhanced Dual Function through Joystick* \$1,550

(HCPCS code E2311)

Additional Switch: Dual Push Button Quad Push Button⁽¹⁹⁾ . . \$97

Dual Toggle Four-Way Toggle

MULTI-FUNCTION SYSTEM - (IE: TILT/RECLINE/POWER LEGS)

TRx Enhanced 2-5 Function Control Box \$1,150

Select: Quad Push Button Four-Way Toggle

Add: Function Display* \$250

Add: Function Display/Step Scanner* \$400

* The display is operated through left/right motion on 4-Way Toggle or buttons 3 and 4 on the Quad Push Button.

Enhanced 2-5 Function through Joystick \$2,350

(HCPCS code E2311)

Includes: Function Display Yes No

Add: Function Display/Step Scanner* \$150

Additional Switch: Quad Push Button Four-Way Toggle . . . \$97

(*Motion Concepts provides the communication interface for Invacare Electronics and includes a function display)

ELECTRONICS OPTIONS - REQUIRES HIGH END ELECTRONICS

TRx M16 Splitter(HCPCS code E2399) \$125

(Allows for multiple switch options with phono plug connection)

TRx Feather Touch Switch - Single function \$97

Dime Size Quarter Size

TRx Feather Touch Switch Multi function \$135

Dual Dime Size Quad Dime Size

Stealth Green Egg (ESGRN) \$77

(For use with Single function electronics only)

TRx Power Shut Off Switch⁽²⁰⁾ - Toggle Switch Only . . . \$175

(HCPCS code E2399)

TRx Lazarus Jump Starter \$150

(HCPCS code A9270)

Accessory Port Power Supply \$195

(3 amp capacity)

TRx Lighting Effects - Through Switch \$350

TRx Front Safety Lights (pair)

Operated through switches \$495

Operated through Joystick \$845

TRx Rear Safety Lights (pair)

Operated through switches \$495

Operated through Joystick \$845

TRX Front & Rear Light Combo Package Credit* . . . (\$150)

*Purchase both sets and receive a credit

MEMORY SEATING OPTIONS

Single Function Memory Seating* Suggested Retail \$495

.RFQ Required

(For use with Enhanced electronics only - Not available to control Elevating seat modules)

*Please call for specification sheet 1-888-433-6818

Multiple Function Memory Seating* . . Suggested Retail \$995

.RFQ Required

(For use with Enhanced electronics only - Not available to control Elevating seat modules)

*Please call for specification sheet 1-888-433-6818

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

O2 HOLDERS

- Non-Articulating O² Gas Holder (Cylinder Style) \$895
 - Articulating O² Holder⁽²¹⁾ \$1,250
- Please Select One:
- O² Gas (Cylinder Style) Liquid O² (6.5" x 4.5")

BATTERIES

- MK U1 Battery (MU-1 SLD G) Each \$140
Quantity () (HCPCS code E2365)
- MK GROUP 22NF Battery (M22NF SLD G) Each
Quantity () (HCPCS code E2361) \$173
- MK GROUP 24 Battery (M24 SLD G) Each
Quantity () (HCPCS code E2363) \$215
- MK GROUP 34 Battery (M34 SLD G) Each
Quantity () \$215

INSTALLATION

- Installation by Motion Concepts⁽²²⁾ N/C
- Dealer Installation (Subject to approval) (\$100)

Total Retail Value \$ _____

NOTES

- 1 Seat-to-floor height is measured at the front of the seat pan based on a standard seat size with a 0° pre-tilt. Seat to floor heights may vary +/- .25" dependent on conditions of the power base. Suspension castors may affect seat-to-floor height up to an additional .5"-1".
- 2 Systems are designed for use with a 2"-3" cushion.
- 3 Tilt range is 45° up to 400 lbs and 30° for 400 lbs or greater.
- 4 Back angle is adjustable from -76° to 116° degrees in increments.
- 5 Total degrees of recline is 174 ie: if 30 degrees of precline is used, then 144 degrees of recline will be available.
- 6 This option allows the back angle to adjust from 81° to 121° in increments.
- 7 This back height is not recommended with Recline Systems, please call for further details. Power Sliding Back (PSB) is not available on these back heights.
- 8 Included with all Recline Systems when no back is selected. May also be used as an interface for after market backs.
- 9 Not available with Recline/ESR or Recline/PSB.
- 10 Laterals are not available with the MaTRx PB Backs.
- 11 Standard front rigging will only be provided when specified. Please provide or order separately.
- 12 Includes a flip up foot platform. Increases seat-to-floor height.
- 13 This item is only warrantied for client weights up to 375 lbs.
- 14 90° front rigging will collide with front casters on the Storm 3G, S626/S646.
- 15 Under normal usage Power Legrest actuators tend to be louder than Tilt/Recline Actuators.
- 16 Individual legs count as two functions.
- 17 Heavy duty durable construction (4.5lbs ea.). Lift and swing away setup.
- 18 90° front rigging will collide with front casters on the Storm 3G.
- 19 Non-reversing relays.
- 20 Must order an Enhanced control box.
- 21 Articulating O2 holders are not standard configuration on Precline Systems. Please call for quotation.
- 22 Batteries may be required at time of installation, please call to confirm.