



X3

"F" Series CG Tilt System

US Price List - Effective August 2017



REQUIRED INFORMATION

Date of Order: _____
 Dealer Name: _____
 Dealer Account #: _____
 PO #: _____
 Tag: _____
 Purchasing Contact: _____
 Phone/Fax: _____
 RTS/Therapist: _____
 Ship to Address: _____

 City: _____ Zip: _____
 Client Info: Gender: _____
 Special Client Conditions: _____

REQUIRED MEASUREMENTS

To ensure system is accurately configured please fill in all required measurements.

Client Height:	
Client Weight (lbs):	
A. Seat to Shoulder:	
B. Trunk Depth:	
C. Chest Width:	
D. Knee to Back:	
E. Seat to Top of Head:	
F. Elbow to Hand:	
G. Seat to Elbow:	
H. Hip Width:	
I. Knee to Heel:	
Cushion Thickness:	

BASE COLOR OPTIONS

Select Color: Red Delicious Blizzard White Original Orange
 Party Pink Matte Black Motion Green
 Bright Blue Solar Yellow

WHEELCHAIR TRANSPORT BRACKETS

Wheelchair Transport Brackets⁽¹⁾ Part # XBKTS \$250

SINGLE POWER OPTION VR2 - K0856

ROVI X3⁽²⁾ HCPCS Code: K0856 \$6,795
 VR2 Non-Expandable Joystick - 2 Actuator Control . . .N/C

SINGLE POWER OPTION R-NET - K0856

ROVI X3⁽²⁾ HCPCS Code: K0856 \$7,395
 R-Net Expandable Controller⁽³⁾ HCPCS Code: E2377 \$615
 R-Net Harness for Expandable Controller⁽³⁾ \$595
 HCPCS Code: E2313

OPTIONAL R-NET JOYSTICKS

R-Net LED Remote Joystick⁽⁴⁾ \$400
 R-Net LED Remote Joystick for Lights⁽⁴⁾ \$500
 R-Net Color Remote Joystick with Monojacks⁽⁴⁾ \$600
 R-Net Color Remote Joystick with Monojacks for Lights⁽⁴⁾ \$650
 R-Net Omni Display for Specialty Drive Input⁽⁴⁾ \$950

JOYSTICK MOUNT Select Mounting Position: Left Right

Fixed N/C
 Swing-Away HCPCS Code: E1028 \$235
 Retractable Swing-Away Joystick HCPCS Code: E1028 \$325
 Adjustable, Swing-Away Mount for Omni HCPCS Code: E1028 \$295
 Select Mounting Position: Left Right

JOYSTICK HANDLES

"T" Handle for Joystick Knob HCPCS Code: E2323 \$115
 Ball Handle for Joystick Knob HCPCS Code: E2323 \$115
 BodyPoint Joystick Handle 3" U-Shaped \$120
 HCPCS Code: E2323
 BodyPoint Joystick Handle 4" U-Shaped \$120
 HCPCS Code: E2323
 BodyPoint Joystick Handle, Domed Rubber \$120
 HCPCS Code: E2323

ELECTRONIC OPTIONS

R-Net Output Module (ECM)⁽⁴⁾⁽⁵⁾ \$995
 R-Net Bluetooth Module⁽⁴⁾⁽⁵⁾ \$895
 Select OS: iOS Apple Android/PC
 Light Package (Fully integrated LED lighting kit: headlights, taillights and turn signals)⁽⁶⁾ \$1,195

ELECTRONIC PROGRAMMERS

R-Net PC Programmer with Dongle \$450
 Diagnostic Test Tool (DTT) with RNet & VR2 Cables \$1,400

SEAT TO FLOOR HEIGHT

Select Seat-to-Floor height:⁽⁷⁾ 17" 18" 19" 20"
 Part # STF17 STF18 STF19 STF20

- NOTES:**
- For unoccupied use only.
 - Price of the powerbase includes: Active Ride Control suspension, 8 amp offboard charger, flat free inserts and positioning belt, Weight capacity is 300 lbs.
 - Required for R-Net Remote Joysticks and Specialty Drive Controls.
 - Requires an R-Net Controller and Harness - order under Single Power Option R-Net.
 - Must order R-Net Color Remote Joystick with these options.
 - Must select Joystick for Lights feature or Omni for operation.
 - Seat to floor height is measured at the front of the seat pan based on a standard seat size. Seat to floor heights may vary +/- .25" dependent on conditions of the power base and on seat depths greater than 18".

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

"RFQ" = Request For Quotation - Please call customer service at 1-888-433-6818 and/or fax 1-888-433-6834

SINGLE FUNCTION SYSTEM

NOTE: Systems Include Seat Pan, Back Pan, Standard Armrests and Front Rigging. They are designed for use with a 2"-3" cushion.

F-Series CG Tilt Only System (50°Tilt)
(HCPCS code E1002) Part # FTR\$4,700

5° Pre-Tilt⁽¹⁾ Part # PRETILTN/C

NOTES:

1. This option increases seat-to-floor height by approximately 1.5"-2" minimum. Please call for confirmation.

ELEVATING SEAT MODULE

Elevating Scissor Lift Module (Scissor Mechanism)
12" Elevation - 250lbs max. Part # ESM12\$2,995

NOTES:

1. Uses an integrated Tilt/Lift mechanism, tilt is 45° and the seat to floor height is not affected by the Elevating Lift Seat.
2. Not available on seat sizes greater than 20" wide by 20" deep, call for RFQ.
3. For weights greater than 175lbs tilt while elevated is not available, call Customer Service for details.

STANDARD SEAT SPECIFICATIONS

Seat Width - Width adjustable range is 15"-21"

15" 16" 17" 18" 19" 20" 21"
Part# SWF15 SWF16 SWF17 SWF18 SWF19 SWF20 SWF21

22"
Part# SWF22\$595

Seat Depth

15" 16" 17" 18" 19" 20" 21"
Part# SDF15 SD16 SD17 SD18 SD19 SD20 SD21

22"\$595
Part# SD22

Select depth adjustability range: From 15" to 19" Part# SDA1519
 From 18" to 21" Part# SDA1821

SEAT CUSHION OPTIONS

- Matrx-Libra Seat Cushion* 16" w-20" w**(HCPCS Code E2624) \$485
- Matrx-Libra Seat Cushion* 21" w**(HCPCS Code E2624)\$525
- Matrx-Libra Seat Cushion* 22" w**(HCPCS Code E2625)\$595
**comes standard with Startex fabric, optional Spacetex fabric available below*
- Matrx-PS Seat Cushion 15" w-20" w** (HCPCS Code E2605) ..\$345
- Matrx-PS Seat Cushion 21" w-22" w** (HCPCS Code E2606) ..\$695
- Matrx-Vi Seat Cushion 15" w-20" w** (HCPCS Code E2607) ..\$425
- Matrx-Vi Seat Cushion 21" w-22" w** (HCPCS Code E2608) ..\$475
- Matrx-Flo-Tech Cushion 14" w-22" w**(HCPCS Code K0734) ..\$475
- "ON-CHAIR" Cushion Credit**(\$20)

FABRIC OPTIONS FOR CUSHIONS

- Spacetex Fabric Option for Libra Cushion**N/C
- Onyx Fabric Upgrade for Matrx-Vi or PS Cushion**\$100

LAP BELTS AND CHEST STRAP

- Padded Lap Belt 2 point⁽¹⁾** Part # PLB2 (HCPCS code E0978) \$195
 - Padded Lap Belt 4 point⁽¹⁾** Part # PLB4 (HCPCS code E0978) \$245
 - Push Button Style Seat Belt**STD
- Select: 48" Part # BELT48 60" Part # BELT60 71" Part # BELT71

NOTE:

1. 67" length only option.

BACK OPTIONS

- Standard Rehab Back and Back cushion** Part # SRB
Finished height from 18" to 27"N/C
- Additional Super Soft/HR Foam⁽¹⁾** Part # SS\$100
- Matrx Fabric Upcharge** Part # MFU\$100
- OMIT: Rigid Back Pan & Cushion⁽⁴⁾** Part # ORBC(\$75)
- OMIT: Back Cushion Only** Part # OBC(\$25)

FINISHED BACK HEIGHT⁽¹⁾

Finished Back height for Tilt Systems⁽¹⁾

18" 19" 20" 21" 22" 23" 24" 25"
Part# FBT18 FBT19 FBT20 FBT21 FBT22 FBT23 FBT24 FBT25

NOTES:

1. Comes with standard cane height.

UPGRADEABLE BACK OPTIONS⁽¹⁾

- Contour Back 16"-21" wide** (HCPCS Code E2613)\$495
Specify actual back height dimension: 12" 15" 17" 21"
- Contour Back 22" wide** (HCPCS Code E2614)\$595
Specify actual back height dimension: 15" 17" 21"
- High Back 16"-21" wide** (HCPCS Code E2620)\$595
Specify actual back height dimension: 16" 18" 20"
- Matrx PB Elite Back 15"-20" wide (3" Contour)**\$625
(HCPCS Code E2620) Specify actual back height dimension:
 10" 12" 14" 16" 18" 20"
- Matrx PB Elite HD Back 21"-22" wide**\$795
(HCPCS Code E2620 for 21" w E2621 for 22" w)
Actual back height dimension: 16" 20"
- Matrx PB Elite TR Back 16"-20" wide** (HCPCS Code E2620) ..\$650
Actual back height dimension: 16" 18" 20"
- Matrx PB Elite Deep Back 15"-21" wide (6" contour)** ..\$760
(HCPCS Code E2620) Specify actual back height dimension:
 12" 16" 18" 20"
- Matrx PB Back 16"-20" wide (5" Contour)**\$625
(HCPCS Code E2620) Specify actual back height dimension:
 12" 16" 20"
- Matrx PB Deep Back 16"-20"(7" Contour)**\$645
(HCPCS Code E2615) Specify actual back height dimension:
 12" 16" 20"
- Matrx PB Heavy Duty Back (5" Contour)**\$835
(For Systems 20"-22" wide and over 300lbs)(HCPCS Code E2615 for 20" w E2616 for 22" w)
Specify actual back height dimension: 16" 20"
- Matrx PB Heavy Duty Deep Back (7" Contour)**\$935
(For Systems 20"-22" wide and over 300lbs)(HCPCS Code E2615 for 20" w E2616 for 22" w)
Specify actual back pan dimensions: 16" 20"
- "ON-CHAIR" Matrx Back Credit**(\$50)
- Super Soft Foam for Elite Back Upgrade**\$100
- Onyx Fabric Upgrade** Contour Back High Back ..\$100

NOTES:

1. 4" gap is recommended for Tilt systems to allow for the seat cushion.

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

"RFQ" = Request For Quotation - Please call customer service at 1-888-433-6818 and/or fax 1-888-433-6834

TILT ONLY BACK CANE SPECIFICATIONS

Standard Tilt Back Canes⁽¹⁾ - Measured from Seat pan to top of back cane

- 22" 24"N/C
 Part # SBC22 SBC24
 Select Cane Angle 75° 79° 85° 90° 95° 101° 106° 112° 116°
 Part # BA75 BA79 BA85 BA90 BA95 BA101 BA106 BA112 BA116

- Angled (5°) Tilt Back Canes⁽¹⁾ 22" 24"N/C
 Part # ABC22 ABC24
 Select Cane Angle 81° 84° 91° 96° 101° 106° 111° 117° 121°
 Part # BA81 BA84 BA91 BA96 BA101 BA106 BA111 BA117 BA121

NOTES:
 1. Back angle is adjustable from in 5 degree increments.

MATRIX ELAN HEADREST PADS (HCPCS code E0955)

- SELECT ONE:
 Matrix Elan Standard Pad\$305
 Part # 6" Pad ESP6 10" Pad ESP10 14" Pad ESP14
 Matrix Elan Standard Pad with Infection Control\$305
 Part # 6" Pad ESP6-IC 10" Pad ESP10-IC 14" Pad ESP14-IC
 Matrix Elan Occipital Pad\$305
 Part # 9" Pad E0P9 12" Pad E0P12
 Matrix Elan Occipital Pad with Infection Control\$305
 Part # 9" Pad E0P9-IC 12" Pad E0P12-IC
 Matrix Elan 4-Point Pad 11"W x 10"H\$305
 Part # Standard Cover E4POINT Infection Control E4POINT-IC

MOTION CONCEPTS HEADREST PADS

- Motion Concepts Standard Part # MCSH\$345
 Motion Concepts Auto Style⁽¹⁾ Part # MCOH\$435

STEALTH HEADRESTS

- Stealth All Positioning -\$315
 Large Part # SAPL X-Large Part # SAPXL
 Stealth Nino Comfort Plus 6" Part # STSHN\$315
 Stealth Comfort Plus\$315
 8" Part # STSHXS 10" Part # STSHS 14" Part # STSHL
 Stealth's World's Best Removable Mounting
 (HCPCS code E1028) Part # WBHU\$285

HEADREST MOUNTING HARDWARE

- SELECT ONE:
 Fixed Mounting (non-removable)⁽¹⁾⁽²⁾ Part # FMN/C
 Elan Headrest Hardware - Multi-Axis Removable
 (HCPCS code E1028) Part # MEHW\$280
 Motion Concepts Multi-Axis Removable Mounting
 (HCPCS code E1028) Part # MAHU\$250
 OMIT: Headrest Part # OHRN/C

NOTES:
 1. Fixed Mounting hardware is not available on Matrix PB Backs.
 2. Unless selecting a "Removable" hwr option, this Headrest will come standard w/ Fixed Straight Mounting Hardware (non-removable).

ARMRESTS

- Adult Dual Post Adjustable Height With Quick Height
 Adjustment Lever (pair)STD
 Select Size: (9.75"-12.75") (12.75"-15.75") (15.75"-18.75")
 Part # DP1 DP2 DP3
 Select: Flat pouch Glove Boxes
 Part # FP GB
 Add: OutBack Arm Option⁽¹⁾ (each) Specify:
 Select: Right Part # OAR Left Part # OAL\$285
 1" Offset Armrest Spacer⁽²⁾ (each)\$125
 Left Part # OAS2L Right Part # OAS2R
 1/2" Offset Armrest Spacer⁽²⁾ (each)\$125
 Left Part # OAS1L Right Part # OAS1R

NOTES:
 1. Allows outward rotation while flipping back the armrest.
 2. Allows seat width to increase by 1" or 2" depending on offset.

ARMPADS

- Modular Arm Pads⁽¹⁾ - (each)\$75
 Left Full(14") Part # MAPTFL Desk(10") Part # MAPTDL
 Must select Pad style
 Iskin - Integrated soft skin foam Part # ISSFL
 Startex Covered Visco foam Part # SCVFL
 Right Full(14") Part # MAPTFR Desk(10") Part # MAPTDR
 Must select Pad style
 Iskin - Integrated soft skin foam Part # ISSFR
 Startex Covered Visco foam Part # SCVFR

- Standard Pads-(each)STD
 Select: Full (2.25"x14")⁽²⁾ Left Part # SAPFL Right Part # SAPFR
 Desk (2.25"x10") Left Part # SAPDL Right Part # SAPDR

- Waterfall Arm Pads-(each)\$45
 Select: Full (2.25"x14")⁽²⁾ Left Part # WAPFL Right Part # WAPFR
 Desk (2.25"x10") Left Part # WAPDL Right Part # WAPDR

- Flat Pads-(each)\$75
 Select: Full (3"x14") Left Part # MPAFL Right Part # MPAFR
 Desk (3"x10") Left Part # MPADL Right Part # MPADR

- Gel Pads-(each) Select:
 Narrow (2"x12") Left Part # GAPNL Right Part # GAPNR ..\$100
 Wide (3.5"x12") Left Part # GAPWL Right Part # GAPWR ..\$110
 Long (3.5"x14") Left Part # GAPLL Right Part # GAPLR ..\$120

- Ergonomic Arm Trough-LEFT(4.5"x13") Part # EATL\$110
 (HCPCS Code E2209)
 Ergonomic Arm Trough-RIGHT(4.5"x13") Part # EATR ..\$110
 (HCPCS Code E2209)
 Flat Hand Pad-LEFT(4.5"x5") (HCPCS Code K0108) Part # FHPL ..\$75
 Flat Hand Pad-RIGHT(4.5"x5") (HCPCS Code K0108) Part # FHPR ..\$75

NOTES:
 1. Requires interchangeable pad insert.

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

"RFQ" = Request For Quotation - Please call customer service at 1-888-433-6818 and/or fax 1-888-433-6834

ARMPADS - ACCESSORIES

- Elbow Block LEFT⁽¹⁾**\$260
 Select Pad Size:
XS(3.25"x4.25") S(3.5"x5.25") M(3.75"x5.75") L(4.25"x6.5")
 Part # EBXL EBSL EBML EBLL
 Vertical Extension: S(3") L(5") Horizontal: S(5") L(7")
 Part # V3 V5 H5 H7
Elbow Block RIGHT⁽¹⁾\$260
 Select Pad Size:
XS(3.25"x4.25") S(3.5"x5.25") M(3.75"x5.75") L(4.25"x6.5")
 Part # EBXR EBSR EBMR EBLR
 Vertical Extension: S(3") L(5") Horizontal: S(5") L(7")
 Part # V3 V5 H5 H7
Multi Axis Upper Extremity Support⁽²⁾ (each)\$260
 With Quick Adjustment Lever
Left Part # MACESL Right Part # MACESR

- NOTES:**
 1. Prevents arm from slipping.
 2. Allows rotation of the armpad, only compatible with Gel pad, Flat pad or Ergonomic Arm Troughs.

MATRIX FIXED LATERALS - FOR ELITE AND ELITE TR FULL BACKS

- Matrix Standard Fixed LEFT**(HCPCS code E0956)\$120
 Select pad : Standard (7.5" x 5") Part # PBLSL
 Small (5.5" x 3") Part # PBSSL
 Matrix Standard Fixed RIGHT(HCPCS code E0956)\$120
 Select pad : Standard (7.5" x 5") Part # PBLSR
 Small (5.5" x 3") Part # PBSSR
 Matrix Offset Fixed LEFT(HCPCS code E0956)Part # PBOLL . . . \$130
 Matrix Offset Fixed RIGHT(HCPCS code E0956)Part # PBOLR . \$130

MATRIX SWING AWAY LATERALS - FOR ELITE AND ELITE TR FULL BACKS

- Matrix Lateral Trunk Support LEFT pad (each)**\$115
 (HCPCS code E0956)
 Select Size:
XS(3.25"x4.25") S(3.5"x5.25") M(3.75"x5.75") L(4.25"x6.5")
 Part # EPADXL EPADSL EPADML EPADLL
 Matrix Lateral Trunk Support RIGHT pad (each)\$115
 (HCPCS code E0956)
 Select Size:
XS(3.25"x4.25") S(3.5"x5.25") M(3.75"x5.75") L(4.25"x6.5")
 Part # EPADXR EPADSR EPADMR EPADLR

SELECT SWING AWAY HARDWARE TO CONVERT ABOVE:

- Matrix Elite Back Swing Away Hardware**
 LEFT (HCPCS code E1028) Part # SAHEBL\$215
 RIGHT (HCPCS code E1028) Part # SAHEBR\$215
 Matrix Offset Elite Swing Away Hardware
 LEFT (HCPCS code E1028) Part # OSAHL\$215
 RIGHT (HCPCS code E1028) Part # OSAHR\$215

LATERALS - FOR ALL OTHER BACKS (EXCLUDING MATRIX PB ELITE AND PB BACKS)

- Fixed Mounted Lateral Trunk Support LEFT**
 Pad Size: XS(3.25"x4.25") S(3.5"x5.25") M(3.75"x5.75") L(4.25"x6.5")
 Part # LATXL LATSL LATML LATLL
 (HCPCS code E0956).....\$115
 Fixed Mounted Lateral Trunk Support RIGHT
 Pad Size: XS(3.25"x4.25") S(3.5"x5.25") M(3.75"x5.75") L(4.25"x6.5")
 Part # LATXR LATSR LATMR LATLR
 (HCPCS code E0956)\$115

SELECT SWING AWAY HARDWARE TO CONVERT ABOVE SUPPORTS:

- Maxx Style Swingaway Multi-Adjustable Mounting Hardware⁽¹⁾**
 (Standard with 1" telescoping link gives 0-2.75" medial offset⁽²⁾)
 LEFT (HCPCS code E1028) Part # SALATL\$250
 RIGHT (HCPCS code E1028) Part # SALATR\$250
 Additional Telescoping Link(ea) Quantity () . . . \$100
 (Adding an additional 1" telescoping link gives up to 5.5" medial offset⁽²⁾)
LEFT Part # ATLL RIGHT Part # ATLR

- NOTES:**
 1. These laterals are not available with Matrix Elite Deep and PB Deep Backs.
 2. The amount of offset available may be impacted by back thickness.

FIXED HIP/KNEE SUPPORTS

- LEFT Lateral Hip Support Pad with Fixed Mounting** \$145
 (HCPCS code E0956)
 Select Size:
 Mounted at Hip (5.5"W x 4"H) (7"W x 6"H) (9"W x 4"H) (13"W x 4"H)
 Part # HSPXLH HSPSLH HSPMLH HSPLLH
 Mounted at Knee (5.5"W x 4"H) (7"W x 6"H) (9"W x 4"H) (13"W x 4"H)
 Part # HSPXLK HSPSLK HSPMLK HSPLLK
RIGHT Lateral Hip Support Pad with Fixed Mounting . . . \$145
 (HCPCS code E0956)
 Select Size:
 Mounted at Hip (5.5"W x 4"H) (7"W x 6"H) (9"W x 4"H) (13"W x 4"H)
 Part # HSPXRH HSPSRH HSPMRH HSPLRH
 Mounted at Knee (5.5"W x 4"H) (7"W x 6"H) (9"W x 4"H) (13"W x 4"H)
 Part # HSPXRK HSPSRK HSPMRK HSPLRK

REMOVEABLE AND SWING AWAY HARDWARE UPGRADE

HARDWARE TO CONVERT ABOVE HIP SUPPORT TO UPGRADE TO REMOVEABLE/SWINGAWAY HIP OR KNEE SUPPORTS

- Lift Off Removable Hardware for Hip Support (Hip position only)**
 (HCPCS code E1028)
 LEFT Part # LOHL\$205
 RIGHT Part # LOHR\$205
 Swingaway Removable Hardware for Hip Supports
 (HCPCS code E1028)
 LEFT Mounted at Hip Part # SHSHLH\$265
 RIGHT Mounted at Hip Part # SHSHRH\$265
 LEFT Mounted at Knee Part # SHSHLK\$265
 RIGHT Mounted at Knee Part # SHSHRK\$265

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

"RFQ" = Request For Quotation - Please call customer service at 1-888-433-6818 and/or fax 1-888-433-6834

STANDARD FRONT RIGGING

CALL TO CONFIRM ON SEAT SIZES BELOW 16"W X 16"D

Legrest Receivers OnlySTD

<p>**Please Specify Front rigging to be used on TRx System⁽¹⁾</p>

Fixed Center Mount Foot Platform 70° 90°N/C
 Seat pan to footplate: S (9"-15") L (16"-21") Set at _____

OR

Invacare 70° Swingaway Footrests (Pair)⁽³⁾N/C

NOTES:

- Standard receivers will be provided. Please provide front rigging or order separately from legrest section.
- Comes with rubber coated foot platform.
- Includes composite footplates and are not covered by Motion Concepts Lifetime warranty.

POWER CENTER MOUNT FOOT PLATFORM OPTIONS

MUST BE COMPLETED	Measurement in Inches
Knee to Heel:	
Cushion Thickness: (Subtract)	
Cushion Compression: (Add)	
Seat Pan to Footplate: (Total)	

CENTER MOUNTS ARE RESTRICTED TO SEAT DEPTHS 16" AND GREATER**

LNX Power Center Mount Foot Platform⁽¹⁾ Part # LNX \$2,795
 (HCPCS code E1012) 13"-19" unlimited adjustabilityN/C

Add: One Piece Malleable Calf Panel Part # 0PMCP\$250

NOTES:

- Includes as standard individual calf pads and rubber coated foot platform. There is 8" of anti-shear extension.
- 2.5" of ground clearance is required.

MANUAL CENTER MOUNT FRONT RIGGING

Maxx Style Fixed Center Mount Foot Platform⁽¹⁾⁽²⁾ Part # MSFCMFP \$595

With Quick Release Height Adjustment

Set at 70° 90° 97°

Part # SA70 SA90 SA97

Seat pan to footplate: 13"-19"

NOTES:

- Depth and height adjustable. Includes individual flat calf pads and a rubber coated foot platform.


FOOT PLATFORM OPTIONS FOR MAXX/LNX CMFP



Small Flip Foot Platform⁽¹⁾
 (12"W x 8.5"D) Part # SFFP \$235



Medium Flip-up Foot Platform⁽¹⁾
 (12"W x 10.75"D) Part # MFFP \$240



Large Flip-up Foot Platform⁽¹⁾
 (14.5"W x 13"D) Part # LFFP \$250



Extra Large Flip-up Foot Platform⁽¹⁾
 (17"W x 13"D) Part # XLFFP \$325



Extra Extra Large Flip-up Foot Platform⁽¹⁾
 (21"W x 13"D) Part # XXLFFP \$325




NOTE:

- May cause interference with certain bases and seat widths call to confirm configuration.

FOOT PLATE OPTIONS FOR CENTER MOUNTS

Individual Footplates for LNX and Maxx Style Fixed Center Mount

In cast aluminum with built in heel cups and rubber mat



Left: S(5"W x 7.5"D) M(5.5"Wx9.5"D) L(6"Wx11.5"D) .. \$175
 Part # IFSL IFML IFLL

Right: S(5"W x 7.5"D) M(5.5"Wx9.5"D) L(6"Wx11.5"D) .. \$175
 Part # IFSR IFMR IFLR

POWER FRONT RIGGING OPTIONS

Maxx Style Swing-Away Power Elevating Pivot Legrests⁽¹⁾ Part # MSPPL \$2,750
 (HCPCS code E1010)

Select Receivers: Straight Flared outward 1"
 Part # RS Part # RF

Seat pan to footplate: S (10"-14") M (13"-17") L (16"-20")
 Part # SPTF5 SPTF2 SPTF7

Select Footplates on page 6 under Foot plate options:

Select Function: Individual⁽²⁾ Combined Synchronized W/Recline
 Part # INDA COMA SYNA

NOTES:

- Includes Swing Away receivers, Pivot Plus Legrests and actuators. Legs swing away manually
- Individual legs count as two functions.

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

"RFQ" = Request For Quotation - Please call customer service at 1-888-433-6818 and/or fax 1-888-433-6834

MANUAL ELEVATING LEGRESTS

- Maxx Style Swing-Away Manual Elevating Pivot Legrests** Part # MMPPSA\$720
- Select Receivers: Straight Part # RS Flared outward 1" Part # RF
- Seat pan to footplate: S (10"-14") Part # SPTF5 M (13"-17") SPTF2 L (16"-20") SPTF7

Select Footplates on next page under Foot plate options:

SWING-AWAY FRONT RIGGING OPTIONS

- Heavy Duty 70° Swing Away Footrests⁽¹⁾ (Pair)**
- Part # HD70\$495
- Select Receivers: Straight Part # RS Flared outward 1" Part # RF
- Seat pan to footplate: XS (6"-11") Part # SPTF11 S (11"-15") SPTF8 M (13"-17") SPTF2 L (16"-20") SPTF7
- Calf pads for HD Swing Away Footrests (pair)**\$275
- Select: Flat Pads Part # HD70FCP Curved Calf pads (inside width is 6") Part # HD70CCP


Select Footplates on page 8 under Foot plate options:

- Residual Limb Support (6.75"W x 9.75"L x 2"H)⁽²⁾**\$495 (HCPCS code E1020)
- Select: Left Part # RLSL Right Part # RLSR
- Dynamic Coil Springs (Pair) for Heavy Duty 70° Swingaway Footrests** Part # DCS\$190

- NOTE:**
1. Heavy duty durable construction (4.5lbs ea.).
 2. Swing away and removable, 50° of angle adjustment, 3.5" height adjustment and 2" of for/aft adjustment.

FOOT PLATE OPTIONS FOR ELEVATING AND SWING AWAY FRONT RIGGING⁽¹⁾

SELECT FOOTPLATE SIZE FOR BELOW HARDWARE:




- Adjustable Angle Foot Plates(Each)**(HCPCS code K0040)\$175
- Left S(4.25"W x 8"D) M(5.25"Wx8"D) L(6.25"Wx8"D) XL(9"Wx11"D)
- Part # AFSL AFML AFLL AFXL
- Right S(4.25"W x 8"D) M(5.25"Wx8"D) L(6.25"Wx8"D) XL(9"Wx11"D)
- Part # AFSR AFMR AFLR AFXR

- Multi-Axis Adjustable Angle Foot Plates(Each)**(HCPCS code K0040)\$225
- Left S(4.25"W x 8"D) M(5.25"Wx8"D) L(6.25"Wx8"D) XL(9"Wx11"D)
- Part # MFSL MFML MFLM MFXL
- Right S(4.25"W x 8"D) M(5.25"Wx8"D) L(6.25"Wx8"D) XL(9"Wx11"D)
- Part # MFSR MFMR MFLR MFXR

- Heel Loops (Each)**(HCPCS code E0951)\$20
- Left Part # HLL Right Part # HLR

FOOT PLATE OPTINS FOR ELEVATING AND SWING AWAY FRONT RIGGING⁽¹⁾



- Adjustable Angle Foot Plates(Each)**(HCPCS code K0040)\$175
- Left S(7.25"W x 9"D) M(7.5"W x 11"D) L(10"W x 14"D)
- Part # AFCSL AFCSM AFCSL AFCLL
- Right S(7.25"W x 9"D) M(7.5"W x 11"D) L(10"W x 14"D)
- Part # AFCSR AFCSR AFCLR

- Multi-Axis Adjustable Angle Foot Plates(Each)**(HCPCS code K0040) ...\$225
- Left S(7.25"W x 9"D) M(7.5"W x 11"D) L(10"W x 14"D)
- Part # MFCSL MFCSM MFCSL MFCLL
- Right S(7.25"W x 9"D) M(7.5"W x 11"D) L(10"W x 14"D)
- Part # MFCSR MFCSR MFCLR

- Heel Loops (Each)**(HCPCS code E0951)\$20
- Left Part # HLL Right Part # HLR

- NOTE:**
1. Certain footplate sizes may not fit all seat widths depending on legrests selected.

CALF PANEL AND STRAP OPTIONS⁽¹⁾

- Airmesh Calf Strap for chairs 15"-18" wide** - 3.5" pad height
- Part # CSM\$135
- Airmesh Calf Strap for chairs 18"-22" wide**-3.5" pad height
- Part # CSL\$135
- Airmesh Calf Panel for chairs 15"-18" wide** 35.5"w x 9"h pad
- Part # CPM\$195
- Airmesh Calf Panel for chairs 18"-22" wide** 38.5"w x 9"h pad
- Part # CPL\$210

ADDITIONAL OPTIONS AND ACCESSORIES

- Concepts Back Pack -** Single Part # SBP\$160
- Double Part # DBP\$255
- Transfer Handles (Height Adjustable)**\$455
- Select: 4"-6.5" Part # THS 9"-11.5" Part # THL
- Flip & Fold Away 1/2 Lap Tray⁽²⁾** (14"W x 10"D)\$995
- Select Mounting: Left Part # FFLTL Right Part # FFLTR

- NOTES:**
1. May increase seat to floor heights on the certain bases.
 2. Not available on joystick side and only available with a flat arm pad or modular arm pad.

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

"RFQ" = Request For Quotation - Please call customer service at 1-888-433-6818 and/or fax 1-888-433-6834

ELECTRONICS SWITCH MOUNTING

Select: Left Right Do Not Mount Specialty Controls*

Part # ESML ESMR ESMO ESMS

SINGLE-FUNCTION ELECTRONICS

Single/Dual Function Through Switch/Toggle⁽¹⁾STD

- Please Select One: Single Push Button Part # SFCBSPB
- Single Toggle Part # SFCBST
- Single Toggle/Bat extension Part # SFCBSTB
- Dual Push Button Part # SFCBDPB
- Omit Switch Part # SFCBOS

Add: Reduced Drive for Elevating Seat Part # RDE \$400

Single/Dual Function Through Joystick⁽¹⁾⁽²⁾ Part # SFCBJ \$1,200
(HCPCS code E2310)

Add: Reduced Drive for Elevating Seat Part # RDE \$400

- Additional Switch: Single Push Button Part # SPBS \$97
- Single Toggle Part # ST
- Single Toggle/Bat extension Part # STB
- Dual Push Button Part # DPB

- NOTE:**
1. These electronics include drive lock-out and tilt limit.
2. Toggle action only.

MULTI-FUNCTION ELECTRONICS

Multi Function Through Switches/Toggles⁽¹⁾ \$1,550

- Select: Quad Push Button Part # MFCBQPB
- Four-Way Toggle Part # MFCBFWT
- 8-Way Rocker⁽²⁾ Part # MFCBEWR

Add: Reduced Drive for Elevating Seat Part # RDE \$400

Multi Function Through Drive Control⁽¹⁾ \$2,600
(HCPCS code E2311)

- Through Joystick Part # MFCBJJ
- Through Specialty Control* Part # MFCBJS

*Please contact Customer Service for a specialty control configuration sheet.

Includes: Reduced Drive for Elevating Seat Yes Part # RDY
 No Part # RDN

- Additional Switch: Quad Push Button Part # QPB \$97
- Four-Way Toggle Part # FWT \$97
- 8-Way Rocker Part # EWR \$150

- NOTE:**
1. These electronics include attendant control, drive lock-out, tilt/recline limit and actuator speed control.
2. Dual function is not available with 8-Way Rocker.

ELECTRONICS OPTIONS

M16 Splitter⁽¹⁾ (HCPCS code E2399) Part # M16S \$165
(Allows for multiple switch options with phono plug connection)

Stealth Egg Switch-Single function \$105

Green Part # SESG Black Part # SESB

Power Shut Off Switch -Toggle Switch Only Part # PSOS . . . \$230
(HCPCS code E2399)

Lazarus Jump Starter Part # LJS \$150
(HCPCS code A9270)

Accessory Port Power Supply (3 amp capacity) Part # APPS . \$195

- NOTE:**
1. Only available with R-Net.

BATTERIES

MK GROUP 34 Battery (M34 SLD G) Each
Quantity () (HCPCS code E2359) \$330

INSTALLATION

Installation by Motion Concepts N/C

Total Retail Value \$ _____

OTHER
