


**REQUIRED INFORMATION**

Date of Order: \_\_\_\_\_ Dealer Account #: \_\_\_\_\_  
 Dealer Name: \_\_\_\_\_  
 PO #: \_\_\_\_\_ Tag \_\_\_\_\_  
 Base PO #: \_\_\_\_\_  
 After Market Items shipping to Motion Concepts:  Yes  No  
 List Items: \_\_\_\_\_  
 PO #: \_\_\_\_\_  
 Name of the Manufacturer: \_\_\_\_\_  
 Purchasing Contact: \_\_\_\_\_  
 Phone: \_\_\_\_\_ Fax: \_\_\_\_\_  
 E-mail: \_\_\_\_\_  
 RTS/Therapist: \_\_\_\_\_  
 Ship to Address: \_\_\_\_\_  
 City: \_\_\_\_\_ State: \_\_\_\_\_  
 Zip: \_\_\_\_\_ Client Gender:  M  F  
 Special Client Conditions: \_\_\_\_\_

**REQUIRED MEASUREMENTS**

To ensure system is accurately configured please fill in all required measurements.

Client Height:	
Client Weight (lbs):	
A. Seat to Shoulder:	
B. Trunk Depth:	
C. Chest Width:	
D. Knee to Back:	
E. Seat to Top of Head:	
F. Elbow to Hand:	
G. Seat to Elbow:	
H. Hip Width:	
I. Knee to Heel:	
Cushion Thickness:	



*Weight Capacity is up to 300lbs unless specified otherwise*



<b>HCPCS CODES</b>	K0848/, K0849
<b>SUSPENSION SYSTEM</b>	Active Ride Control (ARC)
<b>BASE WIDTH</b>	23.25"
<b>BASE LENGTH</b>	36.5"
<b>DRIVE WHEEL SIZE <sup>(1)</sup></b>	14"
<b>CASTER WHEEL SIZE (front/rear)</b>	6"
<b>TURNING RADIUS</b>	20"
<b>DRIVETRAIN</b>	4-pole EAD Motors

<b>BATTERIES</b>	2 x M34
<b>BATTERY CHARGER</b>	8A, Off Board
<b>BATTERIES WEIGHT (each)</b>	38.5 lbs (M34)
<b>RANGE <sup>(2)</sup></b>	17.7 miles
<b>WEIGHT CAPACITY</b>	300 lbs
<b>BASE WEIGHT (w/o batteries)</b>	157 lbs
<b>GROUND CLEARANCE</b>	3"
<b>MAXIMUM SPEED</b>	6.2 mph

**BASE VITALS**

<sup>(1)</sup> All X3 wheels come with solid inserts  
<sup>(2)</sup> Actual driving range and speed may vary due to factors such as user weight, type and grade of terrain, battery condition & charge level; type and condition of the drive wheels, and variations in the drive and control systems.

**SOLID SEAT PAN - K0848**



**CAPTAIN SEAT - K0849**



**Contoured Static Seating on ROVI X3 Base**

Specific Client Measurements Required - Please call customer service at 1-888-433-6818 and/or fax 1-888-433-6834 to ensure that we have all the client details we require for this type of option.

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**BASE COLOR**

Select Color: Standard With Lights<sup>(1)</sup> . . . . . **\$1,195**

- Blizzard White  Part # 800914-00200-03  Part # 800914-00201-03
- Matte Black  Part # 800914-00200-50  Part # 800914-00201-50
- Solar Yellow  Part # 800914-00200-69  Part # 800914-00201-69
- Motion Green  Part # 800914-00200-36  Part # 800914-00201-36
- Original Orange  Part # 800914-00200-23  Part # 800914-00201-23
- Bright Blue  Part # 800914-00200-40  Part # 800914-00201-40
- Red Delicious  Part # 800914-00200-21  Part # 800914-00201-21
- Party Pink  Part # 800914-00200-2C  Part # 800914-00201-2C
- Purple Reign  Part # 800914-00200-80  Part # 800914-00201-80

Matching Battery Doors (select fender color first) . . . **\$495**

**NOTES:**

1. Fully integrated LED lighting kit. Includes headlights, taillights and turn signals. Must select Joystick for Lights feature or Omni for operation.



**CUSTOM FENDER OPTIONS**

Custom Fender Design by Bob Rice . . . . . **\$500**

Please refer to Page 4 for Custom fender options by Bob Rice

- Standard Fender Configuration . . . . . N/C
- Fenders w/Lights Configuration . . . \$1195



Dipped Fender Designs . . . . . **\$395**

Please refer to Page 5 for Custom Dipped fender options

- Standard Fender Configuration . . . . . N/C
- Fenders w/Lights Configuration . . . \$1195



- Front/Rear Battery Doors . . . . . **\$750**
- Front Battery Doors . . . . . **\$395**
- Rear Battery Doors . . . . . **\$395**
- Fenders & Front/Rear Battery Doors . . . . . **\$1100**

**TIRE OPTIONS**

- Standard Solid Tires Part # STIRES **STD**
- Pneumatic Tires Part # PTIRES **\$199**

**MILITARY DECAL**

**3" Diameter Decal (Uninstalled) . . . . . \$50**



- Army  Navy
- Air Force  Marines

**FULL CAPTAIN-STYLE SEAT - K0849**

All bases below come standard with Full Captain-Style Seat in Grey or Black Vinyl and 22" high backrest. Includes a Semi-Recline hinge of 120° with left side only adjustable lever. Standard adjustable height, flip-back armrests and a foot platform that measures 12" wide by 9" deep.

ROVI X3<sup>(1)</sup> HCPCS Code: K0849 Part # 888WNLEX3-SP-SC. . . . . **\$6,695**  
Select One:

- Standard 5.0 mph High Torque Motors . . . . . **N/C**  
Part # HTM5
- Optional 5.8 mph Quick Speed Motors<sup>(2)</sup> . . . . . **\$200**  
Part # GSM58
- Optional 6.2 mph High Speed Motors<sup>(2)</sup> . . . . . **\$1,550**  
with OnTraxx Enhanced Driving Module Part # HSM62

- Select:  Grey Seat Size:  18"x18" Part # V18171/222RG  
 20"x20" Part # SV2121211/2RG
- Black Seat Size:  18"x18" Part # SV18171/222RB  
 20"x20" Part # SV2121211/2RB
- Select STF<sup>(2)</sup>:  20"  21"  22"  23"

**NOTES:**

1. Price of the powerbase includes: Active Ride Control suspension, tie-down hardware for 4-point strap system, 8 amp off-board charger, flat free inserts and positioning belt, Weight capacity is 300 lbs.
2. Seat to floor height is measured at the front of the seat pan based on a standard seat size. Seat to floor heights may vary +/- .25" dependent on conditions of the power base and on seat depths greater than 18".

**SOLID SEAT PAN & CAPTAIN-STYLE SEAT - K0848**

All bases above come standard with a Solid Seat Pan style Captain Seat in Grey or Black Vinyl and 22" high backrest. Includes a Semi-Recline hinge of 120° with left side only adjustable lever. Standard adjustable height, flip-back armrests and a foot platform that measures 12" wide by 9" deep.

ROVI X3<sup>(1)</sup> HCPCS Code: K0848 Part # 888WNLEX3-SP-SS. . . . . **\$6,550**  
Select One:

- Standard 5.0 mph High Torque Motors . . . . . **N/C**  
Part # HTM5
- Optional 5.8 mph Quick Speed Motors<sup>(2)</sup> . . . . . **\$200**  
Part # GSM58
- Optional 6.2 mph High Speed Motors<sup>(2)</sup> . . . . . **\$1,550**  
with OnTraxx Enhanced Driving Module Part # HSM62

- Select:  Grey Seat Size:  18"x18" Part # SP181821RG  
 20"x20" Part # SP202023RG
- Black Seat Size:  18"x18" Part # SP181821RB  
 20"x20" Part # SP202023RB
- Select STF<sup>(2)</sup>:  16.5"  17.5"  18.5"  19.5"

**NOTES:**

1. Price of the powerbase includes: Active Ride Control suspension, tie-down hardware for 4-point strap system, 8 amp off-board charger, flat free inserts and positioning belt, Weight capacity is 300 lbs.
2. Seat to floor height is measured at the front of the seat pan based on a standard seat size. Seat to floor heights may vary +/- .25" dependent on conditions of the power base and on seat depths greater than 18".

**Contoured Static Seating on ROVI X3 Base**

Specific Client Measurements Required - Please call customer service at 1-888-433-6818 and/or fax 1-888-433-6834 to ensure that we have all the client details we require for this type of option.







Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**SEAT CUSHION REFERENCE GUIDE**

						
						
Cushion	Meshtex	Startex	Infection	Spacetex	Fusion	Onyx
Libra	N/A	STD	\$50	N/C	N/A	N/A
PS	N/A	STD - Reversible	STD - Reversible	N/A	N/A	N/A
Vi	STD	\$50	\$50	N/A	\$150	\$100

**SEAT CUSHION OPTIONS - FOR SOLID SEAT**

- Matrx Libra Seat Cushion 18" w-20" w** HCPCS Code E2624 **.\$485**
- Matrx PS Seat Cushion 18" w - 20" w** HCPCS Code E2605 **..\$345**
- Matrx Vi Seat Cushion 18" w - 20" w** HCPCS Code E2607 **..\$425**

**FABRIC OPTIONS FOR CUSHIONS**

- Spacetex Fabric Option for Libra Cushion Part # SFLF . . . . . N/C
- Infection Control Fabric Option for Libra Cushion Part # ICFL . . . . . \$50
- Fusion Fabric Upgrade for Matrx Vi Cushion Part # FFV . . . \$150
- Onyx Fabric Upgrade for Matrx Vi Cushion Part # OCUMA . . . \$100
- Startex Fabric Upgrade for Matrx Vi Cushion
- Smooth Side up Part # SFVS . . . . . \$50
- Fabric Side up Part # SFVF . . . . . \$50

**SINGLE FUNCTION CONTROL OPTIONS**

- VR2 Non-Expandable Joystick** Part # D50680 . . . . . **STD**
- R-Net LED Remote Joystick<sup>(1)(2)</sup>** Part # D51122SC . . . . . **.\$295**
- R-Net LED Remote Joystick for Lights<sup>(1)</sup>** Part # D51108 . . . . . **.\$500**
- R-Net Color Remote Joystick with Monojacks<sup>(1)(2)</sup>** Part # D51165 . . . . . **.\$600**
- R-Net Color Remote Joystick with Monojacks for Lights<sup>(1)</sup>** Part # D50982 . . . . . **.\$650**

 **R-Net Advanced Joystick<sup>(1)</sup>** Part # D51646. . . . . **.\$1295**  
 3.5" high-resolution VGA color display | IR Control as standard | Optional, Integrated Bluetooth module available. **NOW \$600**

- R-Net Expandable Controller<sup>(3)</sup>** . . . . . **.\$615**  
HCPCS Code: E2377 Part # D51109
- R-Net Harness for Expandable Controller<sup>(3)</sup>** . . . . . **.\$415**  
HCPCS Code: E2313 Part # REH

- NOTES:**
1. Requires an R-Net Controller and Harness
  2. Light package not available with the joystick.
  3. Required for R-Net Remote Joysticks.

**JOYSTICK MOUNT**

- Square Joystick Mount Receiver** Part # 1343-00701-00. **STD**
  - Square Joystick Mount** Part # 1343-00707-00. . . . . **STD**
  - Swing Away Mount Bracket** Part # 1343-00709-A3 . . . . **STD**
- Select Mounting Position:  Left  Right  
 Part # ML MR

**JOYSTICK HANDLES**

- "T" Handle for Joystick Knob** . . . . . **.\$115**  
HCPCS Code: E2323 Part # SA80941
- Ball Handle for Joystick Knob** . . . . . **.\$115**  
HCPCS Code: E2323 Part # SA72441
- BodyPoint Joystick Handle** HCPCS Code: E2323 . . . . . **.\$120**
  - 3" U-Shaped Part # PC101B-20
  - 4" U-Shaped Part # PC102B-20
  - Domed Rubber Part # PC107B-20

**ELECTRONIC OPTIONS**

- R-Net Output Module(EMC)<sup>(1)(2)</sup>** . . . . . **.\$995**  
Part # D50883
- R-Net Bluetooth Module<sup>(1)(2)</sup>** . . . . . **.\$895**  
Select OS:  iOS Apple Part # D51580  Android/PC Part # D51111
- Smart Phone Holder** Part # SPH . . . . . **N/C**

- NOTES:**
1. Requires an R-Net Controller and Harness order above
  2. Must order R-Net Color Remote joystick with these options

**ARMRESTS (MUST BE SELECTED)**

- Adjustable Height Flip Back Arm Left** Part # 104300-00112. **.\$125**
- Adjustable Height Flip Back Arm Right** Part # 104400-00112 **.\$125**

**OTHER OPTIONS**

- Wheelchair Transport Brackets<sup>(1)</sup>** Part # XBKTS . . . . . **.\$250**

- NOTES:**
1. For unoccupied use only

**BATTERIES**

- 2 x MK GROUP 34 Batteries** Part # M34 SLD G (Two) . . . **.\$650**  
HCPCS Code: E2359

**Total Retail Value \$**

**Contoured Static Seating on ROVI X3 Base**

Specific Client Measurements Required - Please call customer service at 1-888-433-6818 and/or fax 1-888-433-6834 to ensure that we have all the client details we require for this type of option.

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**CUSTOM FENDER DESIGN BY BOB RICE**



Flames



Personalized



Camouflage



Silver Concrete



U.S. Navy



U.S. Coast Guard



U.S. Air Force



U.S. Marines



U.S. Army

**DESIGN DESCRIPTION**

(PLEASE USE IF YOU SELECTED PERSONALIZED OPTION)

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

**DEALER SIGN OFF**

*This concept is agreed to and accepted by:*

Dealer Name: \_\_\_\_\_

PO #: \_\_\_\_\_

Approved with no changes

Approved with changes

Review again

Date: \_\_\_\_\_

Signature: \_\_\_\_\_



**Contoured Static Seating on ROVI X3 Base**

Specific Client Measurements Required - Please call customer service at 1-888-433-6818 and/or fax 1-888-433-6834 to ensure that we have all the client details we require for this type of option.

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**CUSTOM DIPPED FENDER DESIGNS**



U.S. Flag



Camouflage



Electric Circuit



Flames



80's Green Rock



80's Black Rock



Pink Camo



Flower Power



Graphity



Animal Paw Print



**CUSTOM DIPPED BATTERY DOOR DESIGN**

Please select Fender option first, if battery doors are to match the fenders. If not make the design selection below

- Matching Front Battery Door ..... \$395
  - Matching Rear Battery Door ..... \$395
  - Front Battery Door ..... \$395
  - Rear Battery Door ..... \$395
- |   |   |
|---|---|
| <input type="checkbox"/> U.S. Flag        | <input type="checkbox"/> 80's Black Rock  |
| <input type="checkbox"/> Camouflage       | <input type="checkbox"/> Pink Camo        |
| <input type="checkbox"/> Electric Circuit | <input type="checkbox"/> Flower Power     |
| <input type="checkbox"/> Flames           | <input type="checkbox"/> Graphity         |
| <input type="checkbox"/> 80's Green Rock  | <input type="checkbox"/> Animal Paw Print |