

Dealer Name: \_\_\_\_\_  
 Dealer Account Number: \_\_\_\_\_  
 PO #: \_\_\_\_\_  
 Tag: \_\_\_\_\_

**PLEASE NOTE:**

These Accessory Additions are only available with the purchase of a new power positioning system and completed order for **TRD0335 Rev U** or **TRD0360 Rev R** with all Client Specifications completed. Subject to approval.

**LATERAL TILT MODULE<sup>(1)(2)(3)(4)(5)(6)</sup>**

**Power Lateral Tilt Module (up to 15°)** Part # LATM . . . \$3,995  
 All options must be selected on the system order form prior to selecting this option.

Reduced Drive while in Lateral Tilt?  Yes Part # LATDY  
 No Part # LATDM

**NOTES:**

1. Weight Capacity is 250 lbs.
2. Only available on 3G Storm, TDX-SP, 6000ZHD and Edge HD bases.
3. When ordering an lateral tilt module on a Storm short base frame, standard fixed casters must be ordered on the base.
4. Only available on 16v to 20" seat widths and 16" to 22" seat depths.
5. Only available with Multi-Function electronics.
6. Must order the Group 24 battery tray for the TDX-SP, which will utilize Group 22 batteries to accomodate the lateral tilt module.

**POWER ANTERIOR TILT MODULE 45°<sup>(1)(2)(3)(4)(5)(6)(7)</sup>**

**Power Anterior Tilt Module up to 45°** Part # PAT45 . . . \$3,995  
 All options must be selected on the system order form prior to selecting this option.

**Synchronized Option for Center Mount Front**  
**Rigging with Anterior Tilt<sup>(8)</sup>** Part # SOCM . . . . \$1350

For LNX, Maxx Style Fixed and Power Elevating Calf Panel

**Synchronized Option Recline with**  
**Anterior Tilt<sup>(8)</sup>** Part # SOR . . . . . \$500  
 for use on Tilt/Recline Systems

**NOTES:**

1. This module may limit front rigging options.
2. Weight capacity 250 lbs.
3. Only available on 16"-20" seat widths and 16"-22" seat depths.
4. Only available with Multi-Function electronics.
5. All front rigging will be set a 70° unless using the Synchronized Center Mount option.
6. Minimum 18" seat to floor height at zero degrees.
7. Minimum finished knee to heel to touch floor is 13".
8. This is an upgradeable option only.

**POWER SWING AWAY JOYSTICKS**

**Power Swing Away Joystick<sup>(1)(2)</sup>** . . . . . \$2,995  
 (Joystick Power "on" operation)

Select Movement:

**Horizontal** (Quad link style) Part # PSAJH  
 **Vertical Flip up** (Mid line style) Part # PSAJH

Select Mounting Position:  Left Side  Right Side  
 Part # ML Part # MR

**Operate Power Joystick/Chin Boom with**  
**Joystick Power "off" option<sup>(3)</sup>** Part # JPO . . . . . \$750

**NOTE:**

1. Comes with an egg switch.
2. Counts as an additional electronic function.
3. Allows power joystick to operate when chair is turned off.

**POWER SWING AWAY CHIN BOOM**

**Power Swing Away Chin Boom<sup>(1)(2)</sup>** Part # PSACB . . . \$2,995  
 (Joystick Power "on" operation)

Select Mounting Position:  LEFT Side  RIGHT Side

**Operate Power Joystick/Chin Boom with**  
**Joystick Power "off" option<sup>(3)</sup>** Part # JPO . . . . . \$750

**NOTE:**

1. Comes with an egg switch.
2. Counts as an additional electronic function.
3. Allows Power Chin Boom to operate when chair is turned off.

**POWER COMMUNICATION DEVICE<sup>(1)(2)(3)(4)(5)(6)</sup>**

**Power Communication Device - Single Actuator...** \$2,995  
 (Joystick Power "on" operation)

Select Mounting Position:  LEFT Side  RIGHT Side  
 Part # PCDSL Part # PCDSR

**Power Communication Device - Dual Actuator ....** \$3,995  
 (Joystick Power "on" operation)

Select Mounting Position:  LEFT Side  RIGHT Side  
 Part # PGDDL Part # PGDDR

**Operate Power Communication Device with**  
**Joystick Power "off" option<sup>(7)</sup>** Part # JPO . . . . . \$750

**Additional Electronics<sup>(3)</sup>** Part # AE . . . . . \$1,550

**NOTES:**

1. Considered an additional electronic function and requires Multi-function electronics.
2. Considered an additional two functions.
3. May require additional electronics if system has more than 6 actuators (functions).
4. Minimum seat depth available is 16".
5. Weight capacity 20 lbs.
6. Comes with an egg switch.
7. Allows Power Communication Device to operate when chair is turned off.

**TRANSFER ARMS<sup>(1)(2)(3)(4)(5)(6)</sup>**

**Power Transfer Armrest** . . . . . \$1,995  
 Select Mounting Position:  LEFT Side  RIGHT Side  
 Part # PTAL Part # PTAR

**Power Flip Back Armrest** . . . . . \$1,995  
 Select Mounting Position:  LEFT Side  RIGHT Side  
 Part # PFBAL Part # PFBAR

**Manual Transfer Armrest** . . . . . \$995  
 Select Mounting Position:  LEFT Side  RIGHT Side  
 Part # MTAL Part # MTAR

**NOTES:**

1. Weight capacity is 200 lbs.
2. Only with Adult Dual Post Style Armrest arm for Tilt/Recline, arm height 9"-13".
3. Requires Multi-Function electronics ordered on chair.
4. Minimum seat depth is 16" and 16"-20" seat width
5. Minimum seat to floor height 18".
6. Transfer only available in the chair home position.

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**POWER FOOT PLATFORM OPTIONS FOR LNX**

For use with LNX Foot platform ordered separately.

Power Telescoping Foot Platform for LNX<sup>(1)(3)(4)(5)</sup> .. \$2,250

Powers lower leg length vertically up to a maximum of 4" extension  
Part # LNXPTF

Power Flip up Foot Platform for LNX<sup>(1)(2)(3)(4)(5)</sup> ... \$1,995

Part # LNXPFUF

**NOTES:**

1. Available with one piece calf panel only and not available with a finished knee-to-heel of 13" or less. Requires 2-6 function electronics. Not available with Split Footplates.
2. Comes standard with a 12" x 10.75" footplate made of durable aluminum 1/4" thick with heel loops.
3. The front rigging will be set at approximately 80° on the Storm series.
4. Minimum finished knee to heel is 13".
5. Minimum seat to floor height is 17".

**POWER FOOT PLATFORM OPTIONS**

Fixed Centre Mount footplatform is included

Power Telescoping Foot Platform with Fixed Center Mount<sup>(1)(3)(4)(5)</sup>

Part # FCMFPTF .....\$2,845  
Powers lower leg length vertically up to a maximum of 4" extension

Power Flip up Foot Platform with Fixed Center Mount<sup>(1)(2)(4)(5)</sup>

Part # FCMFPFUF .....\$2,590

**NOTES:**

1. Includes the Fixed Center Mount Foot platform which does not elevate.
2. Comes standard with a 12" x 10.75" footplate made of durable aluminum 1/4" thick with heel loops.
3. The front rigging will be set at approximately 80° on the Storm series.
4. Minimum finished knee to heel is 13".
5. Minimum seat to floor height is 17".

**STORM FLARED CASTER OPTION**

Flared Front Casters for Storm Short Base<sup>(1)</sup> Part # FFCS \$950

**NOTE:**

1. Allows 90° front rigging on Short Storm bases comes with 10" x 10" footplate.

**LATITUDE STYLE FOOTPLATFORM**

Latitude Style Foot Platform Part # LSFP .....\$2,995

For Finished Knee to Heel Measurements of 6"-12"

Specify Finished knee to heel\* \_\_\_\_\_  
\*(including cushion immersion)

Latitude Style Padded Foot Plate Part # LSPFP .....\$250

**NOTE:**

1. Minimum seat width is 16" wide.

**POWER CALF PANEL<sup>(1)(2)(3)(4)</sup>**


Power Elevating Calf Panel Part # PECP .....\$3,995

**NOTES:**

1. No footplate available.
2. Calf Panel measures up to 15" in length and measures 13" wide at the bottom and tapers up to within 2" narrower than the seat width.
3. Minimum seat depth 16" deep.
4. Not available for use with Dual Front Rigging adaptor.

**CUSTOMIZED FOOT PLATFORM OPTIONS**

Custom Foot Platform<sup>(1)</sup> Part # CFP .....\$550



Onyx Pad for Any Footplatform<sup>(2)(3)</sup> Part # OCUFP .....\$450

**NOTES:**

1. Available for Maxx/LNX Center Mount
2. Available for all Maxx/LNX Center Mount Foot Platform options (Page 8 of QRD0335)
3. Pad thickness is 1"

**POWER SWING AWAY FOOTREST MECHANISM<sup>(1)(2)(3)(4)</sup>**

Power Assisted Legrest Swing Away Mechanism ... \$1,995

Specify :  Switch or  Joystick

Part # PALSAMS Part # PALSAMJ

**NOTES:**

1. Only available with Multi-Function electronics.
2. Compatible with Motion HD 70° Swing Away footrest, Maxx Style Swing Away Legrest Power or Manual.
3. Legrests are not removeable.
4. Minimum 16" seat width and minimum 16" seat depth.

**SPECIALTY PAD OPTIONS<sup>(1)(2)</sup>**

Optional Sized Flat One Piece Calf Panel for LNX ... \$750

Maximum pad size is seat width less 2" and 10" in length. Part # OPCC

With minimum 13" finished knee to heel.

Specify preferred size each Height: \_\_\_\_\_ Width: \_\_\_\_\_<sup>(1)</sup>

Lateral Leg Supports for LNX One Piece Calf Panel ... \$650

Select Pad Size: Part # LLSLNX

(5.5"W x 4"H)  (7"W x 6"H)  (9"W x 4"H)  (13"W x 4"H)

Lateral Leg Supports for Power Elevating Calf Panel ... \$650

Select Pad Size: Part # LLSPCP

(5.5"W x 4"H)  (7"W x 6"H)  (9"W x 4"H)  (13"W x 4"H)

Optional Sized Individual Calf Pads for Front Rigging<sup>(3)</sup> \$650

Part # OSICP

Specify preferred size each Height: \_\_\_\_\_ Width: \_\_\_\_\_

Optional Sized Flat Arm Pads .....\$650

Part # OSFAP

Specify preferred size each Length: \_\_\_\_\_ Width: \_\_\_\_\_

Optional Sized Lateral Pad or Elbow Block .....\$650

Part # OSLPEB

Specify preferred size each Length: \_\_\_\_\_ Width: \_\_\_\_\_

**NOTES:**

1. Subject to clearance available.
2. Startex fabric over polyurethane foam is standard.
3. Available on Motion Concepts front rigging options.

**Ultra Low Maxx CG Power Positioning Accessories**

"RFQ" = Request For Quotation - Please call customer service at 1-888-433-6818 and/or fax 1-888-433-6834

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**VENT TRAYS<sup>(1)(2)(3)</sup>**

- Vent Tray - Articulating (Without Battery Tray) . . . . . \$1,350
  - Horizontal (14W X 15D X 10H) Part # ARTVTH
  - Vertical (14.625W X 10.125D X 13H) Part # ARTVTV
- Vent Tray - Trilogy 100/200 (14W x 6D) . . . . . \$1,350  
(14W X 15D X 10H) - (Without Battery Tray) Part # TRILVT

**NOTE:**  
 1. The user weight limit is 225 lbs.  
 2. The maximum weight capacity of the vent trays is 35 lbs.  
 3. Vents are not standard configuration on Elevating or Precline Systems or with Single Posts Backs.

**COMPACT-LAPTOP VENT TRAYS<sup>(1)</sup>**

- Vent Tray (Without Battery Tray) Part # CVT . . . . . \$1,350
- Vent Bag Holder<sup>(2)</sup> - Articulating -without battery tray Part # VBHA \$990
- Vent Bag Holder<sup>(2)</sup> - Non Articulating  
(Without Battery Tray) Part # VBHNA . . . . . \$525

**NOTE:**  
 1. Vents are not standard configuration on Elevating or Precline Systems or with Single Posts Backs.  
 2. Motion Concepts does not supply the Vent Bag.

**ADDITIONAL BATTERY TRAYS**

- Storm Long Base Additional Front Mounted Battery  
Tray Group 40 Battery Part # SLBT40 . . . . . \$500
- Storm Short Base Additional Front Mounted Battery  
Tray for U1 Battery Part # SLBTU1 . . . . . \$500
- Storm Base - 3 x Group 22 Battery Tray Conversion<sup>(3)</sup>  
Part # S322 . . . . . \$950
- TDX-SP Base - 3 x Group 22 Battery Tray<sup>(1)</sup>  
Part # TDX322 . . . . . \$500
- TDX-SP Base Rear Caster Mounted U1 Battery Tray<sup>(2)</sup>  
Part # TDXU1 . . . . . \$500
- Rovi Base - 3 M40 Battery Tray Part # R340 . . . . . \$450
- Edge HD - 3 x Group 22 Battery Tray Conversion<sup>(1)</sup>  
Part # EDGE322 . . . . . \$500

**NOTE:**  
 1. Comes with steel battery door cover.  
 2. Must have factory Transport brackets TBKTS  
 3. Deletes factory battery boxes.

**ADDITIONAL ELECTRONIC OPTIONS<sup>(1)</sup>**

- Step Down Converter<sup>(1)</sup> Part # SDC . . . . . \$500  
-converts 24V to 12V
- Richardson Powertech Vent Power Center<sup>(1)</sup>  
Part # RPVPC . . . . . \$1,950

**NOTE:**  
 1. Should this product be used in life support applications Motion Concepts will not assume any risk and shall be held harmless against all damages.

**O2 HOLDERS<sup>(1)</sup>**

- Non-Articulating O<sup>2</sup> Gas Holder Part # NA02H . . . . . \$395  
(Cylinder Style)
- Articulating O<sup>2</sup> Holder . . . . . \$1,250  
Please Select One:  O<sup>2</sup> Gas (Cylinder Style) Part # A0HG  
 Liquid O<sup>2</sup> (6.5" x 4.5") Part # A0HL
- Non-Articulating IV Pole . . . . . \$395  
 Users Left Part # NA1PL  Users Right Part # NA1PR
- Articulating IV Pole Part # AIVP . . . . . \$1,250

**NOTE:**  
 1. Articulating O<sup>2</sup> holders are not standard configurations on Precline Systems, please call for quotation.

**MEMORY SEATING OPTIONS<sup>(1)(2)(3)</sup>**

- Single Function Memory Seating Part # SFMS . . . . . \$495  
\*Please call for specification sheet 1-888-433-6818
- Multiple Function Memory Seating Part # MFMS . . . . . \$995  
\*Please call for specification sheet 1-888-433-6818

**NOTES:**  
 1. Does not operate elevation function, anterior tilt function or lateral tilt function.  
 2. Comes with a quad button switch with LED indicator.  
 3. Only available with Multi-Function electronics.

**SAFETY LIGHTS**

- Front Safety Lights (pair)
  - Operated through switches<sup>(1)</sup> Part # FSLS . . . . . \$495
  - Operated through Joystick<sup>(2)</sup> Part # FSLJ . . . . . \$845
- Rear Safety Lights (pair)
  - Operated through switches<sup>(1)</sup> Part # RSLS . . . . . \$495
  - Operated through Joystick<sup>(2)</sup> Part # RSLJ . . . . . \$845
- Front & Rear Light Combo Package Credit\* Part # LPC (\$150)  
\*Purchase both sets and receive a credit
- Ground Effects Lighting Through Switch Part #GEL . . . \$350

**NOTE:**  
 1. Operated through a separate power supply and single push button switch.  
 2. Only available with Multi-Function electronics.

**OTHER ACCESSORIES<sup>(1)</sup>**

- Power Back Pack Holder with Single Back Pack . . \$2,995  
Select Mounting Position:  Left Side  Right Side  
Part # PBPHL Part # PBPHR
- Concepts Back Pack -  Single Part # SBP . . . . . \$125  
 Double Part # DBP . . . . . \$195

**NOTE:**  
 1. Weight limit 20 lbs.

**Total Retail Value \$** \_\_\_\_\_