



A3

Contour Static Seating on ROVI A3 Base US PRICE LIST AND ORDER FORM

Price Effective November 8, 2022



PRINT FORM

SAVE

CLEAR FORM

Contact Motion Concepts Customer Service: 1.888.433.6818 | quotes@motionconcepts.com

REQUIRED INFORMATION

☐ Quote ☐ Order

Date of Order: _____ Dealer Account #: _____

Dealer Name: _____

PO #: _____ Tag: _____

Purchasing Contact: _____

Phone: _____ Fax: _____

E-mail: _____

ATP/Therapist: _____

Ship to Address: _____

City: _____ State: _____


Zip: _____ Client Gender: ☐ M ☐ F

Special Client Conditions: _____

REQUIRED MEASUREMENTS

To ensure system is accurately configured
please fill in all required measurements

Client Height:	
Client Weight (lbs):	
A. Seat to Shoulder:	
B. Trunk Depth:	
C. Chest Width:	
D. Knee to Back:	
E. Seat to Top of Head:	
F. Elbow to Hand:	
G. Seat to Elbow:	
H. Hip Width:	
I. Knee to Heel:	
Cushion Thickness:	



ROVI

A3

HCPCS CODES	K0848, K0849
SUSPENSION SYSTEM	Active Ride Control (ARC)
BASE WIDTH	25.5"
BASE LENGTH	36.5"
DRIVE WHEEL SIZE ⁽¹⁾	14"
CASTER WHEEL SIZE (front/rear)	6"
BASE ONLY TURNING RADIUS	20.5"
BASE VITALS	DRIVETRAIN 4-pole EAD Motors

⁽¹⁾ All A3 wheels come with solid inserts

⁽²⁾ Actual driving range and speed may vary due to factors such as user weight, type and grade of terrain, battery condition & charge level; type and condition of the drive wheels, and variations in the drive and control systems.

BATTERIES	2 x M34 3 x M50 option for Vents
BATTERY CHARGER	8A, Off Board
BATTERIES WEIGHT (each)	38.5 lbs
RANGE ⁽²⁾	17.7 miles
WEIGHT CAPACITY	300 lbs
BASE WEIGHT (w/o batteries)	157 lbs
GROUND CLEARANCE	3"
MAXIMUM SPEED	6.2 mph

SOLID SEAT PAN - K0848



CAPTAIN SEAT - K0849



FULL CAPTAIN-STYLE SEAT - K0849

All bases below come standard with Full Captain-Style Seat in Grey or Black Vinyl and 22" high backrest. Includes a Semi-Recline hinge of 120° with left side only adjustable lever. Standard adjustable height, flip-back armrests and a foot platform that measures 12" wide by 9" deep.

☐ 888WNLEA3-MP6-CS ROVI A3 5.8 MPH Quick Speed Motors⁽¹⁾
HCPCS code K0849 \$7,095

☐ HSM62 Optional 6.2 MPH HIGH SPEED MOTORS \$950
with OnTraxx Enhanced Driving Module

☐ OEDMC OnTraxx Enhanced Driving Module HCPCS code K0108 \$1,550

SELECT Seat Size:

☐ SV18171/222RG Grey 18" x 18"

☐ SV18171/222RB Black 18" x 18"

SELECT Seat to floor for 18" x 18" seat size⁽²⁾:

☐ STF19 19" Seat-To-Floor N/C

☐ STF20 20" Seat-To-Floor N/C

☐ STF21 21" Seat-To-Floor N/C

☐ STF22 22" Seat-To-Floor N/C

SELECT Seat Size:

☐ SV2121211/2RG Grey 20" x 20"

☐ SV2121211/2RB Black 20" x 20"

SELECT Seat to floor for 20" x 20" seat size⁽²⁾:

☐ STF20 20" Seat-To-Floor N/C

☐ STF21 21" Seat-To-Floor N/C

☐ STF22 22" Seat-To-Floor N/C

☐ STF23 23" Seat-To-Floor N/C

SOLID SEAT PAN & CAPTAIN-STYLE SEAT - K0848

All bases below come standard with a Solid Seat Pan style Captain Seat in Grey or Black Vinyl and 22" high backrest. Includes a Semi-Recline hinge of 120° with left side only adjustable lever. Standard adjustable height, flip-back armrests and a foot platform that measures 12" wide by 9" deep.

☐ 888WNLEA3-MP6-SS ROVI A3 5.8 MPH Quick Speed Motors⁽¹⁾
HCPCS code K0849 \$6,950

☐ HSM62 Optional 6.2 MPH HIGH SPEED MOTORS \$950
with OnTraxx Enhanced Driving Module

☐ OEDMC OnTraxx Enhanced Driving Module HCPCS code K0108 \$1,550

SELECT Seat Size:

☐ SP181821RG Grey 18" x 18"

☐ SP202023RG Grey 20" x 20"

☐ SP181821RB Black 18" x 18"

☐ SP202023RB Black 20" x 20"

SELECT Seat to floor⁽³⁾:

☐ STF17 17" Seat-To-Floor N/C

☐ STF18 18" Seat-To-Floor N/C

☐ STF19 19" Seat-To-Floor N/C

☐ STF20 20" Seat-To-Floor N/C

1. Price of the powerbase includes: Active Ride Control suspension, tie-down hard ware for 4-point strap system, 8 amp off-board charger, flat free inserts and positioning belt. Weight capacity is 300 lbs.

2. Seat to floor height is measured at the seam of cushion on 20" widths and at the top of the leg trough on 18" widths. Seat to floor heights may vary +/- .25" dependent on base conditions.

3. Seat to floor height is measured at the front of the seat pan based on a stand seat size. Seat to floor heights may vary +/- .25" dependent on conditions of the power base and on seat depths greater than 18"

FENDER COLOR OPTIONS

- | | |
|--|--|
| <input type="checkbox"/> 800914-00200-03 | <input type="checkbox"/> Blizzard White |
| <input type="checkbox"/> 800914-00200-69 | <input type="checkbox"/> Solar Yellow |
| <input type="checkbox"/> 800914-00200-23 | <input type="checkbox"/> Original Orange |
| <input type="checkbox"/> 800914-00200-21 | <input type="checkbox"/> Red Delicious |
| <input type="checkbox"/> 800914-00200-2C | <input type="checkbox"/> Party Pink |



NEW Style Fender Colors⁽¹⁾

- | | |
|------------------------------|--|
| <input type="checkbox"/> BCF | <input type="checkbox"/> Black Carbon (Matte) |
| <input type="checkbox"/> BPF | <input type="checkbox"/> Black Panther (Gloss) |
| <input type="checkbox"/> PRF | <input type="checkbox"/> Purple Reign |
| <input type="checkbox"/> IBF | <input type="checkbox"/> Indigo Blue |
| <input type="checkbox"/> MGF | <input type="checkbox"/> Motion Green |



1. The new style fender design

FENDER COLOR WITH LIGHTS⁽¹⁾⁽²⁾

SPECIAL OFFER

- | | | |
|--|--|-----------|
| <input type="checkbox"/> 800914-00201-03 | <input type="checkbox"/> Blizzard White | Now \$700 |
| <input type="checkbox"/> BCFL | <input type="checkbox"/> Black Carbon (Matte) | Now \$700 |
| <input type="checkbox"/> BPFL | <input type="checkbox"/> Black Panther (Gloss) | Now \$700 |
| <input type="checkbox"/> 800914-00201-69 | <input type="checkbox"/> Solar Yellow | Now \$700 |
| <input type="checkbox"/> 800914-00201-36 | <input type="checkbox"/> Motion Green | Now \$700 |
| <input type="checkbox"/> 800914-00201-23 | <input type="checkbox"/> Original Orange | Now \$700 |
| <input type="checkbox"/> 800914-00201-40 | <input type="checkbox"/> Bright Blue | Now \$700 |
| <input type="checkbox"/> 800914-00201-21 | <input type="checkbox"/> Red Delicious | Now \$700 |
| <input type="checkbox"/> 800914-00201-2C | <input type="checkbox"/> Party Pink | Now \$700 |
| <input type="checkbox"/> 800914-00201-80 | <input type="checkbox"/> Purple Reign | Now \$700 |

1. Fully integrated LED lighting kit. Includes headlights, taillights and turn signals. Must select Joystick for Lights feature or Omni for operation (page 3).

2. These fenders are not available on VR2 Bases.



MATCHING BATTERY DOORS

- ☐ MBD Matching Battery Doors \$495
(Select fender color first)

CUSTOM DESIGNS

Please see form number TRD0607 - ROVI Custom Designs Price Sheet, for Custom/Dipped Fender colors, battery door and other personalized options.

TIRE OPTIONS

- | | | |
|---------------------------------|----------------------|-------|
| <input type="checkbox"/> STIRES | Standard Solid Tires | STD |
| <input type="checkbox"/> PTIRES | Pneumatic Tires | \$199 |

MILITARY DECAL

- | | | |
|---------------------------------|---|------|
| <input type="checkbox"/> DECALA | Army Decal 3" Diameter (Uninstalled) | \$50 |
| <input type="checkbox"/> DECALF | Air Force Decal 3" Diameter (Uninstalled) | \$50 |
| <input type="checkbox"/> DECALN | Navy Decal 3" Diameter (Uninstalled) | \$50 |
| <input type="checkbox"/> DECALM | Marines Decal 3" Diameter (Uninstalled) | \$50 |



Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

SINGLE FUNCTION CONTROLLER OPTIONS

- ☐ D51122SC R-Net LED Joystick⁽¹⁾⁽²⁾ STD
- ☐ D51108 R-Net LED Joystick for Lights⁽¹⁾ N/C
- ☐ D51165 R-Net Color Remote Joystick with Monojacks⁽¹⁾⁽³⁾ \$600
- ☐ D50982 R-Net Color Remote Joystick with Monojacks for Lights⁽¹⁾ \$650
- ☐ D51711 R-Net Advanced Joystick⁽¹⁾⁽³⁾ \$1,295
- ☐ D51109 R-Net Expandable Controller HCPCS Code: E2377 \$615
- ☐ REH R-Net Harness for Expandable Controller HCPCS Code: E2313 \$415

1. Requires an R-Net Controller and Harness of Page 2

2. Light package not available with this joystick

3. 3.5" high-resolution VGA colour display, IR Control as standard, includes integrated Bluetooth Module.

MATRIX SEAT CUSHION OPTIONS

- Matrix PS Seat Cushion wide/reversable cover**
- HCPCS code E2605 \$345
- ☐ PS1818 ☐ PS2020

- Matrix PSP Seat Cushion wide/reversable cover**
- HCPCS code E2607 \$395
- ☐ PSP1818 ☐ PSP2020







- Matrix Multi Cushion wide/reversable cover**
- HCPCS code E2607 \$445
- ☐ MM1818 ☐ MM2020

- Matrix Vi Seat Cushion** HCPCS code E2607 \$425
- ☐ MA1818-VI ☐ MA2020-VI

- Matrix Libra Seat Cushion** HCPCS code E2624 \$485
- ☐ LC1818 ☐ LC2020

- Matrix Libra Adjust Seat Cushion** HCPCS code E2624 \$625
- ☐ LC1818-ADJ ☐ LC2020-ADJ

SEAT CUSHION REFERENCE GUIDE

						
						
Cushion	Meshtex	Startex	Infection	Spacetex	Onyx	Fusion
Libra/Adjust	N/A	STD	\$50	N/C	N/A	N/A
PS/PSP	N/A	STD-Reversible	STD-Reversible	N/A	N/A	N/A
Vi	STD	\$50	\$50	N/A	\$100	\$150

FABRIC OPTIONS FOR CUSHIONS

- ☐ SFLF Spacetex Fabric Option for Libra Cushion. . . N/C
- ☐ ICFL Infection Control Fabric Upgrade for Libra Cushion \$50
- ☐ FFV Fusion Fabric Upgrade for Matrix Vi Cushion.. \$150
- ☐ OCUMA Onyx Fabric Upgrade for Matrix Vi Cushion... \$100
- ☐ SFVS Smooth side Startex Fabric Upgrade for Matrix Vi Cushion \$50
- ☐ SFVF Fabric side Startex Fabric Upgrade for Matrix Vi Cushion \$50

JOYSTICK MOUNT

- ☐ 1343-00707-00L Fixed for Solid Pan/Captain Seat (Left Hand) .. N/C
- ☐ 1343-00707-00R Fixed for Solid Pan/Captain Seat (Right Side) .. N/C
- ☐ SAQR Motion Height Adjustable Swing-Away Quad Link (Right Hand)⁽¹⁾ HCPCS code E1028 .. \$235
- ☐ 1343-00709-00-A3 Swing-Away Mount Bracket (Left Hand) .. \$235

JOYSTICK KNOB - OPTIONAL

- ☐ SA80941 "T" Handle for Joystick Knob HCPCS code E2323 \$115
- ☐ SA72441 Ball Handle for Joystick Knob HCPCS code E2323 \$115
- ☐ PC101B-20 BodyPoint Joystick Handle 3" U-Shaped HCPCS code E2323 .. \$120
- ☐ PC102B-20 BodyPoint Joystick Handle 4" U-Shaped HCPCS code E2323... \$120
- ☐ PC107B-20 BodyPoint Joystick Handle Domed Rubber HCPCS code E2323 \$120

ELECTRONIC OPTIONS

- ☐ D50883 R-Net Output Module (EMC)⁽¹⁾⁽²⁾ \$995
- ☐ D51580 R-Net Bluetooth Module IOS Apple⁽¹⁾⁽²⁾ \$895
- ☐ D51111 R-Net Bluetooth Module Android/PC⁽¹⁾⁽²⁾ \$895
- ☐ SPH Smart Phone Holder \$50
- ☐ MUSB Motion Dual USB Charger N/C

1. Requires an R-Net Controller

2. Must order R-Net Colour Remote joystick with these options

ARMRESTS (MUST BE SELECTED)

- ☐ 104300-00112 Adjustable Height Flip Back Arm Left \$125
- ☐ 104400-00112 Adjustable Height Flip Back Arm Right \$125

WHEELCHAIR TRANSPORT BRACKETS

- ☐ XBKTS Wheelchair Transport Brackets⁽¹⁾ \$250

1. For unoccupied use only

BATTERIES

- ☐ M34 SLD G MK Group 34 Battery (ea) Quantity⁽¹⁾ (2)
- HCPCS code E2359 \$305

1. 2 batteries required.

Total Retail Value \$ _____

SUBMIT**PRINT FORM****SAVE****CLEAR FORM**