



PRINT FORM

SAVE

CLEAR FORM

Contact Motion Concepts Customer Service: 1.888.433.6818 | quotes@motionconcepts.com



Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

### CELL PHONE HOLDER - ARMREST MOUNT <sup>(1)(2)</sup>

- ☐ RCPAL Left .....\$85  
☐ RCPAR Right .....\$85



#### RAM® X-Grip® Phone Mount

The specialized mounting brackets allow the Phone to be seamlessly integrated via BALL MOUNT to the Arm Tube of your Motion Concepts Power Positioning System. The mounts can be reconfigured to accommodate a variety of mounting position.

1. Not available with TRRO, TRROE, XRO, XROE  
2. Not available on the joystick side

### CELL PHONE HOLDER - RESOLVE SWING AWAY MOUNT <sup>(1)(2)</sup>

- ☐ RCPRL Left .....\$350  
☐ RCPRR Right .....\$350



#### RAM® X-Grip® Phone Mount

To fully personalize positioning of your Phone, the Resolve Swing-Away Mount offers up to 4" (2" up or 2" down) of vertical travel relative to the arm pad as well as rotation and angle adjustability and an opportunity to position the accessory in any orientation that best accommodates your personal needs.

1. Not available with TRRO, TRROE, XRO, XROE  
2. Not available on the joystick side

### TABLET HOLDER - SEAT RAIL MOUNT <sup>(1)</sup>

- ☐ RTL Left .....\$270  
☐ RTR Right .....\$270



#### RAM® X-Grip® Wheelchair Seat Track Mount for 9"-10" Tablets

The specialized mounting brackets allow the Tablet to be seamlessly integrated via BALL MOUNT to the Ultra Rail of your Motion Concepts Power Positioning System. The mounts can be reconfigured to accommodate a variety of mounting positions.

1. Not available with TRRO, TRROE, XRO, XROE

### CUP HOLDER - ARMREST MOUNT <sup>(1)(2)</sup>

- ☐ RCHAL Left .....\$85  
☐ RCHAR Right .....\$85



#### RAM® Level Cup™ XL 32oz Drink Holder

The specialized mounting brackets allow the Cup Holder to be seamlessly integrated via BALL MOUNT to the Arm Tube of your Motion Concepts Power Positioning System. The mounts can be reconfigured to accommodate a variety of mounting position. Designed for bigger insulated cups, plastic water bottles or even your reliable steel coffee cup, the Level Cup XL keeps your drink upright, regardless of the conditions

1. Not available with TRRO, TRROE, XRO, XROE  
2. Not available on the joystick side

### CUP HOLDER - RESOLVE SWING AWAY MOUNT <sup>(1)(2)</sup>

- ☐ RCHRL Left .....\$350  
☐ RCHRR Right .....\$350

#### RAM® Level Cup™ XL 32oz Drink Holder

To fully personalize positioning of your Cup Holder, the Resolve Swing-Away Mount offers up to 4" (2" up or 2" down) of vertical travel relative to the arm pad as well as rotation and angle adjustability and an opportunity to position the accessory in any orientation that best accommodates your personal needs. Designed for bigger insulated cups, plastic water bottles or your steel coffee cup, the Level Cup XL keeps your drink upright, regardless of the conditions.

1. Not available with TRRO, TRROE, XRO, XROE  
2. Not available on the joystick side



Total Retail Value \$ \_\_\_\_\_

**SUBMIT**

PRINT FORM

SAVE

CLEAR FORM