

INVACARE® MATRX®

# PSP (Posture Seat Polymer) CUSHION

---



The **Invacare Matrix PSP** (Posture Seat Polymer) cushion combines the stability and support of the PS foam base with a polymer overlay for extra protection and shear reduction. An ultra-lightweight cushion providing zero maintenance, worry free skin protection and positioning.

## Call to order

**1.800.333.6900**  
invacare.com

**1.888.433.6818**  
motionconcepts.com



**matrx®**  
SEATING SERIES




### Features and Functional Benefits

- Startex™ outer cover – Reversible can be used smooth or fabric side out
- Moisture resistant and breathable to keep cool and dry
- Includes hook and loop strips
- Startex™ inner cover for extra moisture protection
- Molded high resilient polyurethane foam made with Ultra-fresh™ anti-microbial
- Rear cushion radius and waterfall front edge
- Gently contoured to provide neutral leg position
- Molded polymer insert for extra protection and shear reduction
- Unique waffled relief zone allows immersion in ischial area
- Designed to provide off-loading to hips and thighs for support and stability

**HCPCS Code: Pediatric/Standard - E2607**  
**Heavy-Duty - E2607 (up to 21" W)/E2608( 22" W and greater)**

### Optional Accessories



**POK – Pelvic Obliquity Kit**  
 (2 pieces) accommodates up to one inch of obliquity



**Cushion Rigidizer –** lightweight and low profile

**Invacare Corporation**  
[www.invacare.com](http://www.invacare.com)

**USA**  
 One Invacare Way  
 Elyria, Ohio  
 44035-4190  
 (800) 333-6900

**Canada**  
 570 Matheson Blvd. E.  
 Unit 8  
 Mississauga, Ontario  
 L4Z 4G4 Canada  
 (800) 668-5324

**Motion Concepts**  
[www.motionconcepts.com](http://www.motionconcepts.com)  
[info@motionconcepts.com](mailto:info@motionconcepts.com)

**USA**  
 700 Ensminger Rd., Suite 112  
 Tonawanda, NY 14150  
 (888) 433-6818  
 Fax: (888) 433-6834

**Canada**  
 84 Citation Dr., Unit 1-7  
 Concord, ON L4K 3C1  
 (866) 748-7943  
 (905)-695-0134  
 Fax: 905-695-0138

© 2021  
 Invacare Corporation.  
 All rights reserved.  
 Trademarks are identified by the symbols™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Matrix is a registered trademark of Motion Concepts. Otto Bock is a trademark of Otto Bock HealthCare. Specifications are subject to change without notification.  
 Form No. 18-296 rev. 1/12  
 160959



| Model #          | Width              | Depth              | Height    | Height at Leg Abductor/Adductor | Product Trough | Weight  | Limited Capacity |
|------------------|--------------------|--------------------|-----------|---------------------------------|----------------|---------|------------------|
| Warranty         |                    |                    |           |                                 |                |         |                  |
| PSP (Pediatric)  | 10" to 14"         | 10", 12", 14", 16" | 3.5"/3.5" | 2.75"                           | 1.2 lb. and up | 150 lb. | 2 Years          |
| PSP (Standard)   | 16", 17", 18", 20" | 16", 18", 20"      | 3.5"/3.5" | 2.75"                           | 2.2 lb. and up | 300 lb. | 2 Years          |
| PSP (Heavy-Duty) | 22", 24", 26"      | 18", 20", 22"      | 3.5"/3.5" | 3.00"                           | 4.6 lb. and up | 600 lb. | 2 Years          |

