

8.0 - Accessories:

V. Headrests

Standard Headrest

The standard headrest is made from sealed cell urethane foam. The headrest assembly allows for easy adjustments for both height and depth, and it can be locked at virtually any desired angle due to its pivoting design.



Onyx Headrest

The Onyx headrest is an automotive style headrest with a leatherette fabric cover. **(Note:** a full Onyx seating system is available in the leatherette finish. The Onyx line consists of a seat cushion, a back cushion and headrest). The headrest is height adjustable and can be adjusted to different angles via the pivot block design. The Onyx headrest is designed to fit most standard multi offset headrest hardware.



Stealth Headrest

The Stealth headrest utilizes Stealth Products "World's Best" technology to produce a fully adjustable headrest bracket. The design of the bracket allows the headrest to be extended, rotated and angled to achieve any desired position. These durable joints are also designed to remain locked following each adjustment.

NOTES