

AVAILABLE
NOW!

UPfront
maxx

be **UPFRONT** with
INDEPENDENCE
& **ACCESSIBILITY**
EXCLUSIVELY ON

ROVI



The all-NEW **UPfront Maxx Module** from Motion Concepts was exclusively designed for **X3 and A3 Power Bases** from ROVI. It allows people to get a bit closer to their everyday! By thoughtfully combining a fully functioning power positioning system with 5.5" of forward shift, consumers are able to explore new ways of accessing previously difficult to reach surfaces and objects. The **UPfront module** also provides confidence and stability by automatically tucking the user's feet rearward, under the front of the seat pan and locking out the front casters for peace of mind and greater flexibility. All this access is available with the durability, function and modularity you have come to expect from Motion Concepts products.

Call to order

888.433.6818
motionconcepts.com



ROVI A3 | X3



UPfront maxX



BASE VITALS

	ROVI A3	ROVI X3
HCPCS Codes	K0856, K0861	K0856, K0861
Suspension System	Active Ride Control (ARC)	Active Ride Control (ARC)
Base Width	25.5"	23.25"
Base Length	36.5"	36.5"
Drive Wheel Size ⁽¹⁾	14"	14"
Caster Wheel Size (front/rear)	6"	6"
Base only Turning Radius	20.5"	19.5"
Drivetrain	4-pole EAD Motors	4-pole EAD Motors
Batteries	2 x M34 3 x M34*	2x M34 3 x M34*
	(*option for Vent application)	(*option for Vent application)
Battery Charger	8A, Off Board	8A, Off Board
Batteries Weight (each)	38.5 lbs (M34)	38.5 lbs (M34)
Range ⁽²⁾	17.7 miles	17.7 miles
Weight Capacity	300 lbs	300 lbs
Base Weight (w/o batteries)	159 lbs	157 lbs
Ground Clearance	3"	3"
Maximum Speed	6.2 mph	6.2 mph

⁽¹⁾ All A3 wheels come with solid inserts

⁽²⁾ Actual driving range and speed may vary due to factors such as user weight, type and grade of terrain, battery condition & charge level; type and condition of the drive wheels, and variations in the drive and control systems.

SYSTEM VITALS

- Fully modular, featuring up to 45° of CG Tilt, 170° of Recline and up to 45° of Anterior Tilt
- 16" - 22" W x 16" - 23"D, fully adjustable
- Reduced drive speed in UPfront position
- Recline available from ALL positions
- When system is not in UpFront position the ROVI A/X3's ARC suspension is fully functioning for maximum performance and climbing capability
- Anti-dive locking mechanism offers additional stability and safety by locking front suspension arms as the UpFront function is engaged
- Available with the wide range of Motion Concepts' accessories and Invacare® Matrx® Seating products, mounted via our Ultra Rail and Ultra Track slotted interfaces
- Seat-to-Floor heights as low as 17.5"
- Up to 250 lb weight capacity



Motion Concepts USA
700 Ensminger Rd., Suite 112
Tonawanda, New York 14150
888.433.6818

Motion Concepts CANADA
84 Citation Dr., Unit 1
Concord, ON L4K 1C3
866.748.7943

www.motionconcepts.com
info@motionconcepts.com

© 2019 Motion Concepts. All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Motion Concepts unless otherwise noted. Specifications are subject to change without notification.

CONTACT YOUR SALES PROFESSIONAL

