

# X3

## Oh Captain, My Captain . . .

The innovative and exciting ROVI X3 boasts the narrowest wheelbase in the industry and all of the function and durability you have come to expect from Motion Concepts. Now the ROVI X3 is available with a stylish and functional Captain's Seat.

- Offered in captain's seat and solid seat pan configurations
- Available in two-tone gray vinyl or black
- Semi-recline angle of 120°
- Low-profile design
- Flip back, adjustable height, adjustable width arms
- High back with headrest
- Built in seat pouch on backrest
- 300 lb. weight capacity
- Seat to floor heights as low as 18"
- 18" x 18" and 20" x 20" seat sizes available
- USB charger included



Follow us on . . .



Toll Free Tel: 1.888.433.6818 ► Toll Free Fax: 1.888.433.6834  
► [www.motionconcepts.com](http://www.motionconcepts.com)



CAPTAIN-STYLE SEAT WITH SOLID SEAT PAN



FULL CAPTAIN-STYLE SEAT-Grey



FULL CAPTAIN-STYLE SEAT - Black

## CUSTOM DIPPED SHROUD DESIGNS



\* Other options are also available

## CUSTOM SHROUD DESIGNS BY Veteran Bob Rice



\* Personalized options are also available

## X3 BASE VITALS

<b>HCPCS CODES</b>	K0848, K0849
<b>SUSPENSION SYSTEM</b>	Active Ride Control (ARC)
<b>BASE WIDTH</b>	23.25"
<b>BASE LENGTH</b>	36.5"
<b>DRIVE WHEEL SIZE <sup>(1)</sup></b>	14"
<b>CASTER WHEEL SIZE (front/rear)</b>	6"
<b>TURNING RADIUS</b>	20"
<b>DRIVETRAIN</b>	4-pole EAD Motors
<b>BATTERIES</b>	2 x M34 3 x M50 option for Vents
<b>BATTERY CHARGER</b>	8A, Off Board
<b>BATTERIES WEIGHT (each)</b>	42 lbs
<b>RANGE <sup>(2)</sup></b>	17.7 miles
<b>WEIGHT CAPACITY</b>	300 lbs
<b>BASE WEIGHT (w/o batteries)</b>	157 lbs
<b>GROUND CLEARANCE</b>	3"
<b>MAXIMUM SPEED</b>	5 mph

<sup>(1)</sup> All X3 wheels come with solid inserts

<sup>(2)</sup> Actual driving range and speed may vary due to factors such as user weight, type and grade of terrain, battery condition & charge level; type and condition of the drive wheels, and variations in the drive and control systems.

## X3 BASE SHROUD COLORS\*



\* Wheel rims are always black