

**10" SEAT  
ELEVATION  
NOW  
AVAILABLE!**

# RISE TO THE OCCASION

**ULTRA LOW MAXX BARIATRIC  
POWER POSITIONING SYSTEM  
ON INVACARE® TDX® SP2 HD SERIES**

- From straightforward to the most complex, Ultra Low Maxx Bariatric Power Positioning Systems (BPPS) are offered in numerous configurations to best accommodate the varied individual needs of today's larger clients
- Up to 50° of CG Tilt for efficient pressure relief
- Up to 150° of Recline with Extended Shear Reduction (ESR) to help reduce damaging shear forces and maintain posture
- **NEW** 10" Power Adjustable Seat Height Module with up to 350 lb. weight capacity
- Easy, 'in-clinic' dimensional adjustments to help ensure optimal system set-up and performance
- Seat-to-Floor heights as low as 17"
- Fully adjustable: 20" - 26" wide x 16" - 23" deep
- Optional LNX Power Center Mount Foot Platform or Maxx Style Fixed Center Mount Legrest for smooth and effective leg elevation and positioning
- Available with the full range of Motion Concepts' accessories and Invacare® Matrx® Bariatric Seating products, mounted via our Ultra Rail and Ultra Track slotted interfaces
- Comes with SureStep® Suspension and Stability Lock that work in synergy to help maintain balance and provide stability when transitioning over obstacles
- Up to 400 lb. weight capacity



Modern Design Meets Timeless Performance in the Next Generation of TDX (Total Driving eXperience)

The Ultra Low Maxx Bariatric Power Positioning System for Invacare® TDX® SP2 HD Series and new LiNX® electronic platform is designed for people seeking exceptional driving performance, elegant design, comfort and remarkable stability.

**Call to order**

**888.433.6818**  
motionconcepts.com



**REHAB SEATING SYSTEM**

**TDXSP2HD**

K0850



**SINGLE ACTUATOR SYSTEM**

**TDXSP2HD-CG**

K0858



**MULTIPLE ACTUATOR SYSTEM**

**TDXSP2HD-MCG**

K0862



**TDX® SP2 HD Series Base Vitals**

<b>Suspension System</b>	SureStep® Suspension and Stability Lock
<b>Base Width</b>	25.5"
<b>Base Length</b>	35.6"
<b>Drive Wheel Size</b>	14"
<b>Caster Wheel Size (front/rear)</b>	6"
<b>Base only Turning Radius</b>	20"
<b>Motor</b>	4-pole Sealed Housing
<b>Batteries</b>	Group 24 x 2
<b>Batteries Weight (each)</b>	51.5 lbs (GP24)
<b>Battery Range <sup>(1)</sup></b>	Up to 20 miles
<b>Weight Capacity (base only)</b>	Up to 400 lbs
<b>Ground Clearance</b>	> 2.5"
<b>Maximum Speed</b>	5.8 mph Standard Motor
<b>Warranty</b>	1 Year Electronics/Motors 1 Year Base & Seat Frame 6 Months Batteries

<sup>(2)</sup> Actual driving range and speed may vary due to factors such as user weight, type and grade of terrain, battery condition & charge level; type and condition of the drive wheels, and variations in the drive and control systems.

**Ultra Low Maxx Bariatric Power Positioning System Vitals**

<b>Seat Width</b>	20" - 26" (Fully adjustable)
<b>Seat Depth</b>	16" - 23" (Fully adjustable)
<b>Seat-to-Floor</b>	as low as 17"
<b>Weight Capacity</b>	Up to 400 lbs

**REHAB SEATING SYSTEM: TDXSP2HD K0850**

- Bariatric Rehab Seat •

**SINGLE ACTUATOR SYSTEM: TDXSP2HD-CG K0858**

- Bariatric 50° Tilt Only Up to 400 lbs

**MULTIPLE ACTUATOR SYSTEM: TDXSP2HD-MCG K0862**

- Bariatric 50° Tilt & 150° Power Recline with ESR\* Up to 400 lbs
- Bariatric 150° Power Recline with ESR \* Up to 400 lbs
- 10" Power Adjustable Seat Height Module Up to 350 lbs

**ADD-ON MODULES**

- 8° Precline •
- 30° Precline •

**MATRIX SEATING OPTIONS**

- Bariatric Cushions •
- Bariatric Backs •

<sup>(\*)</sup> Extended Shear Reduction

**CUSTOMIZE FOR YOUR PERSONAL STYLE\***

- All black frame with complementary black tires and wheels for a modern, elegant design
- Customized, contemporary colored shrouds with reflectors and wheel accent colors
- New rear swing arm design for better controller access
- Standard transport brackets for secure, unoccupied wheelchair transportation



SpaceStation Silver



Monster Green



RockStar Red



Grape JellyBean



Lights Out Black - Matte



Island Blue



Tangy Orange



Invacare® Blue



WhiteOut



Black Ice - Glossy

**Motion Concepts USA**  
700 Ensinger Rd., Suite 112  
Tonawanda, New York 14150  
888.433.6818

**Motion Concepts CANADA**  
84 Citation Dr., Unit 1  
Concord, ON L4K 1C3  
866.748.7943

[www.motionconcepts.com](http://www.motionconcepts.com)  
[info@motionconcepts.com](mailto:info@motionconcepts.com)



TRD0433 Rev B

© 2019 Motion Concepts. All rights reserved. Trademarks are identified by the symbols ™ and ©. All trademarks are owned by or licensed to Motion Concepts unless otherwise noted. Specifications are subject to change without notification.