



SYSTEM WEIGHT CALCULATOR



Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

BASE with Mounting Brackets (no batteries) WEIGHT

- ROVI X3 167 lbs
- ROVI A3 169 lbs
- Invacare TDX SP2 NARROW Group 22 158 lbs
- Invacare TDX SP2 Group 24 160 lbs
- Invacare TDX SP2 HD Group 24 160 lbs
- Invacare AVIVA FX Group 24 170 lbs
- Invacare AVIVA FX Group 34 170 lbs
- Invacare AVIVA Storm RX NARROW Group 50 . 153 LBS
- Invacare AVIVA Storm RX Group 24 155 LBS
- Invacare AVIVA FX Group 24 170 LBS

BATTERY WEIGHT

- Group 22 (Pair) 75 lbs
- Group 24 (Pair) 102 lbs
- Group 24 Flat Top (Pair) 101 lbs
- Group 34 (Pair) 77 lbs
- M 40 (Each) QTY: 1 2 3 22 lbs
- Group 50 (Each) QTY: 1 2 3 33 lbs
- U1 (Each) QTY: 1 2 3 23 lbs

JOYSTICK WEIGHT

- Joystick 1 lb
- Maxx Resolve Swing-Away Joystick Mount 2.2 lb
- Height Adjustable Swing-Away Quad Link 1.5 lb
- Attendant Control & Mount 1.5 lb

SYSTEMS WEIGHT

- Rehab Seat System 50 lbs
- Ultra Low Maxx CG Tilt (50°) Only System 41 lbs
- Ultra Low Maxx CG Tilt (50°) & Power Recline (170°) 55 lbs
- Ultra Low Maxx Recline (170°) Only System 55 lbs

MODULES WEIGHT

- Ultra Low Maxx CG Tilt (50°) Module 27 lbs
- Ultra Low Maxx CG Tilt (45°) & 12" Elevate Lift Seat .55 lbs
- Ultra Low Maxx CG Tilt (50°) & Power Recline (170°) 27 lbs
- Ultra Low Maxx CG Tilt (50°) Power Recline (170°) & 12" Elevate Lift 55 lbs
- Ultra Low Maxx 12" Elevate Lift Seat 35 lbs
- Ultra Low Maxx 10" HD Elevate Lift Seat 45 lbs
- Lateral Tilt 33 lbs

CUSHION OPTIONS WEIGHT

- Matrx Libra Seat Cushion 14"w - 20"w 8 lbs
- Matrx Libra Adjust Seat Cushion 14"w - 20"w .. 8 lbs
- Matrx PS Seat Cushion 14"w - 20"w 5 lbs
- Matrx PSP Seat Cushion 14"w - 20"w 5 lbs
- Matrx Vi Seat Cushion 14"w - 20"w 6 lbs

HEADREST HARDWARE WEIGHT

- Fixed Mounting 2 lbs
- Elan Headrest Hardware - Multi-Axis Removable . 3 lbs
- Loxx Headrest Hardware 1.3 lbs
- Motion Concepts Multi - Axis Removable MTG ... 3 lbs

HEADREST PAD OPTIONS WEIGHT

- Matrx Elan Standard Pad 1 lb
- Matrx Elan Occipital Pad 1 lb
- Matrx Elan 4-Point Pad 1.5 lbs
- Motion Concepts Standard 1 lb
- Motion Concepts Auto Style 2 lbs

BACK OPTIONS WEIGHT

- Standard Rehab Back and Back cushion 7 lbs
- Planar interface plate 8 lbs
- Contour Back 15"w 5 lbs
- Contour Back 16" - 20"w 5.5 lbs
- Contour Back 21" - 22" w 6.5 lbs
- High Back 16" - 22"w 8.5 lbs
- Fixed Space Saver Hardware must be selected for ALL Matrx Backs used on Tilt Only Systems 5 lbs
- Matrx Elite/E2 Back 14" - 20" w 5 lbs
- Matrx Elite HD Back 21" - 24"w 12.5 lbs
- Matrx Elite TR Back 15" - 20" w 8.5 lbs
- Matrx Elite TR HD Back 21" - 24"w 11.5 lbs
- Matrx Elite/E2 Deep Back 14" - 20" w 5.5 lbs
- Matrx Elite Deep HD Back 21" - 24" w 8.5 lbs
- Matrx PB Back 14" - 20"w 6 lbs
- Matrx PB HD Back 21" - 24"w 8.5 lbs
- Matrx PB Deep Back 14" - 20"w 8.5 lbs
- Matrx PB HD Deep Back 21" - 24"w 10 lbs

ARMRESTS WEIGHT

- Adult Dual Post Adjustable Height With Quick Height Adjustment Lever (Pair) 9 lbs
- Ultra Rail Mounted Flip Back Cantilever Maxx Tilt Arm . 9 lbs
- Reclining Adjustable Height 10 lbs
- Maxx Style Cane Mounted Height & Width & Angle Adjustable Cantilever Flip Back Armrest for Recline Systems (Pair) 9 lbs
- Cane Mounted Height & Angle Adjustable Cantilever Flip Back Armrest with/Quick Lock Mechanism for Recline Systems (Pair) 9 lbs
- Outback Cantilever Recline Arms (Pair) 9 lbs
- Outback Cantilever Tilt Arms (Pair) 10.5 lbs

ARM PADS WEIGHT

- Modular Arm Pads (Pair) 2 lb
- Standard Pads (Pair) 2 lb
- Waterfall Arm Pads (Pair) 2 lb
- Flat Multi Position Pads (Pair) 2 lb
- Gel Pads (Pair) 3 lbs
- Ergonomic Arm Trough (Pair) 3 lbs
- Flat Hand Pad (Pair) 1 lbs
- Arm Pad Storage/Cup Holder (Each) 2 lbs

LATERALS WITH HARDWARE **WEIGHT**

- Fixed Mounted Lateral Trunk Support w/ Pad (Pair) 2 lbs
- Maxx Style Swing-away Multi-Adjustable Mounting Hardware w/Pad (Pair) 3 lbs
- Matrx Standard Fixed (Pair) 2 lbs
- Matrx Offset Fixed (Pair) 2 lbs
- Matrx Elite Back Swing Away Hardware w/Pad (Pair) 3 lbs

FIXED HIP/KNEE SUPPORTS WITH PAD **WEIGHT**

- Lateral Hip Support Pad with Fixed Mounting (Pair) 5 lbs
- Swing Away Removable Hardware for Hip Supports w/pad (Pair) 7 lbs
- Maxx Style Quick Release, Removable, Multi-Axis Mounting Hardware for Hip Supports w/Pad (Pair) 7 lbs
- Extended Maxx Style Quick Release, Removable, Multi-Axis Mounting Hardware for Hip Supports (Pair) 2 lbs

FRONT RIGGING **WEIGHT**

- Fixed Center Mount Foot Platform 9 lbs
- Invacare 70° Swingaway Footrests (Pair) 7 lbs
- LNX Power Center Mount Foot Platform 24 lbs
- Latitude Style Foot Platform 12 lbs
- Latitude Style Padded Foot Plate 1 lbs
- Maxx Style Fixed Center Mount Foot Platform 11 lbs
- Invacare Center Mount Adjustable Knee Angle 10 lbs
- Maxx Style Swing-Away Power Elevating Pivot Legrests (Pair) 27 lbs
- Maxx Style Swing-Away Manual Elevating Pivot Legrests (Pair) 27 lbs
- Heavy Duty 70° Swing away Footrests (Pair) 13 lbs
- Vertical Lifting Platform 40 lbs
- Residual Limb Support 3.5 lbs

FOOT PLATFORM/FOOT PLATES **WEIGHT**

- Multi-Axis Adjustable Angle Foot Plates(each) . . 2 lbs

VENT AND BATTERY TRAYS **WEIGHT**

- Vent Tray - Articulating⁽⁴⁾ (Without battery tray) . . 17 lbs
- Vent Tray - Trilogy 100/200 (14W x 6D) 15 lbs
- Compact-Laptop Vent Tray 16 lbs
- Compact-LaptopVent Bag Holder - Articulating . . 15 lbs
- Compact-LaptopVent Bag Holder - Non-Articulating 12 lbs
- M50 Battery Tray for 3 M50 Batteries. 12 lbs
- Battery Tray for 2 GR34 and 1 M50 Battery 12 lbs

O² HOLDERS **WEIGHT**

- Non-Articulating O² Gas Holder 8 lbs
- Articulating O² Gas Holder. 12 lbs
- Non-Articulating IV Pole 5 lbs
- Articulating IV Pole 9 lbs



NOTE: Due to system variance weight might vary ± 10 lbs