



A FLARE FOR FUNCTION

- Flared Front Casters position tight lower extremities and maintain ideal balance of Storm Series® Power bases, no front loading
- Only available on Invacare® Storm Series® Power bases
- Casters flared outward by 3.25" maintaining base width and increasing length by 3"; turning short base in to a long base
- Motion Concepts' 90° Manual and Power Center Mount Foot Platforms can tuck under seat for clients with limited knee extension
- Provides clearance for larger, wider foot platforms on standard center mounts
- Ideal for Anterior Tilt and Uplift Systems

