

... ANOTHER SOLUTION

matrx[®] GENERA

The **matrx**[®] Genera back provides support and comfort at a price so affordable that no client needs to make do with sling upholstery.

NEW – USER FRIENDLY HARDWARE —
installed and adjusted in minutes!!



Features:

- ✓ Adjustable depth (1") and recline (20°)
- ✓ Contoured shape provides comfort and midline support
- ✓ Hardware adjusts to fit chairs 1" wider than back model width
- ✓ Portable – quick release hardware
- ✓ Durable design with 300lb. weight capacity

BACK	Width	Height
Sizes	16", 18", 20"	16", 20"



MEDICARE CODED E2611 – GENERAL USE BACK

www.motionconcepts.com

usa

700 Ensminger Rd. Suite 112
Tonawanda, NY 14150
Toll Free Tel: 1.888.433.6818
Toll Free Fax: 1.888.433.6834

canada

84 Citation Dr.,
Concord, ON L4K 3C1
Toll Free Tel: 1.866.748.7943
Tel: 905.695.0134
Fax: 905.695.0138

Motion Concepts