



The complete solution, designed with the highest level of skin protection, positioning and adjustability for unique user needs.

HCPCS code: E2624/E2625

**Skin Protection and Positioning
Wheelchair Seat Cushion,
Adjustable**

INNER LINER

Highly flexible Platilon liner provides additional layer of moisture protection

FOAM BASE

High resilient foam made with Ultra-Fresh™ anti-microbial and odor protection. Anatomically contoured to provide balance and stability to the pelvis, optimizing body function



OUTER COVER


Moisture resistant, high stretch, breathable Startex™ cover includes hook and loop fastener strips, lifting strap, convenient zippered side pocket with easy pull zipper tab. Optional Spacer Mesh cover available

FLUID SAC

Dual layer, sealed compartments for maximum protection around bony prominences. Non-temperature sensitive, no kneading or maintenance required

LIBRA FIT KIT ACCESSORIES

Accessories to enhance immersion, correct and accommodate orthopedic issues also available



- Recommended for all users including those at highest risk of skin breakdown
- Ab/adductor contour promotes comfortable midline leg position
- Pre-ischial contour prevents sliding to reduce shear
- Trochanteric shelf promotes level pelvis, lateral stability
- Extended rear cushion radius helps prevent posterior pelvic tilt
- Waterfall front edge accommodates tight hamstrings, foot propulsion
- Standard Sizes: 12"- 20" wide; 300 lb. weight capacity
- Heavy-Duty: 21"- 24" wide; 500 lb. weight capacity



www.motionconcepts.com

Follow us on . . .



USA

700 Ensminger Rd. Suite 112
Tonawanda, NY 14150
Toll Free Tel: 1.888.433.6818
Toll Free Fax: 1.888.433.6834

CANADA

84 Citation Dr.,
Concord, ON L4K 3C1
Toll Free Tel: 1.866.748.7943
Tel: 905.695.0134
Fax: 905.695.0138