


... ANOTHER SOLUTION

matrx™ PS POSTURE SEAT



P.S. - I love it!



Features:

- ✓ Positioning, stability and amazing comfort!
- ✓ Unique waffled ischial recess stabilizes pelvis to prevent sliding
- ✓ Reversible cover! Dartex™ outer cover can be used fabric or smooth side out.
- ✓ Foam contains Ultra-Fresh™ for control of odors and bacteria.
- ✓ Sizes 16", 18" and 20" wide x 16", 18" and 20" deep.
- ✓ 300 lb. weight capacity.
- ✓ Medicare coded E2605
- ✓ ADP approved SES1M0060

www.motionconcepts.com

usa
700 Ensminger Rd. Suite 112
Tonawanda, NY 14150
Toll Free Tel: 1.888.433.6818
Toll Free Fax: 1.888.433.6834

canada
84 Citation Dr.,
Concord, ON L4K 3C1
Toll Free Tel: 1.866.748.7943
Tel: 905.695.0134
Fax: 905.695.0138

