

# ... ANOTHER SOLUTION

## ELITE LATERALS - OFFSET SWING-AWAY



**Features:**

- ✓ Provides an additional 3" of vertical position (higher or lower) relative to mounting slots
- ✓ Mount inside back shell without interfering with depth/angle adjustment
- ✓ Easy push button release
- ✓ 4 pad sizes: XS, S, M & L

## ELITE LATERALS - SWING-AWAY



**Features:**

- ✓ 2½" of vertical adjustment on back pan
- ✓ Mount inside back shell without interfering with depth/angle adjustment
- ✓ Easy push button release
- ✓ 4 pad sizes: XS, S, M & L

## ELITE LATERALS - FIXED



**Features:**

- ✓ 2½" of vertical adjustment on back pan
- ✓ Mount inside back shell without interfering with depth/angle adjustment
- ✓ Available in Standard and Pediatric pad sizes



## ELITE LATERALS - OFFSET



**Features:**

- ✓ Provides an additional 3" of vertical position (higher or lower) relative to mounting slots

[www.motionconcepts.com](http://www.motionconcepts.com)

*usa*  
 700 Ensminger Rd. Suite 112  
 Tonawanda, NY 14150  
 Toll Free Tel: 1.888.433.6818  
 Toll Free Fax: 1.888.433.6834

*canada*  
 84 Citation Dr.,  
 Concord, ON L4K 3C1  
 Toll Free Tel: 1.866.748.7943  
 Tel: 905.695.0134  
 Fax: 905.695.0138

