



ADDITIONAL MEASUREMENTS FOR BARIATRIC CLIENTS

Dealer Name: _____ Dealer Account #: _____

PO #: _____ Tag: _____

Client's Name: _____ Client's current weight: _____

Client's Weight History: _____



1. **Buttock/thigh depth:** distance from the back of the buttocks to behind the knee/calf (seat pan depth)
Left side: _____ Right side: _____
2. **Sitting depth:** distance from the most posterior aspect of the trunk to behind the knee/calf (for seat depth)
Left side: _____ Right side: _____
3. **Elbow height:** vertical distance from the top of the sitting surface to the inferior aspect of the elbow in the sitting position
Left side: _____ Right side: _____
4. **Buttocks height:** vertical distance from the top of the sitting surface to the top of gluteal tissue
Left side: _____ Right side: _____

5. **External foot width:** distance between the outside borders of the feet when in desired seated position _____
6. **External lower leg width:** distance between the outside borders of the lower legs at their widest point, when in desired seated position _____
7. **Hip/thigh width:** width at widest point used to determine seat width _____
8. **External elbow width:** _____
9. **Occiput depth:** horizontal distance between the back of the head and the back of the trunk in the sitting position _____