


ADDITIONAL MEASUREMENTS FOR CLIENTS WITH PANNICULUS

Dealer Name: _____ Dealer Account #: _____

PO #: _____ Tag: _____

Client's Name: _____ Client's current weight: _____

Client's Weight History: _____



1. **Trunk width:** overall width of trunk, at lower ribs or widest point _____

2. **External thigh width:** width between lateral aspect of thighs, at widest point _____

3. **Pannus width:** Overall width of Pannus, at widest point _____

4. **External lower leg width:** width between lateral aspect of lower legs, at widest point _____

5. **Pannus to floor distance:** distance from lowest interior surface of unsupported pannus to floor _____

6. **Pannus depth:** dimension of pannus from most anterior to most posterior aspect
Left side: _____ Right side: _____

7. **Buttock depth to pannus:** distance from most posterior aspect of buttocks to most posterior aspect of pannus
Left side: _____ Right side: _____

8. **Supported pannus to floor distance:** distance from interior surface of supported "lifted" pannus to floor _____

9. **Buttock/thigh depth:** dimension from most posterior aspect of buttocks to popliteal fossa
Left side: _____ Right side: _____