



## Thank you for using Motion Concepts eForms.

Motion Concepts is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email\* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.



### Adobe Acrobat Reader DC

\* Interactive functions of our new forms work best with the latest version of Adobe Acrobat Reader DC visit <https://get.adobe.com/reader/> to download and install on your PC or Mac or visit Google Play/ iTunes to download the Adobe Acrobat Reader DC app for your device.

#### Save

Adobe Acrobat Reader DC allows you to save this form with your content - to complete later or use as a starting point for your next form. Please note that content must be added and saved in Acrobat - saving content from completed forms in the browser may not be possible.

**SAVE**

#### Submit

Adobe Acrobat Reader DC allows you to submit this form electronically via your email client.

Simply click the submit button below and step through the simple process.

**SUBMIT**

**CLEAR FORM**

#### Print

If you do not have access to Adobe Acrobat Reader DC simply print this form and complete it by hand and fax it to our Customer Service Department at:

**905-695-0583**

**PRINT FORM**

**PRINT FORM**      **SAVE**      **CLEAR FORM**

**REQUIRED INFORMATION**

Quote       Order

Date of Order: \_\_\_\_\_ Dealer Account #: \_\_\_\_\_

Dealer Name: \_\_\_\_\_

PO #: \_\_\_\_\_ Tag \_\_\_\_\_

Purchasing Contact: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

E-mail: \_\_\_\_\_

RTS/Therapist: \_\_\_\_\_

Ship to Address: \_\_\_\_\_

City: \_\_\_\_\_ Province: \_\_\_\_\_


Postal Code: \_\_\_\_\_ Client Gender:  M  F

Special Client Conditions: \_\_\_\_\_

**REQUIRED MEASUREMENTS**

To ensure system is accurately configured please fill in all **REQUIRED\*** measurements

<b>Client Height*</b>	
<b>Client Weight (lbs)*</b>	
<b>A. Seat to Shoulder</b>	
<b>B. Trunk Depth</b>	
<b>C. Chest Width</b>	
<b>D. Knee to Back</b>	
<b>E. Seat to Top of Head</b>	
<b>F. Back to Joystick</b>	
<b>G. Seat to Elbow</b>	
<b>H. Hip Width</b>	
<b>I. Knee to Heel*</b>	
Cushion Thickness	



**Customer Service: 1.866.748.7943 | Fax: 1.905.695.0583 | orderscdn@motionconcepts.com**

**ROVI**  
**BASE STATISTICS**

ADP CODES	WP0001287
SUSPENSION SYSTEM	Active Ride Control (ARC)
BASE WIDTH X <sub>3</sub>	23.25"
BASE WIDTH A <sub>3</sub>	25.5"
BASE LENGTH	36.5"
DRIVE WHEEL SIZE <sup>(1)</sup>	14"
CASTER WHEEL SIZE (front/rear)	6"
TURNING RADIUS X <sub>3</sub>	20"
TURNING RADIUS A <sub>3</sub>	20.5"
DRIVETRAIN	4-pole EAD Motors

BATTERIES	2 x M34
BATTERY CHARGER	8A, Off Board
BATTERIES WEIGHT (each)	38.5 lbs (M34) / 17.5 kg
RANGE <sup>(2)</sup>	17.7 miles / 28.5 km/h
WEIGHT CAPACITY	300 lbs / 136 kg
BASE WEIGHT X <sub>3</sub> (w/o batteries)	157 lbs / 71 kg
BASE WEIGHT A <sub>3</sub> (w/o batteries)	159 lbs / 72 kg
GROUND CLEARANCE	3"
MAXIMUM SPEED	6.2 mph

**X3**



**A3** **NEW**



<sup>(1)</sup> All X3 and A3 wheels come with solid inserts

<sup>(2)</sup> Actual driving range and speed may vary due to factors such as user weight, type and grade of terrain, battery condition & charge level; type and condition of the drive wheels, and variations in the drive and control systems.

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**X3/A3 BASE OPTIONS**

- 888WNLEX3-RNET58  
**ROVI X<sub>3</sub> 5.8 MPH QUICK SPEED MOTORS<sup>(1)(2)</sup> . \$7,395**
- 888WNLEX3-RNET62  
**ROVI X<sub>3</sub> 6.2 MPH HIGH SPEED MOTORS<sup>(1)(2)</sup> . \$7,395**  
*must order OnTraxx Enhanced Driving Module*
- 888WNLEA3-RNET58  
**ROVI A<sub>3</sub> 5.8 MPH QUICK SPEED MOTORS<sup>(1)(2)</sup> \$6,800**
- 888WNLEA3-RNET62  
**ROVI A<sub>3</sub> 6.2 MPH HIGH SPEED MOTORS<sup>(1)(2)</sup> . \$6,800**  
*must order OnTraxx Enhanced Driving Module*

- OEDM OnTraxx Enhanced Driving Module . . . . . \$1,175
- D51109 R-Net Expandable Controller . . . . . N/C

1. Price of the powerbase includes: Active Ride Control suspension, tie-down hardware for 4-point strap system, 8 amp off-board charger, flat free inserts and positioning belt. Weight capacity is 300 lbs.  
2. When selecting an increased speed motor package there is a reduction in torque output by roughly 15%.

**FENDER COLOUR**

- 800914-00200-03 **Blizzard White**
- 800914-00200-69 **Solar Yellow**
- 800914-00200-23 **Original Orange**
- 800914-00200-21 **Red Delicious**
- 800914-00200-2C **Party Pink**



**NEW Style Fender Colors<sup>(1)</sup>**

- BCF **Black Carbon (Matte)**
- BPF **Black Panther (Gloss)**
- PRF **Purple Reign**
- IBF **Indigo Blue**
- MGF **Motion Green**



1. The new style fender design for VR2/R-Net Bases. For R-Net bases with operational lights please select Fender Color with Lights option below.

**FENDER COLOR WITH LIGHTS<sup>(1)</sup>**

- 800914-00201-03 **Blizzard White** . . . . . \$995
- BCFL **Black Carbon (Matte)** . . . . . \$995
- BPFL **Black Panther (Gloss)** . . . . . \$995
- 800914-00201-69 **Solar Yellow** . . . . . \$995
- 800914-00201-36 **Motion Green** . . . . . \$995
- 800914-00201-23 **Original Orange** . . . . . \$995
- IBFL **Indigo Blue** . . . . . \$995
- 800914-00201-21 **Red Delicious** . . . . . \$995
- 800914-00201-2C **Party Pink** . . . . . \$995
- PRFL **Purple Reign** . . . . . \$995

1. Fully integrated LED lighting kit. Includes headlights, taillights and turn signals. Must select Joystick for Lights feature or Omni for operation (page 3).



**WHEELCHAIR TRANSPORT BRACKETS**

- XBKTS Wheelchair Transport Brackets<sup>(1)</sup> . . . . . STD
1. For unoccupied use only

**TIRE OPTIONS**

- STIRES Standard Solid Tires . . . . . STD
- PTIRES Pneumatic Tires . . . . . N/C

**BATTERIES**

- M34 SLD G MK Group Battery x 2 . . . . . \$650

**JOYSTICK SELECTION**

- D51122 R-Net LED Joystick<sup>(1)(2)</sup> . . . . . N/C
- D51108 R-Net LED Joystick for Lights<sup>(1)</sup> . . . . . \$500
- D51165 R-Net Color Remote Joystick with Monojacks<sup>(1)(2)</sup> . . . . \$875
- D50982 R-Net Color Remote Joystick with Monojacks for Lights<sup>(1)</sup> . . . . . \$975
- D51711 R-Net Advanced Joystick<sup>(1)(3)</sup> . . . . . \$1,595
- D50882 R-Net Attendant Control . . . . . \$995

1. Requires an R-Net Controller  
2. Light package not available with this joystick  
3. 3.5" high-resolution VGA colour display, IR Control as standard, includes integrated Bluetooth Module.

**OMNI 2 CONTROL MODULE FOR R-NET**

- OMNI2 OMNI 2 Control Module for R-Net<sup>(1)</sup> . . . . . \$2,295
- Includes:  D51709 R-Net Omni 2 Display Section Bluetooth  
 D51710 R-Net Omni2 Input Section USB

Select Cable Assembly Length:

- SA81774MK1500-1500 mm (4' 11")
- SA81774MK2000- 2000 mm (6' 7")
- SA81774MK2500- 2500 mm (8' 2")

- ARMS260-ROVI-ML Adjustable Swing-Away Mount - Omni Left . . \$339
- ARMS260-ROVI-MR Adjustable Swing-Away Mount - Omni Right . \$339

1. Requires an R-Net Controller and Harness

**JOYSTICK MOUNT**

- 800914-001-16-ML Fixed Joystick Mounting Left . . . . . N/C
- 800914-001-16-MR Fixed Joystick Mounting Right . . . . . N/C
- MRR Maxx Resolve Swing-Away Joystick Mount (Right Hand) . . . . . \$450
- MRL Maxx Resolve Swing-Away Joystick Mount (Left Hand) . . . . . \$450
- SAQR Height Adjustable Swing-Away Quad Link (Right Hand) . . . . . \$244
- SAQL Height Adjustable Swing-Away Quad Link (Left Hand) . . . . . \$244

**JOYSTICK KNOB - OPTIONAL**

- PC101A Bodypoint U Shaped Handle 3"<sup>(1)</sup> . . . . . \$120
- PC102A Bodypoint U Shaped Handle 4"<sup>(1)</sup> . . . . . \$120
- PC104A BodyPoint Straight Handle . . . . . \$120
- PC107A Bodypoint Rubber Dome<sup>(1)</sup> . . . . . \$120
- PC110A Grooved Mushroom Joystick Handle 3/16" Stem<sup>(1)</sup> . . . . . \$120

**Ultra Low Maxx Power Positioning**  
on ROVI X3/A3 Base

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**ELECTRONIC OPTIONS**

- D50883 R-Net Output Module (EMC)<sup>(1)(2)</sup> . . . . . \$995
- D51580 R-Net Bluetooth Module IOS Apple<sup>(1)(2)</sup> . . . . . \$895
- D51111 R-Net Bluetooth Module Android/PC<sup>(1)(2)</sup> . . . . . \$895
- SPH Smart Phone Holder<sup>(3)</sup> . . . . . \$50
- MUSB Motion Dual USB Charger . . . . . N/C

1. Requires an R-Net Controller
2. Must order R-Net Colour Remote joystick with these options
3. Not available with OMNI2

**ELECTRONIC PROGRAMMER**

- D51812 R-Net Wireless Programmer . . . . . \$750

**SEAT-TO-FLOOR HEIGHT<sup>(1)</sup>**

- STF17 17" Seat-To-Floor . . . . . N/C
- STF18 18" Seat-To-Floor . . . . . N/C
- STF19 19" Seat-To-Floor . . . . . N/C
- STF20 20" Seat-To-Floor . . . . . N/C

1. Seat-to-floor height is measured at the front of the seat pan based on a standard seat size. Seat-to-floor height may vary +/- .25" dependent on condition of the power base and seat depths greater than 18" .

**REHAB SEAT**

- D0195 Ultra Low Maxx Rehab Seat<sup>(1)</sup> . . . . . \$2,995

1. Not available with 12" Elevating Seat Lift Module, Order Ultra Low Maxx Elevating Lift Seat Only.

**POWER POSITIONING SYSTEMS**

- MDCGT-50 Ultra Low Maxx CG Tilt (50°) Only . . . . . \$6,969
- MDCGT-45 Ultra Low Maxx CG Tilt (45°) Only<sup>(1)</sup> . . . . . \$6,969
- MDCGTRESR-50 Ultra Low Maxx CG Tilt (50°) & Power Recline (170°) w/Extended Shear Reduction(ESR) . . \$11,569
- MDCGTRESR-45<sup>(1)</sup> Ultra Low Maxx CG Tilt (45°) & Power Recline (170°) w/Extended Shear Reduction(ESR) . \$11,569
- MDCGTRPSB-50 Ultra Low Maxx CG Tilt (50°) & Power Recline (170°) w/Power Sliding Back(PSB) . . . . . \$12,069
- MDCGTRPSB-45<sup>(1)</sup> Ultra Low Maxx CG Tilt (45°) & Power Recline (170°) w/Power Sliding Back(PSB)<sup>(1)</sup> . . . . . \$12,069
- MDRESR Ultra Low Maxx Power Recline (170°) w/Extended Shear Reduction(ESR) . . . . \$10,350
- MDRPSB Ultra Low Maxx Power Recline (170°) w/Power Sliding Back(PSB)<sup>(2)</sup> . . . . . \$10,850
- MDES12 ROVI X<sub>3</sub> Ultra Low Maxx Elevate Lift Seat Only<sup>(3)(4)</sup> . . . . . \$4,995  
Weight capacity up to 250lbs

1. Tilt will be reduced to 45° if adding an Elevating Lift Seat Module
2. Only available on seat sized 17" wide and greater
3. Not available on seat sizes greater than 20" wide by 20" deep.
4. Only available on Rovi X3 Base.

**ADD: ELEVATING LIFT SEAT MODULE**

- ESM12 12" Elevating Lift Seat Module<sup>(1)(2)(3)</sup> . . . \$2,495  
Weight capacity up to 250lbs

1. Uses an integrated Tilt/Lift mechanism, tilt is 45° and the seat to floor height is not affected by the Elevating Lift Seat.
2. Not available on seat sizes greater than 20" wide by 20" deep
3. For weights greater than 175lbs tilt while elevated is not available.

**ADD: HD ELEVATING LIFT SEAT MODULE**

- ESM10 10" Elevating Lift Seat Module<sup>(1)(2)</sup> . . . \$2,995  
Weight capacity up to 300bs

1. Only Available with Tilt. Uses an integrated Tilt/Lift mechanism, tilt is 45°
2. Not available on seat sizes greater than 22" wide by 22" deep

**MODULES - USE ON SYSTEM**

- 5PT 5° Fixed Pre-Tilt<sup>(1)</sup> . . . . . N/C
- 8PR 8° Precline Module<sup>(2)</sup> only available with Recline . . . . N/C
- 30PR 30° Precline Module<sup>(2)</sup> only available with Recline . . . \$500

1. This Fixed Pre-tilt option adds approximately 2" to the front seat to floor height
2. Total degrees of recline is 170° ie. if 30° of precline is used, then 140° of recline will be available.

**ADULT SEAT WIDTH SETTINGS**

- SW116 16" Wide (Adjustable from 16" to 20") . . . . . N/C
- SW117 17" Wide (Adjustable from 16" to 20") . . . . . N/C
- SW118 18" Wide (Adjustable from 16" to 20") . . . . . N/C
- SW119 19" Wide (Adjustable from 16" to 20") . . . . . N/C
- SW120 20" Wide (Adjustable from 16" to 20") . . . . . N/C
- W219 19" Wide (Adjustable from 19" to 22") . . . . . \$460
- W220 20" Wide (Adjustable from 19" to 22") . . . . . \$460
- W221 21" Wide (Adjustable from 19" to 22") . . . . . \$460
- W222 22" Wide (Adjustable from 19" to 22") . . . . . \$460
- BW123 23" Wide (Adjustable from 20" to 24"<sup>(1)</sup>) . . . . . \$460
- BW124 24" Wide (Adjustable from 20" to 24"<sup>(1)</sup>) . . . . . \$460

1. This option is only available on the Rovi A3 Base.

**ADULT ULTRA RAIL SEAT DEPTH SETTINGS**

- SD116 16" Deep (Adjustable from 16" to 21") . . . . . N/C
- SD117 17" Deep (Adjustable from 16" to 21") . . . . . N/C
- SD118 18" Deep (Adjustable from 16" to 21") . . . . . N/C
- SD119 19" Deep (Adjustable from 16" to 21") . . . . . N/C
- SD120 20" Deep (Adjustable from 16" to 21") . . . . . N/C
- D219 19" Deep (Adjustable from 19" to 23") . . . . . \$460
- D220 20" Deep (Adjustable from 19" to 23") . . . . . \$460
- D221 21" Deep (Adjustable from 19" to 23") . . . . . \$460
- D222 22" Deep (Adjustable from 19" to 23") . . . . . \$460
- D223 23" Deep (Adjustable from 19" to 23") . . . . . \$460

**Ultra Low Maxx Power Positioning**  
on ROVI X3/A3 Base

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**SMALL ADULT SEAT WIDTH SETTINGS<sup>(1)</sup>**

- W314 14" Wide (Adjustable from 14" to 15") ..... \$250
- W315 15" Wide (Adjustable from 14" to 15") ..... \$250

**SMALL ADULT SEAT DEPTH SETTINGS<sup>(1)</sup>**

- D314 14" Deep (Adjustable from 14" to 18") ..... \$250
- D315 15" Deep (Adjustable from 14" to 18") ..... \$250
- D316 16" Deep (Adjustable from 14" to 18") ..... \$250
- D317 17" Deep (Adjustable from 14" to 18") ..... \$250
- D318 18" Deep (Adjustable from 14" to 18") ..... \$250

1. Small Adult Seat Size - weight limit is 175 lbs on these seat sizes and if recline is selected it may be limited to 165° of recline. You must select both small adult seat width and depth together.

**MATRIX SEATING OPTIONS**

- Matrix Cushions available on page 10
- Matrix Backs available on page 11
- Matrix Headrests available on page 12



**CHEST STRAPS**

- MFUS BP Monoflex, Underarm RLS S(3" x 15.5") . . \$195
- MFUM BP Monoflex, Underarm RLS M(4" x 18.25") \$195
- MFUL BP Monoflex, Underarm RLS L(4.5" x 21.5") . \$195
- MFCS BP Monoflex, Center RLS S(2.5" x 15.5") . . \$195
- MFCM BP Monoflex, Center RLS M(3.25" x 18.25") \$195
- MFCL BP Monoflex, Center RLS L(3.5" x 21.5") . . . \$195

**LAP BELTS**

- BELT48 Push Button Style Lap Belt, Length 48" . . . . . STD
- BELT60 Push Button Style Lap Belt, Length 60" . . . . . N/C
- BELT71 Push Button Style Lap Belt, Length 71" . . . . . N/C
- PLB2 Bodypoint Padded Lap Belt 2 Point<sup>(1)</sup> . . . . . \$250
- PLB4 Bodypoint Padded Lap Belt 4 Point<sup>(1)</sup> . . . . . \$315

1. 67" length only option

**MOTION CONCEPTS HEADREST OPTIONS**

- MCSH Motion Concepts Standard<sup>(1)</sup> ..... N/C
- MCOH Motion Concepts Auto Style<sup>(1)</sup> ..... \$235
- FM Fixed Mounting (non-removable)<sup>(2)</sup>..... N/C
- OHR Omit: Headrest.....N/C

1. Must select Fixed Mounting (FM) or Matrix Headrest Hardware on page 12.  
2. Only compatible with Standard Rehab, Contour or High Back Models. All other upgradeable backs require removable hardware.

**MATRIX HEADREST AND HARDWARE OPTIONS**

- Matrix Elan Headrests and Hardware available on page 12



**FINISHED BACK HEIGHT<sup>(1)(2)</sup> FOR TILT ONLY/ELEVATE ONLY**

- 18" FBT18
- 19" FBT19
- 20" FBT20
- 21" FBT21
- 22" FBT22
- 23" FBT23
- 24" FBT24
- 25" FBT25

1. 4" Gap from the seat pan to the bottom of back pan is recommended.  
2. Finished back height is measured from seat pan to top of backrest.

**FINISHED BACK HEIGHT<sup>(1)(2)</sup> FOR RECLINE OR TILT/RECLINE**

- 20" FBR20
- 21" FBR21
- 22" FBR22
- 23" FBR23
- 24" FBR24
- 25" FBR25
- 26" FBR26
- 27" FBR27

1. 6" Gap from the seat pan to the bottom of back pan is recommended.  
2. Finished back height is measured from seat pan to top of backrest.

**BACK CANES (NON RECLINE SYSTEMS)**

- SBC18 18" Straight Back Canes . . . . . STD
- SBC20 20" Straight Back Canes . . . . . STD
- SBC22 22" Straight Back Canes . . . . . STD
- SBC24 24" Straight Back Canes . . . . . STD

Select Cane Angle  75° BA75  79° BA79  85° BA85  90° BA90  95° BA95  101° BA101  106° BA106  112° BA112  116° BA116

- ABC20 20" Angled (5°) Back Canes . . . . . N/C
- ABC22 22" Angled (5°) Back Canes . . . . . N/C
- ABC24 24" Angled (5°) Back Canes . . . . . N/C

Select Cane Angle  81° BA81  84° BA84  91° BA91  96° BA96  101° BA101  106° BA106  111° BA111  117° BA117  121° BA121

- EBC22 22" Endomorph Back Canes. . . . . \$250
- EBC24 24" Endomorph Back Canes. . . . . \$250
- EBCN1 1" Narrower . . . . . N/C
- EBCN2 2" Narrower . . . . . N/C
- EBCN4 4" Narrower . . . . . N/C

**BACK OPTIONS**

- SRB Standard Rehab Back & Cushion . . . . . N/C
- PIPF Flat Planar Interface Plate<sup>(1)</sup> . . . . . N/C
- PIPR Recessed Planar Interface Plate<sup>(1)(2)</sup> . . . . . N/C
- ORBC Omit Rigid Back & Cushion<sup>(3)</sup> . . . . . N/C

1. Not available on tilt only or tilt/elevate systems, for use with Recline Module.  
2. Included with all Recline Systems when no back is selected. May also be used as an interface for after market backs.  
3. Not available with Recline/ESR or Recline/PSB

**MATRIX BACKREST OPTIONS**

- Matrix Backs available on page 11



Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**ARMRESTS (NON RECLINE SYSTEMS)**

- DP1 Adult Dual Post Adj. 9.75"-12.75" . . . . . STD
- DP2 Adult Dual Post Adj. 12.75"-15.75" . . . . . STD
- DP3 Adult Dual Post Adj. 15.75"-18.75" . . . . . STD
- Select:  FP Flat Pouch . . . . . STD
- GB Glove Boxes . . . . . N/C
- OAL ADD: OutBack Arm Option Left<sup>(1)</sup> each . . . . \$195
- OAR ADD: OutBack Arm Option Right<sup>(1)</sup> each . . . . \$195
- RMCAN Ultra Rail Seat Mnt Flip Back Cantilever  
MAXX Tilt Arm (9-12.5")<sup>(2)(3)</sup> . . . . . \$488
- RMCAN2 Ultra Rail Seat Mnt Flip Back Cantilever  
MAXX Tilt Arm (11-14")<sup>(2)(3)</sup> . . . . . \$488
- DPM ADD: Dual Post Module to RMCAN (pair) . . . \$250
- OCTA1 Outback Cantilever Tilt Armrest 9"-12" . . . . \$488
- OCTA2 Outback Cantilever Titl Armrest 12"-15" . . . \$488
- OCTA3 Outback Cantilever Titl Armrest 6"-12" . . . . \$488
- CAFEB Cantilver Arms for Endomorph  
Back Canes (pair)<sup>(3)</sup> . . . . . N/C

1. Only available with Dual Post Armrests.  
2. Mounts to Seat Rail. Slotted upper tube for infinite fore/aft adjustment of arm pad.  
3. Not available with Small Adult Sizes.

**RECLINE ARMRESTS**

- Reclining Adjustable Height . . . . . STD
- RC1 9.5"-13" . . . . . STD
- RC2 13"-16" . . . . . STD
- OAL Outback Arm Option Left . . . . . \$195
- OAR Outback Arm Option Right . . . . . \$195
- Maxx Style Cane Mounted Height, Width, Angle Adj.  
Cantilever Flip Back Armrests
- MCAN1 9.5"-13" . . . . . \$488
- MCAN2 12.5"-16" . . . . . \$488
- Outback Cantilever Recline Armrest
- OCRA1 9"-12"<sup>(1)</sup> . . . . . \$488
- OCRA2 12"-15"<sup>(1)</sup> . . . . . \$488
- OCRA3 7" - 9"<sup>(1)</sup> . . . . . \$488
- Cane Mounted Height and Angle Adj.  
Cantilver Flip Back Armrest
- CANQ1 9.5"-12"<sup>(1)</sup> . . . . . \$488
- CANQ2 10"-13"<sup>(1)</sup> . . . . . \$488
- CANQ3 13"-16"<sup>(1)</sup> . . . . . \$488
- Adult Dual Post On Recline
- DPOR1 9.75" - 12.75" . . . . . \$250
- DPOR2 12.75" - 15.75" . . . . . \$250
- DPOR3 15.75" - 18.75" . . . . . \$250
- FP Flat Pouch . . . . . STD
- GB Glove Boxes . . . . . N/C

1. Not available with Small Adult Sizes.

**ARMREST TUBE LENGTH - PLEASE SPECIFY**

- ATFL Left Full . . . . . STD
- ATDL Left Desk . . . . . STD
- ATRF Right Full . . . . . STD
- ATDR Right Desk . . . . . STD

**ARMPADS**

**Standard**

- SAPFL Left Full (2.25" x 14") . . . . . STD
- SAPDL Left Desk (2.25" x 10") . . . . . STD
- SAPFR Right Full (2.25" x 14") . . . . . STD
- SAPDR Right Desk (2.25" x 10") . . . . . STD

**Waterfall Arm Pads<sup>(1)</sup>**

- WAPFL Left Full (13") . . . . . N/C
- WAPDL Left Desk (9") . . . . . N/C
- WAPFR Right Full (13") . . . . . N/C
- WAPDR Right Desk (10") . . . . . N/C

**Modular Arm Pad - Left<sup>(2)</sup>**

- MAPTFL Left Tray Full (14") . . . . . \$95
- MAPTDL Left Tray Desk (10") . . . . . \$95

Select Pad Insert below:

- ISSFL Iskin - Integrated soft skin foam
- SCVFL Startex Covered Visco foam

**Modular Arm Pad - Right<sup>(2)</sup>**

- MAPTFR Right Tray Full (14") . . . . . \$95
- MAPTDR Right Tray Desk (10") . . . . . \$95

Select Pad Insert below:

- ISSFR Iskin - Integrated soft skin foam
- SCVFR Startex Covered Visco foam

**Flat Multi Position Arm Pads**

- MPAFL Left Full (14") . . . . . \$150
- MPADL Left Desk (10") . . . . . \$150
- MPAFR Right Full (14") . . . . . \$150
- MPADR Right Desk (10") . . . . . \$150

**Gel Pads**

- GAPNL Narrow Left (2" x 12") . . . . . \$88
- GAPNR Narrow Right (2" x 12") . . . . . \$88
- GAPWL Wide Left (3.5" x 12") . . . . . \$96
- GAPWR Wide Right (3.5" x 12") . . . . . \$96
- GAPLL Long Left (3.5" x 14") . . . . . \$103
- GAPLR Long Right (3.5" x 14") . . . . . \$103

**Arm Troughs**

- EATL Ergonomic Arm Trough - Left . . . . . \$121
- EATR Ergonomic Arm Trough - Right . . . . . \$121
- FHPL Flat Hand Pad - Left<sup>(3)</sup> . . . . . \$91
- FHPR Flat Hand Pad - Right<sup>(3)</sup> . . . . . \$91

**Arm Pad Storage/Cup Holder<sup>(4)</sup>**

- APSL Left . . . . . \$295
- APSR Right . . . . . \$295
- OMIT Omit Armpads(7/8" diameter tube) . . . . . N/C

1. Full length pad not compatible with the desk length armrest tube if using Ultra Rail mounted Tilt Arms.  
2. Requires interchangeable pad insert.  
3. Must order with ergonomic arm trough.  
4. Not available on the joystick side.

**ELBOW BLOCKS**

- Elbow Blocks available on page 14

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**ARM ADDUCTOR PAD<sup>(1)(2)</sup>**

- AADSL Left Small (6" x 2.25") . . . . . \$135
- AADSR Right Small (6" x 2.25") . . . . . \$135
- AADLL Left Large (8.5" x 2.25") . . . . . \$135
- AADLR Right Large (8.5" x 2.25") . . . . . \$135

1. Only available with Modular or Flat Multi Position Arm Pads.  
2. 1" thick pad, unlimited fore/aft adjustment.

**MANUAL CENTER MOUNT FRONT RIGGING**


- MSFCMFP Maxx Style Fixed Center Mount Foot Platform<sup>(1)(2)</sup> . . . . . \$275
- MIFPN Rubber Coated Footplate 10.5" W x 10" . . . . . N/C
- ICP Individual calf pads . . . . . \$100

Select seat pan to footplate:

- SPTF1 9" - 13" . . . . . N/C
- SPTF14 13" - 19" . . . . . N/C

Set at:

- SA83 Set at 83° . . . . . N/C
- SA90 Set at 90° . . . . . N/C
- SA97 Set at 97° . . . . . N/C



1. Depth adjustable, and easy slide adjustment for seat pan to foot plate extension.  
2. 3" ground clearance required.

**POWER CENTER MOUNT FOOT PLATFORM**

MUST BE COMPLETED	Measurement in Inches
Knee to Heel:	
Cushion Thickness: (Subtract)	
Cushion Compression: (Add)	
Seat Pan to Footplate: (Total)	

**LNx Power Center Mount Foot Platform**


- LN<sub>X</sub> 9" Articulation<sup>(1)(2)</sup> . . . . . \$3,000
    - LSPTF1 10" - 14" Seat pan to footplate
    - LSPTF2 13" - 19" Seat pan to footplate
  - LN<sub>XA</sub> 6" Articulation<sup>(1)(2)(5)</sup> . . . . . \$3,000
    - LASPTF2 10" - 14" Seat pan to footplate
- Options:
- LIFPN Rubber coated footplate 10.5" W x 10" D . . . . . STD
  - LICP Individual calf pads . . . . . STD
  - OPCP One Piece Malleable Calf Panel<sup>(3)</sup> . . . . . \$250
  - LNXPF Plantar Flex Option for LN<sub>X</sub> Foot platform<sup>(4)</sup>. \$250

**Latitude Style Foot Platform**

- LSFP Latitude Style Foot Platform. . . . . \$2,995  
For finished knee to heel measurements of 6.5"-12"
- LSPFP Latitude Style Padded Foot Plate . . . . . \$250

1. 3" ground clearance required.  
2. Only available on adult range seat systems 16" deep or greater.  
For shorter seat depths call Customer Service for other options.  
3. Not available with Plantar Flex Option for LN<sub>X</sub> Foot platform, or with LN<sub>XA</sub>.  
4. Only available with Foot Platform Options.  
5. For clients under 4' only.


**FOOT PLATFORM OPTIONS FOR FIXED/LN<sub>X</sub> CMFP**



- SFFP Small Foot Platform<sup>(1)</sup>(12"w x 8.5"d) . . . . . \$180
- MFFP Medium Platform<sup>(1)</sup> (12"w x 10.75"d) . . . . . \$180
- LFFP Large Platform<sup>(1)</sup> (14.5"w x 13"d) . . . . . \$250
- XLFFP Extra Large Platform<sup>(1)</sup> (17"w x 13"d) . . . . . \$250
- XXLFFP Extra Extra Large Platform<sup>(1)</sup> (21"w x 13"d). \$250

1. May need to adjust leg position to avoid interference with front castors

**INDIVIDUAL FOOT PLATES FOR FIXED/LN<sub>X</sub> CMFP**



- IFSL Small<sup>(1)</sup> - (5" x 7.5" deep) - Left . . . . . \$165
- IFSR Small<sup>(1)</sup> - (5" x 7.5" deep) - Right . . . . . \$165
- IFML Medium<sup>(1)</sup> - (5.5" x 9.5" deep) - Left . . . . . \$165
- IFMR Medium<sup>(1)</sup> - (5.5" x 9.5" deep) - Right . . . . . \$165
- IFLL Large<sup>(1)</sup> - (6" x 11.5" deep) - Left . . . . . \$165
- IFLR Large<sup>(1)</sup> - (6" x 11.5" deep) - Right . . . . . \$165
- MFOL Left Optional Footplate (6"W x 10"D) . . . . . \$225
- MFOR Right Optional Footplate (6"W x 10"D) . . . . . \$225

1. Foot plates may interfere with casters in certain configurations

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**MAXX STYLE POWER ELEVATING LEGRESTS**

MSPPL MAXX Swing Away Power Elevating<sup>(1)(2)</sup> ..... \$2,490

Select Legrest Function

INDA Individual Movement ..... N/C  
 COMA Combined Movement ..... N/C  
 SYNA Synchronized/Recline Movement ..... N/C

Select Receivers:

RS Straight Receivers ..... N/C  
 RF Flared 1" Outward Receivers ..... N/C

Select Seat Pan to Footplate:

SPTF5 Length: 10"-14" ..... N/C  
 SPTF2 Length: 13"-17" ..... N/C  
 SPTF7 Length: 16"-20" ..... N/C

FCPMS Flat Calf Pad Option (pair) ..... \$275

CCPMS Curved Calf Pad Option (inside with 6")<sup>(3)</sup> (pair) \$275

1. Must choose footplate on page 8
2. Comes standard with moulded i-skin calf pads.
3. Not available on 14 and 15" seat widths. For 16" and 17" seat widths you must select flared receivers.

**MAXX STYLE ELEVATING LEGRESTS**

MMPPSA MAXX Swing Away Manual Elevating<sup>(1)(2)</sup> ..... \$680

Select Receivers:

RS Straight Receivers ..... N/C  
 RF Flared 1" Outward Receivers ..... N/C

Select Seat Pan to Footplate:

SPTF5 Length: 10"-14" ..... N/C  
 SPTF2 Length: 13"-17" ..... N/C  
 SPTF7 Length: 16"-20" ..... N/C

FCPMS Flat Calf Pad Option (pair) ..... \$275

CCPMS Curved Calf Pad Option (inside with 6")<sup>(3)</sup> (pair) \$275

1. Must choose footplate on page 8
2. Comes standard with moulded i-skin calf pads
3. Not available on 14 and 15" seat widths. For 16" and 17" seat widths you must select flared receivers.

**SWING AWAY LEGREST**

HD70 Heavy Duty 70° Swing Away Footrest<sup>(1)(2)</sup> ..... \$495

Select Receivers:

RS Straight Receivers ..... N/C  
 RF Flared 1" Outward Receivers ..... N/C

Select Seat Pan to Footplate:

SPTF11 Length: 6"-10" ..... N/C  
 SPTF8 Length: 11"-15" ..... N/C  
 SPTF2 Length: 13"-17" ..... N/C  
 SPTF7 Length: 16"-20" ..... N/C

HDFCP Flat Calf Pad Option (pair) ..... \$275

HD70CCP Curved Calf Pad Option<sup>(3)</sup> (inside with 6") (pair) \$275

DCS Dynamic Coil Springs (Pair) for Heavy Duty 70° Swingaway Footrests ..... \$190

I70SAF Invacare 70° Swingaway Footrests (Pair) ..... \$325  
 I70-CF Comes with composite foot plates ..... N/C

70HDAIH Invacare 70° Non-Tapered HD Swingaway Footrests with AT5543 Adjustable Angle Foot plates, Heel Loops and Impact Guards<sup>(1)</sup> (Pair) ..... \$457

70HTAPAIH Invacare Tapered 70° Swingaway Footrests with AT5543 Adjustable Angle Footplates, Heel Loops and Impact Guards(Pair) ..... \$497

ALPT Longer Pivot and Slide Tube for 70° Swingaways (Pair) ..... \$132

RLST Residual Limb Support - Left<sup>(4)</sup> (6.75"W x 9.75"L x 2"H) (each) ..... \$495

RLST Residual Limb Support - Right<sup>(4)</sup> (6.75"W x 9.75"L x 2"H) (each) ..... \$495

1. Must choose footplate on page 8
2. Heavy duty durable construction (4.5lbs ea.).
3. Not available on 14 and 15" seat widths. For 16" and 17" seat widths you must select flared receivers.
4. Swing away and removable, 50° of angle adjustment, 3.5" height adjustment and 2" of for/aft adjustment.




Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**FOOTPLATE FOR SWING AWAY/ELEVATING LEGS**



1. Certain footplate sizes may not fit all seat widths depending on legrests selected.

**ADJUSTABLE ANGLE FOOTPLATES**

<input type="checkbox"/> AFSL	Left Small Footplate (4.25"W x 8"D) .....	\$165
<input type="checkbox"/> AFSR	Right Small Footplate (4.25"W x 8"D) .....	\$165
<input type="checkbox"/> AFML	Left Medium Footplate (5.25"W x 8"D) .....	\$165
<input type="checkbox"/> AFMR	Right Medium Footplate (5.25"W x 8"D) .....	\$165
<input type="checkbox"/> AFLL	Left Large Footplate (6.25"W x 8"D) .....	\$165
<input type="checkbox"/> AFLR	Right Large Footplate (6.25"W x 8"D) .....	\$165
<input type="checkbox"/> AFXL	Left Extra Large Footplate (9"W x 11"D) .....	\$165
<input type="checkbox"/> AFXR	Right Extra Large Footplate (9"W x 11"D) .....	\$165


**MULTI-AXIS ADJUSTABLE ANGLE FOOTPLATES**

<input type="checkbox"/> MFSL	Left Small Footplate (4.25"W x 8"D) .....	\$225
<input type="checkbox"/> MFSR	Right Small Footplate (4.25"W x 8"D) .....	\$225
<input type="checkbox"/> MFML	Left Medium Footplate (5.25"W x 8"D) .....	\$225
<input type="checkbox"/> MFMR	Right Medium Footplate (5.25"W x 8"D) .....	\$225
<input type="checkbox"/> MFLL	Left Large Footplate (6.25"W x 8"D) .....	\$225
<input type="checkbox"/> MFLR	Right Large Footplate (6.25"W x 8"D) .....	\$225
<input type="checkbox"/> MFXL	Left Extra Large Footplate (9"W x 11"D) .....	\$225
<input type="checkbox"/> MFXR	Right Extra Large Footplate (9"W x D) .....	\$225

**HEEL LOOPS, CUPS AND TOE STRAPS**

<input type="checkbox"/> HLL	Left Heel Loop (each) .....	\$20
<input type="checkbox"/> HLR	Right Heel Loop (each) .....	\$20
<input type="checkbox"/> HCL	Left Cup Loop (each) .....	\$20
<input type="checkbox"/> HCR	Right Cup Loop (each) .....	\$20
<input type="checkbox"/> TSL	Left Toe Strap (each) .....	\$20
<input type="checkbox"/> TSR	Right Toe Strap (each) .....	\$20

**FOOTPLATE/CUT OUTS SIZE SELECTION**



1. Certain footplate sizes may not fit all seat widths depending on legrests selected

**ADJUSTABLE ANGLE FOOTPLATES/CUT OUTS**

<input type="checkbox"/> AFCSL	Left Small Footplate (7.25"W x 9"D) .....	\$165
<input type="checkbox"/> AFCSR	Right Small Footplate (7.25"W x 9"D) .....	\$165
<input type="checkbox"/> AFCML	Left Medium Footplate (7.5"W x 11"D) .....	\$165
<input type="checkbox"/> AFCMR	Right Medium Footplate (7.5"W x 11"D) .....	\$165
<input type="checkbox"/> AFCLL	Left Large Footplate (10"W x 14"D) .....	\$165
<input type="checkbox"/> AFCLR	Right Large Footplate (10"W x 14"D) .....	\$165

**MULTI-AXIS ADJUSTABLE ANGLE FOOTPLATES/ CUT OUTS**

<input type="checkbox"/> MFCSL	Left Small Footplate (7.25"W x 9"D) .....	\$225
<input type="checkbox"/> MFCSR	Right Small Footplate (7.25"W x 9"D) .....	\$225
<input type="checkbox"/> MFCML	Left Medium Footplate (7.5"W x 11"D) .....	\$225
<input type="checkbox"/> MFCMR	Right Medium Footplate (7.5"W x 11"D) .....	\$225
<input type="checkbox"/> MFCLL	Left Large Footplate (10"W x 14"D) .....	\$225
<input type="checkbox"/> MFCLR	Right Large Footplate (10"W x 14"D) .....	\$225

**CALF STRAPS AND CALF PANELS**

<input type="checkbox"/> CSM	Body Point Airmesh Calf Strap for Chairs 15"-18" wide with 3.5" H Pad <sup>(1)</sup> .....	\$200
<input type="checkbox"/> CSL	Body Point Airmesh Calf Strap for Chairs 18"-22" wide with 3.5" H Pad <sup>(1)</sup> .....	\$200
<input type="checkbox"/> CPM	Body Point Airmesh Calf Panel for Chairs 15"-18" wide with 35.5" W x 9" H Pad <sup>(1)</sup> . . .	\$295
<input type="checkbox"/> CPL	Body Point Airmesh Calf Panel for Chairs 18"-22" wide with 38.5" W x 9" H Pad <sup>(1)</sup> . . .	\$315

**THORACIC LATERALS AND HIP/KNEE SUPPORTS**

- Matrx Laterals please see page 12*
- Maxx Laterals and Hip/Knee supports please see page 13*

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**ELECTRONIC SWITCH MOUNTING**

- ESML Left Side (Armpad)
- ESMR Right Side (Armpad)
- ESMJL Joystick Left (Mounts to Joystick Bracket)
- ESMJR Joystick Right (Mounts to Joystick Bracket)
- ESMD Do not mount

**SINGLE-DUAL FUNCTION ELECTRONICS**

- SFCB Single/Dual Function through Switch/Toggle<sup>(1)</sup> . . . . . STD

Select Switch:

- SPBO Single Push Button
- STO Single Toggle
- DPBO Dual Push Button
- SOS Omit Switch

- SFCBA Single/Dual Function through Drive Control<sup>(1)</sup> . . . . . \$750

Select Control:

- JOY Through Joystick
- SC Through Specialty Control

Optional:

- RDE Reduced Drive for Elevating Seat<sup>(2)</sup> . . . . . \$400

- SPBSA Additional Single Push Button . . . . . \$97
- STA Additional Single Toggle Switch . . . . . \$97
- DPBA Additional Dual Push Button . . . . . \$97

1. These electronics include drive lock-out and tilt limit.  
2. If reduced drive is not selected the system will be in drive lock-out when elevated.

**MULTI- FUNCTION ELECTRONICS**

- MFCB Multi-Function through Switch/Toggle<sup>(1)</sup> . . . . . \$1,950

Select Switch:

- QPBO Quad Push Button
- FWTO Four-Way Toggle
- EWRO 8-Way Rocker<sup>(2)</sup>

- MFCBD Multi-Function through Drive Control<sup>(1)</sup> . . . . . \$1,950

Select Control:

- JOY Through Joystick
- SC Through Specialty Control

Optional:

- RDEY Include Reduced Drive for Elevating Seat<sup>(3)</sup>
- RDEN Do not Induce Reduced Drive for Elevating Seat<sup>(3)</sup>

- QPBA Additional Quad Button. . . . . \$150
- FWTA Additional Four-Way Toggle Switch . . . \$150
- EWRA Additional 8-Way Rocker . . . . . \$255

1. These electronics include attendant control, drive lock-out, tilt/recline limit and actuator speed control  
2. Dual function is not available with 8-Way Rocker.  
3. If reduced drive is not selected system will be in drive lock-out when elevated.

**ELECTRONICS OPTIONS**

- M616S M616 Converter<sup>(1)</sup> . . . . . \$125
- FTSQ Feather Touch Switch Single Function Quarter Size . . . . . \$295
- FTMD Feather Touch Switch Multi-Function Dual Dime Size . . . . . \$295
- FTMQ Feather Touch Switch Multi-Function Quad Dime Size . . . . . \$295
- ASL320R ASL Pulse Mechanical Switch Red. . . . . \$200
- ASL320G ASL Pulse Mechanical Switch Green. . . . \$200
- ASL320BK ASL Pulse Mechanical Switch Black. . . . \$200
- PSOS Power Shut Off Switch - Toggle<sup>(1)</sup> . . . . . \$175
- LJS Lazarus Jump Starter<sup>(1)</sup> . . . . . \$150
- APPS Accessory Port Power Supply<sup>(1)</sup> (3 amp capacity) . . . . . \$195
- TPSMH 12 Piece Switch Mounting Hardware Kit (additional kits). . . . . \$39

1. Not available with Single-Dual Function Electronics.

**VENT TRAYS, O2 HOLDERS, ACCESSORIES**

- For Vent Trays, O2 Holders, IV Pole and Additional Accessories  
Please see page 14

**INSTALLATION**

- MCI Motion Concepts Installation . . . . . Included

**PLEASE SPECIFY PAGES INCLUDED WITH ORDER**

- Pages 1-9 Excludes Matrix and Maxx Seating Accessories Thoraic Supports, Hip Supports, Vent Trays and Accessories
- Pages 1-14 Complete Order Form

**Total Retail Value \$**

\_\_\_\_\_

**SUBMIT**

**PRINT FORM**

**SAVE**

**CLEAR FORM**



Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**CUSHION COVER REFERENCE GUIDE**

Meshtex	Startex	Infection Control			
Spacetex	Onyx				
Cushion	Meshtex	Startex	Infection	Spacetex	Onyx
Libra/Adjust	N/A	STD	\$50	N/C	N/A
PS/PSP	N/A	STD-Reversible	STD-Reversible	N/A	N/A
Vi	STD	\$103	\$103	N/A	\$103

**MATRIX SEAT CUSHION OPTIONS**

**Matrix Libra Adjust Seat Cushion 14" -20" wide . . . . . \$985**

- LC1416-ADJ    LC1516-ADJ    LC1518-ADJ    LC1520-ADJ
- LC1616-ADJ    LC1618-ADJ    LC1620-ADJ    LC1716-ADJ
- LC1718-ADJ    LC1720-ADJ    LC1816-ADJ    LC1818-ADJ
- LC1820-ADJ    LC1916-ADJ    LC1918-ADJ    LC1920-ADJ
- LC2016-ADJ    LC2018-ADJ    LC2020-ADJ

**Matrix Libra Adjust Seat Cushion 22" depth . . . . . \$1,265**

- LC1622-ADJ    LC1722-ADJ    LC1822-ADJ
- LC1922-ADJ    LC2022-ADJ

**Matrix Libra Adjust Seat Cushion 21" - 24" wide . . . . . \$1,265**

- LC2118-ADJ    LC2120-ADJ    LC2122-ADJ
- LC2218-ADJ    LC2220-ADJ    LC2222-ADJ
- LC2318-ADJ    LC2320-ADJ    LC2322-ADJ
- LC2418-ADJ    LC2420-ADJ    LC2422-ADJ

**Matrix Libra Seat Cushion 14" - 20" wide . . . . . \$745**

- LC1416    LC1516    LC1518    LC1520    LC1616
- LC1618    LC1620    LC1716    LC1718    LC1720
- LC1816    LC1818    LC1820    LC1916    LC1918
- LC1920    LC2016    LC2018    LC2020

**Matrix Libra Seat Cushion 22" depth . . . . . \$1,025**

- LC1622    LC1722    LC1822    LC1922    LC2022

**Matrix Libra Seat Cushion 21" - 24" wide . . . . . \$1,265**

- LC2118    LC2120-ADJ    LC2122-ADJ
- LC2218-ADJ    LC2220-ADJ    LC2222-ADJ
- LC2318-ADJ    LC2320-ADJ    LC2322-ADJ
- LC2418-ADJ    LC2420-ADJ    LC2422-ADJ

**Matrix PS Seat Cushion 14" - 20" wide . . . . . \$569**

- PS1416    PS1616    PS1618    PS1620    PS1816
- PS1818    PS1820    PS2016    PS2018    PS2020

**Matrix PS HD Seat Cushion 22" - 24" wide . . . . . \$975**

- PS2218    PS2220    PS2222    PS2418    PS2420
- PS2422

**MATRIX SEAT CUSHION OPTIONS**

**Matrix PSP Seat Cushion 14" - 20" wide . . . . . \$645**

- PSP1416    PSP1616    PSP1618    PSP1620
- PSP1816    PSP1818    PSP1820    PSP2016
- PSP2018    PSP2020

**Matrix PSP HD Seat Cushion 22" - 22" wide . . . . . \$995**

- PSP2218    PSP2220    PSP2222

**Matrix Vi Seat Cushion 14" - 20" wide . . . . . \$624**

- MA1415-VI    MA1417-VI    MA1616-VI    MA1618-VI
- MA1620-VI    MA1816-VI    MA1818-VI    MA1820-VI

**Matrix Vi Seat Cushion 21" - 22" wide . . . . . \$681**

- MA2016-VI    MA2018-VI    MA2020-VI    MA2218-VI
- MA2220-VI

**Matrix Vi HD Seat Cushion 22" - 22" wide . . . . . \$840**

- MA2222-VI    MA2418-VI    MA2420-VI    MA2422-VI

**Matrix Multi Cushion 14" - 20" wide/reversible cover**

HCPCS code E2607 . . . . . \$745

- MM1616    MM1618    MM1620    MM1716    MM1718
- MM1720    MM1816    MM1818    MM1820    MM2016
- MM2018    MM2020

**FABRIC OPTIONS FOR CUSHIONS**

- SFLF   Spacetex Fabric Option for  
Libra Cushion . . . . . N/C
- ICFL   Infection Control Fabric Upgrade for  
Libra Cushion . . . . . \$50
- OCUMA   Onyx Fabric Upgrade for  
Matrix Vi Cushion. . . . . \$100
- SFVS   Smooth side Startex Fabric Upgrade for  
Matrix Vi Cushion . . . . . \$50
- SFVF   Fabric side Startex Fabric Upgrade for  
Matrix Vi Cushion . . . . . \$50

**OTHER UPGRADABLE BACK OPTIONS**

- EBSS   Super Soft Foam for Elite Back Upgrade . . . . \$103
- OCUHG   Onyx Fabric Upgrade for High Backs. . . . . \$103
- WSCU   Infection Control Cover for Elite/PB Backs . . \$103  
*purchased with back*
- WSC   Infection Control Cover for Elite/PB Backs . . \$150  
*purchased with without back*
- CHIEB   Chest Harness Interface Kit for Elite Back<sup>(1)(2)</sup> . \$103
- CHIEB   Chest Harness Interface Kit for PB Back<sup>(1)(2)</sup> . . \$52

1. Does not include the chest strap.  
2. This Interface Kit is only Required for 4 point harnesses.






Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**BACK COVER REFERENCE GUIDE**

			
Back	Meshtex	Infection	Onyx
High	STD	N/A	\$103
Elite/ Elite HD/Elite Deep	STD	\$103	N/A
Elite TR/ Elite TR HD	STD	\$103	N/A
PB/PB HD/PB Deep/PB HD Deep	STD	\$103	N/A

**UPGRADABLE BACK OPTIONS**

- NEW** **Matrix Elite E2 Back 14" - 20" wide/14"/16"H(3" Contour) . \$965**
- E2S1414     E2S1416     E2S1514     E2S1516
  - E2S1614     E2S1616     E2S1714     E2S1716
  - E2S1814     E2S1816     E2S1914     E2S1916
  - E2S2014     E2S2016

- NEW** **Matrix Elite E2 Back 14" - 20" wide/18"/20"H(3" Contour) \$995**
- E2S1418     E2S1420     E2S1518     E2S1520
  - E2S1618     E2S1620     E2S1718     E2S1720
  - E2S1818     E2S1820     E2S1918     E2S1920
  - E2S2018     E2S2020

**UPGRADABLE BACK OPTIONS**

**High Back 16" - 22" wide . . . . . \$595**

- HB1616     HB1618     HB1620     HB1716
- HB1718     HB1720     HB1816     HB1818
- HB1820     HB1916     HB1918     HB1920
- HB2016     HB2018     HB2020     HB2116
- HB2118     HB2120     HB2216     HB2218
- HB2220     HB2316     HB2318     HB2320
- HB2416     HB2418     HB2420

**NEW** **Matrix Elite E2 Deep Back 14" - 20" wide/14"/16"H (6" Contour) \$1,125**

- E2D1414     E2D1416     E2D1514     E2D1516
- E2D1614     E2D1616     E2D1714     E2D1716
- E2D1814     E2D1816     E2D1914     E2D1916
- E2D2014     E2D2016

**NEW** **Matrix Elite E2 Deep Back 14" - 20" wide/18"/20"H (6" Contour) \$1,155**

- E2D1418     E2D1420     E2D1518     E2D1520
- E2D1618     E2D1620     E2D1718     E2D1720
- E2D1818     E2D1820     E2D1918     E2D1920
- E2D2018     E2D2020

**Matrix Elite Back 14" - 20" wide (3" Contour) . . . . . \$766**

- PBE1414     PBE1416     PBE1418 . . .  PBE1420
- PBE1514     PBE1516     PBE1518 . . .  PBE1520
- PBE1614     PBE1616     PBE1618 . . .  PBE1620
- PBE1714     PBE1716     PBE1718 . . .  PBE1720
- PBE1814     PBE1816     PBE1818 . . .  PBE1820
- PBE1916     PBE1918     PBE1920 . . .  PBE2014
- PBE2016     PBE2018     PBE2020

**Matrix Elite TR Back 15" - 20" wide (3" Contour) . . . . . \$866**

- PBE1516-TR     PBE1518-TR     PBE1520-TR
- PBE1616-TR     PBE1618-TR     PBE1620-TR
- PBE1716-TR     PBE1718-TR     PBE1720-TR
- PBE1816-TR     PBE1818-TR     PBE1820-TR
- PBE1916-TR     PBE1918-TR     PBE1920-TR
- PBE2016-TR     PBE2018-TR     PBE2020-TR

**Matrix Elite HD Back 21" - 22" wide (3" Contour) . . . . \$1,095**

- PBE2116-HD     PBE2120-HD     PBE2216-HD
- PBE2220-HD     PBE2316-HD     PBE2320-HD
- PBE2416-HD     PBE2420-HD

**Matrix Elite TR HD Back 21" - 24" wide (3" Contour) . . \$1,151**

- PBE2120-TRHD     PBE2220-TRHD     PBE2320-TRHD
- PBE2420-TRHD

**Matrix Elite Deep Back 14" - 20" wide (6" Contour) . . . . \$918**

- EDB1416     EDB1418     EDB1420     EDB1514
- EDB1516     EDB1518     EDB1520     EDB1614
- EDB1616     EDB1618     EDB1620     EDB1714
- EDB1716     EDB1718     EDB1720     EDB1814
- EDB1816     EDB1818     EDB1820     EDB1914
- EDB1916     EDB1918     EDB1920     EDB2014
- EDB2016     EDB2018     EDB2020

**Matrix PB Posture Back 14" - 20" wide (5" Contour) . . . . . \$734**

- MPB1412     MPB1612     MPB1616     MPB1620
- MPB1712     MPB1716     MPB1720     MPB1812
- MPB1816     MPB1820     MPB2012     MPB2016
- MPB2020

**Matrix PB Posture Back Deep 14" - 20" wide (7" Contour) . \$892**

- PBD1412     PBD1612     PBD1616     PBD1620
- PBD1712     PBD1716     PBD1720     PBD1812
- PBD1816     PBD1820     PBD2012     PBD2016
- PBD2020

**Matrix Elite Deep HD Back 21" - 24" wide (6" Contour) . \$1,095**

- EDB2116-HD     EDB2120-HD     EDB2216-HD
- EDB2220-HD     EDB2316-HD     EDB2320-HD
- EDB2416-HD     EDB2420-HD

**Matrix PB Heavy Duty Back 20" - 22" wide (5" Contour) . \$1,051**

- HDB2016     HDB2020     HDB2216     HDB2220

**Matrix PB Heavy Duty Deep Back 20"-22" wide (5" Contour) \$1,209**

- DHDB2016     DHDB2020     DHDB2216     DHDB2220

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**MATRX ELAN HEADREST PAD<sup>(1)</sup>**

Standard Pads

- ESP6      Matrx Elan Standard Pad 6" ..... \$310
- ESP10     Matrx Elan Standard Pad 10" ..... \$310
- ESP14     Matrx Elan Standard Pad 14" ..... \$310
- ESP6-HUG   Matrx Elan Hug Standard Pad 6" ..... \$310
- ESP10-HUG   Matrx Elan Hug Standard Pad 10" ..... \$310
- ESP14-HUG   Matrx Elan Hug Standard Pad 14" ..... \$310
- SNUG14    Matrx SNUG Standard Pad 14" ..... \$625
- SNUG14-A   Matrx SNUG Asymmetrical Pad 14" ..... \$645
- SNUG18    Matrx SNUG Standard Pad 18" ..... \$625
- SNUG18-A   Matrx SNUG Asymmetrical Pad 18" ..... \$645

Standard Pads with Infection Control Cover

- ESP6-IC     Matrx Elan Infection Control Pad 6" ..... \$335
- ESP10-IC    Matrx Elan Infection Control Pad 10" ..... \$335
- ESP14-IC    Matrx Elan Infection Control Pad 14" ..... \$335
- ESP6-IC-HUG   Matrx Elan Hug Infection Control Pad 6" ..... \$335
- ESP10-IC-HUG   Matrx Elan Hug Infection Control Pad 10" ..... \$335
- ESP14-IC-HUG   Matrx Elan Hug Infection Control Pad 14" ..... \$335
- SNUG14-IC   Matrx SNUG Standard Infection Control Pad 14" ..... \$665
- SNUG14-A-IC   Matrx SNUG Asymmetrical Infection Control Pad 14" ..... \$685
- SNUG18-IC   Matrx SNUG Standard Infection Control Pad 18" ..... \$665
- SNUG18-A-IC   Matrx SNUG Asymmetrical Infection Control Pad 18" ..... \$685

NEW

NEW

Occipital /4-Point Pads

- EOP9      Matrx Elan Occipital Pad 9" ..... \$341
- EOP12     Matrx Elan Occipital Pad 12" ..... \$341
- E4PTSM    Matrx Elan 4-Point Pad 9"W x 7"H ..... \$448
- E4POINT    Matrx Elan 4-Point Pad 11"W x 10"H ..... \$448

Occipital / 4-Point Pads with Infection Control Cover

- EOP9-IC    Matrx Elan Occipital Pad 9" with Infection Control ..... \$366
- EOP12-IC   Matrx Elan Occipital Pad 12" with Infection Control ..... \$366
- E4PTSM-IC   Matrx Elan 4-Point Pad 9"W x 7"H with Infection Control ..... \$473
- E4POINT-IC   Matrx Elan 4-Point Pad 11"W x 10"H with Infection Control ..... \$473

1. If mounting on PB Elite must order removable hardware

**HEADREST MOUNTING HARDWARE OPTIONS**

- LXSHW      Loxx Standard Hardware 15" ..... \$275
- LXMMH      Loxx Mini Hardware 11" ..... \$275
- MEHWRE     Elan Hardware 15" ..... \$285
- MEMH       Elan Mini Hardware 11" ..... \$285
- MEHW-EXT   Elan Extended Hardware ..... \$320
- FPLXSHW    FLIPP Loxx Standard Hardware 15" ..... \$705
- FPLXSHW-EX   FLIPP Loxx Extended Hardware ..... \$705
- FPLXMMH    FLIPP Loxx Mini Hardware 11" ..... \$705
- FPLXMMH-EX   FLIPP Loxx Mini Extended Hardware ..... \$705
- FPMEHW     FLIPP Elan Hardware 15" ..... \$715
- FPMEMH     FLIPP Elan Mini Hardware 11" ..... \$715
- FPMEHW-EXT   FLIPP Elan Extended Hardware ..... \$780

**MATRX FIXED LATERALS - FOR ELITE AND ELITE TR BACKS**

- PBSSL      Matrx Small Lateral Fixed-Left(3x5.5") ..... \$113
- PBSSR      Matrx Small Lateral Fixed-Right(3x5.5") ..... \$113
- PBLSL      Matrx Standard Lateral Fixed-Left(7.5x5") ..... \$116
- PBLSR      Matrx Standard Lateral Fixed-Right (7.5x5") ..... \$116
- PBOLL      Matrx Offset Lateral Fixed Left ..... \$125
- PBOLR      Matrx Offset Lateral Fixed Right ..... \$125

**MATRX SWING AWAY LATERALS - FOR ELITE AND ELITE TR BACKS**

Matrx Lateral Trunk Support Left<sup>(1)</sup>

- EPADXL     Size XS (3.25x4.25") ..... \$115
- EPADSL     Size S (3.5x5.25") ..... \$115
- EPADML     Size M (3.75x5.75") ..... \$115
- EPADLL     Size L (4.25x6.5") ..... \$115

Matrx Lateral Trunk Support Right<sup>(1)</sup>

- EPADXR     Size XS (3.25x4.25") ..... \$115
- EPADSR     Size S (3.5x5.25") ..... \$115
- EPADMR     Size M (3.75x5.75") ..... \$115
- EPADLR     Size L (4.25x6.5") ..... \$115

1. Must choose swing away hardware below

**MATRX SWING AWAY HARDWARE - FOR ELITE AND ELITE TR BACKS**

- SAEHBL     Matrx Elite Swing Away Hardware Left ..... \$250
- SAHEBR     Matrx Elite Swing Away Hardware Right ..... \$250
- OSAHL      Matrx Offset Elite Swing Away Hardware Left ..... \$250
- OSAHR      Matrx Offset Elite Swing Away Hardware Right ..... \$250

**OTHER MATRX ACCESSORIES**

Multi Axis Upper Extreimity Support

- MACESL     Multi Axis Arm Support - Left<sup>(1)</sup> ..... \$225
- MACESR     Multi Axis Arm Support - Right<sup>(1)</sup> ..... \$225

1. Allows rotation of the armpad, only compatible with Gel pad, Flat pad or Ergonomic Arm Troughs

**MAXX LATERALS TRUNK SUPPORTS<sup>(1)</sup>**

**Left Fixed Mount Lateral Trunk Support<sup>(1)</sup>**

Select Left Trunk Support Size

- LATXL Size XS (3.25" x 4.25") . . . . . \$100
- LATSL Size S (3.5" x 5.25") . . . . . \$100
- LATML Size M (3.75" x 5.75") . . . . . \$100
- LATLL Size L (4.25" x 6.5") . . . . . \$100

**Right Fixed Mount Lateral Trunk Support<sup>(1)</sup>**

Select Right Trunk Support Size

- LATXR Size XS (3.25" x 4.25") . . . . . \$100
- LATSR Size S (3.5" x 5.25") . . . . . \$100
- LATMR Size M (3.75" x 5.75") . . . . . \$100
- LATLR Size L (4.25" x 6.5") . . . . . \$100

1. These laterals are not available with Matrix PB, PB Elite Deep and PB Deep Backs.

**MAXX STYLE SWING-AWAY TRUNK SUPPORT HARDWARE<sup>(1)</sup>**

- SALATL Left MAXX Style Swing Away HDW . . . . . \$195
- Select Left Trunk Support Size
- LATXL Size XS (3.25" x 4.25") . . . . . \$100
  - LATSL Size S (3.5" x 5.25") . . . . . \$100
  - LATML Size M (3.75" x 5.75") . . . . . \$100
  - LATLL Size L (4.25" x 6.5") . . . . . \$100

- SALATR Right MAXX Style Swing Away HDW<sup>(1)</sup> . . . \$195
- Select Right Trunk Support Size
- LATXR Size XS (3.25" x 4.25") . . . . . \$100
  - LATSR Size S (3.5" x 5.25") . . . . . \$100
  - LATMR Size M (3.75" x 5.75") . . . . . \$100
  - LATLR Size L (4.25" x 6.5") . . . . . \$100

- ATLL Additional Telescoping Link - Left<sup>(2)</sup> . . . \$75
- ATLR Additional Telescoping Link - Right<sup>(2)</sup> . . \$75

1. These laterals are not available with Matrix PB, PB Elite Deep and PB Deep Backs  
2. Up to 2.75" medial offset

**FIXED HIP/KNEE SUPPORTS**

**Lateral Hip Support with Fixed Mounting Hardware Left<sup>(1)</sup>**

Select Pad Size and mounting position at Hip:

- HSPXLH Size XS (5.5" W x 4" H) mounted at Hip . \$120
- HSPSMLH Size S/M (7" W x 4" H) mounted at Hip . \$120
- HSPSLH Size S (7"W x 6"H) mounted at Hip . . . . \$120
- HSPMLH Size M (9" W x 4"H) mounted at Hip . . . \$120
- HSPLLH Size L (13" W x 4"H) mounted at Hip . . . \$120

Select Pad Size and mounting position at Knee:

- HSPXLK Size XS (5.5"W x 4" H) mounted at Knee . . . \$120
- HSPSMLK Size S/M (7"W x 4" H) mounted at Knee . . . \$120
- HSPSLK Size S (7"W x 6" H) mounted at Knee . . . . . \$120
- HSPMLK Size M (9"W x 4"H) mounted at Knee . . . . . \$120
- HSPLLK Size L (13"W x 4"H) mounted at Knee . . . . . \$120

**FIXED HIP/KNEE SUPPORTS - CONTINUED**

**Lateral Hip Support with Fixed Mounting Hardware Right<sup>(1)</sup>**

Select Pad Size and mounting position:

- HSPXRH Size XS (5.5"W x 4"H) mounted at Hip . . . . \$120
- HSPSMRH Size S/M (7"W x 4"H) mounted at Hip . . . . \$120
- HSPSRH Size S (7"W x 6"H) mounted at Hip . . . . . \$120
- HSPMRH Size M (9"W x 4"H) mounted at Hip . . . . . \$120
- HSPLRH Size L (13"W x 4"H) mounted at Hip . . . . . \$120

Select Pad Size and mounting position at Knee:

- HSPXRK Size XS (5.5"W x 4"H) mounted at Knee . . \$120
- HSPSMRK Size S/M (7"W x 4"H) mounted at Knee . . \$120
- HSPSRK Size S (7"W x 6"H) mounted at Knee . . . . . \$120
- HSPMRK Size M (9"W x 4"H) mounted at Knee . . . . . \$120
- HSPLRK Size L (13"W x 4"H) mounted at Knee . . . . . \$120

1. Must choose swing away hardware below if required.

**REMOVABLE HIP/KNEE SUPPORT HARDWARE**

**Lift Off Removable Hardware for Hip Support**

- LOHL Left . . . . . \$195
- LOHR Right . . . . . \$195

**MAXX Style Quick Release, Removable, Multi-Axis Mounting Hardware for Hip Supports**

Select mounting position:

- MHSHLH Left mounted at Hip . . . . . \$215
- MHSHLK Left mounted at Knee . . . . . \$215
- MHSHRH Right mounted at Hip . . . . . \$215
- MHSHRK Right mounted at Knee . . . . . \$215

**Swing Away Removable Hardware for Hip Supports**

Select mounting position:

- SHSHLH Left mounted at Hip . . . . . \$215
- SHSHLK Left mounted at Knee . . . . . \$215
- SHSHRH Right mounted at Hip . . . . . \$215
- SHSHRK Right mounted at Knee . . . . . \$215

**Extended Maxx Style Quick Release, Removable, Multi-Axis Mounting Hardware for Hip Supports<sup>(1)</sup>**

Select mounting position:

- EMQHLH Left mounted at Hip . . . . . \$215
- EMQHLLK Left mounted at Knee . . . . . \$215
- EMQHRH Right mounted at Hip . . . . . \$215
- EMQHRK Right mounted at Knee . . . . . \$215

1. Recommended when finished seat cushion thickness will be greater than 3.5" high.

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**ELBOW BLOCKS**

**Left Elbow Block<sup>(1)</sup>**

Select Pad Size and Both Extensions:

- EBXL Size XS (3.25" x 4.25") . . . . . \$195
- EBSL Size S (3.5" x 5.25") . . . . . \$195
- EBML Size M (3.75" x 5.75") . . . . . \$195
- EBLL Size L (4.25" x 6.5") . . . . . \$195
  
- V3 Left Verticle Extension 3"
- V5 Left Verticle Extension 5"
- H5 Left Horizontal Extension 5"
- H7 Left Horizontal Extension 7"

**Right Elbow Block<sup>(1)</sup>**

Select Pad Size and Both Extensions:

- EBXL Size XS (3.25" x 4.25") . . . . . \$195
- EBSL Size S (3.5" x 5.25") . . . . . \$195
- EBML Size M (3.75" x 5.75") . . . . . \$195
- EBLL Size L (4.25" x 6.5") . . . . . \$195
  
- V3 Right Verticle Extension 3"
- V5 Right Verticle Extension 5"
- H5 Right Horizontal Extension 5"
- H7 Right Horizontal Extension 7"

1. Prevents arm from slipping.

**VENT TRAYS<sup>(1)(2)(3)</sup>**

- ARTVTH Vent Tray - Articulating Horizontal<sup>(4)</sup> . . \$1,350  
(14"W x 15"D x 10"H)
- ARTVTV Vent Tray - Articulating Vertical<sup>(4)</sup> . . . . \$1,350  
(14.625"W x 10.125"D x 13"H)
- TRILVT Vent Tray - Trilogy 100/200 . . . . . \$1,350  
(14"W x 6"D)

1. The user weight limit is 225 lbs.
2. The maximum weight capacity of the vent tray is 35 lbs.
3. Vents are not standard configuration on Elevating or Precline Systems or with Single Post Backs.
4. Required for use on Recline Systems

**COMPACT-LAPTOP VENT TRAYS<sup>(1)(2)(3)</sup>**

- CVT Vent Tray . . . . . \$1,350
- VBHA Vent Bag Holder Articulating<sup>(4)(5)</sup> . . . . . \$990
- VBHNA Vent Bag Holder Non-Articulating . . . . . \$990

1. The user weight limit is 225 lbs.
2. The maximum weight capacity of the vent tray is 20 lbs.
3. Vents are not standard configuration on Elevating or Precline Systems or with Single Post Backs.
4. Required for use on Recline Systems.
5. Motion Concepts does not supply the Vent Bag

**SUBMIT**

**TRAY OPTION<sup>(1)(2)(3)(4)</sup>**

- USAT Universal Style Articulating Tray . . . . . \$1,350  
(14"W x Adjustable 6"-10"D)

1. The user weight limit is 225 lbs.
2. The maximum weight capacity of the vent tray is 20 lbs.
3. Trays are not standard configuration on Elevating or Precline Systems or with Single Post Backs.
4. Required for use on Recline Systems.

**ADDITIONAL BATTERY TRAYS**

- R350 M50 Battery Tray for 3 M50 Batteries<sup>(1)</sup> . . \$450
- R234150 Battery Tray for 3 G34 Batteries<sup>(2)</sup> . . . . . \$450

1. Only available on the Rovi X3 Base.
2. Only available on the Rovi A3 Base.

**O2 HOLDERS<sup>(1)</sup>**

- NAO2H Non-Articulating O<sup>2</sup> Gas Holder . . . . . \$395
- AOHG Articulating O<sup>2</sup> Gas Holder  
Cylinder Style<sup>(2)(3)</sup> . . . . . \$1,250
- AOHL Articulating Liquid O<sup>2</sup> Holder  
(6.5" x 4.5")<sup>(2)(3)</sup> . . . . . \$1,250
- NAIPL Non-Articulating IV Pole User's Left . . . . \$395
- NAIPR Non-Articulating IV Pole User's Right . . . \$395
- AIVP Articulating IV Pole<sup>(2)(3)</sup> . . . . . \$1,250

1. The user weight limit is 225 lbs. The maximum weight capacity of the O2 holder is 20 lbs.
2. Required for use on Recline Systems
3. Articulating O2 Holders are not standard configuration on Elevating or Precline Systems or with Single Post Backs.

**ADDITIONAL OPTIONS & ACCESSORIES**

- UCH Universal Cup Holder<sup>(1)</sup> . . . . . \$50
- THS Transfer Handles 4" - 6.5" . . . . . \$350
- THL Transfer Handles 9" - 11.5" . . . . . \$350
- RPFS Recline Push Handles Straight<sup>(1)</sup> . . . . . \$350
- RPHA Recline Push Handles Angled<sup>(1)</sup> . . . . . \$350
- FFLTL Flip and Fold Away 1/2 Lap Tray  
Left Mounting(14"W x 10"D)<sup>(2)(3)</sup> . . . . . \$795
- FFLTR Flip and Fold Away 1/2 Lap Tray  
Right Mounting(14"W x 10"D)<sup>(2)(3)</sup> . . . . . \$795

1. 4" of total vertical mounting adjustment
2. Not available on joystick side and only available with a flat arm pad or modular arm pad.
3. The maximum weight capacity of the tray is 10 lbs.

PRINT FORM

SAVE

CLEAR FORM