



Thank you for using Motion Concepts eForms.

Motion Concepts is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.



Adobe Acrobat Reader DC

* Interactive functions of our new forms work best with the latest version of Adobe Acrobat Reader DC visit <https://get.adobe.com/reader/> to download and install on your PC or Mac or visit Google Play/ iTunes to download the Adobe Acrobat Reader DC app for your device.

Save

Adobe Acrobat Reader DC allows you to save this form with your content - to complete later or use as a starting point for your next form. Please note that content must be added and saved in Acrobat - saving content from completed forms in the browser may not be possible.

SAVE

Submit

Adobe Acrobat Reader DC allows you to submit this form electronically via your email client.

Simply click the submit button below and step through the simple process.

SUBMIT

CLEAR FORM

Print

If you do not have access to Adobe Acrobat Reader DC simply print this form and complete it by hand and fax it to our Customer Service Department at:

888-433-6818

PRINT FORM



ULTRA LOW MAXX CG POWER POSITIONING SYSTEM ON TDX-SP2 BASE

CDN PRICE LIST AND ORDER FORM - Price Effective February 1, 2026



PRINT FORM

SAVE

CLEAR FORM

Motion Concepts Customer Service: 1-866-748-7943 quotescdn@motionconcepts.com

REQUIRED INFORMATION

Quote Order

Date of Order: _____ Dealer Account #: _____

Dealer Name: _____

PO #: _____ Tag _____

Purchasing Contact: _____

Phone: _____ Fax: _____

E-mail: _____

RTS/Therapist: _____

Ship to Address: _____

City: _____ Province: _____


Postal Code: _____ Client Gender: M F

Special Client Conditions: _____

REQUIRED MEASUREMENTS

To ensure system is accurately configured please fill in all REQUIRED* measurements

Client Height*	
Client Weight (lbs)*	
A. Seat to Shoulder	
B. Trunk Depth	
C. Chest Width	
D. Knee to Back	
E. Seat to Top of Head	
F. Back to Joystick	
G. Seat to Elbow	
H. Hip Width	
I. Knee to Heel*	
Cushion Thickness*	



POWER BASE OPTIONS

TDXSP2-CG TDX-SP2 Base for Single Actuator System . . . \$7,595

BASE WIDTH

NARROWFRM Narrow Base 24" (MK 22 NF Batteries)N/C

STDFRM Standard Base 25.5" (MK 24 NF Batteries)N/C

BATTERIES

22NFBATTERY 22 NF BATTERY (ON CHAIR) PAIR \$700

24BATTERY 24 GEL BATTERY (ON CHAIR) PAIR \$840

SEAT-TO-FLOOR HEIGHT⁽¹⁾

SLOW 16.75" Seat-To-Floor⁽²⁾ \$250

LOW 17.25" Seat-To-FloorN/C

MED 18.25" Seat-To-FloorN/C

TALL 19.25" Seat-To-FloorN/C

1. Seat-to-floor height is measured at the front of the seat pan based on a standard seat size. Seat-to-floor height may vary +/- .25" dependent on condition of the power base and seat depths greater than 18"

2. This seat to floor height is only available with the Group 24 Battery Tray (24TRY) and on 17" seat depths and greater.

USER WEIGHT LIMIT

U300 Weight Capacity up to 300lbs STD

WHEELCHAIR TRANSPORT BRACKETS

TRBKTS Wheelchair Transport Brackets⁽¹⁾ STD

TRRO Transport Ready Option⁽²⁾⁽⁴⁾ \$499

TRRO-E Transport Ready Option for Elevate⁽²⁾⁽³⁾⁽⁴⁾ \$499

1. For unoccupied use only

2. TRRO meets ISO-7176 Part 19

3. Required for Elevate Systems to meet ISO-7176 Part 19, includes two additional tie-down loops on rail.

4. Must select a Motion Concepts or Matrix Elan headrest with hardware and Motion Concepts Standard Rehab, High Back or Matrix Upgradeable back. Comes with pelvic safety belt standard, please select TROBELT on page 5.

TIRE OPTIONS

B1431-3 14" x 3" Black Tire with Gel Foam Inserts STD

B1430-3 14" x 3" Black Tire - Pneumatic TireN/C

FORK OPTIONS

DSFK Double-Sided Fork Package STD

BATTERY CHARGER OPTION

I10CHARGER 110 Volt Battery ChargerN/C

L8900 Omit Battery ChargerN/C

FRAME FINISH

170P Obsidian BlackN/C

171P Ruby RedN/C

172P Peacock BlueN/C

174P Purple PassionN/C

175P Denim NavyN/C

176P Lime TwistN/C

177P Arctic BlueN/C

RIM INSERT COLORS - MUST PICK ONE

170PR Obsidian BlackN/C

171PR Ruby RedN/C

172PR Peacock BlueN/C

173PR Pink FlamingoN/C

174PR Purple PassionN/C

175PR Denim NavyN/C

176PR Lime TwistN/C

177PR Arctic BlueN/C

178PR Banana SplitN/C

179PR Orange MarmaladeN/C

180PR Morning MistN/C

Ultra Low Maxx CG Power Positioning System TDX-SP2

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

HIGHSPEED MOTORS AND ENHANCED DRIVING

- HSMTR Performance Motors 7.5 mph\$995
- GTRAC-LX LINX G-TRAC Module⁽¹⁾ \$1,200

1. Recommended for use with Specialty Drive Controls

LIGHTS AND INDICATORS

- LIGHTS LED Lights and Indicators⁽¹⁾\$995

1. Must order REM215, REM216, REM400 or REM500

NON-EXPANDABLE SINGLE/DUAL JOYSTICKS

- REM110 LED Drive Only Remote/Joystick⁽¹⁾N/C
Single drive profile
- REM210 LED Non-Expandable Remote/Joystick⁽¹⁾N/C
Three drive profiles, control of two seat functions

1. Not available with lights and indicators.

EXPANDABLE JOYSTICKS⁽¹⁾

- REM211 LED Expandable Remote/Joystick⁽²⁾\$875
Three drive profiles, control of up to 6 actuators
Not for use with specialty controls
- REM216 LED Expandable Remote/Joystick
with Lights⁽³⁾\$875
Three drive profiles, control of up to 6 actuators
and light operation. Not for use with specialty controls
- REM400 Color 3.5" Touch Screen
Remote/Joystick⁽⁴⁾\$1,100
Four drive profiles, control of up to 6 actuators through the
drive control and light operation.
- SW400TGL Toggle Kit for REM400⁽⁵⁾ \$395

1. Must choose EXPC Expandable Controller and PWH Harness.
2. Not available with lights and indicators
3. Must choose lighting package separately.
4. REM400 or REM500 is required when using specialty controls, choose specialty controls purchased separately.
5. Adds 2 toggle switches to REM400.

EXPANDABLE SPECIALTY DISPLAY

- REM500 LINX Color 3.5" Touch Screen
Display⁽¹⁾⁽²⁾⁽³⁾\$1,100

Select mounting position for REM500

- R5MR RightN/C
- R5ML LeftN/C
- R5RPB Rear Push Bar\$140

1. Must choose EXPC Expandable Controller and PWH Harness.
2. Driver Control not included.
3. REM400 or REM500 is required when using specialty controls, choose specialty controls purchased separately.

SPECIALTY PROPORTIONAL CONTROLS⁽¹⁾⁽²⁾

- CREM LINX Compact Remote HCPCS code \$1,097
- CREM-LF LINX Compact Remote Low Force \$1,097

1. Requires REM400 or REM500
2. Select Mount in Mounting For Driver Control section

ATTENDANT DRIVE CONTROL CHOICES

- ACU Proportional Attendant Control⁽¹⁾ \$1,024

Select mounting position for ACU

- RPBMR Right
- RPBML Left
- RPBMC Center

JOYSTICK MOUNT

Arm Mounted

- FR Fixed Joystick Mounting - Right STD
- FL Fixed Joystick Mounting - Left STD
- MRR Maxx Resolve Swing-Away Joystick Mount
(Right Hand) \$450
- MRL Maxx Resolve Swing-Away Joystick Mount
(Left Hand) \$450
- SAQR Motion Height Adjustable Swing-Away
Quad Link (Right Hand) \$244
- SAQL Motion Height Adjustable Swing-Away
Quad Link (Left Hand) \$244

Rear Push Bar Mounted

- RPBMR Right \$190
- RPBML Left \$190
- RPBMC Center \$190

JOYSTICK KNOB - OPTIONAL

- PC101A BodyPoint 3" U-Shaped \$120
- PC102A BodyPoint 4" U-Shaped \$120
- PC104A BodyPoint Straight Handle \$120
- PC107A BodyPoint Domed Rubber \$120
- PC110A BodyPoint Grooved Mushroom \$120

1. To remove the joystick knob from the REM400, pull straight up on the joystick (DO NOT TWIST) otherwise damage may occur and may void the warranty.

ELECTRONICS - EXPANDABLE MODULES

- OUT500 Output Module \$645
- IN200 LINX Switch Input Module with
Power ON/OFF for REM200 Series⁽¹⁾⁽²⁾ \$645
- IN500 LINX 9 Pin Input Module with
Sip-N-Puff Nozzle⁽²⁾⁽³⁾⁽⁴⁾⁽⁶⁾ \$1,800
- PKG32666 Therafin Sip-N-Puff Breath Tube Kit \$440
HCPCS code E2326
- COVER-L Large Recline Shroud Cover⁽⁵⁾ N/C

1. Not for use with Alternative Drive Controls. (Uses switch ports).
2. Only one input module per chair.
3. Input module also required for any 9-pin specialty control.
4. 9 pin specialty control and Sip-N-Puff cannot both be used together.
5. May be required when choosing OUT500, IN200 or IN500 options.
6. Requires REM400 or REM500 for compatibility.

LINX ELECTRONIC ACCESSORIES

- LAK LINX Access Key⁽¹⁾ \$50

1. Required for programming

Ultra Low Maxx CG Power Positioning System
TDX-SP2

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

POWER POSITIONING SYSTEMS

- ULMT Ultra Low Maxx CG Tilt (50°) Only System⁽¹⁾
.....\$6,969
- ULMTRE Ultra Low Maxx CG Tilt (50°) & Power Recline
(170°) w/Extended Shear Reduction(ESR)⁽¹⁾
.....\$11,569
- ULMTRP Ultra Low Maxx CG Tilt (50°) & Power Recline
(170°) w/Power Sliding Back(PSB)⁽¹⁾⁽²⁾
.....\$12,069
- ULMRE Ultra Low Maxx Power Recline (170°)
w/Extended Shear Reduction(ESR)
.....\$10,350
- ULMRP Ultra Low Maxx Power Recline (170°)
w/Power Sliding Back(PSB)⁽²⁾
.....\$10,850
- ULME Ultra Low Maxx 12" Elevating
Lift Seat Only⁽³⁾\$4,995
Weight capacity up to 250lbs
- ULME10 Ultra Low Maxx 10" Elevating
Lift Seat Only⁽⁴⁾\$5,995
Weight capacity up to 300lbs

1. Tilt will be reduced to 45° if adding an Elevating Lift Seat Module
2. Requires Multi-Function M610 through Drive Control (MFCBD).
3. Not available on seat sizes greater than 20" wide by 20" deep.
4. Not available on seat sizes greater than 22" wide by 22" deep.

ADD: POWER ADJUSTABLE SEAT HEIGHT MODULE

- ESM12 12" Power Adjustable Seat Height Module⁽¹⁾⁽²⁾⁽³⁾
.....\$2,495
Weight capacity up to 250lbs
- ESM10 10" HD Power Adjustable Seat Height Module⁽⁴⁾⁽⁵⁾
.....\$2,995
Weight capacity up to 300lbs

1. Uses an integrated Tilt/Lift mechanism, tilt is 45° and the seat to floor height is not affected by the Elevating Lift Seat.
2. Not available on seat sizes greater than 20" wide by 20" deep
3. For weights greater than 175lbs tilt while elevated is not available.
4. Only available with tilt, uses an integrated Tilt/Lift mechanism, tilt is 45°.
5. Not available on seat sizes greater than 22" wide by 22" deep. Only available on Standard Base (STDFRM).

MODULES - MUST SELECT SYSTEM TO ADD MODULES

- PRETILT 5° Fixed Pre-Tilt⁽¹⁾N/C
- PRM8 8° Precline Module⁽²⁾N/C
- PRM30 30° Precline Module⁽²⁾\$500
- ERSSP Programable Electronic Recline Function
Start & Stop⁽³⁾\$400

Specify: Start: _____degrees Stop: _____degrees

- MRSTD 30° Manual Recline⁽⁴⁾ For Users under 200lbs **NEW!**
.....\$590
- MRHD 30° Manual Recline⁽⁴⁾⁽⁵⁾ For Users 200-250lbs
.....\$590

1. This Fixed Pre-tilt option adds approximately 2" to the front seat to floor height.
2. Total degrees of recline is 170° ie. if 30° of precline is used, then 140° of recline will be available. Must order Recline System
3. Requires Multi-Function M610 through Drive Control (MFCBD). The range is dependent on the Recline module selected. Will be set to factory settings unless specified.
4. For Attendant use. Not available with Recline Arms/ Cane Mounted Cantilever arms, or Vent Trays. Uses Dual-Post Tilt Arms or RMCAN Arms. Comes with Standard Back Canes with adustable back angle range from 95°-125°.
5. Not available with Elevate Module

ADULT SEAT WIDTH SETTINGS

- SW116 16" Wide (Adjustable from 16" to 20")N/C
- SW117 17" Wide (Adjustable from 16" to 20")N/C
- SW118 18" Wide (Adjustable from 16" to 20")N/C
- SW119 19" Wide (Adjustable from 16" to 20")N/C
- SW120 20" Wide (Adjustable from 16" to 20")N/C
- SW219 19" Wide (Adjustable from 19" to 22")\$460
- SW220 20" Wide (Adjustable from 19" to 22")\$460
- SW221 21" Wide (Adjustable from 19" to 22")\$460
- SW222 22" Wide (Adjustable from 19" to 22")\$460

ADULT ULTRA RAIL SEAT DEPTH SETTINGS

- SD116 16" Deep (Adjustable from 16" to 21")N/C
- SD117 17" Deep (Adjustable from 16" to 21")N/C
- SD118 18" Deep (Adjustable from 16" to 21")N/C
- SD119 19" Deep (Adjustable from 16" to 21")N/C
- SD120 20" Deep (Adjustable from 16" to 21")N/C
- SD121 21" Deep (Adjustable from 16" to 21")N/C
- SD219 19" Deep (Adjustable from 19" to 23")\$460
- SD221 21" Deep (Adjustable from 19" to 23")\$460
- SD222 22" Deep (Adjustable from 19" to 23")\$460
- SD223 23" Deep (Adjustable from 19" to 23")\$460

SMALL ADULT SEAT WIDTH SETTINGS⁽¹⁾

- SW314 14" Wide (Adjustable from 14" to 15")\$250
 - SW315 15" Wide (Adjustable from 14" to 15")\$250
1. Weight limit limited to 175 lbs. If recline is selected it may be limited to 165° of recline. You must select both small adult seat width and depth together.

SMALL ADULT SEAT DEPTH SETTINGS⁽¹⁾

- SD314 14" Deep (Adjustable from 14" to 18")\$250
 - SD315 15" Deep (Adjustable from 14" to 18")\$250
 - SD316 16" Deep (Adjustable from 14" to 18")\$250
 - SD317 17" Deep (Adjustable from 14" to 18")\$250
 - SD318 18" Deep (Adjustable from 14" to 18")\$250
1. Weight limit limited to 175 lbs. If recline is selected it may be limited to 165° of recline. You must select both small adult seat width and depth together.

Ultra Low Maxx CG Power Positioning System
TDX-SP2

NOTE: Systems Include Seat Pan, Back Pan, and Standard Armrests. They are designed for use with a 2" - 3" cushion

MATRIX SEAT CUSHION OPTIONS

Matrix PS Seat Cushion 14" - 20" wide/reversable cover \$569
 PS1414 PS1416 PS1616 PS1618 PS1620
 PS1716 PS1718 PS1720 PS1816 PS1818
 PS1820 PS2016 PS2018 PS2020

Matrix PS HD Seat Cushion 21"-22" wide \$975
 PS2218 PS2220 PS2222

Matrix PSP Seat Cushion 14" - 20" wide/reversable cover \$645
 PSP1414 PSP1416 PSP1616 PSP1618 PSP1620
 PSP1716 PSP1718 PSP1720 PSP1816 PSP1818
 PSP1820 PSP2016 PSP2018 PSP2020

Matrix PSP HD Seat Cushion 21"-22" wide \$995
 PSP2218 PSP2220 PSP2222

Matrix Vi Seat Cushion 14" - 20" wide \$624
 KA1414-VI MA1415-VI MA1417-VI MA1616-VI MA1618-VI
 MA1620-VI MA1716-VI MA1718-VI MA1720-VI MA1816-VI
 MA1818-VI MA1820-VI MA2016-VI MA2018-VI
 MA2020-VI

Matrix Vi Seat Cushion 21" - 22" wide \$681
 MA2218-VI MA2220-VI

Matrix Vi HD Seat Cushion 22" wide \$840
 MA2222-VI

Matrix Libra Seat Cushion 14" - 20" wide \$745
 LC1414 LC1416 LC1516 LC1518 LC1520
 LC1616 LC1618 LC1620 LC1716 LC1718
 LC1720 LC1816 LC1818 LC1820 LC1916
 LC1918 LC1920 LC2016 LC2018 LC2020

Matrix Libra HD Seat Cushion 21" wide \$1,025
 LC2118 LC2120

Matrix Libra Seat Cushion 22" deep \$1,025
 LC1622 LC1722 LC1822 LC1922 LC2022
 LC2122

Matrix Libra HD Seat Cushion 22" wide \$1,025
 LC2218 LC2220 LC2222

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

Matrix Libra Adjust Seat Cushion 14" -20" wide \$985
 LC1414-ADJ LC1416-ADJ LC1516-ADJ LC1518-ADJ
 LC1520-ADJ LC1616-ADJ LC1618-ADJ LC1620-ADJ
 LC1716-ADJ LC1718-ADJ LC1720-ADJ LC1816-ADJ
 LC1818-ADJ LC1820-ADJ LC1916-ADJ LC1918-ADJ
 LC1920-ADJ LC2016-ADJ LC2018-ADJ LC2020-ADJ

Matrix Libra Adjust Seat Cushion 22" deep \$1,265
 LC1622-ADJ LC1722-ADJ LC1822-ADJ LC1922-ADJ
 LC2022-ADJ LC2122-ADJ

Matrix Libra Adjust HD Seat Cushion 22" wide \$1,265
 LC2118-ADJ LC2120-ADJ

Matrix Multi Cushion 14" - 20" wide w/reversable cover \$745
 MM1414 MM1416 MM1616 MM1618 MM1620
 MM1716 MM1718 MM1720 MM1816 MM1818
 MM1820 MM2016 MM2018 MM2020

Matrix Multi Cushion 22" wide w/reversable cover \$995
 MM2216 MM2218 MM2220 MM2222

FABRIC OPTIONS FOR CUSHIONS

Cushion	Meshtex	Startex	Infection Control	Spacetex	Onyx	Fusion
Libra/Adjust	N/A	STD	\$50	N/A	N/A	N/A
PS/PSP	N/A	STD Reversible	STD Reversible	N/A	N/A	N/A
Vi	STD	\$50	\$50	N/A	\$100	\$150

LC-SP Spacetex Fabric Option for Libra Cushion N/C
 LC-IC Infection Control Fabric Upgrade for Libra Cushion \$50
 FSU-MA Fusion Fabric Upgrade for Matrix Vi Cushion... \$150
 OCU-MA Onyx Fabric Upgrade for Matrix Vi Cushion \$100
 STARTEX-S Smooth side Startex Fabric Upgrade for Matrix Vi Cushion \$50
 STARTEX-F Fabric side Startex Fabric Upgrade for Matrix Vi Cushion \$50

LAP BELTS

BELT48 Push Button Style Lap Belt, Length 48"⁽¹⁾ STD
 BELT60 Push Button Style Lap Belt, Length 60"⁽¹⁾ STD
 BELT71 Push Button Style Lap Belt, Length 71"⁽¹⁾ STD
 PLB2 Bodypoint Padded Lap Belt 2 Point⁽²⁾ \$250
 PLB4 Bodypoint Padded Lap Belt 4 Point⁽²⁾ \$315
TROBELT Safety Belt for use only with TRRO & TRRO-E... N/C

1. Not available with TRRO and TRRO-E.
2. 67" length only option

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

MATRIX HEADREST PADS⁽¹⁾

Standard Pads

<input type="checkbox"/>	ESP6	Matrx Elan Standard Pad 6"	\$310
<input type="checkbox"/>	ESP10	Matrx Elan Standard Pad 10"	\$310
<input type="checkbox"/>	ESP14	Matrx Elan Standard Pad 14"	\$310
<input type="checkbox"/>	ESP6-HUG	Matrx Elan Hug Standard Pad 6"	\$310
<input type="checkbox"/>	ESP10-HUG	Matrx Elan Hug Standard Pad 10"	\$310
<input type="checkbox"/>	ESP14-HUG	Matrx Elan Hug Standard Pad 14"	\$310
<input type="checkbox"/>	SNUG14	Matrx SNUG Standard Pad 14"	\$625
<input type="checkbox"/>	SNUG14-A	Matrx SNUG Asymmetrical Pad 14"	\$645
<input type="checkbox"/>	SNUG18	Matrx SNUG Standard Pad 18"	\$625
<input type="checkbox"/>	SNUG18-A	Matrx SNUG Asymmetrical Pad 18"	\$645

Standard Pads with Infection Control Cover

<input type="checkbox"/>	ESP6-IC	Matrx Elan Infection Control Pad 6"	\$335
<input type="checkbox"/>	ESP10-IC	Matrx Elan Infection Control Pad 10"	\$335
<input type="checkbox"/>	ESP14-IC	Matrx Elan Infection Control Pad 14"	\$335
<input type="checkbox"/>	ESP6-IC-HUG	Matrx Elan Hug Infection Control Pad 6"	\$335
<input type="checkbox"/>	ESP10-IC-HUG	Matrx Elan Hug Infection Control Pad 10"	\$335
<input type="checkbox"/>	ESP14-IC-HUG	Matrx Elan Hug Infection Control Pad 14"	\$335
<input type="checkbox"/>	SNUG14-IC	Matrx SNUG Standard Infection Control Pad 14"	\$665
<input type="checkbox"/>	SNUG14-A-IC	Matrx SNUG Asymmetrical Infection Control Pad 14"	\$685
<input type="checkbox"/>	SNUG18-IC	Matrx SNUG Standard Infection Control Pad 18"	\$665
<input type="checkbox"/>	SNUG18-A-IC	Matrx SNUG Asymmetrical Infection Control Pad 18"	\$685

Occipital /4-Point Pads

<input type="checkbox"/>	EOP9	Matrx Elan Occipital Pad 9"	\$341
<input type="checkbox"/>	EOP12	Matrx Elan Occipital Pad 12"	\$341
<input type="checkbox"/>	E4PTSM	Matrx Elan 4-Point Pad 9"W x 7"H	\$448
<input type="checkbox"/>	E4POINT	Matrx Elan 4-Point Pad 11"W x 10"H	\$448

Occipital / 4-Point Pads with Infection Control Cover

<input type="checkbox"/>	EOP9-IC	Matrx Elan Occipital Pad 9" with Infection Control	\$366
<input type="checkbox"/>	EOP12-IC	Matrx Elan Occipital Pad 12" with Infection Control	\$366
<input type="checkbox"/>	E4PTSM-IC	Matrx Elan 4-Point Pad 9"W x 7"H with Infection Control	\$473
<input type="checkbox"/>	E4POINT-IC	Matrx Elan 4-Point Pad 11"W x 10"H with Infection Control	\$473

1. If mounting on PB Elite must order removable hardware

MOTION CONCEPTS HEADREST PAD

<input type="checkbox"/>	MCSH	Motion Concepts Standard	N/C
<input type="checkbox"/>	MCOH	Motion Concepts Auto Style ⁽¹⁾	\$235

1. Comes Standard with Fixed Mounting hardware; for use on Matrx PB or Elite Backs you must select the removable hardware upgrade below.

HEADREST MOUNTING HARDWARE

<input type="checkbox"/>	FM	Fixed Mounting	N/C
<input type="checkbox"/>	LXSHW	Loxx Standard Hardware 15"	\$275
<input type="checkbox"/>	LXMMH	Loxx Mini Hardware 11"	\$275
<input type="checkbox"/>	MEHWRE	Elan Hardware 15"	\$285
<input type="checkbox"/>	MEMH	Elan Mini Hardware 11"	\$285
<input type="checkbox"/>	MEHW-EXT	Elan Extended Hardware	\$320
<input type="checkbox"/>	FPLXSHW	FLIPP Loxx Standard Hardware 15"	\$705
<input type="checkbox"/>	FPLXSHW-EX	FLIPP Loxx Extended Hardware	\$705
<input type="checkbox"/>	FPLXMMH	FLIPP Loxx Mini Hardware 11"	\$705
<input type="checkbox"/>	FPLXMMH-EX	FLIPP Loxx Mini Extended Hardware	\$705
<input type="checkbox"/>	FPMEHW	FLIPP Elan Hardware 15"	\$715
<input type="checkbox"/>	FPMEMH	FLIPP Elan Mini Hardware 11"	\$715
<input type="checkbox"/>	FPMEHW-EXT	FLIPP Elan Extended Hardware	\$780
<input type="checkbox"/>	OHR	OMIT: Headrest	N/C

STANDARD BACK CANES - FOR TILT SYSTEMS

<input type="checkbox"/>	SBC18	18" Straight Back Canes	STD
<input type="checkbox"/>	SBC20	20" Straight Back Canes	STD
<input type="checkbox"/>	SBC22	22" Straight Back Canes	STD
<input type="checkbox"/>	SBC24	24" Straight Back Canes ⁽⁷⁾	STD

Select Cane Angle⁽⁷⁾

<input type="checkbox"/>	BA75 75°	<input type="checkbox"/>	BA79 79°	<input type="checkbox"/>	BA85 85°	<input type="checkbox"/>	BA90 90°	<input type="checkbox"/>	BA95 95°
<input type="checkbox"/>	BA101 101°	<input type="checkbox"/>	BA106 106°	<input type="checkbox"/>	BA112 112°	<input type="checkbox"/>	BA116 116°		

<input type="checkbox"/>	ABC20	20" Angled (5°) Back Canes ⁽⁷⁾	N/C
<input type="checkbox"/>	ABC22	22" Angled (5°) Back Canes ⁽⁷⁾	N/C
<input type="checkbox"/>	ABC24	24" Angled (5°) Back Canes ⁽⁷⁾	N/C

Select Cane Angle⁽⁷⁾

<input type="checkbox"/>	BA81 81°	<input type="checkbox"/>	BA84 84°	<input type="checkbox"/>	BA91 91°	<input type="checkbox"/>	BA96 96°	<input type="checkbox"/>	BA101 101°
<input type="checkbox"/>	BA106 106°	<input type="checkbox"/>	BA111 111°	<input type="checkbox"/>	BA117 117°	<input type="checkbox"/>	BA121 121°		

<input type="checkbox"/>	EBC22	22" Endomorph Back Canes ⁽¹⁾⁽⁷⁾	\$250
<input type="checkbox"/>	EBC24	24" Endomorph Back Canes ⁽¹⁾⁽⁷⁾	\$250

Please Select

<input type="checkbox"/>	EBCN1	1" Narrower	N/C
<input type="checkbox"/>	EBCN2	2" Narrower	N/C
<input type="checkbox"/>	EBCN4	4" Narrower	N/C

<input type="checkbox"/>	TBC18	18" 10° Mid-Angle Tilt Back Canes ⁽²⁾⁽⁷⁾	\$250
<input type="checkbox"/>	TBC20	20" 10° Mid-Angle Tilt Back Canes ⁽²⁾⁽⁷⁾	\$250
<input type="checkbox"/>	TBC22	22" 10° Mid-Angle Tilt Back Canes ⁽²⁾⁽⁷⁾	\$250
<input type="checkbox"/>	TBC24	24" 10° Mid-Angle Tilt Back Canes ⁽²⁾⁽⁷⁾	\$250

<input type="checkbox"/>	FDBC	Fold Down Back Canes w/ Pull Pin ⁽⁶⁾⁽⁷⁾	\$250
--------------------------	------	--	-------

1. Comes standard with Cantilever Armrest Cane Mounted.
2. Quick Height Adjustment on Backrest and Maxx Spreader Bars not available.
3. Minimum finished back widths of 16".
4. Handle measures 6.5" high and has .75" height adjustment with angle adjustment.
5. 4" of total vertical mounting adjustment.
6. 250 lbs weight capacity. Must select canes above.
7. Not available with Manual Recline.
6. 250 lbs weight capacity. Must select canes above.

CANE OPTIONS - FOR RECLINE SYSTEMS

Ultra Low Maxx CG Power Positioning System
TDX-SP2

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

- BFS2 Recline Body Fit 2" Narrower⁽¹⁾\$995
- BFS4 Recline Body Fit 4" Narrower⁽¹⁾\$995
- AHH Adjustable Hooking Handles⁽²⁾ (pair) \$450
- RPHS Push Handles for Recline - Straight⁽³⁾ (pair) \$350
- RPHA Push Handles for Recline - Angled⁽³⁾ (pair) \$350

1. Minimum finished back widths of 16".
2. Handle measures 6.5" high and has .75" height adjustment with angle adjustment.
3. 4" of total vertical mounting adjustment.

FINISHED BACK HEIGHT⁽¹⁾⁽²⁾ - FOR TILT ONLY/ELEVATE ONLY

- FBT18 18" FBT19 19" FBT20 20" FBT21 21"
- FBT22 22" FBT23 23" FBT24 24" FBT25 25"

1. 4" Gap from the seat pan to the bottom of back pan is recommended.
2. Finished back height is measured from seat pan to top of backrest.

FINISHED BACK HEIGHT⁽¹⁾⁽²⁾ - FOR RECLINE OR TILT/RECLINE

- FBR20 20" FBR21 21" FBR22 22" FBR23 23"
- FBR24 24" FBR25 25" FBR26 26" FBR27 27"

1. 6" Gap from the seat pan to the bottom of back pan is recommended.
2. Finished back height is measured from seat pan to top of backrest.

STANDARD BACK OPTIONS

- SRB Standard Rehab Back & Cushion N/C
- SS Additional Super Soft/HR Foam⁽¹⁾ \$450
- MFU Matrix Fabric Upcharge \$450
- OBC OMIT: Back Cushion N/C
- PIPF Flat Planar Interface/I-Back Plate⁽²⁾⁽³⁾⁽⁵⁾ N/C
- PIPR Recessed Planar/I-Back Interface Plate⁽²⁾⁽³⁾⁽⁵⁾⁽⁶⁾ N/C
- ORBC OMIT: Rigid Back & Cushion⁽⁴⁾ N/C

1. Additional foam added may impact seat depth
2. Not available on tilt only or tilt/elevate systems, for use with Recline Module.
3. Required with all Recline Systems when no back is selected. May also be used as an interface for after market backs.
4. Not available with Recline/ESR or Recline/PSB.
5. Not available with TRRO and TRRO-E.
6. Fits Comfort Company, Freedom, Pindot and Moulded Backs.

MOTION CONCEPTS SINGLE POST BACKREST⁽¹⁾⁽²⁾⁽³⁾

- SPBMI Single Post Angle Adjustable Fold Down Back interface-Matrix Back.....\$750 (Folds to 25° forward)
- SPBCI Single Post Angle Adjustable Fold Down Back for Contour/High Back⁽⁴⁾\$750 (Folds to 25° forward)

1. Single Post Back is an angle adjustable back between 73° and 109° in 6° increments from 17° forward to 19° back. The back folds down and is only available for 16"-20" seat widths & 21"-24" back height. A 4" gap from bottom of back to seat pan is required.
2. This option is not available with Recline and/or a vent tray and has a 250lb weight capacity.
3. Cane mounted armrest are not available with this option.
4. Not available with TRRO/TRRO-E

MOTION CONCEPTS REHAB SEAT SLING BACK⁽¹⁾

- SLINGB-TT Sling Back with Tension Fabric \$950
- SLINGB-TA Sling Back with Adjustable Fabric \$1,250

1. Only available on Rehab Seat System. Not available with Elevate. Cane height must be 2" taller than Finished Back Height.

UPGRADABLE BACK OPTIONS

Contour Back 15" wide - 20" wide \$495

- CB1512 CB1515 CB1517 CB1612
- CB1615 CB1617 CB1621 CB1712
- CB1715 CB1717 CB1721 CB1812
- CB1815 CB1817 CB1821 CB1912
- CB1915 CB1917 CB1921 CB2012
- CB2015 CB2017 CB2021

Contour Back 21"-22" wide \$595

- CB2115 CB2117 CB2121 CB2215
- CB2217 CB2221

High Back 16" - 22" W \$595

- HB1616 HB1618 HB1620 HB1716
- HB1718 HB1720 HB1816 HB1818
- HB1820 HB1916 HB1918 HB1920
- HB2016 HB2018 HB2020 HB2116
- HB2118 HB2120 HB2216 HB2218
- HB2220

Matrix Elite Back 14" - 20" W (3" Contour)..... \$766

- PBE1414 PBE1416 PBE1418 PBE1420
- PBE1514 PBE1516 PBE1518 PBE1520
- PBE1614 PBE1616 PBE1618 PBE1620
- PBE1714 PBE1716 PBE1718 PBE1720
- PBE1814 PBE1816 PBE1818 PBE1820
- PBE1916 PBE1918 PBE1920 PBE2014
- PBE2016 PBE2018 PBE2020

Matrix Elite HD Back 21" - 22" W (3" Contour) \$1,095

- PBE2116-HD PBE2120-HD PBE2216-HD PBE2220-HD

Matrix Elite TR Back 15" - 20" W (3" Contour) \$866

- PBE1516-TR PBE1518-TR PBE1520-TR PBE1616-TR
- PBE1618-TR PBE1620-TR PBE1716-TR PBE1718-TR
- PBE1720-TR PBE1816-TR PBE1818-TR PBE1820-TR
- PBE1916-TR PBE1918-TR PBE1920-TR PBE2016-TR
- PBE2018-TR PBE2020-TR

Matrix Elite TR HD Back 21" - 22" W (3" Contour) \$1,151

- PBE2120-TRHD PBE2220-TRHD

Matrix Elite Deep Back 14" - 20" W (6" Contour) \$918

- EDB1416 EDB1418 EDB1420 EDB1514
- EDB1516 EDB1518 EDB1520 EDB1614
- EDB1616 EDB1618 EDB1620 EDB1714
- EDB1716 EDB1718 EDB1720 EDB1814
- EDB1816 EDB1818 EDB1820 EDB1914
- EDB1916 EDB1918 EDB1920 EDB2014
- EDB2016 EDB2018 EDB2020

UPGRADABLE BACK OPTIONS

Ultra Low Maxx CG Power Positioning System
TDX-SP2

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

Matrix Elite Deep HD Back 21" - 22" W (6" Contour) \$1,095

- EDB2116-HD EDB2118-HD EDB2120-HD
- EDB2216-HD EDB2218-HD EDB2220-HD

Matrix Elite E2S Back 14" - 20" W x 14"/16"H (3" Contour) \$1,125

- E2S1414 E2S1416 E2S1514 E2S1516
- E2S1614 E2S1616 E2S1714 E2S1716
- E2S1814 E2S1816 E2S1914 E2S1916
- E2S2014 E2S2016

Matrix Elite E2S Back 14" - 20" W x 18"/20"H (3" Contour) \$995

- E2S1418 E2S1420 E2S1518 E2S1520
- E2S1618 E2S1620 E2S1718 E2S1720
- E2S1818 E2S1820 E2S1918 E2S1920
- E2S2018 E2S2020

Matrix Elite E2 Deep Back 14" - 20" W x 14"/16"H (6" Contour) \$1,125

- E2D1414 E2D1416 E2D1514 E2D1516
- E2D1614 E2D1616 E2D1714 E2D1716
- E2D1814 E2D1816 E2D1914 E2D1916
- E2D2014 E2D2016

Matrix Elite E2 Deep Back 14" - 20" W x 18"/20"H (6" Contour) \$1,155

- E2D1418 E2D1420 E2D1518 E2D1520
- E2D1618 E2D1620 E2D1718 E2D1720
- E2D1818 E2D1820 E2D1918 E2D1920
- E2D2018 E2D2020

Matrix PB Posture Back 14" - 20" W (5" Contour) \$734

- MPB1412 MPB1612 MPB1616 MPB1620
- MPB1712 MPB1716 MPB1720 MPB1812
- MPB1816 MPB1820 MPB2012 MPB2016
- MPB2020

Matrix PB Heavy Duty Back⁽¹⁾ 20" - 22" W (5" Contour) \$1,051

- HDB2016 HDB2020 HDB2216 HDB2220

Matrix PB Posture Back Deep⁽¹⁾ 14" - 20" W (7" Contour) \$892

- PBD1412 PBD1612 PBD1616 PBD1620
- PBD1712 PBD1716 PBD1720 PBD1812
- PBD1816 PBD1820 PBD2012 PBD2016
- PBD2020

Matrix PB Heavy Duty Deep Back⁽¹⁾ 20"-22" W (5" Contour) \$1,209

- DHDB2016 DHDB2020 DHDB2216 DHDB2220

Matrix MAC Back⁽²⁾ \$1,495

- MAC1616 MAC1618 MAC1620 MAC1622
- MAC1716 MAC1718 MAC1720 MAC1722
- MAC1816 MAC1818 MAC1820 MAC1822
- MAC1916 MAC1918 MAC1920 MAC1922
- MAC2016 MAC2018 MAC2020 MAC2022

Matrix MAC Back w/ HUG⁽²⁾ \$1,645

- MAC1616-HUG MAC1618-HUG MAC1620-HUG
- MAC1622-HUG MAC1716-HUG MAC1718-HUG
- MAC1720-HUG MAC1722-HUG MAC1816-HUG
- MAC1818-HUG MAC1820-HUG MAC1822-HUG
- MAC1916-HUG MAC1918-HUG MAC1920-HUG
- MAC1922-HUG MAC2016-HUG MAC2018-HUG
- MAC2020-HUG MAC2022-HUG

Matrix MAC Lateral Contour

- MLCSL Left 3.5" MLCSR Right 3.5"
- MLCDL Left 6" MLCDR Right 6"
- MLCXDL Left 9" MLCXDR Right 9"

1. The Heavy Duty PB Back and PB Deep Back are not available on Recline systems with seat widths of 21" and 22"
2. Must select Pad Contour in Matrix MAC Lateral Contour section.

OTHER UPGRADABLE BACK OPTIONS

Cushion	Meshtex	4-Way Stretch Polyester	Infection Control	Onyx	Fusion
Contour	N/A	STD	N/A	\$100	N/A
High	STD	N/A	N/A	\$100	\$100
E2, Elite, Elite HD, Elite Deep	STD	N/A	\$100	N/A	N/A
E2, Elite TR, Elite TR HD	STD	N/A	\$100	N/A	N/A

OTHER UPGRADABLE BACK OPTIONS

- EBSS Super Soft Foam for Elite Back Upgrade..... \$150
- FFUHB Fusion Fabric Upgrade for High Backs..... \$150
- OCUCB Onyx Fabric Upgrade for Contour Backs..... \$150
- OCUHB Onyx Fabric Upgrade for High Backs..... \$150
- CHIEB Chest Harness Interface Kit for Elite Back⁽¹⁾⁽²⁾.... \$103
- CHIPB Chest Harness Interface Kit for PB Back⁽¹⁾⁽²⁾..... \$52
- E2CHI-S Chest Harness Interface Kit for E2 Back for 14"-15" Backs⁽¹⁾⁽²⁾..... \$85
- E2CHI-M Chest Harness Interface Kit for E2 Back for 16"-20" Backs⁽¹⁾⁽²⁾..... \$85
- E2BKB Contoured anatomical PSIS/Sacrolumbar Support..... \$75

1. Does not include the chest strap.
2. This Interface Kit is only Required for 4 point harness

Ultra Low Maxx CG Power Positioning System
TDX-SP2

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

CHEST STRAPS - ABDOMINAL SUPPORT HCPCS code E0960

- MFUS BP Monoflex, Underarm RLS S(3" x 15.5") \$195
- MFUM BP Monoflex, Underarm RLS M(4" x 18.25") \$195
- MFUL BP Monoflex, Underarm RLS L(4.5" x 21.5") \$195
- MFCS BP Monoflex, Center RLS S(2.5" x 15.5") \$195
- MFCM BP Monoflex, Center RLS M(3.25" x 18.25") \$195
- MFCL BP Monoflex, Center RLS L(3.5" x 21.5") \$195
- PFLEX-M PostureFlex Abdominal Support
Height 6.5", length 16" to 23"⁽¹⁾ \$295
- PFLEX-L PostureFlex Abdominal Support
Height 7", length 20" to 32"⁽¹⁾ \$295
- QRTAB Quad release tab for easy independent
application and removal of PostureFlex \$35

1. Only available on the E2 and E2 Deep Backs.

TILT ARMRESTS

Adult Dual Post

- DP1 9.75" - 12.75" STD
- DP2 12.75" - 15.75" STD
- DP3 15.75" - 18.75" STD
- FP Flat Pouch STD
- GB Glove Boxes N/C
- OAL Outback Arm Option Left \$195
- OAR Outback Arm Option Right \$195

Ultra Rail Seat Mount Flip Back Cantilever MAXX Tilt Arm⁽¹⁾

- RMCAN 9" - 12.5" \$488
- RMCAN2 11" - 14" \$488
- DPM Dual Post Module to RMCAN (pair) \$250

Outback Cantilever Tilt Armrest

- OCTA1 9"-12" \$488
- OCTA2 12"-15" \$488
- OCTA2 6"-12" \$488

Cantiliver Arms for Endomorph Back Canes⁽¹⁾

- CAFEB See page 4 STD

1. Not available with Small Adult Sizes.

RECLINE ARMRESTS

Reclining Adjustable Height

- RC1 9.5"-13" STD
- RC2 13"-16" STD
- OAL Outback Arm Option Left \$195
- OAR Outback Arm Option Right \$195

Maxx Style Cane Mounted Height, Width, Angle Adj.

Cantilever Flip Back Armrests⁽¹⁾

- MCAN1 9.5"-13" \$250
- MCAN2 12.5"-16" \$250

Outback Cantilever Recline Armrest

- OCRA1 9"-12" \$350
- OCRA2 12"-15" \$350
- OCRA3 7" - 9" \$350

Cane Mounted Height and Angle Adj.

Cantiliver Flip Back Armrest⁽¹⁾

- CANQ1 9.5"-12" \$350
- CANQ2 10"-13" \$350
- CANQ3 13"-16" \$350

Adult Dual Post On Recline

- DPOR1 9.75" - 12.75" \$250
- DPOR2 12.75" - 15.75" \$250
- DPOR3 15.75" - 18.75" \$250
- FP Flat Pouch STD
- GB Glove Boxes N/C

1. Not available with Small Adult Sizes.

ARMREST TUBE LENGTH - PLEASE SPECIFY

- ATFL Left Full STD
- ATDL Left Desk STD
- ATFR Right Full STD
- ATDR Right Desk STD

ARMPADS⁽⁵⁾

Standard

- SAPFL Left Full (2.25" x 14") STD
- SAPDL Left Desk (2.25" x 10") STD
- SAPML Left Mini (2.25" x 7") STD
- SAPFR Right Full (2.25" x 14") STD
- SAPDR Right Desk (2.25" x 10") STD
- SAPMR Right Mini (2.25" x 7") STD

Waterfall Arm Pads⁽¹⁾

- WAPFL Left Full (13") N/C
- WAPDL Left Desk (9") N/C
- WAPFR Right Full (13") N/C
- WAPDR Right Desk (9") N/C

Modular Arm Pad - Left⁽²⁾

- MAPTFL Left Tray Full (14") \$95
- MAPTDL Left Tray Desk (10") \$95

Select Pad Insert below:

- ISSFL Iskin - Integrated soft skin foam
- SCVFL Startex Covered Visco foam

Modular Arm Pad - Right⁽²⁾

- MAPTFR Right Tray Full (14") \$95
- MAPTDR Right Tray Desk (10") \$95

Select Pad Insert below:

- ISSFR Iskin - Integrated soft skin foam
- SCVFR Startex Covered Visco foam

Flat Multi Position Arm Pads

- MPAFL Left Full (14") \$150
- MPADL Left Desk (10") \$150
- MPAML Left Mini (7") \$150
- MPAFR Right Full (14") \$150
- MPADR Right Desk (10") \$150
- MPAMR Right Mini (7") \$150

Gel Pads

- GAPNL Narrow Left Desk (2" x 12") \$88
- GAPNR Narrow Right Desk (2" x 12") \$88
- GAPWL Wide Left Desk (3.5" x 12") \$96
- GAPWR Wide Right Desk (3.5" x 12") \$96
- GAPLL Long Left Full (3.5" x 14") \$103
- GAPLR Long Right Full (3.5" x 14") \$103

Arm Troughs

- EATL Ergonomic Arm Trough - Left \$121
- EATR Ergonomic Arm Trough - Right \$121
- FHPL Flat Hand Pad - Left⁽³⁾ \$91
- FHPR Flat Hand Pad - Right⁽³⁾ \$91

Ultra Low Maxx CG Power Positioning System
TDX-SP2

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

ARMPADS⁽⁵⁾ - CONTINUED

Arm Pad Storage/Cup Holder⁽⁴⁾

- APSL Left \$295
- APSR Right \$295
- OAP Omit Armpads(7/8" diameter tube) N/C

1. Full length pad not compatible with the desk length armrest tube if using Ultra Rail mounted Tilt Arms.
2. Requires interchangeable pad insert.
3. Must order with ergonomic arm trough.
4. Not available on the joystick side, or with Outback Cantilever arms.
5. Full length arm pads used on the same side as the Joystick result in a "Back to Joystick" measurement of 20" minimum (24" including switch mounting)

ARM ADDUCTOR PAD⁽¹⁾⁽²⁾

- AADSL Left Small (6" x 2.25") \$140
- AADSR Right Small (6" x 2.25") \$140
- AADLL Left Large (8.5" x 2.25") \$140
- AADLR Right Large (8.5" x 2.25") \$140

1. Only available with Modular or Flat Multi Position Arm Pads.
2. 1" thick pad, unlimited fore/aft adjustment.

ELBOW BLOCKS

Left Elbow Block⁽¹⁾

Select Pad Size and Both Extensions:

- EBXL Size XS (3.25" x 4.25") \$195
- EBSL Size S (3.5" x 5.25") \$195
- EBML Size M (3.75" x 5.75") \$195
- EBLL Size L (4.25" x 6.5") \$195

- LV3 Left Vertical Extension 3"
- LV5 Left Vertical Extension 5"
- LH5 Left Horizontal Extension 5"
- LH7 Left Horizontal Extension 7"

Right Elbow Block⁽¹⁾

Select Pad Size and Both Extensions:

- EBXR Size XS (3.25" x 4.25") \$195
- EBSR Size S (3.5" x 5.25") \$195
- EBMR Size M (3.75" x 5.75") \$195
- EBLR Size L (4.25" x 6.5") \$195

- RV3 Right Vertical Extension 3"
- RV5 Right Vertical Extension 5"
- RH5 Right Horizontal Extension 5"
- RH7 Right Horizontal Extension 7"

ADDITIONAL ARM ACCESSORIES

- MACESL Multi Axis Upper Extremity Support⁽¹⁾(each)
Left Side \$225
- MACESR Multi Axis Upper Extremity Support⁽¹⁾ (each)
Right Side \$225
- OAS2L 1" Offset Armrest Spacer Left(each) \$125
- OAS2R 1" Offset Armrest Spacer Right(each) \$125
- OAS1L 1/2" Offset Armrest Spacer Left(each) \$125
- OAS1R 1/2" Offset Armrest Spacer Right(each) \$125
- MPAPA Adapter Block⁽²⁾ (each) Quantity () \$35

1. Allows rotation of the armpad, only compatible with Gel Arm pad, Flat Arm pad or Ergonomic Arm Troughs.
2. Not available with Standard Arm Pads or Waterfall Arm Pads

MAXX LATERAL TRUNK SUPPORTS⁽¹⁾

Lateral Trunk Support Left

- LATXL Size XS (3.25" x 4.25") \$100
- LATSL Size S (3.5" x 5.25") \$100
- LATML Size M (3.75" x 5.75") \$100
- LATLL Size L (4.25" x 6.5") \$100

Add:

- SALATL MAXX Style Swing Away Hardware - Left⁽²⁾ \$195
- ATLL Additional Telescoping Link - Left⁽³⁾ \$75

Lateral Trunk Support Right

- LATXR Size XS (3.25" x 4.25") \$100
- LATSR Size S (3.5" x 5.25") \$100
- LATMR Size M (3.75" x 5.75") \$100
- LATLR Size L (4.25" x 6.5") \$100

Add:

- SALATR MAXX Style Swing Away Hardware - Right⁽²⁾ \$195
- ATLR Additional Telescoping Link - Right⁽³⁾ \$75

1. Comes standard with Fixed Mount Hardware. Not available with Matrix PB, PB Elite Deep, and PB Deep Backs.
2. Comes standard with 1" telescoping link which gives 0-2.75" medial offset.
3. The amount of offset available may be impacted by back thickness.

MATRIX FIXED LATERAL TRUNK SUPPORTS⁽¹⁾

- E2FLL Matrix E2 Fixed Laterals-Left(6.5x5")⁽²⁾ \$195
- E2FLR Matrix E2 Fixed Laterals-Right(6.5x5")⁽²⁾ \$195
- PBSSL Matrix Small Lateral Fixed-Left (3x5.5") \$113
- PBSSR Matrix Small Lateral Fixed-Right (3x5.5") \$113
- PBLSL Matrix Standard Lateral Fixed-Left (7.5x5") \$116
- PBLSR Matrix Standard Lateral Fixed-Right (7.5x5") \$116
- PBOLL Matrix Offset Lateral Fixed Left (7.5x5") \$125
- PBOLR Matrix Offset Lateral Fixed Right (7.5x5") \$125

1. Not available with Matrix PB, PB Elite Deep, and PB Deep Backs.
2. Only available with E2 backs.

MATRIX SWING AWAY LATERAL TRUNK SUPPORTS⁽¹⁾

Matrix Lateral Trunk Support Left

- EPADXL Size XS (3.25" x 4.25") \$115
- EPADSL Size S (3.5" x 5.25") \$115
- EPADML Size M (3.75" x 5.75") \$115
- EPADLL Size L (4.25" x 6.5") \$115

Must select one option:

- SAHEBL Matrix Elite Swing Away Hardware Left \$250
- OSAHL Matrix Offset Elite Swing Away Hardware Left \$250

Matrix Lateral Trunk Support Right

- EPADXR Size XS (3.25" x 4.25") \$115
- EPADSR Size S (3.5" x 5.25") \$115
- EPADMR Size M (3.75" x 5.75") \$115
- EPADLR Size L (4.25" x 6.5") \$115

Must select one option:

- SAHEBR Matrix Elite Swing Away Hardware Right \$250
- OSAHR Matrix Offset Elite Swing Away Hardware Right \$250

1. Not available with Matrix PB, PB Elite Deep, and PB Deep Backs.

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

HIP AND KNEE SUPPORTS

Lateral Hip Support Left

- HSPXLH Size XS (5.5" W x 4" H) \$120
- HSPSMLH Size S/M (7" W x 4" H) \$120
- HSPSLH Size S (7"W x 6"H) \$120
- HSPMLH Size M (9" W x 4"H) \$120
- HSPLLH Size L (13" W x 4"H) \$120

Lateral Hip Support Right⁽¹⁾

- HSPXRH Size XS (5.5"W x 4"H) \$120
- HSPSMRH Size S/M (7"W x 4"H) \$120
- HSPSRH Size S (7"W x 6"H) \$120
- HSPMRH Size M (9"W x 4"H) \$120
- HSPLRH Size L (13"W x 4"H) \$120

Lateral Knee Support Left⁽¹⁾

- HSPXLK Size XS (5.5"W x 4" H) \$120
- HSPSMLK Size S/M (7"W x 4" H) \$120
- HSPSLK Size S (7"W x 6" H) \$120
- HSPMLK Size M (9"W x 4"H) \$120
- HSPLLK Size L (13"W x 4"H) \$120

Lateral Knee Support Right⁽¹⁾

- HSPXRK Size XS (5.5"W x 4"H) \$120
- HSPSMRK Size S/M (7"W x 4"H) \$120
- HSPSRK Size S (7"W x 6"H) \$120
- HSPMRK Size M (9"W x 4"H) \$120
- HSPLRK Size L (13"W x 4"H) \$120

Add Removable Hardware:

MAXX Style Quick Release, Multi-Axis Mounting

Hardware for Hip Supports

- MSHSLH Left mounted at Hip \$215
- MSHSLK Left mounted at Knee \$215
- MSHSRH Right mounted at Hip \$215
- MSHSRK Right mounted at Knee..... \$215

Swing Away Hardware for Hip Supports

- SHSHLH Left mounted at Hip \$215
- SHSHLK Left mounted at Knee \$215
- SHSHRH Right mounted at Hip \$215
- SHSHRK Right mounted at Knee..... \$215

Extended Maxx Style Quick Release, Multi-Axis Mounting

Hardware for Hip Supports⁽²⁾

- EMQHLH Left mounted at Hip \$215
- EMQHLK Left mounted at Knee \$215
- EMQHRH Right mounted at Hip \$215
- EMQHRK Right mounted at Knee..... \$215

Lift Off Hardware for Hip Support

(Hip Position Only)

- LOHL Left \$195
- LOHR Right..... \$195

1. Comes standard with Fixed Mount Hardware.
2. Recommended when finished seat cushion thickness will be greater than 3.5" high.

POWER CENTER MOUNT FOOT PLATFORM

MUST BE COMPLETED	Measurement in Inches
Knee to Heel:	
Cushion Thickness: (Subtract)	
Cushion Compression:(Add)	
Seat Pan to Footplate:(Total)	

LNx Power Center Mount Foot Platform

- LN1 9" Articulation⁽¹⁾⁽²⁾ \$3,000
 - LSPTF1 10" - 14" Seat pan to footplate
 - LSPTF2 13" - 19" Seat pan to footplate
- LN1A 6" Articulation⁽¹⁾⁽²⁾⁽⁶⁾ \$3,000
 - LASPTF2 10" - 14" Seat pan to footplate

Options:

- LIFP Rubber Coated Footplate 11.5" W x 10" D STD
- LICP Individual calf pads STD
- OPCP One Piece Malleable Calf Panel⁽³⁾ \$250
- LNXPFF Plantar Flex Option for LN1 Foot platform⁽⁴⁾ \$250

Latitude Style Foot Platform

- LSFP Latitude Style Foot Platform⁽⁵⁾ \$2,995
For finished knee to heel measurements of 6.5"-12"
- LSPFP Latitude Style Padded Foot Plate \$250

1. 3" ground clearance required.
2. Only available on adult range seat systems 16" deep or greater. For shorter seat depths call Customer Service for other options.
3. Not available with Plantar Flex Option for LN1 Foot platform, or with LN1A.
4. Only available with Foot Platform Options.
5. Not available with XRO and XRO-E.
6. For clients under 4' only.

MANUAL ELEVATING CENTER MOUNT LEGREST

- MECML Manual Elevating Center Mount Legrest
With Plantar Flex \$595

Select legrest set up position:

- 90SU Set at 90° (elevating range 90-40°) STD
- 110TU Tucked at 110° (elevating range 110-60°)⁽¹⁾ N/C

Rubber Coated Footplate:

- MRCP Rubber coated footplate 11.5" W x 10" D N/C
- Select seat pan to footplate for Rubber Coated Footplate
 - CSPTF1 11.5" - 15" N/C
 - CSPTF2 14" - 19" N/C

Individual Footplates:


- IFSL Small - (5" x 7.5" deep) - Left \$165
- IFSR Small - (5" x 7.5" deep) - Right \$165
- IFML Medium - (5.5" x 9.5" deep) - Left \$165
- IFMR Medium - (5.5" x 9.5" deep) - Right \$165
- IFLL Large - (6" x 11.5" deep) - Left \$165
- IFRL Large - (6" x 11.5" deep) - Right \$165

Select seat pan to footplate for Individual Footplates

- CSPTF3 9.5" - 14.5" N/C
- CSPTF4 12" - 19" N/C

Optional:

- CICP Individual calf pads \$200



1. Not available on seat depths 16", 17".

Ultra Low Maxx CG Power Positioning System
TDX-SP2

Dealer Name: _____


Dealer Account Number: _____

PO #: _____

Tag: _____


STANDARD FRONT RIGGING

- MSFCMFP Maxx Style Fixed Center Mount Foot Platform⁽¹⁾⁽²⁾\$275
 - MIFP Comes with a rubber coated footplate 11.5" W x 10" DN/C
 - MICP Individual calf pads\$100
- Select seat pan to footplate:
- MSPTF1 9" - 13"N/C
 - MSPTF14 13" -19"N/C
- Set at:
- MSA70 Set at 70°N/C
 - MSA83 Set at 83°N/C
 - MSA90 Set at 90°N/C
 - MSA97 Set at 97°N/C



1. Depth adjustable, and easy slide adjustment for seat pan to foot plate extension.
2. 3" ground clearance required.

INDIVIDUAL FOOT PLATES FOR FIXED/LNX CMFP



- IFSL Small (5" W x 7.5" D) - Left\$165
 - IFSR Small (5" W x 7.5" D) - Right\$165
 - IFML Medium (5.5" W x 9.5" D) - Left\$165
 - IFMR Medium (5.5" W x 9.5" D) - Right\$165
 - IFLL Large (6" W x 11.5" D) - Left\$165
 - IFLR Large (6" W x 11.5" D) - Right\$165
- Add:
- TSSL Left Toe Strap Small 6.5"\$125
 - TSSR Right Toe Strap Small 6.5"\$125
 - TSML Left Toe Strap Medium 8.5"\$125
 - TSMR Right Toe Strap Medium 8.5"\$125
 - TSLL Left Toe Strap Large 10.5"\$125
 - TSLR Right Toe Strap Large 10.5"\$125
- MFOL Multi-Axis Footplate (6"W x 10"D) - Left\$225
 - MFOR Multi-Axis Footplate (6"W x 10"D) - Right\$225

1. Foot plates may interfere with casters in certain configurations


FOOT PLATFORM OPTIONS FOR FIXED/LNX CMFP



- SFFP Small Foot Platform (12"w x 8.5"d)\$180
- MFFP Medium Platform (12"w x10.75"d)\$180
- LFFP Large Platform (14.5"w x13"d)\$250
- XLFFP Extra Large Platform (17"w x13"d)\$250
- XXLFFP Extra Extra Large Platform (21"w x13"d)\$250

1. May need to adjust leg position to avoid interference with front casters

MAXX INDEPENDENT ADJUSTABLE POSITIONING



- 5.5" - 6" WIDTH ADJUSTMENT
- 9.5" - 11.5" DEPTH ADJUSTMENT
- UP TO 9° ADDUCTION
- UP TO 45° ABDUCTION

- MPFPL Maxx Independent Adjustable Positioning Foot plate⁽¹⁾ - Left\$200
- MPFPR Maxx Independent Adjustable Positioning Foot plate⁽¹⁾ - Right\$200
- MTSL Left Toe Strap (each)\$35
- MTSR Right Toe Strap (each)\$35

1. Foot plates may interfere with casters in certain configurations.

Ultra Low Maxx CG Power Positioning System
TDX-SP2

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

INVACARE CENTER MOUNT⁽¹⁾

Invacare Center Mount with Footboard⁽³⁾

Select seat pan to footplate:

- CTMT-FJ-FB Invacare Center Mount 9.25" - 12" \$760
- CTMT-F-FB Invacare Center Mount 12.25" -15" \$760

Invacare Center Mount with Individual Footplates

Select seat pan to footplate:

- CTMT-FJ Invacare Center Mount 9.25" - 12" \$595
- CTMT-F Invacare Center Mount 12.25" -15" \$595

Select Footplates:

- SM Individual Foot Plates (5.25"W x 7.5"D)..... \$165
- MDM Individual Foot Plates (5.25"W x 10"D)..... \$165
- LRG Individual Foot Plates (5.9"W x 12"D)..... \$165

Optional:

- CALF Invacare Calf Pads for use with CTMT⁽²⁾ \$268

1. 3" ground clearance required.
2. Not available with CTMT-FJ, CTMT-FJ-FB
3. Footboard Included

MAXX STYLE POWER ELEVATING LEGRESTS

- MSPPL MAXX Swing Away Power Elevating⁽¹⁾ \$2,490

Select Power Legrest Function:

- INDA Individual Movement N/C
- COMA Combined Movement N/C
- SYNA Synchronized/Recline Movement..... N/C

Select Receivers:

- MRS Straight Receivers N/C
- MRF Flared 1" Outward Receivers N/C
- MRF2 Flared 2" Outward Receivers..... \$250

Select Seat Pan to Footplate:

- PSPTF5 Length: 10"-14" N/C
- PSPTF2 Length: 13"-17" N/C
- PSPTF7 Length: 16"-20" N/C

Optional:

- PMCP Moulded I-Skin Calf Pads STD
- PFCPMS Flat Calf Pad Option (pair) \$275
- PCCPMS Curved Calf Pad Option⁽²⁾ (pair) \$275

1. Must choose footplate on page 15.
2. Not available on 14 and 15" seat widths. For 16" and 17" seat widths you must select flared receivers.

MAXX STYLE MANUAL ELEVATING LEGRESTS

- MMPPSA MAXX Swing Away Manual Elevating⁽¹⁾ \$680

Select Receivers:

- ERS Straight Receivers N/C
- ERF Flared 1" Outward Receivers N/C
- ERF2 Flared 2" Outward Receivers \$250

Select Seat Pan to Footplate:

- ESPTF5 Length: 10"-14" N/C
- ESPTF2 Length: 13"-17" N/C
- ESPTF7 Length: 16"-20" N/C

Optional:

- EMCP Moulded I-Skin Calf Pads STD
- EFCPMS Flat Calf Pad Option (pair) \$275
- ECCPMS Curved Calf Pad Option⁽²⁾ (pair) \$275

1. Must choose footplate on page 15.
2. Not available on 14 and 15" seat widths. For 16" and 17" seat widths you must select flared receivers.

SWING AWAY LEGREST

- HD70 Heavy Duty 70° Swing Away Footrest(pair)⁽¹⁾⁽²⁾ . \$495

Select Receivers:

- RS Straight Receivers N/C
- RF Flared 1" Outward Receivers N/C
- RF2 Flared 2" Outward Receivers \$250

Select Seat Pan to Footplate:

- HSPTF11 Length: 6"-10" N/C
- HSPTF8 Length: 10"-14" N/C
- HSPTF2 Length: 12"-16" N/C
- HSPTF7 Length: 15"-19" N/C

Options:

- HD70FCP Flat Calf Pad Option (pair) \$275
- HD70CCP Curved Calf Pad Option⁽³⁾ (pair) \$275

- I70SAFT Invacare 70° Swingaway Footrests..... \$325
- I70-CFT Comes with composite foot plates..... N/C

- RLSL Residual Limb Support - Left (8.5" x 6")⁽⁴⁾ \$495
- RLSR Residual Limb Support - Right (8.5" x 6")⁽⁴⁾ \$495

1. Must choose footplate on page 15.
2. Heavy duty durable construction (4.5lbs ea.).
3. Not available on 14 and 15" seat widths. For 16" and 17" seat widths you must select flared receivers.
4. Swing away and removable, 50° of angle adjustment, 4.5" height adjustment, 3" horizontal and 2" of for/aft adjustment.

DUAL FRONT RIGGING ADAPTOR⁽¹⁾

- DFRA Dual Front Rigging Receivers..... \$450

1. Must purchase each set of front rigging separately. Not available on no charge front rigging options.

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

FOOTPLATES - SWING AWAY/ELEVATING LEGRESTS⁽¹⁾




Angle Adjustable Footplates

- AFSL Left Small Footplate (4.25"W x 8"D) \$165
- AFSR Right Small Footplate (4.25"W x 8"D) \$165
- AFML Left Medium Footplate (5.25"W x 8"D) \$165
- AFMR Right Medium Footplate (5.25"W x 8"D) \$165
- AFLL Left Large Footplate (6.25"W x 8"D) \$165
- AFLR Right Large Footplate (6.25"W x 8"D) \$165
- AFXL Left Extra Large Footplate (9"W x 11"D) \$165
- AFXR Right Extra Large Footplate (9"W x 11"D) \$165

Multi-Axis Adjustable Angle Footplates

- MFSL Left Small Footplate (4.25"W x 8"D) \$225
- MFSR Right Small Footplate (4.25"W x 8"D) \$225
- MFML Left Medium Footplate (5.25"W x 8"D) \$225
- MFMR Right Medium Footplate (5.25"W x 8"D) \$225
- MFLL Left Large Footplate (6.25"W x 8"D) \$225
- MFRL Right Large Footplate (6.25"W x 8"D) \$225
- MFXL Left Extra Large Footplate (9"W x 11"D) \$225
- MFXR Right Extra Large Footplate (9"W x 11"D) \$225



Cut Out Angle Adjustable Footplates

- AFCSL Left Small Footplate (7.25"W x 9"D) \$165
- AFCSR Right Small Footplate (7.25"W x 9"D) \$165
- AFCML Left Medium Footplate (7.5"W x 11"D) \$165
- AFCMR Right Medium Footplate (7.5"W x 11"D) \$165
- AFCLL Left Large Footplate (10"W x 14"D) \$165
- AFCLR Right Large Footplate (10"W x 14"D) \$165

Cut Out Multi-Axis Adjustable Angle Footplates

- MFCSL Left Small Footplate (7.25"W x 9"D) \$225
- MFCSR Right Small Footplate (7.25"W x 9"D) \$225
- MFCML Left Medium Footplate (7.5"W x 11"D) \$225
- MFCMR Right Medium Footplate (7.5"W x 11"D) \$225
- MFCLL Left Large Footplate (10"W x 14"D) \$225
- MFCLR Right Large Footplate (10"W x 14"D) \$225

1. Certain footplate sizes may not fit all seat widths depending on legrests selected.

HEEL LOOPS, CUPS AND TOE STRAPS

- HLL Left Heel Loop (each) \$20
- HLR Right Heel Loop (each) \$20
- HCL Left Cup Loop (each) \$20
- HCR Right Cup Loop (each) \$20
- TSL Left Toe Strap (each) \$20
- TSR Right Toe Strap (each) \$20

CALF STRAPS AND CALF PANELS⁽¹⁾

- CSM Body Point Airmesh Calf Strap for Chairs
15"-18" wide with 3.5" H Pad⁽¹⁾ \$200
- CSL Body Point Airmesh Calf Strap for Chairs
18"-22" wide with 3.5" H Pad⁽¹⁾ \$200
- CPM Body Point Airmesh Calf Panel for Chairs
15"-18" wide with 35.5" W x 9" H Pad⁽¹⁾ \$295
- CPL Body Point Airmesh Calf Panel for Chairs
18"-22" wide with 38.5" W x 9" H Pad⁽¹⁾ \$315

1. Not available for Center Mounted Front Rigging.

STANDARD FRONT RECEIVERS

- PHDR Pride HD Drop-in Receivers Only N/C
- PSAR Pride Style #8 Swing Away Receivers Only N/C
- IAR Invacare 70° Swing Away Receivers Only N/C
- Q70R Quickie 70° Swing Away Receivers Only N/C
- Q70HDR Quickie 70HD Swing Away Receivers Only N/C

SINGLE-DUAL FUNCTION ELECTRONICS

- SFCB Single/Dual Function Electronics⁽¹⁾ STD
- Select Switch:
 - SPBSO Single Push Button
 - STO Single Toggle
 - DPBO Dual Push Button
 - OSO Omit Switch
- SFCBA Single/Dual Function Drive Control⁽¹⁾ \$750
- Select Control:
 - JOY Through Joystick
 - SC Through Specialty Controls
- OSO Reduced Drive for Elevating Seat⁽²⁾ \$400

1. These electronics include drive lock-out and tilt limit.
2. If reduced drive is not selected the system will be in drive lock-out when elevated.

MULTI- FUNCTION ELECTRONICS

- MFCB Multi-Function through Switch⁽¹⁾ \$1,950
- Select Switch:
 - QPBO Quad Push Button N/C
 - FWT Four-Way Toggle Switch \$150
 - EWR 8-Way Rocker⁽²⁾ \$205
- MFCBD Multi-Function M610 Through Drive Control⁽¹⁾ \$1,950
- Includes reduced drive for Elevating Seat
- Select Control:
 - JOY Through Joystick
 - SC Through Specialty Controls
- Optional:
 - RDEY Include Reduced Drive for Elevating Seat⁽³⁾
 - RDEN Do Not Include Reduced Drive for Elevating Seat⁽³⁾
- Optional:
 - MFMSR Multi Function Memory Seating N/C

1. These electronics include attendant control, drive lock-out, tilt/recline limit and actuator speed control.
2. Dual function is not available with 8-way Rocker.
3. If reduced drive is not selected system will be in drive lock-out when elevated.

ELECTRONIC SWITCH MOUNTING

- ESML Left Side (Armpad)
- ESMR Right Side (Armpad)
- ESMJL Joystick Left (Mounts to Joystick Bracket)
- ESMJR Joystick Right (Mounts to Joystick Bracket)
- ESMD Do not mount

ELECTRONICS OPTIONS

- Additional Switches⁽⁴⁾
- SPBSA Single Push Button \$150
 - STA Single Toggle Switch⁽³⁾ \$150
 - DPBA Dual Push Button⁽³⁾ \$150
 - QPBA Quad Push Button⁽¹⁾⁽³⁾ \$150
 - 10WA 10-Way Push Button⁽¹⁾⁽³⁾ \$255
 - ASL320R ASL Pulse Mechanical Switch - Red \$200
 - ASL320G ASL Pulse Mechanical Switch - Green \$200
 - ASL320BK ASL Pulse Mechanical Switch - Black \$200
 - SESG Stealth Egg Switch - Green \$135
 - SESB Stealth Egg Switch - Black \$135
 - FTSQ Feather Touch - Single Quarter Size \$295
 - FTMD Feather Touch - Dual Dime Size⁽¹⁾⁽³⁾ \$295
 - FTMQ Feather Touch - Quad Dime Size⁽¹⁾⁽³⁾ \$295
- Additional Switches⁽⁴⁾ (not available with ACT4 electronics)
- FWTA Four-Way Toggle Switch⁽¹⁾⁽³⁾ \$150
 - EWRA 8-Way Rocker⁽¹⁾⁽³⁾ \$255
 - TWRA 10-Way Rocker⁽¹⁾⁽³⁾ \$255
 - TWTA 10-Way Toggle⁽¹⁾⁽³⁾ \$255
 - EWPA 8-Way Push Button⁽¹⁾⁽³⁾ \$255

Specify Use for Additional Switch:

- STG Seating
- JM Joystick Mode⁽²⁾
- JOO Joystick On/Off⁽²⁾

Specify Mounting for Additional Switch:

- ESMLA Left Side (Armpad)
- ESMRA Right Side (Armpad)
- ESMJLA Joystick Left (Mounts to Joystick Bracket)
- ESMJRA Joystick Right (Mounts to Joystick Bracket)
- ESMDA Do not mount

Miscellaneous

- MJBS Motion Junction Box for Switch Control⁽¹⁾⁽⁵⁾
For Mode and/or Emergency Stop. \$175
- S46L Single to Dual Switch Adapter \$265
Allows for dual use of a single switch
- M616S M616 Converter⁽¹⁾ \$125
Allows up to 4 mono phono switches into Multi-Function box
- LJS Lazarus Jump Starter \$150
- APPS Accessory Port Power Supply⁽¹⁾
3 amp capacity \$195
- TPSMH 12 Piece Switch Mounting Hardware Kit (additional kits) \$39

1. Not available with Single-Dual Function Electronics.
2. IN200 is required with REM210, REM211, REM216 joysticks.
3. Not available with REM110 joystick.
4. Not available for Joystick Mode or Joystick On/Off. Must Select Seating (STG)
5. Cannot be used for power On/Off.
6. Not available on Multi Function ACT4-LiNX Electronics

Ultra Low Maxx CG Power Positioning System
TDX-SP2

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

ARTICULATING TRAYS⁽¹⁾

- USAT Universal Style Articulating Tray \$1,350
14"W x Adjustable 6"-10"D
- Add:
- TMA Trilogy Mounting Adapter \$50
- AVT Articulating Vent (LTV) Tray \$1,350
- AVBH Articulating Vent Bag Holder⁽²⁾ \$1,350
- Add:
- ABMC Accessory Bag Mounting and Clips \$100

1. The user weight limit is 225 lbs. The maximum weight capacity of the vent tray is 35 lbs. Vents are not standard configuration on Elevating or Precline Systems or with Single Post Backs. Not available with TRRO and TRRO-E.
2. Motion Concepts does not supply the Vent Bag.

GIMBAL VENT TRAYS⁽¹⁾

- GUT Gimbal Universal Tray \$1,350
14"W x Adjustable 6"-10"D
- Add:
- GTMA Trilogy Mounting Adapter \$50
- GVT Gimbal Vent (LTV) Tray \$1,350
- GVBH Gimbal Vent Bag Holder⁽²⁾ \$1,350
- Add:
- GABMC Accessory Bag Mounting and Clips \$100

1. The Gimbal Vent Trays are not available with Tilt Only and/or Elevate Only Systems. The user weight limit is 200 lbs. Vents are not standard configuration on Elevating or Precline Systems or with Single Post Backs. The maximum weight capacity of the vent tray is 20 lbs. Not available with TRRO and TRRO-E.
2. Motion Concepts does not supply the Vent Bag.

O2 HOLDERS⁽¹⁾

- NAO2H Non-Articulating O2 Gas Holder \$395
- AOHG Articulating O2 Gas Holder Cylinder Style⁽²⁾⁽³⁾ . \$1,250
- AOHL Articulating Liquid O2 Holder (6.5" x 4.5")⁽²⁾⁽³⁾ . \$1,250
- NAIPL Non-Articulating IV Pole User's Left⁽⁴⁾ \$395
- NAIPR Non-Articulating IV Pole User's Right⁽⁴⁾ \$395
- AIVP Articulating IV Pole⁽²⁾⁽³⁾⁽⁴⁾ \$1,250

1. The user weight limit is 225 lbs. The maximum weight capacity of the O2 holder is 20 lbs. Not available with TRRO and TRRO-E.
2. Required for use on Recline Systems
3. Articulating O2 Holders are not standard configuration on Elevating or Precline Systems or with Single Post Backs.
4. The maximum weight capacity of the IV pole is 5 lbs.

VENTILATOR INFORMATION

Required when ordering a Vent Tray

- BIB Built in battery
- DTSVPS Dealer to supply vent power source

REAR STORAGE BASKET⁽¹⁾⁽²⁾

- RSB Removable Rear Storage Basket \$1,395
Exterior dimensions: 18"W x 6"H x 7"D
Interior dimensions: 17.8"W x 6"H x 6.8"D
Interior capacity: 3 gallon
Max Load: 20 lbs

1. Tall items extending beyond the basket's height may be damaged upon tilting.
2. Not available with TRRO, TRRO-E brackets

ACCESSORIES

- MUSB Motion Dual USB Charger N/C
- UCH Universal Cup Holder⁽¹⁾⁽⁵⁾ \$50
- SPH Smart Phone Holder⁽¹⁾ \$50
- PH X-Grip Phone Holder \$95
- AMC Accessory Mount with Clip (QTY ____) ⁽²⁾ \$75
- CRVMR Rear View Mirror Right Mounted \$65
- CRVML Rear View Mirror Left Mounted \$65
- TAG571850 Wheelchair Cover⁽²⁾ \$150
- HK Hydration Kit⁽²⁾ \$350
- MUB Medical Utility Bag w/ Mounting Hardware \$350
- THS Transfer Handles 4" - 6.5" \$350
- THL Transfer Handles 9" - 11.5" \$350
- FFLTL Flip and Fold Away 1/2 Lap Tray
Left Mounting(14"W x 10"D)⁽²⁾⁽³⁾⁽⁴⁾ \$795
- FFLTR Flip and Fold Away 1/2 Lap Tray
Right Mounting(14"W x 10"D)⁽²⁾⁽³⁾⁽⁴⁾ \$795
- CFH Flag Holder ⁽²⁾ \$90
 - Left Mount
 - Right Mount
- CSCC Cane and Crutch Holder \$395
 - Left Mount
 - Right Mount

1. Not available with TRRO, TRRO-E brackets
2. Not available on joystick side and only available with a flat arm pad or modular arm pad. Not available with Outback Cantilever Arms.
3. The maximum weight capacity of the tray is 10 lbs.
4. Not available on the joystick side

INSTALLATION

- MCI Motion Concepts Installation⁽¹⁾ Included

1. Batteries may be required if items on Accessories order form are selected at time of installation, please call to confirm.

RETAIL \$ _____

SUBMIT

PRINT FORM
SAVE
CLEAR FORM