



## Thank you for using Motion Concepts eForms.

Motion Concepts is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email\* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.



## Adobe Acrobat Reader DC

\* Interactive functions of our new forms work best with the latest version of Adobe Acrobat Reader DC visit <https://get.adobe.com/reader/> to download and install on your PC or Mac or visit Google Play/ iTunes to download the Adobe Acrobat Reader DC app for your device.

### Save

Adobe Acrobat Reader DC allows you to save this form with your content - to complete later or use as a starting point for your next form. Please note that content must be added and saved in Acrobat - saving content from completed forms in the browser may not be possible.

SAVE

### Submit

Adobe Acrobat Reader DC allows you to submit this form electronically via your email client.

Simply click the submit button below and step through the simple process.

SUBMIT

CLEAR FORM

### Print

If you do not have access to Adobe Acrobat Reader DC simply print this form and complete it by hand and fax it to our Customer Service Department at:

888-433-6818

PRINT FORM



**TURNING,  
DRIVING  
& RIDE  
COMFORT!**

# ULTRA LOW MAXX CG POWER POSITIONING SYSTEM ON STORM TORQUE RX20 BASE

CDN PRICE LIST AND ORDER FORM - Price Effective October 22, 2025



PRINT FORM

SAVE

CLEAR FORM

Motion Concepts Customer Service: 1-866-748-7943 quotescdn@motionconcepts.com

## REQUIRED INFORMATION

Quote  Order

Date of Order: \_\_\_\_\_ Dealer Account #: \_\_\_\_\_

Dealer Name: \_\_\_\_\_

PO #: \_\_\_\_\_ Tag \_\_\_\_\_

Purchasing Contact: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

E-mail: \_\_\_\_\_

RTS/Therapist: \_\_\_\_\_

Ship to Address: \_\_\_\_\_

City: \_\_\_\_\_ Province: \_\_\_\_\_


Postal Code: \_\_\_\_\_ Client Gender:  M  F

Special Client Conditions: \_\_\_\_\_

## REQUIRED MEASUREMENTS

To ensure system is accurately configured please fill in all REQUIRED\* measurements

Client Height*	
Client Weight (lbs)*	
A. Seat to Shoulder	
B. Trunk Depth	
C. Chest Width	
D. Knee to Back	
E. Seat to Top of Head	
F. Back to Joystick	
G. Seat to Elbow	
H. Hip Width	
I. Knee to Heel*	
Cushion Thickness*	



## BASE OPTIONS

**NEW!**

SRX-20 Storm Torque RX20 Base<sup>(1)</sup> ..... \$7,575



- Enhanced Front Caster Housing and Bearings
- New Front Caster Design
- Improved double-sided Front Forks
- Upgraded Rear Suspension
- Optimized Linx Parameters

1. Comes Standard with Lighting Package

## BASE WIDTH

NARROWFRM Narrow Base 22.4" (MK 50 Batteries) ..... N/C  
 STDFRM Standard Base 24.4" (MK 24 FT Batteries) ..... N/C

## SEAT-TO-FLOOR HEIGHT<sup>(1)</sup>

SLOW 16.75" Seat-To-Floor<sup>(2)</sup> ..... \$250  
 LOW 17.25" Seat-To-Floor ..... N/C  
 MED 18.25" Seat-To-Floor ..... N/C  
 TALL 19.25" Seat-To-Floor ..... N/C

1. Seat-to-floor height is measured at the front of the seat pan based on a standard seat size. Seat-to-floor height may vary +/- .25" dependent on condition of the power base and seat depths greater than 18"
2. This seat to floor height is only available with the Group 24 Battery Tray (24TRY) and on 17" seat depths and greater.

## USER WEIGHT LIMIT

UW250 Weight Capacity up to 250lbs<sup>(1)</sup> ..... STD  
 U300 Weight Capacity up to 300lbs ..... STD

1. Required for Elevate Systems and/or NARROWFRM

## WHEELCHAIR TRANSPORT BRACKETS

TRBKTS Wheelchair Transport Brackets<sup>(1)</sup> ..... STD  
 TRRO Transport Ready Option<sup>(2)(4)</sup> ..... \$499  
 TRRO-E Transport Ready Option for Elevate<sup>(2)(3)(4)</sup> ..... \$499

1. For unoccupied use only
2. TRRO meets ISO-7176 Part 19

## WHEELCHAIR TRANSPORT BRACKETS

3. Required for Elevate Systems to meet ISO-7176 Part 19, includes two additional tie-down loops on rail.
4. Must select a Motion Concepts or Matrix Elan headrest with hardware and Motion Concepts Standard Rehab, High Back or Matrix Upgradeable back. Comes with pelvic safety belt standard, please select TROBELT on page 5.

## TIRE OPTIONS

AB1431-3 14" x 3" Black Tire with Gel Foam Inserts ..... N/C  
 AB1430-3 14" x 3" Black Tire - Pneumatic Tire ..... N/C

## CASTER OPTION

**NEW!**

AB1228-3 8" x 2.5" Black Ribbed Casters with Solid Inlay .. STD

## BATTERY CHARGER OPTION

110CHARGER 110 Volt Battery Charger ..... N/C  
 OMIT Omit Battery Charger ..... N/C

## BATTERIES

M24 SLD G FT MK Group 24 FT Gel Battery (Pair) ..... \$840  
 M50-12 SLD M MK Group 50 Battery (Pair) ..... \$650

## HIGHSPEED MOTORS AND ENHANCED DRIVING

HSMTR Performance Motors 7.5 mph<sup>(1)(2)(3)</sup> ..... \$995  
 GTRAC-LX LINX G-TRAC Module<sup>(4)</sup> ..... \$1,200

1. Not available on the NARROWFRM.
2. Extended high speed use on rugged terrain is not recommended.
3. Performance/High Speed Motors may have less torque and may have less range than standard speed motors.
4. Recommended for use with Specialty Drive Controls

Ultra Low Maxx CG Power Positioning System  
Storm Torque RX20









Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**FRAME FINISH**

<input type="checkbox"/> 170P		Obsidian Black.....	N/C
<input type="checkbox"/> 171P		Ruby Red.....	N/C
<input type="checkbox"/> 172P		Peacock Blue.....	N/C
<input type="checkbox"/> 173P		Pink Flamingo.....	N/C
<input type="checkbox"/> 174P		Purple Passion.....	N/C
<input type="checkbox"/> 175P		Denim Navy.....	N/C
<input type="checkbox"/> 176P		Lime Twist.....	N/C
<input type="checkbox"/> 177P		Arctic Blue.....	N/C

**RIM INSERT COLORS - MUST PICK ONE**

<input type="checkbox"/> 170PR		Obsidian Black.....	N/C
<input type="checkbox"/> 171PR		Ruby Red.....	N/C
<input type="checkbox"/> 172PR		Peacock Blue.....	N/C
<input type="checkbox"/> 173PR		Pink Flamingo.....	N/C
<input type="checkbox"/> 174PR		Purple Passion.....	N/C
<input type="checkbox"/> 175PR		Denim Navy.....	N/C
<input type="checkbox"/> 176PR		Lime Twist.....	N/C
<input type="checkbox"/> 177PR		Arctic Blue.....	N/C
<input type="checkbox"/> 178PR		Banana Split.....	N/C
<input type="checkbox"/> 179PR		Orange Marmalade.....	N/C
<input type="checkbox"/> 180PR		Morning Mist.....	N/C

**NON-EXPANDABLE SINGLE/DUAL JOYSTICK**

<input type="checkbox"/> REM215	LED Non-Expandable Remote with Lights.....	N/C
	Three drive profiles, control of two seat functions and light operation	

**EXPANDABLE JOYSTICKS**

<input type="checkbox"/> REM216	LED Expandable Remote/Joystick with Lights <sup>(2)</sup> .....	\$875
	Three drive profiles, control of up to 6 actuators and light operation. Not for use with specialty controls	
<input type="checkbox"/> REM400	Color 3.5" Touch Screen Remote/Joystick <sup>(1)</sup> .....	\$1,100
	Four drive profiles, control of up to 6 actuators through the drive control and light operation.	
<input type="checkbox"/> SW400TGL	Toggle Kit for REM400 <sup>(2)</sup> .....	\$200

1. REM400 or REM500 is required when using specialty controls, choose specialty controls purchased separately.  
2. Adds 2 toggle switches to REM400.

**EXPANDABLE SPECIALTY DISPLAY**

<input type="checkbox"/> REM500	LiNX Color 3.5" Touch Screen Display <sup>(1)(2)</sup> .....	\$1,100
---------------------------------	--	---------

Select mounting position for REM500

<input type="checkbox"/> R5MR	Right.....	N/C
<input type="checkbox"/> R5ML	Left.....	N/C
<input type="checkbox"/> R5RPB	Rear Push Bar.....	\$140

1. Driver Control not included.  
2. REM400 or REM500 is required when using specialty controls, choose specialty controls purchased separately.

**SPECIALTY PROPORTIONAL CONTROLS<sup>(1)(2)</sup>**

<input type="checkbox"/> CREM	LiNX Compact Remote.....	\$1,097
<input type="checkbox"/> CREM-LF	LiNX Compact Remote Low Force.....	\$1,097

1. Requires REM400 or REM500  
2. Select Mount in Mounting For Driver Control section controls purchased separately.

**ATTENDANT DRIVE CONTROL CHOICES**

<input type="checkbox"/> ACU	Proportional Attendant Control <sup>(1)</sup> .....	\$1,024
------------------------------	---	---------

Select mounting position for ACU

<input type="checkbox"/> RPBMR	Right
<input type="checkbox"/> RPBML	Left
<input type="checkbox"/> RPBMC	Center

**JOYSTICK MOUNT**

Arm Mounted

<input type="checkbox"/> FR	Fixed Joystick Mounting - Right.....	STD
<input type="checkbox"/> FL	Fixed Joystick Mounting - Left.....	STD
<input type="checkbox"/> MRR	Maxx Resolve Swing-Away Joystick Mount (Right Hand).....	\$450
<input type="checkbox"/> MRL	Maxx Resolve Swing-Away Joystick Mount (Left Hand).....	\$450
<input type="checkbox"/> SAQR	Motion Height Adjustable Swing-Away Quad Link (Right Hand).....	\$244
<input type="checkbox"/> SAQL	Motion Height Adjustable Swing-Away Quad Link (Left Hand).....	\$244

Rear Push Bar Mounted

<input type="checkbox"/> RPBMR	Right.....	\$140
<input type="checkbox"/> RPBML	Left.....	\$140
<input type="checkbox"/> RPBMC	Center.....	\$140

**JOYSTICK KNOB - OPTIONAL**

<input type="checkbox"/> PC101A	BodyPoint 3" U-Shaped.....	\$120
<input type="checkbox"/> PC102A	BodyPoint 4" U-Shaped.....	\$120
<input type="checkbox"/> PC104A	BodyPoint Straight Handle.....	\$120
<input type="checkbox"/> PC107A	BodyPoint Domed Rubber.....	\$120
<input type="checkbox"/> PC110A	BodyPoint Grooved Mushroom.....	\$120

1. To remove the joystick knob from the REM400, pull straight up on the joystick (DO NOT TWIST) otherwise damage may occur and may void the warranty.



# Ultra Low Maxx CG Power Positioning System

## Storm Torque RX20

NOTE: Systems Include Seat Pan, Back Pan, and Standard Armrests. They are designed for use with a 2" - 3" cushion

### POWER POSITIONING SYSTEMS<sup>(1)</sup>

- ULMT Ultra Low Maxx CG Tilt (50°) Only System<sup>(1)</sup> ..... \$6,969
- ULMTRE Ultra Low Maxx CG Tilt (50°) & Power Recline (170°) w/Extended Shear Reduction(ESR)<sup>(1)</sup> ..... \$11,569
- ULMTRP Ultra Low Maxx CG Tilt (50°) & Power Recline (170°) w/Power Sliding Back(PSB)<sup>(1)(2)</sup> ..... \$12,069
- ULMRE Ultra Low Maxx Power Recline (170°) w/Extended Shear Reduction(ESR) ..... \$10,350
- ULMRP Ultra Low Maxx Power Recline (170°) w/Power Sliding Back(PSB)<sup>(2)</sup> ..... \$10,850
- ULME Ultra Low Maxx 12 " Elevating Lift Seat Only<sup>(3)</sup> ..... \$4,995  
Weight capacity up to 250lbs

1. Tilt will be reduced to 45° if adding an Elevating Lift Seat Module
2. Requires Multi-Function M610 through Drive Control (MFCBD).
3. Not available on seat sizes greater than 20" wide by 20" deep.

### ADD: POWER ADJUSTABLE SEAT HEIGHT MODULE

- ESM12 12" Power Adjustable Seat Height Module<sup>(1)(2)(3)</sup> ..... \$2,495  
Weight capacity up to 250lbs

1. Uses an integrated Tilt/Lift mechanism, tilt is 45° and the seat to floor height is not affected by the Elevating Lift Seat.
2. Not available on seat sizes greater than 20" wide by 20" deep
3. For weights greater than 175lbs tilt while elevated is not available.

### MODULES - MUST SELECT SYSTEM TO ADD MODULES

- EPRETILT 5° Electronic Pre-Tilt<sup>(1)</sup> ..... N/C
- PRM8 8° Precline Module<sup>(2)</sup> ..... N/C
- PRM30 30° Precline Module<sup>(2)</sup> ..... \$500
- ERSSP Programable Electronic Recline Function Start & Stop<sup>(3)</sup> ..... \$400

Specify:  
Start: \_\_\_\_\_degrees Stop: \_\_\_\_\_degrees

1. This Pre-tilt option can add approximately 2" (5°) to the front seat to floor height. Multi-Function M610 through Drive Control (MFCBD).
2. Total degrees of recline is 170° ie. if 30° of precline is used, then 140° of recline will be available. Must order Recline System
3. Requires Multi-Function M610 through Drive Control (MFCBD). The range is dependent on the Recline module selected. Will be set to factory settings unless specified.

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

### ADULT SEAT WIDTH SETTINGS

- SW116 16" Wide (Adjustable from 16" to 20") ..... N/C
- SW117 17" Wide (Adjustable from 16" to 20") ..... N/C
- SW118 18" Wide (Adjustable from 16" to 20") ..... N/C
- SW119 19" Wide (Adjustable from 16" to 20") ..... N/C
- SW120 20" Wide (Adjustable from 16" to 20") ..... N/C
- SW219 19" Wide (Adjustable from 19" to 22")<sup>(1)(2)</sup> ..... \$460
- SW220 20" Wide (Adjustable from 19" to 22")<sup>(1)(2)</sup> ..... \$460
- SW221 21" Wide (Adjustable from 19" to 22")<sup>(1)(2)</sup> ..... \$460
- SW222 22" Wide (Adjustable from 19" to 22")<sup>(1)(2)</sup> ..... \$460

1. This seat size is not available on the NARROWFRM.
2. This seat range is not available with ESM12

### ADULT ULTRA RAIL SEAT DEPTH SETTINGS

- SD116 16" Deep (Adjustable from 16" to 21") ..... N/C
- SD117 17" Deep (Adjustable from 16" to 21") ..... N/C
- SD118 18" Deep (Adjustable from 16" to 21") ..... N/C
- SD119 19" Deep (Adjustable from 16" to 21") ..... N/C
- SD120 20" Deep (Adjustable from 16" to 21") ..... N/C
- SD121 21" Deep (Adjustable from 16" to 21") ..... N/C
- SD219 19" Deep (Adjustable from 19" to 23")<sup>(1)(2)</sup> ..... \$460
- SD221 21" Deep (Adjustable from 19" to 23")<sup>(1)(2)</sup> ..... \$460
- SD222 22" Deep (Adjustable from 19" to 23")<sup>(1)(2)</sup> ..... \$460
- SD223 23" Deep (Adjustable from 19" to 23")<sup>(1)(2)</sup> ..... \$460

1. This seat size is not available on the NARROWFRM.
2. This seat range is not available with ESM12

### SMALL ADULT SEAT WIDTH SETTINGS<sup>(1)</sup>

- SW314 14" Wide (Adjustable from 14" to 15") ..... \$250
- SW315 15" Wide (Adjustable from 14" to 15") ..... \$250

1. Weight limit limited to 175 lbs. If recline is selected it may be limited to 165° of recline. You must select both small adult seat width and depth together.

### SMALL ADULT SEAT DEPTH SETTINGS<sup>(1)</sup>

- SD314 14" Deep (Adjustable from 14" to 18") ..... \$250
- SD315 15" Deep (Adjustable from 14" to 18") ..... \$250
- SD316 16" Deep (Adjustable from 14" to 18") ..... \$250
- SD317 17" Deep (Adjustable from 14" to 18") ..... \$250
- SD318 18" Deep (Adjustable from 14" to 18") ..... \$250

1. Weight limit limited to 175 lbs. If recline is selected it may be limited to 165° of recline. You must select both small adult seat width and depth together.

Ultra Low Maxx CG Power Positioning System  
Storm Torque RX20

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**MATRIX SEAT CUSHION OPTIONS**

Matrix PS Seat Cushion 14" - 20" wide/reversable cover  
..... \$569

PS1414     PS1416     PS1616     PS1618     PS1620  
 PS1716     PS1718     PS1720     PS1816     PS1818  
 PS1820     PS2016     PS2018     PS2020

Matrix PS HD Seat Cushion 21"-22" wide  
..... \$975

PS2218     PS2220     PS2222

Matrix PSP Seat Cushion 14" - 20" wide/reversable cover  
..... \$645

PSP1414     PSP1416     PSP1616     PSP1618     PSP1620  
 PSP1716     PSP1718     PSP1720     PSP1816     PSP1818  
 PSP1820     PSP2016     PSP2018     PSP2020

Matrix PSP HD Seat Cushion 21"-22" wide  
..... \$995

PSP2218     PSP2220     PSP2222

Matrix Vi Seat Cushion 14" - 20" wide  
..... \$624

KA1414-VI     MA1415-VI     MA1417-VI     MA1616-VI     MA1618-VI  
 MA1620-VI     MA1716-VI     MA1718-VI     MA1720-VI     MA1816-VI  
 MA1818-VI     MA1820-VI     MA2016-VI     MA2018-VI  
 MA2020-VI

Matrix Vi Seat Cushion 21" - 22" wide  
..... \$681

MA2218-VI     MA2220-VI

Matrix Vi HD Seat Cushion 22" wide  
..... \$840

MA2222-VI

Matrix Libra Seat Cushion 14" - 20" wide  
..... \$745

LC1414     LC1416     LC1516     LC1518     LC1520  
 LC1616     LC1618     LC1620     LC1716     LC1718  
 LC1720     LC1816     LC1818     LC1820     LC1916  
 LC1918     LC1920     LC2016     LC2018     LC2020

Matrix Libra HD Seat Cushion 21" wide  
..... \$1,025

LC2118     LC2120

Matrix Libra Seat Cushion 22" deep  
..... \$1,025

LC1622     LC1722     LC1822     LC1922     LC2022  
 LC2122

Matrix Libra HD Seat Cushion 22" wide  
..... \$1,025

LC2218     LC2220     LC2222

Matrix Libra Adjust Seat Cushion 14" -20" wide  
..... \$985

LC1414-ADJ     LC1416-ADJ     LC1516-ADJ     LC1518-ADJ  
 LC1520-ADJ     LC1616-ADJ     LC1618-ADJ     LC1620-ADJ  
 LC1716-ADJ     LC1718-ADJ     LC1720-ADJ     LC1816-ADJ  
 LC1818-ADJ     LC1820-ADJ     LC1916-ADJ     LC1918-ADJ  
 LC1920-ADJ     LC2016-ADJ     LC2018-ADJ     LC2020-ADJ

Matrix Libra Adjust Seat Cushion 22" deep  
..... \$1,265

LC1622-ADJ     LC1722-ADJ     LC1822-ADJ     LC1922-ADJ  
 LC2022-ADJ     LC2122-ADJ

Matrix Libra Adjust HD Seat Cushion 22" wide  
..... \$1,265

LC2118-ADJ     LC2120-ADJ

Matrix Multi Cushion 14" - 20" wide w/reversable cover  
..... \$745

MM1414     MM1416     MM1616     MM1618     MM1620  
 MM1716     MM1718     MM1720     MM1816     MM1818  
 MM1820     MM2016     MM2018     MM2020

Matrix Multi Cushion 22" wide w/reversable cover  
..... \$995

MM2216     MM2218     MM2220     MM2222

**FABRIC OPTIONS FOR CUSHIONS**

Cushion	Meshtex	Startex	Infection Control	Spacetex	Onyx	Fusion
Libra/Adjust	N/A	STD	\$50	N/A	N/A	N/A
PS/PSP	N/A	STD Reversible	STD Reversible	N/A	N/A	N/A
Vi	STD	\$50	\$50	N/A	\$100	\$150

LC-SP    Spacetex Fabric Option for Libra Cushion ..... N/C

LC-IC    Infection Control Fabric Upgrade for Libra Cushion ..... \$50

FSU-MA    Fusion Fabric Upgrade for Matrix Vi Cushion... \$150

OCU-MA    Onyx Fabric Upgrade for Matrix Vi Cushion .... \$100

STARTEX-S    Smooth side Startex Fabric Upgrade for Matrix Vi Cushion ..... \$50

STARTEX-F    Fabric side Startex Fabric Upgrade for Matrix Vi Cushion ..... \$50

**LAP BELTS**

BELT48    Push Button Style Lap Belt, Length 48"<sup>(1)</sup> ..... STD

BELT60    Push Button Style Lap Belt, Length 60"<sup>(1)</sup> ..... STD

BELT71    Push Button Style Lap Belt, Length 71"<sup>(1)</sup> ..... STD

PLB2    Bodypoint Padded Lap Belt 2 Point<sup>(2)</sup> ..... \$250

PLB4    Bodypoint Padded Lap Belt 4 Point<sup>(2)</sup> ..... \$315

TROBELT    Safety Belt for use only with TRRO & TRRO-E... N/C

1. Not available with TRRO and TRRO-E.

2. 67" length only option

**Ultra Low Maxx CG Power Positioning System**  
Storm Torque RX20

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**MATRX HEADREST PADS<sup>(1)</sup>**

**Standard Pads**

<input type="checkbox"/> ESP6	Matrx Elan Standard Pad 6"	\$310
<input type="checkbox"/> ESP10	Matrx Elan Standard Pad 10"	\$310
<input type="checkbox"/> ESP14	Matrx Elan Standard Pad 14"	\$310
<input type="checkbox"/> ESP6-HUG	Matrx Elan Hug Standard Pad 6"	\$310
<input type="checkbox"/> ESP10-HUG	Matrx Elan Hug Standard Pad 10"	\$310
<input type="checkbox"/> ESP14-HUG	Matrx Elan Hug Standard Pad 14"	\$310
<input type="checkbox"/> SNUG14	Matrx SNUG Standard Pad 14"	\$625
<input type="checkbox"/> SNUG14-A	Matrx SNUG Asymmetrical Pad 14"	\$645
<input type="checkbox"/> SNUG18	Matrx SNUG Standard Pad 18"	\$625
<input type="checkbox"/> SNUG18-A	Matrx SNUG Asymmetrical Pad 18"	\$645

**Standard Pads with Infection Control Cover**

<input type="checkbox"/> ESP6-IC	Matrx Elan Infection Control Pad 6"	\$335
<input type="checkbox"/> ESP10-IC	Matrx Elan Infection Control Pad 10"	\$335
<input type="checkbox"/> ESP14-IC	Matrx Elan Infection Control Pad 14"	\$335
<input type="checkbox"/> ESP6-IC-HUG	Matrx Elan Hug Infection Control Pad 6"	\$335
<input type="checkbox"/> ESP10-IC-HUG	Matrx Elan Hug Infection Control Pad 10"	\$335
<input type="checkbox"/> ESP14-IC-HUG	Matrx Elan Hug Infection Control Pad 14"	\$335
<input type="checkbox"/> SNUG14-IC	Matrx SNUG Standard Infection Control Pad 14"	\$665
<input type="checkbox"/> SNUG14-A-IC	Matrx SNUG Asymmetrical Infection Control Pad 14"	\$685
<input type="checkbox"/> SNUG18-IC	Matrx SNUG Standard Infection Control Pad 18"	\$665
<input type="checkbox"/> SNUG18-A-IC	Matrx SNUG Asymmetrical Infection Control Pad 18"	\$685

**Occipital /4-Point Pads**

<input type="checkbox"/> EOP9	Matrx Elan Occipital Pad 9"	\$341
<input type="checkbox"/> EOP12	Matrx Elan Occipital Pad 12"	\$341
<input type="checkbox"/> E4PTSM	Matrx Elan 4-Point Pad 9"W x 7"H	\$448
<input type="checkbox"/> E4POINT	Matrx Elan 4-Point Pad 11"W x 10"H	\$448

**Occipital / 4-Point Pads with Infection Control Cover**

<input type="checkbox"/> EOP9-IC	Matrx Elan Occipital Pad 9" with Infection Control	\$366
<input type="checkbox"/> EOP12-IC	Matrx Elan Occipital Pad 12" with Infection Control	\$366
<input type="checkbox"/> E4PTSM-IC	Matrx Elan 4-Point Pad 9"W x 7"H with Infection Control	\$473
<input type="checkbox"/> E4POINT-IC	Matrx Elan 4-Point Pad 11"W x 10"H with Infection Control	\$473

1. If mounting on PB Elite must order removable hardware

**MOTION CONCEPTS HEADREST PAD**

<input type="checkbox"/> MCSH	Motion Concepts Standard	N/C
<input type="checkbox"/> MCOH	Motion Concepts Auto Style <sup>(1)</sup>	\$235

1. Comes Standard with Fixed Mounting hardware; for use on Matrx PB or Elite Backs you must select the removable hardware upgrade below.

**HEADREST MOUNTING HARDWARE**

<input type="checkbox"/> FM	Fixed Mounting	N/C
<input type="checkbox"/> LXSHW	Loxx Standard Hardware 15"	\$275
<input type="checkbox"/> LXMMH	Loxx Mini Hardware 11"	\$275
<input type="checkbox"/> MEHWRE	Elan Hardware 15"	\$285
<input type="checkbox"/> MEMH	Elan Mini Hardware 11"	\$285
<input type="checkbox"/> MEHW-EXT	Elan Extended Hardware	\$320
<input type="checkbox"/> FPLXSHW	FLIPP Loxx Standard Hardware 15"	\$705
<input type="checkbox"/> FPLXSHW-EX	FLIPP Loxx Extended Hardware	\$705
<input type="checkbox"/> FPLXMMH	FLIPP Loxx Mini Hardware 11"	\$705
<input type="checkbox"/> FPLXMMH-EX	FLIPP Loxx Mini Extended Hardware	\$705
<input type="checkbox"/> FPMEHW	FLIPP Elan Hardware 15"	\$715
<input type="checkbox"/> FPMEMH	FLIPP Elan Mini Hardware 11"	\$715
<input type="checkbox"/> FPMEHW-EXT	FLIPP Elan Extended Hardware	\$780
<input type="checkbox"/> OHR	OMIT: Headrest	N/C

**STANDARD BACK CANES - FOR TILT SYSTEMS**

<input type="checkbox"/> SBC18	18" Straight Back Canes	STD
<input type="checkbox"/> SBC20	20" Straight Back Canes	STD
<input type="checkbox"/> SBC22	22" Straight Back Canes	STD
<input type="checkbox"/> SBC24	24" Straight Back Canes	STD

Select Cane Angle

<input type="checkbox"/> BA75 75°	<input type="checkbox"/> BA79 79°	<input type="checkbox"/> BA85 85°	<input type="checkbox"/> BA90 90°	<input type="checkbox"/> BA95 95°
<input type="checkbox"/> BA101 101°	<input type="checkbox"/> BA106 106°	<input type="checkbox"/> BA112 112°	<input type="checkbox"/> BA116 116°	

<input type="checkbox"/> ABC20	20" Angled (5°) Back Canes	N/C
<input type="checkbox"/> ABC22	22" Angled (5°) Back Canes	N/C
<input type="checkbox"/> ABC24	24" Angled (5°) Back Canes	N/C

Select Cane Angle

<input type="checkbox"/> BA81 81°	<input type="checkbox"/> BA84 84°	<input type="checkbox"/> BA91 91°	<input type="checkbox"/> BA96 96°	<input type="checkbox"/> BA101 101°
<input type="checkbox"/> BA106 106°	<input type="checkbox"/> BA111 111°	<input type="checkbox"/> BA117 117°	<input type="checkbox"/> BA121 121°	

<input type="checkbox"/> EBC22	22" Endomorph Back Canes <sup>(1)</sup>	\$250
<input type="checkbox"/> EBC24	24" Endomorph Back Canes <sup>(1)</sup>	\$250

Please Select

<input type="checkbox"/> EBCN1	1" Narrower	N/C
<input type="checkbox"/> EBCN2	2" Narrower	N/C
<input type="checkbox"/> EBCN4	4" Narrower	N/C

<input type="checkbox"/> TBC18	18" 10° Mid-Angle Tilt Back Canes <sup>(2)</sup>	\$250
<input type="checkbox"/> TBC20	20" 10° Mid-Angle Tilt Back Canes <sup>(2)</sup>	\$250
<input type="checkbox"/> TBC22	22" 10° Mid-Angle Tilt Back Canes <sup>(2)</sup>	\$250
<input type="checkbox"/> TBC24	24" 10° Mid-Angle Tilt Back Canes <sup>(2)</sup>	\$250

<input type="checkbox"/> FDBC	Fold Down Back Canes w/ Pull Pin <sup>(3)</sup>	\$250
-------------------------------	---	-------

1. Comes standard with Cantilever Armrest Cane Mounted.  
2. Quick Height Adjustment on Backrest and Maxx Spreader Bars not available.  
3. 250 lbs weight capacity. Must select canes above.

**Ultra Low Maxx CG Power Positioning System**  
**Storm Torque RX20**

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**CANE OPTIONS - FOR RECLINE SYSTEMS**

- BFS2 Recline Body Fit 2" Narrower<sup>(1)</sup> .....\$995
- BFS4 Recline Body Fit 4" Narrower<sup>(1)</sup> .....\$995
- AHH Adjustable Hooking Handles<sup>(2)</sup> (pair) ..... \$450
- RPBS Push Handles for Recline - Straight<sup>(3)</sup> (pair) ....\$350
- RPBA Push Handles for Recline - Angled<sup>(3)</sup> (pair) ....\$350

1. Minimum finished back widths of 16"
2. Handle measures 6.5" high and has .75" height adjustment with angle adjustment.
3. 4" of total vertical mounting adjustment.

**FINISHED BACK HEIGHT<sup>(1)(2)</sup> - FOR TILT ONLY/ELEVATE ONLY**

- FBT18 18"
- FBT19 19"
- FBT20 20"
- FBT21 21"
- FBT22 22"
- FBT23 23"
- FBT24 24"
- FBT25 25"

1. 4" Gap from the seat pan to the bottom of back pan is recommended.
2. Finished back height is measured from seat pan to top of backrest.

**FINISHED BACK HEIGHT<sup>(1)(2)</sup> - FOR RECLINE OR TILT/RECLINE**

- FBR20 20"
- FBR21 21"
- FBR22 22"
- FBR23 23"
- FBR24 24"
- FBR25 25"
- FBR26 26"
- FBR27 27"

1. 6" Gap from the seat pan to the bottom of back pan is recommended.
2. Finished back height is measured from seat pan to top of backrest.

**STANDARD BACK OPTIONS**

- SRB Standard Rehab Back & Cushion .....N/C
- SS Additional Super Soft/HR Foam<sup>(1)</sup> ..... \$450
- MFU Matrix Fabric Upcharge ..... \$450
- OBC OMIT: Back Cushion .....N/C
- PIPF Flat Planar Interface/I-Back Plate<sup>(2)(3)(5)</sup> ..... \$125
- PIPR Recessed Planar/I-Back Interface Plate<sup>(2)(3)(5)(6)</sup>... \$125
- ORBC OMIT: Rigid Back & Cushion<sup>(4)</sup>.....N/C

1. Additional foam added may impact seat depth
2. Not available on tilt only or tilt/elevate systems, for use with Recline Module.
3. Required with all Recline Systems when no back is selected. May also be used as an interface for after market backs.
4. Not available with Recline/ESR or Recline/PSB.
5. Not available with TRRO and TRRO-E.
6. Fits Comfort Company, Freedom, Pindot and Moulded Backs.

**MOTION CONCEPTS SINGLE POST BACKREST<sup>(1)(2)(3)</sup>**

- SPBMI Single Post Angle Adjustable Fold Down Back interface-Matrix Back.....\$750 (Folds to 25° forward)
- SPBCI Single Post Angle Adjustable Fold Down Back for Contour/High Back<sup>(4)</sup> .....\$750 (Folds to 25° forward)

1. Single Post Back is an angle adjustable back between 73° and 109° in 6° increments from 17° forward to 19° back. The back folds down and is only available for 16"-20" seat widths & 21"-24" back height. A 4" gap from bottom of back to seat pan is required.
2. This option is not available with Recline and/or a vent tray and has a 250lb weight capacity.
3. Cane mounted armrest are not available with this option.
4. Not available with TRRO/TRRO-E

**UPGRADEABLE BACK OPTIONS**

- Contour Back 15" wide - 20" wide ..... \$495
- CB1512
  - CB1515
  - CB1517
  - CB1612
  - CB1615
  - CB1617
  - CB1621
  - CB1712
  - CB1715
  - CB1717
  - CB1721
  - CB1812
  - CB1815
  - CB1817
  - CB1821
  - CB1912
  - CB1915
  - CB1917
  - CB1921
  - CB2012
  - CB2015
  - CB2017
  - CB2021

- Contour Back 21"-22" wide .....\$595
- CB2115
  - CB2117
  - CB2121
  - CB2215
  - CB2217
  - CB2221

- High Back 16" - 22" W .....\$595
- HB1616
  - HB1618
  - HB1620
  - HB1716
  - HB1718
  - HB1720
  - HB1816
  - HB1818
  - HB1820
  - HB1916
  - HB1918
  - HB1920
  - HB2016
  - HB2018
  - HB2020
  - HB2116
  - HB2118
  - HB2120
  - HB2216
  - HB2218
  - HB2220

- Matrix Elite Back 14" - 20" W (3" Contour)..... \$766
- PBE1414
  - PBE1416
  - PBE1418
  - PBE1420
  - PBE1514
  - PBE1516
  - PBE1518
  - PBE1520
  - PBE1614
  - PBE1616
  - PBE1618
  - PBE1620
  - PBE1714
  - PBE1716
  - PBE1718
  - PBE1720
  - PBE1814
  - PBE1816
  - PBE1818
  - PBE1820
  - PBE1916
  - PBE1918
  - PBE1920
  - PBE2014
  - PBE2016
  - PBE2018
  - PBE2020

- Matrix Elite HD Back 21" - 22" W (3" Contour) .....\$1,095
- PBE2116-HD
  - PBE2120-HD
  - PBE2216-HD
  - PBE2220-HD

- Matrix Elite TR Back 15" - 20" W (3" Contour) ..... \$866
- PBE1516-TR
  - PBE1518-TR
  - PBE1520-TR
  - PBE1616-TR
  - PBE1618-TR
  - PBE1620-TR
  - PBE1716-TR
  - PBE1718-TR
  - PBE1720-TR
  - PBE1816-TR
  - PBE1818-TR
  - PBE1820-TR
  - PBE1916-TR
  - PBE1918-TR
  - PBE1920-TR
  - PBE2016-TR
  - PBE2018-TR
  - PBE2020-TR

- Matrix Elite TR HD Back 21" - 22" W (3" Contour) ..... \$1,151
- PBE2120-TRHD
  - PBE2220-TRHD

- Matrix Elite Deep Back 14" - 20" W (6" Contour) .....\$918
- EDB1416
  - EDB1418
  - EDB1420
  - EDB1514
  - EDB1516
  - EDB1518
  - EDB1520
  - EDB1614
  - EDB1616
  - EDB1618
  - EDB1620
  - EDB1714
  - EDB1716
  - EDB1718
  - EDB1720
  - EDB1814
  - EDB1816
  - EDB1818
  - EDB1820
  - EDB1914
  - EDB1916
  - EDB1918
  - EDB1920
  - EDB2014
  - EDB2016
  - EDB2018
  - EDB2020

Ultra Low Maxx CG Power Positioning System  
Storm Torque RX20

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**UPGRADEABLE BACK OPTIONS**

Matrix Elite Deep HD Back 21" - 22" W (6" Contour) ..... \$1,095

- EDB2116-HD     EDB2118-HD     EDB2120-HD
- EDB2216-HD     EDB2218-HD     EDB2220-HD

Matrix Elite E2S Back 14" - 20" W x 14"/16"H (3" Contour) ..... \$1,125

- E2S1414     E2S1416     E2S1514     E2S1516
- E2S1614     E2S1616     E2S1714     E2S1716
- E2S1814     E2S1816     E2S1914     E2S1916
- E2S2014     E2S2016

Matrix Elite E2S Back 14" - 20" W x 18"/20"H (3" Contour) ..... \$995

- E2S1418     E2S1420     E2S1518     E2S1520
- E2S1618     E2S1620     E2S1718     E2S1720
- E2S1818     E2S1820     E2S1918     E2S1920
- E2S2018     E2S2020

Matrix Elite E2 Deep Back 14" - 20" W x 14"/16"H (6" Contour) ..... \$1,125

- E2D1414     E2D1416     E2D1514     E2D1516
- E2D1614     E2D1616     E2D1714     E2D1716
- E2D1814     E2D1816     E2D1914     E2D1916
- E2D2014     E2D2016

Matrix Elite E2 Deep Back 14" - 20" W x 18"/20"H (6" Contour) ..... \$1,155

- E2D1418     E2D1420     E2D1518     E2D1520
- E2D1618     E2D1620     E2D1718     E2D1720
- E2D1818     E2D1820     E2D1918     E2D1920
- E2D2018     E2D2020

Matrix PB Posture Back 14" - 20" W (5" Contour) ..... \$734

- MPB1412     MPB1612     MPB1616     MPB1620
- MPB1712     MPB1716     MPB1720     MPB1812
- MPB1816     MPB1820     MPB2012     MPB2016
- MPB2020

Matrix PB Heavy Duty Back<sup>(1)</sup> 20" - 22" W (5" Contour) ..... \$1,051

- HDB2016     HDB2020     HDB2216     HDB2220

Matrix PB Posture Back Deep<sup>(1)</sup> 14" - 20" W (7" Contour) ..... \$892

- PBD1412     PBD1612     PBD1616     PBD1620
- PBD1712     PBD1716     PBD1720     PBD1812
- PBD1816     PBD1820     PBD2012     PBD2016
- PBD2020

Matrix PB Heavy Duty Deep Back<sup>(1)</sup> 20"-22" W (5" Contour) ..... \$1,209

- DHDB2016     DHDB2020     DHDB2216     DHDB2220

Matrix MAC Back<sup>(2)</sup> ..... \$1,495

- MAC1616     MAC1618     MAC1620     MAC1622
- MAC1716     MAC1718     MAC1720     MAC1722
- MAC1816     MAC1818     MAC1820     MAC1822
- MAC1916     MAC1918     MAC1920     MAC1922
- MAC2016     MAC2018     MAC2020     MAC2022

Matrix MAC Back w/ HUG<sup>(2)</sup> ..... \$1,645

- MAC1616-HUG     MAC1618-HUG     MAC1620-HUG
- MAC1622-HUG     MAC1716-HUG     MAC1718-HUG
- MAC1720-HUG     MAC1722-HUG     MAC1816-HUG
- MAC1818-HUG     MAC1820-HUG     MAC1822-HUG
- MAC1916-HUG     MAC1918-HUG     MAC1920-HUG
- MAC1922-HUG     MAC2016-HUG     MAC2018-HUG
- MAC2020-HUG     MAC2022-HUG

Matrix MAC Lateral Contour

- MLCSL    Left 3.5"     MLCSR    Right 3.5"
- MLCDL    Left 6"     MLCDR    Right 6"
- MLCXDL    Left 9"     MLCXDR    Right 9"

1. The Heavy Duty PB Back and PB Deep Back are not available on Recline systems with seat widths of 21" and 22"
2. Must select Pad Contour in Matrix MAC Lateral Contour section.

**OTHER UPGRADABLE BACK OPTIONS**

Cushion	Meshtex	4-Way Stretch Polyester	Infection Control	Onyx	Fusion
Contour	N/A	STD	N/A	\$100	N/A
High	STD	N/A	N/A	\$100	\$100
E2, Elite, Elite HD, Elite Deep	STD	N/A	\$100	N/A	N/A
E2, Elite TR, Elite TR HD	STD	N/A	\$100	N/A	N/A

**OTHER UPGRADABLE BACK OPTIONS**

- EBSS    Super Soft Foam for Elite Back Upgrade..... \$150
- FFUHB    Fusion Fabric Upgrade for High Backs..... \$150
- OCUCB    Onyx Fabric Upgrade for Contour Backs..... \$150
- OCUHB    Onyx Fabric Upgrade for High Backs..... \$150
- CHIEB    Chest Harness Interface Kit for Elite Back<sup>(1)(2)</sup>..... \$103
- CHIPB    Chest Harness Interface Kit for PB Back<sup>(1)(2)</sup>..... \$52
- E2CHI-S    Chest Harness Interface Kit for E2 Back for 14"-15" Backs<sup>(1)(2)</sup>..... \$85
- E2CHI-M    Chest Harness Interface Kit for E2 Back for 16"-20" Backs<sup>(1)(2)</sup>..... \$85
- E2BKB    Contoured anatomical PSIS/Sacrolumbar Support..... \$75

1. Does not include the chest strap.
2. This Interface Kit is only Required for 4 point harness

Ultra Low Maxx CG Power Positioning System  
Storm Torque RX20

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**CHEST STRAPS - ABDOMINAL SUPPORT**

- MFUS BP Monoflex, Underarm RLS S(3" x 15.5") ..... \$195
- MFUM BP Monoflex, Underarm RLS M(4" x 18.25") ..... \$195
- MFUL BP Monoflex, Underarm RLS L(4.5" x 21.5") ..... \$195
- MFCS BP Monoflex, Center RLS S(2.5" x 15.5") ..... \$195
- MFCM BP Monoflex, Center RLS M(3.25" x 18.25") ..... \$195
- MFCL BP Monoflex, Center RLS L(3.5" x 21.5") ..... \$195
- PFLEX-M PostureFlex Abdominal Support  
Height 6.5", length 16" to 23"<sup>(1)</sup> ..... \$245
- PFLEX-L PostureFlex Abdominal Support  
Height 7", length 20" to 32"<sup>(1)</sup> ..... \$245
- QRTAB Quad release tab for easy independent  
application and removal of PostureFlex ..... \$35

1. Only available on the E2 and E2 Deep Backs.

**TILT ARMRESTS**

Adult Dual Post

- DP1 9.75" - 12.75" ..... STD
- DP2 12.75" - 15.75" ..... STD
- DP3 15.75" - 18.75" ..... STD
- FP Flat Pouch ..... STD
- GB Glove Boxes ..... N/C
- OAL Outback Arm Option Left ..... \$195
- OAR Outback Arm Option Right ..... \$195

Ultra Rail Seat Mount Flip Back Cantilever MAXX Tilt Arm<sup>(1)</sup>

- RMCAN 9" - 12.5" ..... \$488
- RMCAN2 11" - 14" ..... \$488
- DPM Dual Post Module to RMCAN (pair) ..... \$250

Outback Cantilever Tilt Armrest

- OCTA1 9"-12" ..... \$488
- OCTA2 12"-15" ..... \$488
- OCTA2 6"-12" ..... \$488

Cantilever Arms for Endomorph Back Canes<sup>(1)</sup>

- CAFEB See page 4 ..... STD

1. Not available with Small Adult Sizes.

**RECLINE ARMRESTS**

Reclining Adjustable Height

- RC1 9.5"-13" ..... STD
- RC2 13"-16" ..... STD
- OAL Outback Arm Option Left ..... \$195
- OAR Outback Arm Option Right ..... \$195

Maxx Style Cane Mounted Height, Width, Angle Adj.

Cantilever Flip Back Armrests<sup>(1)</sup>

- MCAN1 9.5"-13" ..... \$250
- MCAN2 12.5"-16" ..... \$250

Outback Cantilever Recline Armrest

- OCRA1 9"-12" ..... \$350
- OCRA2 12"-15" ..... \$350
- OCRA3 7" - 9" ..... \$350

Cane Mounted Height and Angle Adj.

Cantilever Flip Back Armrest<sup>(1)</sup>

- CANQ1 9.5"-12" ..... \$350
- CANQ2 10"-13" ..... \$350
- CANQ3 13"-16" ..... \$350

Adult Dual Post On Recline

- DPOR1 9.75" - 12.75" ..... \$250
- DPOR2 12.75" - 15.75" ..... \$250
- DPOR3 15.75" - 18.75" ..... \$250
- FP Flat Pouch ..... STD
- GB Glove Boxes ..... N/C

1. Not available with Small Adult Sizes.

**ARMREST TUBE LENGTH - PLEASE SPECIFY**

- ATFL Left Full ..... STD
- ATDL Left Desk ..... STD
- ATFR Right Full ..... STD
- ATDR Right Desk ..... STD

**ARMPADS<sup>(5)</sup>**

Standard

- SAPFL Left Full (2.25" x 14") ..... STD
- SAPDL Left Desk (2.25" x 10") ..... STD
- SAPML Left Mini (2.25" x 7") ..... STD
- SAPFR Right Full (2.25" x 14") ..... STD
- SAPDR Right Desk (2.25" x 10") ..... STD
- SAPMR Right Mini (2.25" x 7") ..... STD

Waterfall Arm Pads<sup>(1)</sup>

- WAPFL Left Full (13") ..... N/C
- WAPDL Left Desk (9") ..... N/C
- WAPFR Right Full (13") ..... N/C
- WAPDR Right Desk (9") ..... N/C

Modular Arm Pad - Left<sup>(2)</sup>

- MAPTFL Left Tray Full (14") ..... \$95
- MAPTDL Left Tray Desk (10") ..... \$95

Select Pad Insert below:

- ISSFL Iskin - Integrated soft skin foam
- SCVFL Startex Covered Visco foam

Modular Arm Pad - Right<sup>(2)</sup>

- MAPTFR Right Tray Full (14") ..... \$95
- MAPTDR Right Tray Desk (10") ..... \$95

Select Pad Insert below:

- ISSFR Iskin - Integrated soft skin foam
- SCVFR Startex Covered Visco foam

Flat Multi Position Arm Pads

- MPAFL Left Full (14") ..... \$150
- MPADL Left Desk (10") ..... \$150
- MPAML Left Mini (7") ..... \$150
- MPAFR Right Full (14") ..... \$150
- MPADR Right Desk (10") ..... \$150
- MPAMR Right Mini (7") ..... \$150

Gel Pads

- GAPNL Narrow Left Desk (2" x 12") ..... \$88
- GAPNR Narrow Right Desk (2" x 12") ..... \$88
- GAPWL Wide Left Desk (3.5" x 12") ..... \$96
- GAPWR Wide Right Desk (3.5" x 12") ..... \$96
- GAPLL Long Left Full (3.5" x 14") ..... \$103
- GAPLR Long Right Full (3.5" x 14") ..... \$103

Arm Troughs

- EATL Ergonomic Arm Trough - Left ..... \$121
- EATR Ergonomic Arm Trough - Right ..... \$121
- FHPL Flat Hand Pad - Left<sup>(3)</sup> ..... \$91
- FHPR Flat Hand Pad - Right<sup>(3)</sup> ..... \$91

**Ultra Low Maxx CG Power Positioning System**  
**Storm Torque RX20**

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**ARMPADS<sup>(5)</sup> - CONTINUED**

**Arm Pad Storage/Cup Holder<sup>(4)</sup>**

- APSL Left ..... \$295
- APSR Right ..... \$295
- OAP Omit Armpads(7/8" diameter tube) ..... N/C

1. Full length pad not compatible with the desk length armrest tube if using Ultra Rail mounted Tilt Arms.
2. Requires interchangeable pad insert.
3. Must order with ergonomic arm trough.
4. Not available on the joystick side, or with Outback Cantilever arms.
5. Full length arm pads used on the same side as the Joystick result in a "Back to Joystick" measurement of 20" minimum (24" including switch mounting)

**ARM ADDUCTOR PAD<sup>(1)(2)</sup>**

- AADSL Left Small (6" x 2.25") ..... \$140
- AADSR Right Small (6" x 2.25") ..... \$140
- AADLL Left Large (8.5" x 2.25") ..... \$140
- AADLR Right Large (8.5" x 2.25") ..... \$140

1. Only available with Modular or Flat Multi Position Arm Pads.
2. 1" thick pad, unlimited fore/aft adjustment.

**ELBOW BLOCKS**

**Left Elbow Block<sup>(1)</sup>**

Select Pad Size and Both Extensions:

- EBXL Size XS (3.25" x 4.25") ..... \$195
- EBSL Size S (3.5" x 5.25") ..... \$195
- EBML Size M (3.75" x 5.75") ..... \$195
- EBLL Size L (4.25" x 6.5") ..... \$195
  
- LV3 Left Vertical Extension 3"
- LV5 Left Vertical Extension 5"
- LH5 Left Horizontal Extension 5"
- LH7 Left Horizontal Extension 7"

**Right Elbow Block<sup>(1)</sup>**

Select Pad Size and Both Extensions:

- EBXR Size XS (3.25" x 4.25") ..... \$195
- EBSR Size S (3.5" x 5.25") ..... \$195
- EBMR Size M (3.75" x 5.75") ..... \$195
- EBLR Size L (4.25" x 6.5") ..... \$195
  
- RV3 Right Vertical Extension 3"
- RV5 Right Vertical Extension 5"
- RH5 Right Horizontal Extension 5"
- RH7 Right Horizontal Extension 7"

1. Prevents arm from slipping.

**ADDITIONAL ARM ACCESSORIES**

- MACESL Multi Axis Upper Extremity Support<sup>(1)</sup>(each)  
Left Side ..... \$225
- MACESR Multi Axis Upper Extremity Support<sup>(1)</sup> (each)  
Right Side ..... \$225
- OAS2L 1" Offset Armrest Spacer Left(each) ..... \$125
- OAS2R 1" Offset Armrest Spacer Right(each) ..... \$125
- OAS1L 1/2" Offset Armrest Spacer Left(each) ..... \$125
- OAS1R 1/2" Offset Armrest Spacer Right(each) ..... \$125
- MPAPA Adapter Block<sup>(2)</sup> (each) Quantity ( ) ..... \$35

1. Allows rotation of the armpad, only compatible with Gel Arm pad, Flat Arm pad or Ergonomic Arm Troughs.
2. Not available with Standard Arm Pads or Waterfall Arm Pads

**MAXX LATERAL TRUNK SUPPORTS<sup>(1)</sup>**

**Lateral Trunk Support Left**

- LATXL Size XS (3.25" x 4.25") ..... \$100
- LATSL Size S (3.5" x 5.25") ..... \$100
- LATML Size M (3.75" x 5.75") ..... \$100
- LATLL Size L (4.25" x 6.5") ..... \$100

Add:

- SALATL MAXX Style Swing Away Hardware - Left<sup>(2)</sup> ..... \$195
- ATLL Additional Telescoping Link - Left<sup>(3)</sup> ..... \$75

**Lateral Trunk Support Right**

- LATXR Size XS (3.25" x 4.25") ..... \$100
- LATSR Size S (3.5" x 5.25") ..... \$100
- LATMR Size M (3.75" x 5.75") ..... \$100
- LATLR Size L (4.25" x 6.5") ..... \$100

Add:

- SALATR MAXX Style Swing Away Hardware - Right<sup>(2)</sup> ..... \$195
- ATLR Additional Telescoping Link - Right<sup>(3)</sup> ..... \$75

1. Comes standard with Fixed Mount Hardware. Not available with Matrix PB, PB Elite Deep, and PB Deep Backs.
2. Comes standard with 1" telescoping link which gives 0-2.75" medial offset.
3. The amount of offset available may be impacted by back thickness.

**MATRIX FIXED LATERAL TRUNK SUPPORTS<sup>(1)</sup>**

- E2FLL Matrix E2 Fixed Laterals-Left(6.5x5")<sup>(2)</sup> ..... \$195
- E2FLR Matrix E2 Fixed Laterals-Right(6.5x5")<sup>(2)</sup> ..... \$195
- PBSSL Matrix Small Lateral Fixed-Left (3x5.5") ..... \$113
- PBSSR Matrix Small Lateral Fixed-Right (3x5.5") ..... \$113
- PBLSL Matrix Standard Lateral Fixed-Left (7.5x5") ..... \$116
- PBLSR Matrix Standard Lateral Fixed-Right (7.5x5") ..... \$116
- PBOLL Matrix Offset Lateral Fixed Left (7.5x5") ..... \$125
- PBOLR Matrix Offset Lateral Fixed Right (7.5x5") ..... \$125

1. Not available with Matrix PB, PB Elite Deep, and PB Deep Backs.
2. Only available with E2 backs.

**MATRIX SWING AWAY LATERAL TRUNK SUPPORTS<sup>(1)</sup>**

**Matrix Lateral Trunk Support Left**

- EPADXL Size XS (3.25" x 4.25") ..... \$115
- EPADSL Size S (3.5" x 5.25") ..... \$115
- EPADML Size M (3.75" x 5.75") ..... \$115
- EPADLL Size L (4.25" x 6.5") ..... \$115

Must select one option:

- SAHEBL Matrix Elite Swing Away Hardware Left ..... \$250
- OSAHL Matrix Offset Elite Swing Away Hardware Left ..... \$250

**Matrix Lateral Trunk Support Right**

- EPADXR Size XS (3.25" x 4.25") ..... \$115
- EPADSR Size S (3.5" x 5.25") ..... \$115
- EPADMR Size M (3.75" x 5.75") ..... \$115
- EPADLR Size L (4.25" x 6.5") ..... \$115

Must select one option:

- SAHEBR Matrix Elite Swing Away Hardware Right ..... \$250
- OSAHR Matrix Offset Elite Swing Away Hardware Right ..... \$250

1. Not available with Matrix PB, PB Elite Deep, and PB Deep Backs.

**Ultra Low Maxx CG Power Positioning System**  
**Storm Torque RX20**

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**HIP AND KNEE SUPPORTS**

**Lateral Hip Support Left**

- HSPXLH Size XS (5.5" W x 4" H) ..... \$120
- HSPSMLH Size S/M (7" W x 4" H) ..... \$120
- HSPSLH Size S (7"W x 6"H) ..... \$120
- HSPMLH Size M (9" W x 4"H) ..... \$120
- HSPLLH Size L (13" W x 4"H) ..... \$120

**Lateral Hip Support Right<sup>(1)</sup>**

- HSPXRH Size XS (5.5"W x 4"H) ..... \$120
- HSPSMRH Size S/M (7"W x 4"H) ..... \$120
- HSPSRH Size S (7"W x 6"H) ..... \$120
- HSPMRH Size M (9"W x 4"H) ..... \$120
- HSPLRH Size L (13"W x 4"H) ..... \$120

**Lateral Knee Support Left<sup>(1)</sup>**

- HSPXLK Size XS (5.5"W x 4" H) ..... \$120
- HSPSMLK Size S/M (7"W x 4" H) ..... \$120
- HSPSLK Size S (7"W x 6" H) ..... \$120
- HSPMLK Size M (9"W x 4"H) ..... \$120
- HSPLLK Size L (13"W x 4"H) ..... \$120

**Lateral Knee Support Right<sup>(1)</sup>**

- HSPXRK Size XS (5.5"W x 4"H) ..... \$120
- HSPSMRK Size S/M (7"W x 4"H) ..... \$120
- HSPSRK Size S (7"W x 6"H) ..... \$120
- HSPMRK Size M (9"W x 4"H) ..... \$120
- HSPLRK Size L (13"W x 4"H) ..... \$120

Add Removable Hardware:

**MAXX Style Quick Release, Multi-Axis Mounting**

**Hardware for Hip Supports**

- MHSHLH Left mounted at Hip ..... \$215
- MHSHLK Left mounted at Knee ..... \$215
- MHSHRH Right mounted at Hip ..... \$215
- MHSHRK Right mounted at Knee ..... \$215

**Swing Away Hardware for Hip Supports**

- SHSHLH Left mounted at Hip ..... \$215
- SHSHLK Left mounted at Knee ..... \$215
- SHSHRH Right mounted at Hip ..... \$215
- SHSHRK Right mounted at Knee ..... \$215

**Extended Maxx Style Quick Release, Multi-Axis Mounting**

**Hardware for Hip Supports<sup>(2)</sup>**

- EMQHLH Left mounted at Hip ..... \$215
- EMQHLK Left mounted at Knee ..... \$215
- EMQHRH Right mounted at Hip ..... \$215
- EMQHRK Right mounted at Knee ..... \$215

**Lift Off Hardware for Hip Support**

**(Hip Position Only)**

- LOHL Left ..... \$195
- LOHR Right ..... \$195

1. Comes standard with Fixed Mount Hardware.

2. Recommended when finished seat cushion thickness will be greater than 3.5" high.

**POWER CENTER MOUNT FOOT PLATFORM**

MUST BE COMPLETED	Measurement in Inches
Knee to Heel:	
Cushion Thickness: (Subtract)	
Cushion Compression:(Add)	
Seat Pan to Footplate:(Total)	

**LNx Power Center Mount Foot Platform**

- LN<sub>X</sub> 9" Articulation<sup>(1)(2)</sup> ..... \$3,000
  - LSPTF1 10" - 14" Seat pan to footplate
  - LSPTF2 13" - 19" Seat pan to footplate
- LN<sub>XA</sub> 6" Articulation<sup>(1)(2)(6)</sup> ..... \$3,000
  - LASPTF2 10" - 14" Seat pan to footplate

Options:

- LIFPN Rubber coated footplate 10.5" W x 10" D ..... STD
- LICP Individual calf pads ..... STD
- OPCP One Piece Malleable Calf Panel<sup>(3)</sup> ..... \$250
- LN<sub>X</sub>PF Plantar Flex Option for LN<sub>X</sub> Foot platform<sup>(4)</sup> ..... \$250

**Latitude Style Foot Platform**

- LSFP Latitude Style Foot Platform<sup>(5)</sup> ..... \$2,995  
For finished knee to heel measurements of 6.5"-12"
- LSPFP Latitude Style Padded Foot Plate ..... \$250

1. 3" ground clearance required.
2. Only available on adult range seat systems 16" deep or greater.  
For shorter seat depths call Customer Service for other options.
3. Not available with Plantar Flex Option for LN<sub>X</sub> Foot platform, or with LN<sub>XA</sub>.
4. Only available with Foot Platform Options.
5. Not available with TRRO and TRRO-E.
6. For clients under 4' only.

**MANUAL ELEVATING CENTER MOUNT LEGREST**

- MECML Manual Elevating Center Mount Legrest  
With Plantar Flex ..... \$595

Select legrest set up position:

- 90SU Set at 90° (elevating range 90-40°) ..... STD

Rubber Coated Footplate:

- MRCP Rubber Coated Footplate 10.5" W x 10" D ..... N/C

Select seat pan to footplate for One Piece Foot Plate

- CSPTF1 11.5" - 15" ..... N/C
- CSPTF2 14" - 19" ..... N/C

Individual Footplates:


- IFSL Small - (5" x 7.5" deep) - Left ..... \$165
- IFSR Small - (5" x 7.5" deep) - Right ..... \$165
- IFML Medium - (5.5" x 9.5" deep) - Left ..... \$165
- IFMR Medium - (5.5" x 9.5" deep) - Right ..... \$165
- IFLL Large - (6" x 11.5" deep) - Left ..... \$165
- IFLR Large - (6" x 11.5" deep) - Right ..... \$165

Select seat pan to footplate for Individual Footplates

- CSPTF3 9.5" - 14.5" ..... N/C
- CSPTF4 12" - 19" ..... N/C

Optional:

- CICP Individual calf pads ..... \$200



**Ultra Low Maxx CG Power Positioning System**  
**Storm Torque RX20**

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_


Tag: \_\_\_\_\_

**MANUAL CENTER MOUNT FRONT RIGGING**

- MSFCMFP Maxx Style Fixed Center Mount Foot Platform<sup>(1)(2)</sup> .....\$275
- MIFPN Comes with a rubber coated footplate 10.5" W x 10" D. ....N/C
- MICP Individual calf pads .....\$100


Select seat pan to footplate:

- MSPTF1 9" - 13" .....N/C
- MSPTF14 13" -19" .....N/C



1. Depth adjustable, and easy slide adjustment for seat pan to foot plate extension.
2. 3" ground clearance required.


**FOOT PLATFORM OPTIONS FOR FIXED/LNX CMFP**



- SFFP Small Foot Platform (12"w x 8.5"d) .....\$180
- MFFP Medium Platform (12"w x10.75"d) .....\$180
- LFFP Large Platform (14.5"w x13"d) .....\$250
- XLFFP Extra Large Platform (17"w x13"d) .....\$250
- XXLFFP Extra Extra Large Platform (21"w x13"d) .....\$250

1. May need to adjust leg position to avoid interference with front casters

**INDIVIDUAL FOOT PLATES FOR FIXED/LNX CMFP**



- IFSL Small (5" W x 7.5" D) - Left .....\$165
- IFSR Small (5" W x 7.5" D) - Right .....\$165
- IFML Medium (5.5" W x 9.5" D) - Left .....\$165
- IFMR Medium (5.5" W x 9.5" D) - Right .....\$165
- IFLL Large (6" W x 11.5" D) - Left .....\$165
- IFLR Large (6" W x 11.5" D) - Right .....\$165


Add:

- TSSL Left Toe Strap Small 6.5" .....\$125
- TSSR Right Toe Strap Small 6.5" .....\$125
- TSML Left Toe Strap Medium 8.5" .....\$125
- TSMR Right Toe Strap Medium 8.5" .....\$125
- TSLL Left Toe Strap Large 10.5" .....\$125
- TSLR Right Toe Strap Large 10.5" .....\$125

- MFOL Multi-Axis Footplate (6"W x 10"D) - Left .....\$225
- MFOR Multi-Axis Footplate (6"W x 10"D) - Right .....\$225

1. Foot plates may interfere with casters in certain configurations

**MAXX INDEPENDENT ADJUSTABLE POSITIONING**



- 5.5" - 6" WIDTH ADJUSTMENT
- 9.5" - 11.5" DEPTH ADJUSTMENT
- UP TO 9° ADDUCTION
- UP TO 45° ABDUCTION

- MPFPL Maxx Independent Adjustable Positioning Foot plate<sup>(1)</sup> - Left .....\$200
- MPFPR Maxx Independent Adjustable Positioning Foot plate<sup>(1)</sup> - Right .....\$200
- MTSL Left Toe Strap (each) .....\$35
- MTSR Right Toe Strap (each) .....\$35

1. Foot plates may interfere with casters in certain configurations.

**MAXX STYLE POWER ELEVATING LEGRESTS**

- MSPPL MAXX Swing Away Power Elevating<sup>(1)</sup> .....\$2,490

Select Power Legrest Function:

- INDA Individual Movement .....N/C
- COMA Combined Movement .....N/C
- SYNA Synchronized/Recline Movement.....N/C

Select Receivers:

- MRS Straight Receivers .....N/C
- MRF Flared 1" Outward Receivers .....N/C
- MRF2 Flared 2" Outward Receivers.....\$250

Select Seat Pan to Footplate:

- PSPTF5 Length: 10"-14" .....N/C
- PSPTF2 Length: 13"-17" .....N/C
- PSPTF7 Length: 16"-20" .....N/C

Optional:

- PMCP Moulded I-Skin Calf Pads ..... STD
- PFCPMS Flat Calf Pad Option (pair) .....\$275
- PCCPMS Curved Calf Pad Option<sup>(2)</sup> (pair) .....\$275

1. Must choose footplate on page 13.
2. Not available on 14 and 15" seat widths. For 16" and 17" seat widths you must select flared receivers.

**MAXX STYLE MANUAL ELEVATING LEGRESTS**

- MMPPSA MAXX Swing Away Manual Elevating<sup>(1)</sup> .....\$680

Select Receivers:

- ERS Straight Receivers .....N/C
- ERF Flared 1" Outward Receivers .....N/C
- ERF2 Flared 2" Outward Receivers.....\$250

Select Seat Pan to Footplate:

- ESPTF5 Length: 10"-14" .....N/C
- ESPTF2 Length: 13"-17" .....N/C
- ESPTF7 Length: 16"-20" .....N/C

Optional:

- EMCP Moulded I-Skin Calf Pads ..... STD
- EFCPMS Flat Calf Pad Option (pair) .....\$275
- ECCPMS Curved Calf Pad Option<sup>(2)</sup> (pair) .....\$275

1. Must choose footplate on page 13.
2. Not available on 14 and 15" seat widths. For 16" and 17" seat widths you must select flared receivers.

Ultra Low Maxx CG Power Positioning System  
Storm Torque RX20

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**SWING AWAY LEGREST**

HD70 Heavy Duty 70° Swing Away Footrest(pair)<sup>(1)(2)</sup> . \$495

Select Receivers:

- RS Straight Receivers .....N/C
- RF Flared 1" Outward Receivers .....N/C
- RF2 Flared 2" Outward Receivers.....\$250

Select Seat Pan to Footplate:

- HSPTF1 Length: 6"-10" .....N/C
- HSPTF8 Length: 10"-14" .....N/C
- HSPTF2 Length: 12"-16" .....N/C
- HSPTF7 Length: 15"-19" .....N/C

Options:

- HD70FCP Flat Calf Pad Option (pair) .....\$275
- HD70CCP Curved Calf Pad Option<sup>(3)</sup> (pair) .....\$275

I70SAFT Invacare 70° Swingaway Footrests.....\$325  
 I70-CFT Comes with composite foot plates.....N/C

RLSL Residual Limb Support - Left (8.5" x 6")<sup>(4)</sup>  
 .....\$495

RLSR Residual Limb Support - Right (8.5" x 6")<sup>(4)</sup>  
 .....\$495

1. Must choose footplate on page 13.
2. Heavy duty durable construction (4.5lbs ea.).
3. Not available on 14 and 15" seat widths. For 16" and 17" seat widths you must select flared receivers.
4. Swing away and removable, 50° of angle adjustment, 4.5" height adjustment, 3" horizontal and 2" of for/aft adjustment.

**DUAL FRONT RIGGING ADAPTOR<sup>(1)</sup>**

DFRA Dual Front Rigging Receivers..... \$450

1. Must purchase each set of front rigging separately. Not available on no charge front rigging options.

**FOOTPLATES - SWING AWAY/ELEVATING LEGRESTS<sup>(1)</sup>**




Angle Adjustable Footplates

- AFSL Left Small Footplate (4.25"W x 8"D) ..... \$165
- AFSR Right Small Footplate (4.25"W x 8"D) ..... \$165
- AFML Left Medium Footplate (5.25"W x 8"D) ..... \$165
- AFMR Right Medium Footplate (5.25"W x 8"D) ..... \$165
- AFLL Left Large Footplate (6.25"W x 8"D) ..... \$165
- AFLR Right Large Footplate (6.25"W x 8"D) ..... \$165
- AFXL Left Extra Large Footplate (9"W x 11"D) ..... \$165
- AFXR Right Extra Large Footplate (9"W x 11"D) ..... \$165

Multi-Axis Adjustable Angle Footplates

- MFSL Left Small Footplate (4.25"W x 8"D) ..... \$225
- MFSR Right Small Footplate (4.25"W x 8"D) ..... \$225
- MFML Left Medium Footplate (5.25"W x 8"D) ..... \$225
- MFMR Right Medium Footplate (5.25"W x 8"D) ..... \$225
- MFLL Left Large Footplate (6.25"W x 8"D) ..... \$225
- MFLR Right Large Footplate (6.25"W x 8"D) ..... \$225
- MFXL Left Extra Large Footplate (9"W x 11"D) ..... \$225
- MFXR Right Extra Large Footplate (9"W x 11"D) ..... \$225



Cut Out Angle Adjustable Footplates

- AFCSL Left Small Footplate (7.25"W x 9"D) ..... \$165
- AFCSR Right Small Footplate (7.25"W x 9"D) ..... \$165
- AFCML Left Medium Footplate (7.5"W x 11"D) ..... \$165
- AFCMR Right Medium Footplate (7.5"W x 11"D) ..... \$165
- AFCLL Left Large Footplate (10"W x 14"D) ..... \$165
- AFCLR Right Large Footplate (10"W x 14"D) ..... \$165

Cut Out Multi-Axis Adjustable Angle Footplates

- MFCSL Left Small Footplate (7.25"W x 9"D) ..... \$225
- MFCSR Right Small Footplate (7.25"W x 9"D) ..... \$225
- MFCLL Left Medium Footplate (7.5"W x 11"D) ..... \$225
- MFCMR Right Medium Footplate (7.5"W x 11"D) ..... \$225
- MFCLL Left Large Footplate (10"W x 14"D) ..... \$225
- MFCLR Right Large Footplate (10"W x 14"D) ..... \$225

1. Certain footplate sizes may not fit all seat widths depending on legrests selected.

**Ultra Low Maxx CG Power Positioning System**  
**Storm Torque RX20**

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

**HEEL LOOPS, CUPS AND TOE STRAPS**

- HLL Left Heel Loop (each) .....\$20
- HLR Right Heel Loop (each) .....\$20
- HCL Left Cup Loop (each) .....\$20
- HCR Right Cup Loop (each) .....\$20
- TSL Left Toe Strap (each) .....\$20
- TSR Right Toe Strap (each) .....\$20

**CALF STRAPS AND CALF PANELS<sup>(1)</sup>**

- CSM Body Point Airmesh Calf Strap for Chairs  
15"-18" wide with 3.5" H Pad<sup>(1)</sup> ..... \$200
- CSL Body Point Airmesh Calf Strap for Chairs  
18"-22" wide with 3.5" H Pad<sup>(1)</sup> ..... \$200
- CPM Body Point Airmesh Calf Panel for Chairs  
15"-18" wide with 35.5" W x 9" H Pad<sup>(1)</sup> .....\$295
- CPL Body Point Airmesh Calf Panel for Chairs  
18"-22" wide with 38.5" W x 9" H Pad<sup>(1)</sup> ..... \$315

1. Not available for Center Mounted Front Rigging.

**STANDARD FRONT RECEIVERS**

- PHDR Pride HD Drop-in Receivers Only .....N/C
- PSAR Pride Style #8 Swing Away Receivers Only .....N/C
- IAR Invacare 70° Swing Away Receivers Only .....N/C
- Q70R Quickie 70° Swing Away Receivers Only .....N/C
- Q70HDR Quickie 70HD Swing Away Receivers Only .....N/C

**SINGLE-DUAL FUNCTION ELECTRONICS**

- SFCB Single/Dual Function Electronics<sup>(1)</sup> ..... STD  
 Select Switch:  
 SPBSO Single Push Button  
 STO Single Toggle  
 DPBO Dual Push Button  
 OSO Omit Switch
- SFCBA Single/Dual Function Drive Control<sup>(1)</sup> ..... \$750  
 Select Control:  
 JOY Through Joystick  
 SC Through Specialty Controls
- OSO Reduced Drive for Elevating Seat<sup>(2)</sup> ..... \$400

1. These electronics include drive lock-out and tilt limit.  
 2. If reduced drive is not selected the system will be in drive lock-out when elevated.

**MULTI- FUNCTION ELECTRONICS**

- MFCB Multi-Function through Switch<sup>(1)</sup> ..... \$1,950  
 Select Switch:  
 QPBO Quad Push Button .....N/C  
 FWT Four-Way Toggle Switch ..... \$150  
 EWR 8-Way Rocker<sup>(2)</sup> .....\$205
- MFCBD Multi-Function M610 Through  
 Drive Control<sup>(1)</sup> ..... \$1,950  
 Includes reduced drive for Elevating Seat  
 Select Control:  
 JOY Through Joystick  
 SC Through Specialty Controls  
 Optional:  
 RDEY Include Reduced Drive for Elevating Seat<sup>(3)</sup>  
 RDEN Do Not Include Reduced Drive for Elevating Seat<sup>(3)</sup>  
 Optional:  
 MFMSR Multi Function Memory Seating .....N/C

1. These electronics include attendant control, drive lock-out, tilt/recline limit and actuator speed control.  
 2. Dual function is not available with 8-way Rocker.  
 3. If reduced drive is not selected system will be in drive lock-out when elevated.

**ELECTRONIC SWITCH MOUNTING**

- ESML Left Side (Armpad)
- ESMR Right Side (Armpad)
- ESMJL Joystick Left (Mounts to Joystick Bracket)
- ESMJR Joystick Right (Mounts to Joystick Bracket)
- ESMD Do not mount

# Ultra Low Maxx CG Power Positioning System

## Storm Torque RX20

Dealer Name: \_\_\_\_\_

Dealer Account Number: \_\_\_\_\_

PO #: \_\_\_\_\_

Tag: \_\_\_\_\_

### ELECTRONICS OPTIONS

#### Additional Switches<sup>(4)</sup>

<input type="checkbox"/> SPBSA	Single Push Button .....	\$150
<input type="checkbox"/> STA	Single Toggle Switch <sup>(3)</sup> .....	\$150
<input type="checkbox"/> DPBA	Dual Push Button <sup>(3)</sup> .....	\$150
<input type="checkbox"/> QPBA	Quad Push Button <sup>(1)(3)</sup> .....	\$150
<input type="checkbox"/> 10WA	10-Way Push Button <sup>(1)(3)</sup> .....	\$255
<input type="checkbox"/> ASL320R	ASL Pulse Mechanical Switch - Red .....	\$200
<input type="checkbox"/> ASL320G	ASL Pulse Mechanical Switch - Green .....	\$200
<input type="checkbox"/> ASL320BK	ASL Pulse Mechanical Switch - Black .....	\$200
<input type="checkbox"/> SESG	Stealth Egg Switch - Green .....	\$135
<input type="checkbox"/> SESB	Stealth Egg Switch - Black .....	\$135
<input type="checkbox"/> FTSQ	Feather Touch - Single Quarter Size .....	\$295
<input type="checkbox"/> FTMD	Feather Touch - Dual Dime Size <sup>(1)(3)</sup> .....	\$295
<input type="checkbox"/> FTMQ	Feather Touch - Quad Dime Size <sup>(1)(3)</sup> .....	\$295
<input type="checkbox"/> FWTA	Four-Way Toggle Switch <sup>(1)(3)</sup> .....	\$150
<input type="checkbox"/> EWRA	8-Way Rocker <sup>(1)(3)</sup> .....	\$255
<input type="checkbox"/> TWRA	10-Way Rocker <sup>(1)</sup> .....	\$255
<input type="checkbox"/> TWTA	10-Way Toggle <sup>(1)</sup> .....	\$255
<input type="checkbox"/> EWPA	8-Way Push Button <sup>(1)</sup> .....	\$255

#### Specify Use for Additional Switch:

- STG Seating
- JM Joystick Mode<sup>(2)</sup>
- JOO Joystick On/Off<sup>(2)</sup>

#### Specify Mounting for Additional Switch:

- ESMLA Left Side (Armpad)
- ESMRA Right Side (Armpad)
- ESMJLA Joystick Left (Mounts to Joystick Bracket)
- ESMJRA Joystick Right (Mounts to Joystick Bracket)
- ESMDA Do not mount

#### Miscellaneous

<input type="checkbox"/> MJBS	Motion Junction Box for Switch Control <sup>(1)(5)</sup> For Mode and/or Emergency Stop. ....	\$175
<input type="checkbox"/> S46L	Single to Dual Switch Adapter .....	\$265
	Allows for dual use of a single switch	
<input type="checkbox"/> M616S	M616 Converter <sup>(1)</sup> .....	\$125
	Allows up to 4 mono phono switches into Multi-Function box	
<input type="checkbox"/> LJS	Lazarus Jump Starter .....	\$150
<input type="checkbox"/> APPS	Accessory Port Power Supply <sup>(1)</sup> 3 amp capacity .....	\$195
<input type="checkbox"/> TPMSH	12 Piece Switch Mounting Hardware Kit (additional kits) .....	\$39

1. Not available with Single-Dual Function Electronics.
2. IN200 is required with REM210, REM211, REM216 joysticks.
3. Not available with REM110 joystick.
4. Not available for Joystick Mode or Joystick On/Off. Must Select Seating (STG)
5. Cannot be used for power On/Off.

USAT Universal Style Articulating Tray ..... \$1,350  
14"W x Adjustable 6"-10"D

#### Add:

- TMA Trilogy Mounting Adapter ..... \$50
  - AVT Articulating Vent (LTV) Tray ..... \$1,350
  - AVBH Articulating Vent Bag Holder<sup>(2)</sup> ..... \$1,350
- Add:
- ABMC Accessory Bag Mounting and Clips ..... \$100

1. The user weight limit is 225 lbs. The maximum weight capacity of the vent tray is 35 lbs. Vents are not standard configuration on Elevating or Precline Systems or with Single Post Backs. Not available with TRRO and TRRO-E.
2. Motion Concepts does not supply the Vent Bag.

### GIMBAL VENT TRAYS<sup>(1)</sup>

GUT Gimbal Universal Tray ..... \$1,350  
14"W x Adjustable 6"-10"D

#### Add:

- GTMA Trilogy Mounting Adapter ..... \$50
  - GVT Gimbal Vent (LTV) Tray ..... \$1,350
  - GVBH Gimbal Vent Bag Holder<sup>(2)</sup> ..... \$1,350
- Add:
- GABMC Accessory Bag Mounting and Clips ..... \$100

1. The Gimbal Vent Trays are not available with Tilt Only and/or Elevate Only Systems. The user weight limit is 200 lbs. Vents are not standard configuration on Elevating or Precline Systems or with Single Post Backs. The maximum weight capacity of the vent tray is 20 lbs. Not available with TRRO and TRRO-E.
2. Motion Concepts does not supply the Vent Bag.

### O2 HOLDERS<sup>(1)</sup>

<input type="checkbox"/> NAO2H	Non-Articulating O2 Gas Holder .....	\$395
<input type="checkbox"/> AOHG	Articulating O2 Gas Holder Cylinder Style <sup>(2)(3)</sup> ..	\$1,250
<input type="checkbox"/> AOHL	Articulating Liquid O2 Holder (6.5" x 4.5") <sup>(2)(3)</sup> ..	\$1,250
<input type="checkbox"/> NAJPL	Non-Articulating IV Pole User's Left <sup>(4)</sup> .....	\$395
<input type="checkbox"/> NAJPR	Non-Articulating IV Pole User's Right <sup>(4)</sup> .....	\$395
<input type="checkbox"/> AIVP	Articulating IV Pole <sup>(2)(3)(4)</sup> .....	\$1,250

1. The user weight limit is 225 lbs. The maximum weight capacity of the O2 holder is 20 lbs. Not available with TRRO and TRRO-E.
2. Required for use on Recline Systems
3. Articulating O2 Holders are not standard configuration on Elevating or Precline Systems or with Single Post Backs.
4. The maximum weight capacity of the IV pole is 5 lbs.

### VENTILATOR INFORMATION

#### Required when ordering a Vent Tray

- BIB Built in battery
- DTSVPS Dealer to supply vent power source

### ARTICULATING TRAYS<sup>(1)</sup>

### ACCESSORIES

