



Thank you for using Motion Concepts eForms.

Motion Concepts is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.



Adobe Acrobat Reader DC

* Interactive functions of our new forms work best with the latest version of Adobe Acrobat Reader DC visit <https://get.adobe.com/reader/> to download and install on your PC or Mac or visit Google Play/ iTunes to download the Adobe Acrobat Reader DC app for your device.

Save

Adobe Acrobat Reader DC allows you to save this form with your content - to complete later or use as a starting point for your next form. Please note that content must be added and saved in Acrobat - saving content from completed forms in the browser may not be possible.

SAVE

Submit

Adobe Acrobat Reader DC allows you to submit this form electronically via your email client.

Simply click the submit button below and step through the simple process.

SUBMIT

CLEAR FORM

Print

If you do not have access to Adobe Acrobat Reader DC simply print this form and complete it by hand and fax it to our Customer Service Department at:

905-695-0583

PRINT FORM



AVIVA STORM RX-20 BASE WITH ULTRA LOW MINI POWER POSITIONING SYSTEM

CDN PRICE LIST AND ORDER FORM
Price Effective August 12, 2022



PRINT FORM

SAVE

CLEAR FORM

REQUIRED INFORMATION

Quote Order

Date of Order: _____ Dealer Account #: _____

Dealer Name: _____

PO #: _____ Tag _____

Purchasing Contact: _____

Phone: _____ Fax: _____

E-mail: _____

RTS/Therapist: _____

Ship to Address: _____

City: _____ Province: _____


Postal Code: _____ Client Gender: M F

Special Client Conditions: _____

REQUIRED MEASUREMENTS

To ensure system is accurately configured please fill in all required

Client Height:	
Client Weight (lbs):	
A. Seat to Shoulder:	
B. Trunk Depth:	
C. Chest Width:	
D. Knee to Back:	
E. Seat to Top of Head:	
F. Elbow to Hand:	
G. Seat to Elbow:	
H. Hip Width:	
I. Knee to Heel:	
Cushion Thickness:	



Customer Service: 1.800.748.7943 | Fax: 1.905.695.0583 | www.motionconcepts.com

INVACARE POWER BASES

SRX-20 AVIVA Storm RX-20 \$7,575

Select STF⁽¹⁾: LOWP 18.125" MEDP 18.625" TALLP 19.625"
(Base comes standard with LED Lights)

1. Seat-to-floor height is measured at the front of the seat pan based on a standard seat size.

USER WEIGHT LIMITS

U175 Weight Capacity up to 175 lbs STD

BASE WIDTH

NARROWFRM Narrow Base 22.4" (uses 2 M50 Batteries) N/C
STDFRM Standard Base 24.4" (uses 2 Group 24 Batteries) N/C

PERFORMANCE MOTORS

STDMTR Standard Motor - 5.8 mph STD
HSMTR High Speed 4-Pole Motor - 7.5 mph⁽¹⁾⁽²⁾⁽³⁾ \$995

1. Not available on the NARROWFRM.
2. Extended high speed use on rugged terrain is not recommended.
3. Performance/High Speed Motors may have less torque and may have less range than standard speed motors.

TRANSPORT TIE DOWN

TRBKTS Wheelchair Transport Brackets⁽¹⁾ STD
1. For unoccupied use only

CASTER OPTION

AB1228-3 8" x 2.5" Black Casters w/Solid Inlay N/C

TIRE OPTION

AB1430-3 14" x 3" Black Tires Pneumatic N/C
AB1431-3 14" x 3" Black Tires w/Gel Foam Inserts N/C

BATTERY CHARGER OPTION

110CHARGER 110 Volt Battery Charger N/C
OMIT Omit Battery Charger N/C

INSTALLATION OF BATTERIES

M24 SLD G FT MK Group 24 Battery (On Chair) \$840
M50-12 SLD M MK Group 50 Battery (On Chair) \$650

FRAME FINISH

<input type="checkbox"/> 170P		Obsidian Black N/C
<input type="checkbox"/> 171P		Ruby Red N/C
<input type="checkbox"/> 172P		Peacock Blue N/C
<input type="checkbox"/> 173P		Pink Flamingo N/C
<input type="checkbox"/> 174P		Purple Passion N/C
<input type="checkbox"/> 175P		Denim Navy N/C
<input type="checkbox"/> 176P		Lime Twist N/C
<input type="checkbox"/> 177P		Arctic Blue N/C

RIM INSERT COLORS - MUST PICK ONE

<input type="checkbox"/> 170PR		Obsidian Black N/C
<input type="checkbox"/> 171PR		Ruby Red N/C
<input type="checkbox"/> 172PR		Peacock Blue N/C
<input type="checkbox"/> 173PR		Pink Flamingo N/C
<input type="checkbox"/> 174PR		Purple Passion N/C
<input type="checkbox"/> 175PR		Denim Navy N/C
<input type="checkbox"/> 176PR		Lime Twist N/C
<input type="checkbox"/> 177PR		Arctic Blue N/C
<input type="checkbox"/> 178PR		Banana Split N/C
<input type="checkbox"/> 179PR		Orange Marmalade N/C
<input type="checkbox"/> 180PR		Morning Mist N/C

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

CONTROLLERS

- EXPC Expandable Controller STD
- PWH Harness Required for Expandable System N/C
- GTRAC-LX LiNX G-TRAC Module \$1,200

DRIVER CONTROL

- REM215 LED Non -Expandable Remote/Drive Only/Joystick⁽¹⁾ . . . STD

1. Controls up to 2 actuators. Must choose Single/Dual function through drive control and light operation.

- REM216 LED Expandable Remote/Joystick⁽¹⁾ \$875
Three drive functions, control of up to 6 actuators through the drive control and light operation.

- REM400 Color 3.5" Touch Screen Remote/Joystick⁽¹⁾⁽²⁾ \$1,100
 - SW400TGL Toggle Kit for REM400⁽³⁾ \$200

- REM500 LiNX Color 3.5" Touch Screen Display Only - No Driver Control⁽¹⁾⁽²⁾ \$1,100

Select mounting position R5MR Right
 R5ML Left

1. Must choose EXPC Expandable Controller and PWH Harness.
2. REM400 or REM500 is required when using specialty controls, choose specialty controls purchased separately.
3. Adds 2 toggle switches to REM400.

ELECTRONICS - EXPANDABLE MODULES

- OUT500 Output Module \$645
- IN200 LiNX Switch Input Module with Power ON/OFF for REM200 Series⁽¹⁾ \$645
- IN500 LiNX 9 Pin Input Module with Sip-N-Puff Nozzle⁽²⁾⁽³⁾⁽⁴⁾ \$1,800
- PKG32666 Therafin Sip-N-Puff Breath Tube Kit \$440

1. Not for use with Alternative Drive Controls.
2. Only one input module per chair, requires REM400 or REM500.
3. Input module also required for any 9-pin specialty control.
4. 9 pin specialty control and Sip-N-Puff cannot both be used together.

MOUNTING FOR DRIVER CONTROLS

- FR Fixed Joystick Mounting - Right STD
- FL Fixed Joystick Mounting - Left STD
- SAQR Motion Height Adjustable Swing-Away Quad Link (Right Hand)⁽¹⁾ \$344
- SAQL Motion Height Adjustable Swing-Away Quad Link (Left Hand)⁽¹⁾ \$344
- MRR Maxx Resolve Swing-Away Joystick Mount (Right Hand)⁽¹⁾ \$450
- MRL Maxx Resolve Swing-Away Joystick Mount (Left Hand)⁽¹⁾ \$450

1. Not available with the REM500.

ATTENDANT DRIVE CONTROL CHOICES

- ACU Proportional Attendant Control⁽¹⁾ \$1,024
Select mounting position RIGHT-ACU Right
 LEFT-ACU Left

1. Requires Expandable Controller.

SPECIALTY PROPORTIONAL CONTROLS REQUIRES REM400 OR REM500

- CREM LiNX Compact Remote \$1,097
 - CREM-LF LiNX Compact Remote Low Force \$1,097
- Select mounting position
 MR Right
 ML Left

LINX ELECTRONIC ACCESSORIES

- LAK LiNX Access Key⁽¹⁾ \$50
- MUSB Motion Dual USB Charger N/C

1. Required for programming

JOYSTICK TOPS

- PC101A Bodypoint U Shaped Handle 3"⁽¹⁾ \$120
- PC102A Bodypoint U Shaped Handle 4"⁽¹⁾ \$120
- PC107A Bodypoint Rubber Dome⁽¹⁾ \$120
- PC110A Grooved Mushroom Joystick Handle 3/16" Stem⁽¹⁾ \$120
- 1560 T Handle Flexible Joystick Extension⁽¹⁾ \$120
- 1561 Straight Handle Flexible Joystick Extension⁽¹⁾ \$120
- 1826 Chin Cup⁽¹⁾ \$120

1. To remove/disengage the joystick knob from the REM400 Remote pull straight up on the joystick (DO NOT TWIST) otherwise damage may occur and may void the warranty.

NOTE: Systems Include Seat Pan, Back Pan, and Standard Armrests. They are designed for use with a 2"- 3" cushion

REHAB SEAT

PED-RS Mini Pediatric Rehab Seat \$3,995

1. Not available with 12" Elevating Seat Lift Module, Order Ultra Low Maxx Elevating Lift Seat Only.

POWER POSITIONING SYSTEMS

MDCGT-50 Ultra Low Mini Maxx CG Tilt (50°) Only. . . . \$6,969

MDCGT-45 Ultra Low Mini Maxx CG Tilt (45°) Only⁽¹⁾ . . \$6,969

MDCGTRESR-50 Ultra Low Mini Maxx CG Tilt (50°) & Power Recline(163°) w/Extended Shear Reduction(ESR) \$11,569

MDCGTRESR-45 Ultra Low Mini Maxx CG Tilt (45°) & Power Recline(163°) w/Extended Shear Reduction(ESR)⁽¹⁾ \$11,569

MDRESR Ultra Low Mini Maxx Power Recline (163°) w/Extended Shear Reduction(ESR) \$10,350

MDES12 Ultra Low Mini Maxx Elevate Lift Seat Only \$4,995
Weight capacity up to 175lbs

1. Tilt will be reduced to 45° if adding an Elevating Lift Seat Module

ADD: ELEVATING LIFT SEAT MODULE

ESM12 12" Elevating Lift Seat Module⁽¹⁾ \$2,495
Weight capacity up to 175lbs

1. Uses an integrated Tilt/Lift mechanism, tilt is 45° and the seat to floor height is not affected by the Elevating Lift Seat.

MINI MODULES MUST SELECT SYSTEM TO ADD MODULES

EPRETILT 5° Electronic Pre-Tilt⁽¹⁾⁽²⁾ N/C

ERSSP Programable Electronic Recline Function Start & Stop⁽²⁾ Specify Start ___degrees & Stop ___degrees. \$500

PRM8 8° Precline Module⁽³⁾ only available with Recline N/C

PRM30 30° Precline Module⁽³⁾ only available with Recline . . . \$500

1. This Electronic Pre-tilt option adds approximately 2" to the front seat to floor height
2. Requires Multi-Function Electronics through Drive Control (MFCBD). Set range between 90° and 163°. Will be set to factory settings unless specified.
3. Total degrees of recline is 163° ie. if 30° of precline is used, then 133° of recline will be available.

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

SEAT WIDTH SETTINGS

MSW11 11" Wide (Adjustable from 11" to 15") N/C

MSW12 12" Wide (Adjustable from 11" to 15") N/C

MSW13 13" Wide (Adjustable from 11" to 15") N/C

MSW14 14" Wide (Adjustable from 11" to 15") N/C

MSW15 15" Wide (Adjustable from 11" to 15") N/C

ULTRA RAIL SEAT DEPTH SETTINGS

MSD11 11" Deep⁽¹⁾ (Adjustable from 11" to 15") N/C

MSD12 12" Deep (Adjustable from 11" to 15") N/C

MSD13 13" Deep (Adjustable from 11" to 15") N/C

MSD14 14" Deep (Adjustable from 11" to 15") N/C

MSD15 15" Deep (Adjustable from 11" to 15") N/C

1. For seat depths of 11" weight capacity is 150lbs.

MATRIX SEATING OPTIONS

Matrix Cushions available on page 8

Matrix Backs available on page 8

Matrix Headrests available on page 9



LAP BELTS

BELT35 Push Button Style Lap Belt, Length 35" STD

BELT48 Push Button Style Lap Belt, Length 48" . . . N/C

MBPHB Mini BP Padded Hip Belt 2 point⁽¹⁾ \$175

1. Adjustable from 14" to 41" in length, pad is 13" x 1.5" in length

CHEST STRAP

MFUXS BP Monoflex, Underarm RLS - XS(2.75"x12.25") \$195

MFUS BP Monoflex, Underarm RLS - S(3"x15.5") . . \$195

MFCXS BP Monoflex, Centre RLS - XS(2.75"x15.25") . \$195

MFCS BP Monoflex, Centre RLS - S(3"x15.5") . . . \$195

MINIFLEX-XS Miniflex Abdominal Support Band 4.5" H for Matrix Mini backs 11"-14" wide \$245

MINIFLEX-S Miniflex Abdominal Support Band 6" H for Matrix Mini backs 11"-14" wide . . . \$245

MOTION CONCEPTS HEADREST OPTIONS

MCSH Motion Concepts Standard⁽¹⁾ N/C

MCOH Motion Concepts Auto Style⁽¹⁾ \$235

FM Fixed Mounting (non-removable)⁽²⁾ N/C

OHR Omit: Headrest N/C

1. Must select Fixed Mounting (FM) or Matrix Elan Headrest & Hardware on page 9.
2. Only compatible with Standard Rehab, Contour or High Back Models. All other upgradeable backs require removeable hardware.

MATRIX HEADRESTS OPTIONS

Matrix Elan Headrests and Hardware available on page 9



Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

BACK CANES (NON RECLINE SYSTEMS)⁽¹⁾

- SBC14 14" Straight Back Canes STD
- SBC16 16" Straight Back Canes STD
- SBC18 18" Straight Back Canes STD
- SBC20 20" Straight Back Canes STD

Select Cane Angle 75° 79° 85° 90° 95° 101° 106° 112° 116°
BA75 BA79 BA85 BA90 BA95 BA101 BA106 BA112 BA116

1. Canes will be set at 90° angle and are adjustable from 75° to 116°

FINISHED BACK HEIGHT⁽¹⁾⁽²⁾ FOR REHAB, TILT AND/ELEVATE ONLY

- 14" 15" 16" 17" 18" 19" 20" 21" 22"
FBT14 FBT15 FBT16 FBT17 FBT18 FBT19 FBT20 FBT21 FBT22
- 23" 24" 25"
FBT23 FBT24 FBT25

1. 4" Gap from the seat pan to the bottom of back pan is recommended.
2. Finished back height is measured from seat pan to top of backrest.

FINISHED BACK HEIGHT⁽¹⁾⁽²⁾ FOR RECLINE OR TILT/RECLINE

- 16" 17" 18" 19" 20" 21" 22" 23" 24"
FBR16 FBR17 FBR18 FBR19 FBR20 FBR21 FBR22 FBR23 FBR24
- 25"
FBR25

1. 6" Gap from the seat pan to the bottom of back pan is recommended.
2. Finished back height is measured from seat pan to top of backrest.

BACK OPTIONS

- SRB Standard Rehab Back & Cushion \$125
- PIPF Flat Planar/I-Back Interface Plate⁽¹⁾⁽²⁾⁽⁴⁾ \$125
- PIPR Recessed Planar/I-Back Interface Plate⁽¹⁾⁽²⁾⁽⁴⁾ . \$125
- ORBC Omit Rigid Back & Cushion⁽³⁾ N/C

1. Not available on tilt only or tilt/elevate systems, for use with Recline Module.
2. Included with all Recline Systems when no back is selected. May also be used as an interface for after market backs.
3. Not available with Recline/ESR.
4. Fits Comfort Company, Freedom, Pindot and Moulded Backs.

MATRIX BACKREST OPTIONS

Matrix Backs available on page 6/7



ARMRESTS (NON RECLINE SYSTEMS)

- PDP1 Pediatric Dual Post Adj.Height 6.5"-9.5" STD
- PDP2 Pediatric Dual Post Adj.Height 9.75"-12.75" . . . STD
- PDP3 Pediatric Dual Post Adj. Height 12.75"-15.75" . STD

RECLINING ARMRESTS

- PDR1 Pediatric Reclining Adj.Height 7"-8.25" STD
- PDR2 Pediatric Reclining Adj.Height 9.5"-13" STD
- PDR3 Pediatric Reclining Adj. Height 13"-16" STD

ARMREST TUBE LENGTH - PLEASE SPECIFY

- ATFL Left Full STD
- ATDL Left Desk STD
- ATFR Right Full STD
- ATDR Right Desk STD

ARM PADS

Standard

- SAPFL Left Full (2.25" x 14") STD
- SAPDL Left Desk (2.25" x 10") STD
- SAPFR Right Full (2.25" x 14") STD
- SAPDR Right Desk (2.25" x 10") STD

Waterfall Arm Pads

- WAPFL Left Full (13") N/C
- WAPDL Left Desk (9") N/C
- WAPFR Right Full (13") N/C
- WAPDR Right Desk (9") N/C

Modular Arm Pad - Left⁽¹⁾

- MAPTFL Left Tray Full (14") \$95
- MAPTDL Left Tray Desk (10") \$95

Select Pad Insert below:

- ISSFL Iskin - Integrated soft skin foam
- SCVFL Startex Covered Visco foam

Modular Arm Pad - Right⁽²⁾

- MAPTFR Right Tray Full (14") \$95
- MAPTDR Right Tray Desk (10") \$95

Select Pad Insert below:

- ISSFR Iskin - Integrated soft skin foam
- SCVFR Startex Covered Visco foam

Flat Multi Position Arm Pads

- MPAFL Left Full (14") \$150
- MPADL Left Desk (10") \$150
- MPAFR Right Full (14") \$150
- MPADR Right Desk (10") \$150

Gel Pads

- GAPNL Narrow Left (2" x 12") \$88
- GAPNR Narrow Right (2" x 12") \$88
- GAPWL Wide Left (3.5" x 12") \$96
- GAPWR Wide Right (3.5" x 12") \$96
- GAPLL Long Left (3.5" x 14") \$103
- GAPLR Long Right (3.5" x 14") \$103

- OMIT Omit Armpads(7/8" diameter tube) (\$25)

1. Requires interchangeable pad insert.

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

POWER CENTRE MOUNT FOOT PLATFORM

MUST BE COMPLETED	Measurement in Inches
Knee to Heel:	
Cushion Thickness: (Subtract)	
Cushion Compression:(Add)	
Seat Pan to Footplate:(Total)	


LNx Power Center Mount Foot Platform⁽¹⁾⁽²⁾ HCPCS code E1012

- LNxA 10" - 14" Seat pan to footplate \$3,200
- Options
 - LIFPN Comes with a rubber coated footplate (10.5" W x 10") STD
 - LICP Comes standard with individual calf pads. STD
 - OPCP One Piece Calf Panel⁽³⁾ \$250
 - LNXPf Plantar Flex Option for LNx Foot platform⁽²⁾ . \$325
 - LSFP Latitude Style Foot Platform \$2,995
For Finished Knee to heel measurements of 6.5"-12"
 - LSPFP Latitude Style Padded Foot platform \$250

1. 3" ground clearance required.
2. Only available with Foot Platform Options.
3. Not available with TRRO and TRRO-E.

STANDARD FRONT RIGGING


- FCMP Fixed Centre Mount Platform⁽¹⁾ \$395
 - IFP Comes with a rubber coated footplate (10.5" W x 10") N/C
- Select seat pan to footplate:
- SPTF12 6" - 13" N/C



1. 3" ground clearance required.


MANUAL CENTRE MOUNT FRONT RIGGING

- MSFCMFP Maxx Style Fixed Centre Mount Foot Platform⁽¹⁾⁽²⁾ \$395
With Quick Release Height Adjustment
 - IFP Comes with a rubber coated footplate (10.5" W x 10") N/C
 - ICP Individual calf pads (pair) \$200
- Select seat pan to footplate:
- SPTF13 8" - 13" N/C
 - SPTF14 13" -19" N/C



1. Depth and height adjustable.
2. 3" ground clearance required.

FOOT PLATFORM OPTIONS FOR MAXX CMFP



- SFFP Small Foot Platform⁽¹⁾ (12" w x 8.5" d) \$180
 - MFFP Medium Platform⁽¹⁾ (12" w x 10.75" d) \$180
1. May need to adjust leg position to avoid interference with front castors

SWING AWAY FOOTRESTS

- HD70 Heavy Duty 70° Swing Away Footrest^{(pair)(1)(2)} \$495
- Select Receivers:
- RS Straight Receivers N/C
 - RF Flared 1" Outward Receivers N/C
- Select Seat Pan to Footplate:
- SPTF11 Seat pan to footplate 6"-10" N/C
 - SPTF8 Seat pan to footplate 11"-15" N/C
 - SPTF2 Seat pan to footplate 13"-17" N/C
 - HDFCP Flat Calf Pad Option (pair) \$275


1. Must choose footplate on page 6
2. Heavy duty durable construction (4.5lbs ea.).

SWING AWAY FOOTRESTS

- I70SAF Invacare 70° Swingaway Footrests (Pair) \$325
 - I70-CF Comes with composite foot plates N/C
 - 90SWING Invacare 90° Swingaway Footrests with AT5543 Adjustable Angle Footplates, Heel Loops and Impact Guards⁽¹⁾ (Range 4.5"-7.5")(Pair) \$450
- Select Seat Pan to Footplate⁽¹⁾:
- SPTF11 Seat pan to footplate 6"-10" N/C
 - SPTF8 Seat pan to footplate 11"-15" N/C
 - SPTF2 Seat pan to footplate 13"-17" N/C

1. For optional Seat pan to footplate extensions on 90SWING you must order with Motion Concepts Foot plates

FOOTPLATE SIZE SELECTION



1. Certain footplate sizes may not fit all seat widths depending on legrests selected.

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

ADJUSTABLE ANGLE FOOTPLATES

<input type="checkbox"/> AFSL	Left Small Footplate (4.25"W x 8"D)	\$165
<input type="checkbox"/> AFSR	Right Small Footplate (4.25"W x 8"D)	\$165
<input type="checkbox"/> AFML	Left Medium Footplate (5.25"W x 8"D)	\$165
<input type="checkbox"/> AFMR	Right Medium Footplate (5.25"W x 8"D)	\$165
<input type="checkbox"/> AFLL	Left Large Footplate (6.25"W x 8"D)	\$165
<input type="checkbox"/> AFLR	Right Large Footplate (6.25"W x 8"D)	\$165

MULTI-AXIS ADJUSTABLE ANGLE FOOTPLATES

<input type="checkbox"/> MFSL	Left Small Footplate (4.25"W x 8"D)	\$225
<input type="checkbox"/> MFSR	Right Small Footplate (4.25"W x 8"D)	\$225
<input type="checkbox"/> MFML	Left Medium Footplate (5.25"W x 8"D)	\$225
<input type="checkbox"/> MFMR	Right Medium Footplate (5.25"W x 8"D)	\$225
<input type="checkbox"/> MFLL	Left Large Footplate (6.25"W x 8"D)	\$225
<input type="checkbox"/> MFLR	Right Large Footplate (6.25"W x 8"D)	\$225

HEEL LOOPS, CUPS AND TOE STRAPS

<input type="checkbox"/> HLL	Left Heel Loop (each)	\$20
<input type="checkbox"/> HLR	Right Heel Loop (each)	\$20
<input type="checkbox"/> HCL	Left Cup Loop (each)	\$20
<input type="checkbox"/> HCR	Right Cup Loop (each)	\$20
<input type="checkbox"/> TSL	Left Toe Strap (each)	\$20
<input type="checkbox"/> TSR	Right Toe Strap (each)	\$20

SINGLE-DUAL FUNCTION ELECTRONICS

<input type="checkbox"/> SFCB	Single/Dual Function Electronics ⁽¹⁾	STD
-------------------------------	---	-----

Select Switch:

- SPBO Single Push Button
- STO Single Toggle
- DPBO Dual Push Button
- SOS Omit Switch
- JOY Through Joystick
- SC Through Specialty Control

Optional:

<input type="checkbox"/> RDE	Reduced Drive for Elevating Seat ⁽²⁾	\$400
<input type="checkbox"/> SPBSA	Additional Single Push Button	\$150
<input type="checkbox"/> STA	Additional Single Toggle Switch	\$150
<input type="checkbox"/> DPBA	Additional Dual Push Button	\$150

1. These electronics include drive lock-out and tilt limit.
2. If reduced drive is not selected the system will be in drive lock-out when elevated.

MULTI-FUNCTION ELECTRONICS

<input type="checkbox"/> MFCB	Multi-Function through Switch/Toggle ⁽¹⁾	\$1,950
-------------------------------	---	---------

Select Switch:

- QPBO Quad Push Button
- FWTO Four-Way Toggle Switch
- EWRO 8-Way Rocker⁽²⁾

<input type="checkbox"/> MFCBD	Multi-Function through Drive Control ⁽¹⁾	\$1,950
--------------------------------	---	---------

Select Control:

- JOY Through Joystick
- SC Through Specialty Control

Optional:

- RDY
 - RDN
- Include Reduced Drive for Elevating Seat⁽³⁾
Do not Include Reduced Drive for Elevating Seat⁽³⁾

<input type="checkbox"/> QPBA	Additional Quad Button	\$150
<input type="checkbox"/> FWTA	Additional Four-Way Toggle Switch	\$150
<input type="checkbox"/> EWRA	Additional 8-Way Rocker ⁽²⁾	\$205
<input type="checkbox"/> 10WA	Additional 10-Way Switch	\$255

1. These electronics include attendant control, drive lock-out, tilt/recline limit and actuator speed control
2. Dual function is not available with 8-Way Rocker.
3. If reduced drive is not selected system will be in drive lock-out when elevated.

ELECTRONIC SWITCH MOUNTING

<input type="checkbox"/> ESML	Left Side (Armpad)
<input type="checkbox"/> ESMR	Right Side (Armpad)
<input type="checkbox"/> ESMJL	Joystick Left (Mounts to Joystick Bracket)
<input type="checkbox"/> ESMJR	Joystick Right (Mount to Joystick Bracket)
<input type="checkbox"/> ESMDD	Do not mount

ELECTRONICS OPTIONS

<input type="checkbox"/> M616S	M616 Converter ⁽¹⁾	\$125
<input type="checkbox"/> FTSQ	Feather Touch Switch Single Function Quarter Size	\$295
<input type="checkbox"/> FTMD	Feather Touch Switch Multi-Function Dual Dime Size	\$295
<input type="checkbox"/> FTMQ	Feather Touch Switch Multi-Function Quad Dime Size	\$295
<input type="checkbox"/> ASL320R	ASL Pulse Mechanical Switch Red	\$125
<input type="checkbox"/> ASL320G	ASL Pulse Mechanical Switch Green	\$125
<input type="checkbox"/> ASL320BK	ASL Pulse Mechanical Switch Black	\$125
<input type="checkbox"/> SESG	Stealth Egg Switch - Single Function Green	\$135
<input type="checkbox"/> SESB	Stealth Egg Switch - Single Function Black	\$135
<input type="checkbox"/> PSOS	Power Shut Off Switch - Toggle ⁽¹⁾	\$175
<input type="checkbox"/> LJS	Lazarus Jump Starter ⁽¹⁾	\$150

Dealer Name: _____
Dealer Account Number: _____
PO #: _____
Tag: _____

ELECTRONICS OPTIONS CONTINUED

- APPS Accessory Port Power Supply⁽¹⁾
(3 amp capacity) \$195
- TPSMH 12 Piece Switch Mounting Hardware
Kit (additional kits). \$39

1. Not available with Single-Dual Function Electronics.

Total Retail Value \$ _____



GIMBAL VENT TRAYS⁽¹⁾⁽²⁾⁽³⁾⁽⁴⁾

- GUT Gimbal Universal Tray. \$1,350
(14"W x Adjustable 6" - 10"D)
- Add: GTMA Trilogy Mounting Adapter \$50
- GVT Gimbal Vent (LTV)Tray. \$1,350
- GVBH Gimbal Vent Bag Holder⁽⁵⁾ \$1,350
- Add: GABMC Accessory Bag Mounting and Clips \$100

- The Gimbal Vent Trays are not available with Tilt Only Systems.
- The user weight limit is 150 lbs.
- Vents are not standard configuration on Elevating or Precline Systems. Please call Customer Service for RFQ Quotation.
- The maximum weight capacity of the vent tray is 20 lbs.
- Motion Concepts does not supply the Vent Bag.

[PRINT FORM](#) [SAVE](#) [CLEAR FORM](#)

NOTES

ADDITIONAL OPTIONS & ACCESSORIES

- THS Transfer Handles 4" - 6.5" \$350
- THL Transfer Handles 9" - 11.5" \$350
- FFLTL Flip and Fold Away 1/2 Lap Tray
Left Mounting(14"W x 10"D)⁽¹⁾⁽²⁾ \$795
- FFLTR Flip and Fold Away 1/2 Lap Tray
Right Mounting(14"W x 10"D)⁽¹⁾⁽²⁾ \$795

- Not available on joystick side and only available with a flat arm pad or modular arm pad. For seat widths 14" and greater.
- The maximum weight capacity of the tray is 10 lbs.

INSTALLATION

- MCI Installation by Motion Concepts N/C

PLEASE SPECIFY PAGES INCLUDED

- Pages 1-7 Excludes Matrix and Maxx Seating Accessories
- Pages 1-10 Complete Order Form



CUSHION COVER REFERENCE GUIDE

Cushion	Kid*ab*ra	Meshtex	Startex/ Spacetex	Pink Lollipop/ Blue Blizzard	Infection Control
PS/PSP	N/C	N/A	STD	STD	Reversible
Kid*ab*ra	STD	\$103	N/A	\$103	N/A
Vi	\$103	STD	\$103	\$103	N/A
Libra/Libra Adjust	N/A	N/A	STD	N/A	\$50

MATRIX SEAT CUSHION OPTIONS

Matrix Libra Adjust Seat Cushion 14" - 15" wide \$985

- LC1414-ADJ LC1416-ADJ LC1516-ADJ

Matrix Libra Seat Cushion 14" - 15" wide \$745

- LC1414 LC1416 LC1516

Fabric Options for Libra/Libra Adjust Cushions

- SFLF Spacetex N/C
 ICFL Infection Control N/C

Matrix PS Seat Cushion 11" - 14" wide \$569

- PS1112 PS1114 PS1212 PS1214 PS1312 PS1314
 PS1412 PS1414

Matrix PSP Seat Cushion 11" - 14" wide \$645

- PSP1112 PSP1114 PSP1212 PSP1214 PSP1312
 PSP1314 PSP1412 PSP1414

Fabric Options for PS/PSP Cushions

- CV-K Kid*ab*ra (Non Reversible) N/C
 CV-P Pink Lollipop (Non Reversible) N/C
 CV-B Blue Blizzard (Non Reversible) N/C

Matrix Kid*ab*ra Cushion 12" - 14" wide \$576
*(Comes standard with Kid*ab*bra fabric)*

- KA1212-VI KA1214-VI KA1412-VI KA1414-VI

Fabric Options for Kidabra Cushions

- CV-M Meshtex (Non Reversible) \$103
 CV-P Pink Lollipop (Non Reversible) \$103
 CV-B Blue Blizzard (Non Reversible) \$103

Matrix Vi Seat Cushion 14" wide \$624

- MA1415-VI

BACK COVER REFERENCE GUIDE

Back	Kid*ab*ra	Meshtex	Startex	Pink Lollipop	Blue Blizzard	Infection Control
Mini	STD	N/A	N/C	N/C	N/C	N/A
Elite	N/A	N/A	N/A	N/A	N/A	N/A
Elite TR	N/A	N/A	N/A	N/A	N/A	N/A
EliteDeep	N/A	N/A	STD	N/A	N/A	N/A
PB	N/A	N/A	STD	N/A	N/A	N/A
PB Deep	N/A	N/A	STD	N/A	N/A	N/A

UPGRADABLE BACK OPTIONS

Matrix Mini Back 10" - 14" wide \$845
*(Comes standard with Kid*ab*bra fabric)*

- MINI1010 MINI1012 MINI1110 MINI1112
 MINI1210 MINI1212 MINI1214 MINI1310
 MINI1312 MINI1314 MINI1410 MINI1412
 MINI1414

Fabric Options for Mini Back

- MINI-S Startex (Non Reversible) N/C
 MINI-P Pink Lollipop (Non Reversible) N/C
 MINI-B Blue Blizzard (Non Reversible) N/C
 HIMINI Chest Harness Interface Kit for PB Back⁽¹⁾⁽²⁾ ... \$110

1. Does not include the chest strap.
2. This Interface Kit is only Required for 4 point harness

MATRIX MINI FIXED LATERALS - FORMINI BACKS

Matrix MINI Fixed Lateral with Pad

- MINILATL LEFT (4.5" x 3") \$130
 MINILATR RIGHT (4.5" x 3") \$130

MATRIX MINI SWING AWAY LATERALS - FORMINI BACKS

MINI Swing Away Lateral Mounting Hardware

- SAHMINIL LEFT SIDE \$250
 SAHMINIR RIGHT SIDE \$250

Mini Lateral Trunk Support RIGHT(each)

- Select Pad Size:
 EPADXSR XS(3.25"x4.25") \$115
 EPADSR S(3.5"x5.25") \$115

Mini Lateral Trunk Support LEFT(each)

- Select Pad Size:
 EPADXSL XS(3.25"x4.25") \$115
 EPADSL S(3.5"x5.25") \$115

Fabric Options for Mini Lateral Trunk Supports

- MINILAT-K Startex (None Reversible) N/C
 MINILAT-P Pink Lollipop (None Reversible) N/C



Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

UPGRADABLE BACK OPTIONS - CONTINUED

Matrx Elite Back 14" - 15" wide (3" Contour) \$766

- PBE1410 PBE1412 PBE1414 PBE1416 PBE1418
- PBE1420 PBE1510 PBE1512 PBE1514 PBE1516
- PBE1518 PBE1520

Matrx Elite TR Back 15" wide \$866

- PBE1516-TR PBE1518-TR PBE1520-TR

Matrx Elite Deep Back 14" - 15" wide \$918

- EDB1412 EDB1416 EDB1418 EDB1420 EDB1510
- EDB1512 EDB1514 EDB1516 EDB1518 EDB1520

Matrx Elite E2 Back 14" - 15" wide/10"-16"H(3" Contour) . \$965

- E2S1410 E2S1412 E2S1414 E2S1416
- E2S1514 E2S1516

Matrx Elite E2 Back 14" - 15" wide/18"/20"H(3" Contour) ..\$995

- E2S1418 E2S1420 E2S1518 E2S1520

**Matrx Elite E2 Deep Back 14" - 15" wide/10"-16"H (6" Contour)
.....\$1,125**

- E2D1410 E2D1412 E2D1414 E2D1416
- E2D1514 E2D1516

**Matrx Elite E2 Deep Back 14" - 15" wide/18"/20"H (6" Contour)
\$1,155**

- E2D1418 E2D1420 E2D1518 E2D1520

Matrx PB Back 14"(5" Contour) \$734

- MPB1412 MPB1416

Matrx PB Deep Back 14" (7" Contour) \$892

- MPBD1412 MPBD1416

OTHER UPGRADABLE BACK OPTIONS

- EBSS **Super Soft Foam for Elite Back Upgrade . . . \$103**
- WSCU **Infection Control Cover for Elite/PB Backs . . . \$103**
purchased with back
- WSC **Infection Control Cover for Elite/PB Backs . . . \$150**
purchased with without back
- CHIEB **Chest Harness Interface Kit for Elite Back⁽¹⁾⁽²⁾ . \$103**
- CHIEB **Chest Harness Interface Kit for PB Back⁽¹⁾⁽²⁾ . . . \$52**

1. Does not include the chest strap.

2. This Interface Kit is only Required for 4 point harnesses.

MATRIX ELAN HEADREST PAD

- ESP6 **Matrx Elan Standard Pad 6" \$310**
- ESP10 **Matrx Elan Standard Pad 10" \$310**
- ESP14 **Matrx Elan Standard Pad 14" \$310**
- EOP9 **Matrx Elan Occipital Pad 9" \$341**
- EOP12 **Matrx Elan Occipital Pad 12" \$341**
- E4PTSM **Matrx Elan 4-Point Pad 9"W x 7"H \$448**

- E4POINT **Matrx Elan 4-Point Pad 11"W x 10"H \$448**
- ESP6-IC **Matrx Elan Standard Pad 6" with
Infection Control \$335**
- ESP10-IC **Matrx Elan Standard Pad 10" with
Infection Control \$335**
- ESP14-IC **Matrx Elan Standard Pad 14" with
Infection Control Cover \$335**
- EOP9-IC **Matrx Elan Occipital Pad 9" with
Infection Control Cover \$366**
- EOP12-IC **Matrx Elan Occipital Pad 12"with
Infection Control Cover \$366**
- E4PTSM-IC **Matrx Elan 4-Point Pad 9"W x 7"H
with Infection Control Cover \$473**
- E4POINT-IC **Matrx Elan 4-Point Pad 11"W x 10"H
with Infection Control Cover \$473**

HEADREST MOUNTING HARDWARE OPTIONS⁽¹⁾

- LXSHW **Loxx Standard Hardware (15" Vertical Post) . \$275**
- LXMMH **Loxx Mini Hardware (11" Vertical Post) \$275**
- MEMH **Elan Mini Removable HDW \$285**
- MEHWRE **Elan Multi-Axis Removable HDW \$285**

1. Only compatible w/Standard Rehab, Contour & High Back Models. All other head rests REQUIRE an upgrade to Removable Hardware.

MATRIX FIXED LATERALS - FOR ELITE AND ELITE TR BACKS

- PBSSL **Matrx Small Lateral Fixed-Left(3x5.5") \$113**
- PBSSR **Matrx Small Lateral Fixed-Right(3x5.5") . . . \$113**
- PBLSL **Matrx Standard Lateral Fixed-Left(7.5x5") . \$116**
- PBLSR **Matrx Standard Lateral Fixed-Right (7.5x5")\$116**
- PBOLL **Matrx Offset Lateral Fixed Left \$125**
- PBOLR **Matrx Offset Lateral Fixed Right \$125**

MATRIX SWING AWAY LATERALS - FORELITEAND ELITE TRBACKS

- Matrx Lateral Trunk Support Left⁽¹⁾**
- EPADXL **Size XS (3.25x4.25") \$115**
- EPADSL **Size S (3.5x5.25") \$115**
- EPADML **Size M (3.75x5.75") \$115**
- EPADLL **Size L (4.25x6.5") \$115**
- Matrx Lateral Trunk Support Right⁽¹⁾**
- EPADXR **Size XS (3.25x4.25") \$115**
- EPADSR **Size S (3.5x5.25") \$115**
- EPADMR **Size M (3.75x5.75") \$115**
- EPADLR **Size L (4.25x6.5") \$115**

1. Must choose swing away hardware on next page



MATRIX SWING AWAY HARDWARE - FOR ELITE AND ELITE TRUCKS

<input type="checkbox"/> SAEHBL	Matrix Elite Swing Away Hardware Left	\$250
<input type="checkbox"/> SAHEBR	Matrix Elite Swing Away Hardware Right	\$250
<input type="checkbox"/> OSAHL	Matrix Offset Elite Swing Away Hardware Left	\$250
<input type="checkbox"/> OSAHR	Matrix Offset Elite Swing Away Hardware Right	\$250

OTHER MATRIX ACCESSORIES

<input type="checkbox"/> MACESL	Multi Axis Arm Support - Left ⁽¹⁾	\$225
<input type="checkbox"/> MACESR	Multi Axis Arm Support - Right ⁽¹⁾	\$225

1. Allows rotation of the armpad, only compatible with Gel pad, Flat pad or Ergonomic Arm Troughs

MAXX LATERALS TRUNK SUPPORTS⁽¹⁾

Left Fixed Mount Lateral Trunk Support⁽¹⁾
Select Left Trunk Support Size

<input type="checkbox"/> LATXL	Size XS (3.25" x 4.25")	\$100
<input type="checkbox"/> LATSL	Size S (3.5" x 5.25")	\$100
<input type="checkbox"/> LATML	Size M (3.75" x 5.75")	\$100
<input type="checkbox"/> LATLL	Size L (4.25" x 6.5")	\$100

Right Fixed Mount Lateral Trunk Support⁽¹⁾
Select Right Trunk Support Size

<input type="checkbox"/> LATXR	Size XS (3.25" x 4.25")	\$100
<input type="checkbox"/> LATSR	Size S (3.5" x 5.25")	\$100
<input type="checkbox"/> LATMR	Size M (3.75" x 5.75")	\$100
<input type="checkbox"/> LATLR	Size L (4.25" x 6.5")	\$100

1. Not compatible with "deep" matrix backs or sling back.

Swing Away Trunk Support Hardware⁽⁴⁾

<input type="checkbox"/> SALATL	Left MAXX Style Swing Away HDW ⁽¹⁾	\$195
Select Left Trunk Support Size		
<input type="checkbox"/> LATXL	Size XS (3.25" x 4.25")	\$100
<input type="checkbox"/> LATSL	Size S (3.5" x 5.25")	\$100
<input type="checkbox"/> LATML	Size M (3.75" x 5.75")	\$100
<input type="checkbox"/> LATLL	Size L (4.25" x 6.5")	\$100

<input type="checkbox"/> SALATR	Right MAXX Style Swing Away HDW ⁽¹⁾	\$195
Select Right Trunk Support Size		
<input type="checkbox"/> LATXR	Size XS (3.25" x 4.25")	\$100
<input type="checkbox"/> LATSR	Size S (3.5" x 5.25")	\$100
<input type="checkbox"/> LATMR	Size M (3.75" x 5.75")	\$100
<input type="checkbox"/> LATLR	Size L (4.25" x 6.5")	\$100

<input type="checkbox"/> ATLL	Additional Telescoping Link - Left	\$75
<input type="checkbox"/> ATLR	Additional Telescoping Link - Right	\$75

- Up to 2.75" medial offset
- Recommended for all backs except rehab back, PB and elite deep backs
- Medial offset depends on back thickness
- Must choose fixed mount and pads

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

HIP SUPPORTS

Lateral Hip Support with Fixed Mounting Hardware Left⁽¹⁾

Select Pad Size and mounting position at Hip:

<input type="checkbox"/> HSPXLH	Size XS (5.5" W x 4" H) mounted at Hip . . .	\$120
<input type="checkbox"/> HSPSMLH	Size S/M (7" W x 4" H) mounted at Hip . . .	\$120
<input type="checkbox"/> HSPSLH	Size S (7" W x 6" H) mounted at Hip	\$120
<input type="checkbox"/> HSPMLH	Size M (9" W x 4" H) mounted at Hip	\$120
<input type="checkbox"/> HSPLLH	Size L (13" W x 4" H) mounted at Hip	\$120

Select Pad Size and mounting position at Knee:

<input type="checkbox"/> HSPXLK	Size XS (5.5" W x 4" H) mounted at Knee . . .	\$120
<input type="checkbox"/> HSPSMLK	Size S/M (7" W x 4" H) mounted at Knee . . .	\$120
<input type="checkbox"/> HSPSLK	Size S (7" W x 6" H) mounted at Knee	\$120
<input type="checkbox"/> HSPMLK	Size M (9" W x 4" H) mounted at Knee	\$120
<input type="checkbox"/> HSPLLK	Size L (13" W x 4" H) mounted at Knee	\$120

Lateral Hip Support with Fixed Mounting Hardware Right⁽¹⁾

Select Pad Size and mounting position:

<input type="checkbox"/> HSPXRH	Size XS (5.5" W x 4" H) mounted at Hip	\$120
<input type="checkbox"/> HSPSMRH	Size S/M (7" W x 4" H) mounted at Hip	\$120
<input type="checkbox"/> HSPSRH	Size S (7" W x 6" H) mounted at Hip	\$120
<input type="checkbox"/> HSPMRH	Size M (9" W x 4" H) mounted at Hip	\$120
<input type="checkbox"/> HSPLRH	Size L (13" W x 4" H) mounted at Hip	\$120

Select Pad Size and mounting position at Knee:

<input type="checkbox"/> HSPXRK	Size XS (5.5" W x 4" H) mounted at Knee . . .	\$120
<input type="checkbox"/> HSPSMRK	Size S/M (7" W x 4" H) mounted at Knee . . .	\$120
<input type="checkbox"/> HSPSRK	Size S (7" W x 6" H) mounted at Knee	\$120
<input type="checkbox"/> HSPMRK	Size M (9" W x 4" H) mounted at Knee	\$120
<input type="checkbox"/> HSPLRK	Size L (13" W x 4" H) mounted at Knee	\$120

1. Must choose swing away hardware below if required.

REMOVABLE HIP SUPPORT HARDWARE

Lift Off Removable Hardware for Hip Support

<input type="checkbox"/> LOHL	Left	\$195
<input type="checkbox"/> LOHR	Right	\$195

MAXX Style Quick Release, Removable, Multi-Axis Mounting Hardware for Hip Supports

Select mounting position:

<input type="checkbox"/> MHSHLH	Left mounted at Hip	\$215
<input type="checkbox"/> MHSHLK	Left mounted at Knee	\$215
<input type="checkbox"/> MHSHRH	Right mounted at Hip	\$215
<input type="checkbox"/> MHSHRK	Right mounted at Knee	\$215

Swing Away Removable Hardware for Hip Supports

Select mounting position:

<input type="checkbox"/> SHSHLH	Left mounted at Hip	\$215
<input type="checkbox"/> SHSHLK	Left mounted at Knee	\$215
<input type="checkbox"/> SHSHRH	Right mounted at Hip	\$215
<input type="checkbox"/> SHSHRK	Right mounted at Knee	\$215

Extended Maxx Style Quick Release, Removable, Multi-Axis Mounting Hardware for Hip Supports⁽¹⁾

Select mounting position: