



Invacare Storm Torque RX20 Power Wheelchair with Captains and Solid Seat

US PRICE LIST AND ORDER FORM - Price Effective March 2, 2026



CONTACT MOTION CONCEPTS CUSTOMER SERVICE: 1.888.433.6818 / QUOTES@MOTIONCONCEPTS.COM

REQUIRED INFORMATION

Quote Order

Date of Order: _____ Dealer Account #: _____

Dealer Name: _____

PO #: _____ Tag _____

Purchasing Contact: _____

Phone: _____ Fax: _____

E-mail: _____

ATP/Therapist: _____

Ship to Address: _____

City: _____ State: _____

Zip: _____ Client Gender: M F

Special Client Conditions: _____

REQUIRED MEASUREMENTS

To ensure system is accurately configured please fill in all required measurements.

Client Height:		
Client Weight:		
Use the table below to calculate your Seat Pan to Footplate Measurement. Then select the appropriate range in the legrest section.	Solid Seat	Captain Seat
Knee to Heel:		
Cushion Thickness: (Subtract)		4"
Cushion Compression:(Add)		

Notes:

HCPCS CODES	K0848, K0849
BASE WIDTH NARROW	22.4"
BASE WIDTH STANDARD	24.4"
BASE LENGTH	33.5"
DRIVE WHEEL SIZE	14"
CASTER WHEEL SIZE	8"
BASE ONLY TURNING RADIUS	24.8"
DRIVETRAIN	4-pole Sealed Housing
RANGE ⁽¹⁾	19.7 miles

BATTERIES NARROW	M50
BATTERIES STANDARD	Group 24
BATTERY CHARGER	8A, Off Board
WEIGHT CAPACITY NARROW	300 lbs
WEIGHT CAPACITY STANDARD	300 lbs
GROUND CLEARANCE	3.2"
MAXIMUM SPEED NARROW	5.8 mph
MAXIMUM SPEED STANDARD	7.5 mph

⁽¹⁾ Actual driving range and speed may vary due to factors such as user weight, type and grade of terrain, battery condition & charge level; type and condition of the drive wheels, and variations in the drive and control systems.

SOLID SEAT



CAPTAIN SEAT



Invacare Storm Torque RX20
Captain and Solid Seat

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

INVACARE POWER BASES

NEW!

- SRX-20C Storm **Torque** RX20 Base for Captain's Seat
HCPCS code K0849 \$8,705
- SRX-20S Storm **Torque** RX20 Base with Solid Seat Pan
HCPCS code K0848 \$8,705



- Enhanced Front Caster Housing and Bearings
- New Front Caster Design
- Improved double-sided Front Forks
- Upgraded Rear Suspension
- Optimized Linx Parameters

SELECT SEAT TO FLOOR HEIGHT (inches)	STF Height	Captain's Seat		Solid Seat	
		No Elevate	With Elevate	No Elevate	With Elevate
<input type="checkbox"/> RSLW <input type="checkbox"/> RSLW-S	Super Low ⁽¹⁾	20 1/4	22	16 1/2	18 1/4
<input type="checkbox"/> RLOW <input type="checkbox"/> RLOW-S	Low	20 3/4	22 1/2	17	18 3/4
<input type="checkbox"/> RMED <input type="checkbox"/> RMED-S	Medium	21 3/4	23 1/2	18	19 3/4
<input type="checkbox"/> RTALL <input type="checkbox"/> RTALL-S	Tall	22 3/4	24 1/2	19	20 3/4

1. RSLW / RSLW-S \$363 up charge. Not available on Narrow Base

SEAT STYLE

- SEMIR Semi-Recline (35°) Captain's Seat STD
- SEMIRB Semi-Recline (35°) with Solid Seat \$145

SEAT COVER

- BSC Black Slip Cover \$573

1. Includes seat, back, and headrest cover

USER WEIGHT LIMITS

- UW250 Weight Capacity up to 250 lbs⁽¹⁾ STD
- U300 Weight Capacity up to 300 lbs STD

1. Required for Narrow Frame Base.

BASE WIDTH

- NARROWFRM Narrow Base 22.4"⁽¹⁾ N/C
- STDFRM Standard Base 24.4"⁽²⁾ N/C

1. Uses (2) M50 Batteries. User weight limited to 250 lbs.

2. Uses (2) 24 Gel Batteries.

PERFORMANCE MOTORS

- HSMTR Performance Motors - 7.5 mph⁽¹⁾⁽²⁾⁽³⁾ \$1,470

1. Must order STDFRM Battery Tray Type. Not available on Narrow base.

2. Extended high speed use on rugged terrain is not recommended.

3. Performance/High Speed Motors may have less torque and may have less range than standard speed motors.

TRANSPORT TIE DOWN

- TRBKTS Wheelchair Transport Brackets⁽¹⁾ STD

1. For unoccupied use only

TIRE OPTIONS

- AB1431-3 14" x 3" Black Tire with Gel Foam Inserts STD
- AB1430-3 14" x 3" Black Tire - Pneumatic Tire N/C

CASTER OPTION

NEW!

- AB1228-3 8" x 2.5" Black Ribbed Casters w/ Solid Inlay N/C

INSTALLATION OF BATTERIES

- M50-12 SLD M MK Group 50 Battery (ea) Quantity()
HCPCS code E2361 \$355
- M24 SLD G FT MK Group 24 Gel Battery (ea) Quantity()
HCPCS code E2363 \$529

BATTERY CHARGER OPTION

- 110CHARGER 110 Volt Battery Charger N/C
- L8900 Omit Battery Charger (\$84)

FRAME FINISH

- 170P Obsidian Black N/C
- 171P Ruby Red N/C
- 172P Peacock Blue N/C
- 173P Pink Flamingo N/C
- 174P Purple Passion N/C
- 175P Denim Navy N/C
- 176P Lime Twist N/C
- 177P Arctic Blue N/C

RIM INSERT COLORS - MUST PICK ONE

- 170PR Obsidian Black N/C
- 171PR Ruby Red N/C
- 172PR Peacock Blue N/C
- 173PR Pink Flamingo N/C
- 174PR Purple Passion N/C
- 175PR Denim Navy N/C
- 176PR Lime Twist N/C
- 177PR Arctic Blue N/C
- 178PR Banana Split N/C
- 179PR Orange Marmalade N/C
- 180PR Morning Mist N/C

ELEVATING SEAT

- ESM12C Elevate lift 12" \$3,504
Weight capacity up to 300lbs

- RDE Reduced Drive for Elevating Seat⁽¹⁾ \$580
1. If reduced drive is not selected the system will be in drive lock-out when elevated.

ELECTRONICS FOR ELEVATING SEAT

- SFCB Single Function Electronics⁽¹⁾ STD
- Select Switch:
- SPBSO Single Push Button
 - STO Single Toggle Switch
- SFCBA Single/Dual Function Drive Control⁽¹⁾
HCPCS code E2310 \$2,190

**Invacare Storm Torque RX20
Captain and Solid Seat**

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

ADDITIONAL SWITCH FOR POWER SEATING

- SPBSA Single Push Button \$218
- STA Single Toggle Switch..... \$218
- ASL320R ASL Pulse Mechanical Switch - Red..... \$181
- ASL320G ASL Pulse Mechanical Switch - Green \$181
- ASL320BK ASL Pulse Mechanical Switch - Black..... \$181
- SESG Stealth Egg Switch - Green \$196
- SESB Stealth Egg Switch - Black..... \$196

MOUNTING FOR POWER SEATING SWITCH

- ESML Left Side (Armpad)
- ESMR Right Side (Armpad)
- ESMJL Joystick Left (Mounts to Joystick Bracket)
- ESMJR Joystick Right (Mounts to Joystick Bracket)
- ESMD Do not mount

SEAT WIDTH SETTING

- W18 18" Wide (Actual Width 17.7" Solid / 19" Captain) N/C
- W20 20" Wide (Actual Width 18.5" Solid / 20" Captain) N/C
- W22 22" Wide (Actual Width 20.5" Solid / 22" Captain) N/C

SEAT DEPTH SETTING

- D16 16" Deep N/C
- D17 17" Deep N/C
- D18 18" Deep N/C
- D19 19" Deep N/C
- D20 20" Deep..... N/C
- D21 21" Deep N/C

Seat Width and Depth Compatibility

	16D	17D	18D	19D	20D	21
18W	✓	✓	✓	✓	N/A	N/A
20W	N/A	N/A	N/A	✓	✓	✓
22W	N/A	N/A	N/A	✓	✓	✓

NON-EXPANDABLE SINGLE DRIVE JOYSTICKS

- REM215N LED Non-Expandable Remote..... N/C

EXPANDABLE JOYSTICKS AND CONTROLLER

- EXPC Expandable Controller
HCPCS code E2377 \$775
- PWH Harness Required for Expandable System
HCPCS code E2313 \$595
- GTRAC-LX LiNX G-TRAC Module
HCPCS code K0108 \$1,045
- REM216N LED Expandable Remote with Lights⁽¹⁾ N/C
- REM400 Color 3.5" Touch Screen Remote⁽¹⁾ \$1,050
- SW400TGL Toggle Kit for REM400⁽²⁾ \$415

1. Must choose EXPC Expandable Controller and PWH Harness.
2. Adds 2 toggle switches to REM400.

MOUNTING FOR JOYSTICKS

- FR Fixed Joystick Mounting - Right STD
- FL Fixed Joystick Mounting - Left..... STD
- MRR Maxx Resolve Swing-Away Joystick Mount
(Right Hand) HCPCS code E1032 \$573
- MRL Maxx Resolve Swing-Away Joystick Mount
(Left Hand) HCPCS code E1032 \$573
- SAQR Motion Height Adjustable Swing-Away
Quad Link (Right Hand) HCPCS code E1032..... \$341
- SAQL Motion Height Adjustable Swing-Away
Quad Link (Left Hand) HCPCS code E1032 \$341

LINX ELECTRONIC ACCESSORIES

- LAK LiNX Access Key⁽¹⁾ \$105
- MUSB Motion Dual USB Charger N/C

1. Required for programming

JOYSTICK TOPS⁽¹⁾

HCPCS code E2323

- PC101A Bodypoint U Shaped Handle 3" \$126
- PC102A Bodypoint U Shaped Handle 4" \$126
- PC107A Bodypoint Rubber Dome..... \$126
- PC110A Grooved Mushroom Joystick Handle \$126

1. To remove/disengage the joystick knob from the REM400 Remote pull straight up on the joystick (DO NOT TWIST) otherwise damage may occur and may void the warranty.

LAP BELTS

- BELT48 Push Button Style Lap Belt, Length 48" STD
- BELT60 Push Button Style Lap Belt, Length 60"..... N/C
- BELT71 Push Button Style Lap Belt, Length 71" N/C
- PLB2 Bodypoint Padded Lap Belt 2 Point⁽¹⁾ \$222
- PLB4 Bodypoint Padded Lap Belt 4 Point⁽¹⁾ \$275

1. 67" length only option

ARM RESTS

HCPCS code E0973

- LVAHF Left Full Length, Adjustable Height..... \$167
- RVAHF Right Full Length, Adjustable Height..... \$167

**Invacare Storm Torque RX20
Captain and Solid Seat**

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

MATRIX SEAT CUSHIONS - FOR SOLID SEAT SYSTEMS

Matrix PS Seat Cushion 14" - 20" wide/reversible cover
 HCPCS code E2605. \$383
 PS1816 PS1818 PS2018 PS2020

Matrix PS Seat Cushion 21"-22" wide
 HCPCS code E2606. \$605
 PS2218 PS2220

Matrix PSP Seat Cushion 14" - 20" wide/reversible cover
 HCPCS code E2607. \$438
 PSP1816 PSP1818 PSP2018 PSP2020

Matrix PSP Seat Cushion 21"-22" wide
 HCPCS code E2608. \$660
 PSP2218

Matrix Vi Seat Cushion 14" - 20" wide
 HCPCS code E2607. \$497
 MA1816-VI MA1818-VI MA2018-VI MA2020-VI

Matrix Vi Seat Cushion 21" - 22" wide
 HCPCS code E2608. \$556
 MA2218-VI MA2220-VI

Matrix Libra Seat Cushion 14" - 20" wide
 HCPCS code E2624. \$567
 LC1816 LC1818 LC2018 LC2020

Matrix Libra HD Seat Cushion 21" wide
 HCPCS code E2624. \$696
 LC2218 LC2220

Matrix Libra Adjust Seat Cushion 14" -20" wide
 HCPCS code E2624. \$731
 LC1816-ADJ LC1818-ADJ LC2018-ADJ LC2020-ADJ

Matrix Libra Adjust Seat Cushion 22" deep
 HCPCS code E2624. \$860
 LC2218-ADJ LC2220-ADJ

Matrix Multi Cushion 14" - 20" wide w/reversible cover
 HCPCS code E2607. \$521
 MM1816 MM1818 MM2018 MM2020

Matrix Multi Cushion 22" wide w/reversible cover
 HCPCS code E2608. \$696
 MM2218 MM2220 MM2222

OCCC "ON-CHAIR" Cushion Credit⁽¹⁾(\$105)


1. Not available on Matrix PS or PSP Cushions.

STANDARD FRONT RIGGING

- FCMP Fixed Center Mount Platform(1)(2)N/C
- FIFPN Rubber Coated Footplate (10.5" W x 10")N/C

Select seat pan to footplate:

- FSPTF3 10.25" - 14.25"N/C
- FSPTF4 14.25" - 20.25"N/C



1. 3" ground clearance required.
 2. Fixed Center Mount will be set at factory seating to prevent caster interference.

MANUAL ELEVATING CENTER MOUNT LEGREST

- MECML Manual Elevating Center Mount Legrest
 With Plantar Flex\$863

Select legrest set up position:

- 90SU Set at 90° (elevating range 90-40°) STD

Rubber Coated Footplate:

- MRCPN Rubber Coated Footplate 10.5" W x 10" D \$165

Select seat pan to footplate for One Piece Foot Plate

- CSPTF1 12.75" - 16.25"N/C
- CSPTF2 15.25" - 20.25"N/C

Individual Footplates:


- IFSL Small - (5" x 7.5" deep) - Left \$239
- IFSR Small - (5" x 7.5" deep) - Right \$239
- IFML Medium - (5.5" x 9.5" deep) - Left \$239
- IFMR Medium - (5.5" x 9.5" deep) - Right \$239
- IFLL Large - (6" x 11.5" deep) - Left \$239
- IFLR Large - (6" x 11.5" deep) - Right \$239

Select seat pan to footplate for Individual Footplates

- CSPTF3 10" - 15"N/C
- CSPTF4 12.5" - 19.5"N/C

Optional:

- CICP Individual calf pads \$290



1. Seat pan to footplate measurements shown are for Solid Seat systems.
 Captain Seat measurements are an additional 4"

