



ROVI X3 / A3 Power Wheelchair with Captains Seat

US PRICE LIST AND ORDER FORM
Price Effective January 1, 2025



CONTACT MOTION CONCEPTS CUSTOMER SERVICE: 1.888.433.6818 / QUOTES@MOTIONCONCEPTS.COM

REQUIRED INFORMATION

Quote Order

Date of Order: _____ Dealer Account #: _____

Dealer Name: _____

PO #: _____ Tag _____

Purchasing Contact: _____

Phone: _____ Fax: _____

E-mail: _____

ATP/Therapist: _____

Ship to Address: _____

City: _____ State: _____

Zip: _____ Client Gender: M F

Special Client Conditions: _____

REQUIRED MEASUREMENTS

To ensure system is accurately configured please fill in all required measurements.

Client Height:		
Client Weight:		
Use the table below to calculate your Seat Pan to Footplate Measurement. Then select the appropriate range in the legrest section.	Solid Seat	Captain Seat
Knee to Heel:		
Cushion Thickness: (Subtract)		4"
Cushion Compression:(Add)		
Seat Pan to Footplate:(Total)		

Notes:

BASE	X3	A3
HCPCS CODES	K0856, K0861	K0848, K0849
BASE WIDTH	23.25"	25.5"
BASE LENGTH	36.5"	36.5"
DRIVE WHEEL SIZE ⁽¹⁾	14"	14"
CASTER WHEEL SIZE	6" (front/rear)	6" (front/rear)
BASE ONLY TURNING RADIUS	19.5"	20.5"
DRIVETRAIN	4-pole Motors	4-pole Motors
BATTERIES	2 x M34	2 x M34

BASE	X3	A3
BATTERY CHARGER	8A, 110V	8A, 110V
BATTERIES WEIGHT	38.5 lbs (each)	38.5 lbs (each)
RANGE ⁽²⁾	17.7 miles	17.7 miles
SUSPENSION SYSTEM	Active Ride Control (ARC)	
WEIGHT CAPACITY	300 lbs	300lbs
BASE WEIGHT	157 lbs	159 lbs
GROUND CLEARANCE	3"	3"
MAXIMUM SPEED	6.2 mph	6.2 mph

⁽¹⁾ All A3/X3 wheels come with solid inserts

⁽²⁾ Actual driving range and speed may vary due to factors such as user weight, type and grade of terrain, battery condition & charge level; type and condition of the drive wheels, and variations in the drive and control systems.

SOLID SEAT



CAPTAIN SEAT



ROVI POWER BASE

888WNLEX3-SEMIR
 ROVI X3 Base for Captain Seat
 HCPCS code K0849. \$7,395

888WNLEX3-SEMIRB
 ROVI X3 Base for Solid Seat
 HCPCS code K0848. \$7,395

888WNLEA3-SEMIR
 ROVI A3 Base for Captain Seat
 HCPCS code K0849. \$7,595

888WNLEA3-SEMIRB
 ROVI A3 Base for Solid Seat
 HCPCS code K0848. \$7,595

Select seat to floor height: (inches)		STF Height	Captain's Seat		Solid Seat	
			No Elevate	With Elevate	No Elevate	With Elevate
<input type="checkbox"/> TLOW	<input type="checkbox"/> TLOW-S	Low	20 ^{3/4}	22 ^{1/2}	17	18 ^{3/4}
<input type="checkbox"/> TMED	<input type="checkbox"/> TMED-S	Medium	21 ^{3/4}	23 ^{1/2}	18	19 ^{3/4}
<input type="checkbox"/> TTALL	<input type="checkbox"/> TTALL-S	Tall	22 ^{3/4}	24 ^{1/2}	19	20 ^{3/4}
<input type="checkbox"/> TXTALL	<input type="checkbox"/> TXTALL-S	X-Tall	23 ^{3/4}	N/A	20	N/A

SEAT STYLE

SEMIR Semi-Recline (35°) Captain's Seat STD
 SEMIRB Semi-Recline (35°) with Solid Seat \$100

SEAT COVER

BSC Black Slip Cover \$395
 1. Includes seat, back, and headrest cover

PERFORMANCE MOTORS

HSM62 Optional 6.2 mph High Speed Motors. \$950
 Must order OnTraxx Enhanced Driving Module
 OEDM OnTraxx Enhanced Driving Module
 HCPCS code K0108 \$750

TRANSPORT TIE DOWN

XBKTS Wheelchair Transport Brackets⁽¹⁾ N/C
 1. For unoccupied use only

TIRE OPTIONS

STIRES Standard Solid Tires STD
 PTIRES Pneumatic Tires N/C

INSTALLATION OF BATTERIES

HCPCS code E2359
 M34 SLD G MK Group 34 Battery (ea) Quantity () . . . \$305

LIGHTED FENDER

Lights operate through single push button

- BWFL  Blizzard White (Metallic)
- OCFL  Orange Crush
- MGFL  Motion Green
- BCFL  Black Carbon (Matte)
- BPFL  Black Panther (Gloss)
- CARFL  Candy Apple Red (Metallic)
- GGFL  Graphite Grey (Metallic)
- IBFL  Indigo Blue
- MYFL  Mellow Yellow
- NBFL  Nordic Blue
- RDPFL  Razzle Dazzle Pink (Metallic)
- PRFL  Purple Reign

MATCHING BATTERY DOORS

MBD Matching Battery Doors. \$495
 (Select fender color first)

CUSTOM DESIGNS

Please see form number TRD0607 - ROVI Custom Designs Price Sheet, for Custom/Dipped Fender colors, battery door and other personalized options.

ELEVATING SEAT

ESM12C Elevate lift 12" HCPCS code E2298. \$2,995
 Weight capacity up to 300lbs
 RDE Reduced Drive for Elevating Seat⁽¹⁾ \$400

1. If reduced drive is not selected the system will be in drive lock-out when elevated.

ELECTRONICS FOR ELEVATING SEAT

SFCB Single Function Electronics STD
 Select Switch:
 SPBSO Single Push Button
 STO Single Toggle Switch
 SFCBA Single/Dual Function Drive Control⁽¹⁾
 HCPCS code E2310 \$1,795

ADDITIONAL SWITCH FOR POWERED SEATING

- ASL320R ASL Pulse Mechanical Switch - Red \$125
- ASL320G ASL Pulse Mechanical Switch - Green \$125
- ASL320BK ASL Pulse Mechanical Switch - Black \$125
- SESG Stealth Egg Switch - Green \$135
- SESB Stealth Egg Switch - Black \$135

MOUNTING FOR POWERED SEATING SWITCH

- ESML Left Side (Armpad)
- ESMR Right Side (Armpad)
- ESMJL Joystick Left (Mounts to Joystick Bracket)
- ESMJR Joystick Right (Mounts to Joystick Bracket)
- ESMD Do not mount

SEAT WIDTH SETTING

- W18 18" Wide (Actual Width 17.7" Solid / 19" Captain) . . N/C
- W20 20" Wide (Actual Width 18.5" Solid / 20" Captain) . . N/C
- W22 22" Wide (Actual Width 20.5" Solid / 22" Captain) . . N/C

SEAT DEPTH SETTING

- D16 16" Deep N/C
- D17 17" Deep N/C
- D18 18" Deep N/C
- D19 19" Deep N/C
- D20 20" Deep N/C
- D21 21" Deep N/C

Seat Width and Depth Compatibility						
	16D	17D	18D	19D	20D	21D
18W	✓	✓	✓	✓	N/A	N/A
20W	N/A	N/A	N/A	✓	✓	✓
22W	N/A	N/A	N/A	✓	✓	✓

JOYSTICK SELECTION

- D50677 VR2 Non Expandable Joystick N/C
- D51122 R-Net LED Joystick⁽¹⁾ N/C
- D51165 R-Net Color Remote Joystick with Monojacks⁽¹⁾ \$600
- D51711 R-Net Advanced Joystick⁽¹⁾⁽²⁾ \$995

1. Requires Expandable Controller and Harness (D51109 and REH)
 2. 3.5" high-resolution VGA colour display, IR Control as standard, includes integrated Bluetooth Module.

EXPANDABLE CONTROLLER⁽¹⁾

- D51109 R-Net Expandable Controller
HCPCS code E2377 \$615
- REH R-Net Harness for Expandable Controller
HCPCS code E2313 \$595

1. Required with R-Net Joysticks

MOUNTING FOR DRIVER CONTROL

- FR Fixed Joystick Mounting - Right STD
- FL Fixed Joystick Mounting - Left STD
- MRR Maxx Resolve Swing-Away Joystick Mount (Right Hand) \$395
- MRL Maxx Resolve Swing-Away Joystick Mount (Left Hand) \$395
- SAQR Motion Height Adjustable Swing-Away Quad Link (Right Hand) \$235
- SAQL Motion Height Adjustable Swing-Away Quad Link (Left Hand) \$235

ELECTRONIC OPTIONS

- D50883 R-Net Output Module (EMC)⁽¹⁾ \$995
- D51580 R-Net Bluetooth Module IOS Apple⁽¹⁾ \$895
- D51111 R-Net Bluetooth Module Android/PC⁽¹⁾ \$895
- SPH Smart Phone Holder \$50
- MUSB Motion Dual USB Charger N/C

1. Requires D51165 or D51711 Joystick.

OPTIONAL JOYSTICK KNOB

HCPCS code E2323

- SA80941 "T" Handle for Joystick Knob \$115
- SA72441 Ball Handle for Joystick Knob \$115
- PC101B-20 BodyPoint Joystick Handle 3" U-Shaped \$120
- PC102B-20 BodyPoint Joystick Handle 4" U-Shaped \$120
- PC104B-20 BodyPoint Straight Handle \$120
- PC107B-20 BodyPoint Joystick Handle Domed Rubber \$120

LAP BELTS

- BELT48 Push Button Style Lap Belt, Length 48" STD
- BELT60 Push Button Style Lap Belt, Length 60" N/C
- BELT71 Push Button Style Lap Belt, Length 71" N/C
- PLB2 Bodypoint Padded Lap Belt 2 Point⁽¹⁾ \$190
- PLB4 Bodypoint Padded Lap Belt 4 Point⁽¹⁾ \$235

1. 67" length only option

ARMRESTS

HCPCS code E0973

- LVAHF Left Full Length, Adjustable Height \$115
- RVAHF Right Full Length, Adjustable Height \$115

MATRIX SEAT CUSHIONS - FOR SOLID SEAT SYSTEMS

Matrix PS Seat Cushion 14" - 20" wide/reversible cover
 HCPCS code E2605 \$345
 PS1816 PS1818 PS2018 PS2020

Matrix PS Seat Cushion 21"-22" wide
 HCPCS code E2606 \$545
 PS2218 PS2220

Matrix PSP Seat Cushion 14" - 20" wide/reversible cover
 HCPCS code E2607 \$395
 PSP1816 PSP1818 PSP2018 PSP2020

Matrix PSP Seat Cushion 21"-22" wide
 HCPCS code E2608 \$595
 PSP2218

Matrix Vi Seat Cushion 14" - 20" wide
 HCPCS code E2607 \$425
 MA1816-VI MA1818-VI MA2018-VI MA2020-VI

Matrix Vi Seat Cushion 21" - 22" wide
 HCPCS code E2608 \$475
 MA2218-VI MA2220-VI

Matrix Libra Seat Cushion 14" - 20" wide
 HCPCS code E2624 \$485
 LC1816 LC1818 LC2018 LC2020

Matrix Libra HD Seat Cushion 22" wide
 HCPCS code E2625 \$595
 LC2218 LC2220

Matrix Libra Adjust Seat Cushion 14" -20" wide
 HCPCS code E2624 \$625
 LC1816-ADJ LC1818-ADJ LC2018-ADJ LC2020-ADJ

Matrix Libra Adjust HD Seat Cushion 22" wide
 HCPCS code E2625 \$735
 LC2218-ADJ LC2220-ADJ

Matrix Multi Cushion 14" - 20" wide/reversible cover
 HCPCS code E2607 \$445
 MM1816 MM1818 MM2018 MM2020

Matrix Multi Cushion 22" wide/reversible cover
 HCPCS code E2608 \$595
 MM2218 MM2220 MM2222

OCCC "ON-CHAIR" Cushion Credit⁽¹⁾ (\$90)

1. Not available on Matrix PS or PSP Cushions.

STANDARD FRONT RIGGING


FCMP Fixed Center Mount Platform⁽¹⁾ N/C
 FIFP Rubber Coated Footplate (11.5" W x 10") N/C
 FBIFP Large Rubber Coated Footplate (17" w x 13"d) . . \$250

Select seat pan to footplate:

FSPTF3 10.25" - 14.25" N/C
 FSPTF4 14.25" - 20.25" N/C

Set at:

FSA83 Set at 83° N/C
 FSA90 Set at 90°⁽²⁾ N/C
 FSA97 Set at 97°⁽²⁾ N/C



1. 3" ground clearance required.
 2. Not available on 18" W seats.

MANUAL ELEVATING CENTER MOUNT LEG REST

MECML Manual Elevating Center Mount Legrest
 With Plantar Flex \$595

Select legrest set up position:

90SU Set at 90° (elevating range 90-40°) N/C
 110TU Tucked at 110° (elevating range 110-60°)⁽¹⁾ . . N/C

Rubber Coated Footplate:

MRCP Rubber Coated Footplate (11.5" W x 10") \$165

Select seat pan to footplate for Rubber Coated Footplate

CSPTF1 12.75" - 16.25" N/C
 CSPTF2 15.25" - 20.25" N/C

Individual Footplates:


IFSL Small - (5" x 7.5" deep) - Left \$165
 IFSR Small - (5" x 7.5" deep) - Right \$165
 IFML Medium - (5.5" x 9.5" deep) - Left \$165
 IFMR Medium - (5.5" x 9.5" deep) - Right \$165
 IFLL Large - (6" x 11.5" deep) - Left \$165
 IFLR Large - (6" x 11.5" deep) - Right \$165

Select seat pan to footplate for Individual Footplates

CSPTF3 10" - 15" N/C
 CSPTF4 12.5" - 19.5" N/C

Optional:

CICP Individual calf pads \$200



1. Not available on 18" W seats.

MAXX STYLE FRONT RIGGING


- MSFCMFP Maxx Style Fixed Center Mount Platform⁽¹⁾⁽²⁾ . \$395
 - MIFFP Rubber Coated Footplate (11.5" W x 10") N/C
 - MBIFP Large Rubber Coated Footplate (17" w x 13"d) . . \$250
- MICP Individual calf pads. N/C

Select seat pan to footplate:

- MSPTF3 10.25" - 14.25" N/C
- MSPTF4 14.25" - 20.25" N/C

Set at:

- MSA83 Set at 83° N/C
- MSA90 Set at 90°⁽³⁾ N/C
- MSA97 Set at 97°⁽³⁾ N/C



1. Depth adjustable, and easy slide adjustment for seat pan to foot plate extension.
 2. 3" ground clearance required.
 3. Not available on 18" W seats.

FOOT PLATFORM OPTIONS ⁽¹⁾



- SFFP Small Foot Platform (12" x 8.5" deep) \$180
- MFFP Medium Platform (12" x 10.75" deep) \$180
- LFFP Large Platform (14.5" x 13" deep) \$250
- XLFFP Extra Large Platform (17" x 13" deep) \$250
- XXLFFP Extra Extra Large Platform (21" x 13" deep) . \$250

1. May need to adjust leg position to avoid interference with front casters

ADDITIONAL ACCESSORIES

- PH X-Grip Phone Holder \$95
- UCH Universal Cup Holder ⁽¹⁾ \$35
- CSCC Cane and Crutch Holder \$395
- CSO2 Non-Articulating O² Gas Holder \$250
- CFH Flag Holder \$65
- CRVMR Rear View Mirror Right Mounted \$45
- CRVML Rear View Mirror Left Mounted \$45
- TAG571850 Captain Seat Wheelchair Cover \$105
- AMC Accessory Mount with Clip \$50
- CSWH Captain Seat Walker Holder ⁽²⁾ \$200

1. Not available on the joystick side.
 2. Maximum load is 20 lbs

MILITARY DECAL⁽¹⁾

- DECALA Army Decal 3" Diameter \$50
- DECALF Air Force Decal 3" Diameter \$50
- DECALN Navy Decal 3" Diameter \$50
- DECALM Marines Decal 3" Diameter \$50



1. Comes uninstalled off chair.

NOTES

Total Retail Value \$ _____

SUBMIT

PRINT FORM

SAVE

CLEAR FORM