

Date of Order: _____	Ship to Address: _____
Dealer Name: _____	_____
Dealer Account #: _____	City: _____
PO #: _____	State: _____ Zip: _____
Tag: _____	Comments: _____
Purchasing Contact: _____	_____
Phone: _____ Fax: _____	_____

Pediatric Backs

► Page

Invacare [®] Matrx [®] Mini Back/Extra Outer Cover	2
--	---

Standard Backs

► Page

Invacare [®] Matrx [®] PB Back/Extra Outer Cover	3
Invacare [®] Matrx [®] PB Deep Back/Extra Outer Cover	3
Invacare [®] Matrx [®] Elite Back/Extra Outer Cover	6/7
Invacare [®] Matrx [®] Elite Full Back/Extra Outer Cover	7
Invacare [®] Matrx [®] Elite Tilt/Recline Back/Extra Outer Cover	8
Invacare [®] Matrx [®] Elite Deep Back/Extra Outer Cover	9
Hardware Upgrades for Matrx Elite and PB Backs	12/13
Invacare [®] Matrx [®] MX1 Back/Extra Outer Cover	16/17
Hardware Upgrades for Matrx MX1 Back	18
Invacare [®] Matrx [®] Genera Back/Extra Outer Cover	18

Bariatric Backs

► Page

Invacare [®] Matrx [®] PB Heavy-Duty Back/Extra Outer Cover	4
Invacare [®] Matrx [®] PB Extra Wide Heavy-Duty Back/Extra Outer Cover	4
Invacare [®] Matrx [®] PB Deep Heavy-Duty Back/Extra Outer Cover	5
Invacare [®] Matrx [®] PB Extra Wide Deep Heavy-Duty Back/Extra Outer Cover	5
Invacare [®] Matrx [®] Elite Heavy-Duty Back/Extra Outer Cover	10
Invacare [®] Matrx [®] Elite Extra Wide Heavy-Duty Back/Extra Outer Cover	10
Hardware Upgrades for Matrx Elite and PB Backs	12/13

Modifications & Accessories

► Page

Invacare [®] Matrx [®] Elite Back Laterals	11
Invacare [®] Matrx [®] Mini Back Laterals	11
Matrx [®] PB & Elite Back Options and Modifications	14/15
Headrest Supports & Hardware	19
Other Accessories	20

INVACARE® MATRX® BACKS – INTERNATIONAL PRICE LIST/ORDER FORM

MATRX® MINI BACK: Pediatric – Positioning Back

- Mini Set Hardware fits ¾", 7/8", 1" diameter back canes
- 3" Contour Depth
- Weight capacity up to 150 lbs
- Backs are drilled for headrest attachment
- Extra outer cover includes foam

- Standard Kid*ab*ra Fabric upholstery or no charge cover options of Pink or Blue



Mini Back



QUANTITY	MODEL	FITS WIDTH	FITS HEIGHT	MSRP
<input type="checkbox"/> QTY ____	MINI1007	10", 11"	7"	\$625.00
<input type="checkbox"/> QTY ____	MINI1010	10", 11"	10"	\$625.00
<input type="checkbox"/> QTY ____	MINI1012	10", 11"	12"	\$625.00
<input type="checkbox"/> QTY ____	MINI1107	11", 12"	7"	\$625.00
<input type="checkbox"/> QTY ____	MINI1110	11", 12"	10"	\$625.00
<input type="checkbox"/> QTY ____	MINI1112	11", 12"	12"	\$625.00
<input type="checkbox"/> QTY ____	MINI1207	12", 13"	7"	\$625.00
<input type="checkbox"/> QTY ____	MINI1210	12", 13"	10"	\$625.00
<input type="checkbox"/> QTY ____	MINI1212	12", 13"	12"	\$625.00
<input type="checkbox"/> QTY ____	MINI1214	12", 13"	14"	\$625.00
<input type="checkbox"/> QTY ____	MINI1307	13", 14"	7"	\$625.00
<input type="checkbox"/> QTY ____	MINI1310	13", 14"	10"	\$625.00
<input type="checkbox"/> QTY ____	MINI1312	13", 14"	12"	\$625.00
<input type="checkbox"/> QTY ____	MINI1314	13", 14"	14"	\$625.00
<input type="checkbox"/> QTY ____	MINI1407	14", 15"	7"	\$625.00
<input type="checkbox"/> QTY ____	MINI1410	14", 15"	10"	\$625.00
<input type="checkbox"/> QTY ____	MINI1412	14", 15"	12"	\$625.00
<input type="checkbox"/> QTY ____	MINI1414	14", 15"	14"	\$625.00

COVER OPTIONS:

SPECIFY: Kid*ab*ra Pink Blue

Mini Extra Outer Cover (Comes standard with foam)



MODEL	FITS WIDTH	FITS HEIGHT	MSRP
MC	REFER TO BACK SIZES ABOVE		\$135.00

SPECIFY: WIDTH ____ HEIGHT ____ QTY ____

COVER OPTIONS:

SPECIFY: Kid*ab*ra Pink Blue

INVACARE® MATRX® BACKS – INTERNATIONAL PRICE LIST/ORDER FORM

MATRX® PB BACK: Standard – Positioning Back

- Easy Set Hardware fits ¾", 7/8", 1" diameter back canes
- See Page 12 for other no-charge hardware options
- 5" Contour Depth
- Weight capacity up to 300 lbs
- Backs include reinforcement plate for headrest

- Standard Meshtex Fabric upholstery



PB Back



QUANTITY	MODEL	FITS WIDTH	FITS HEIGHT	MSRP
<input type="checkbox"/> QTY ___	MPB1412	13"-15"	12"	\$625.00
<input type="checkbox"/> QTY ___	MPB1416	13"-15"	16"	\$625.00
<input type="checkbox"/> QTY ___	MPB1612	15"-17"	12"	\$625.00
<input type="checkbox"/> QTY ___	MPB1616	15"-17"	16"	\$625.00
<input type="checkbox"/> QTY ___	MPB1620	15"-17"	20"	\$625.00
<input type="checkbox"/> QTY ___	MPB1712	16"-18"	12"	\$625.00
<input type="checkbox"/> QTY ___	MPB1716	16"-18"	16"	\$625.00
<input type="checkbox"/> QTY ___	MPB1720	16"-18"	20"	\$625.00
<input type="checkbox"/> QTY ___	MPB1812	17"-19"	12"	\$625.00
<input type="checkbox"/> QTY ___	MPB1816	17"-19"	16"	\$625.00
<input type="checkbox"/> QTY ___	MPB1820	17"-19"	20"	\$625.00
<input type="checkbox"/> QTY ___	MPB2012	19"-21"	12"	\$625.00
<input type="checkbox"/> QTY ___	MPB2016	19"-21"	16"	\$625.00
<input type="checkbox"/> QTY ___	MPB2020	19"-21"	20"	\$625.00

PB Extra Outer Cover

MODEL	FITS WIDTH	FITS HEIGHT	MSRP
PBC	REFER TO BACK SIZES ABOVE		\$100.00
SPECIFY:	WIDTH _____	HEIGHT _____	<input type="checkbox"/> QTY _____

MATRX® PB DEEP BACK: Standard – Positioning Back

- 7" Contour Depth

PB Deep Back



QUANTITY	MODEL	FITS WIDTH	FITS HEIGHT	MSRP
<input type="checkbox"/> QTY ___	PBD1412	13"-15"	12"	\$645.00
<input type="checkbox"/> QTY ___	PBD1416	13"-15"	16"	\$645.00
<input type="checkbox"/> QTY ___	PBD1612	15"-17"	12"	\$645.00
<input type="checkbox"/> QTY ___	PBD1616	15"-17"	16"	\$645.00
<input type="checkbox"/> QTY ___	PBD1620	15"-17"	20"	\$645.00
<input type="checkbox"/> QTY ___	PBD1712	16"-18"	12"	\$645.00
<input type="checkbox"/> QTY ___	PBD1716	16"-18"	16"	\$645.00
<input type="checkbox"/> QTY ___	PBD1720	16"-18"	20"	\$645.00
<input type="checkbox"/> QTY ___	PBD1812	17"-19"	12"	\$645.00
<input type="checkbox"/> QTY ___	PBD1816	17"-19"	16"	\$645.00
<input type="checkbox"/> QTY ___	PBD1820	17"-19"	20"	\$645.00
<input type="checkbox"/> QTY ___	PBD2012	19"-21"	12"	\$645.00
<input type="checkbox"/> QTY ___	PBD2016	19"-21"	16"	\$645.00
<input type="checkbox"/> QTY ___	PBD2020	19"-21"	20"	\$645.00

PB Deep Extra Outer Cover

MODEL	FITS WIDTH	FITS HEIGHT	MSRP
PBDC	REFER TO BACK SIZES ABOVE		\$110.00
SPECIFY:	WIDTH _____	HEIGHT _____	<input type="checkbox"/> QTY _____

MATRX® PB HEAVY-DUTY BACK: Bariatric – Positioning Back

- Easy Set Hardware fits 3/4", 7/8", 1" diameter back canes
- See Page 13 for other no-charge hardware options
- 5" Contour Depth
- Weight capacity up to 500 lbs
- Backs include reinforcement plate for headrest

- Standard Meshtex Fabric upholstery



PB Heavy-Duty Back



QUANTITY	MODEL	FITS WIDTH	FITS HEIGHT	MSRP
<input type="checkbox"/> QTY ____	HDB2016	20", 21"	16"	\$835.00
<input type="checkbox"/> QTY ____	HDB2020	20"-21"	20"	\$835.00
<input type="checkbox"/> QTY ____	HDB2216	22", 23"	16"	\$835.00
<input type="checkbox"/> QTY ____	HDB2220	22", 23"	20"	\$835.00
<input type="checkbox"/> QTY ____	HDB2416	24", 25"	16"	\$835.00
<input type="checkbox"/> QTY ____	HDB2420	24", 25"	20"	\$835.00

PB Heavy-Duty Extra Outer Back Cover

MODEL	FITS WIDTH	FITS HEIGHT	MSRP
HDC	REFER TO BACK SIZES ABOVE		\$100.00
SPECIFY:	WIDTH ____	HEIGHT ____	<input type="checkbox"/> QTY ____

PB Extra Wide Heavy-Duty Back (Additional one week lead time)*



MODEL	FITS WIDTH	FITS HEIGHT	MSRP
HDBWW16	REFERENCE CHART BELOW		\$1435.00
SPECIFY:	WIDTH ____	<input type="checkbox"/> QTY ____	
HDBWW20	REFERENCE CHART BELOW		\$1435.00
SPECIFY:	WIDTH ____	<input type="checkbox"/> QTY ____	

W I D T H (WW)	
	25" 26" 27" 28" 29" 30"
HEIGHT	
16"	X X X X X X
20"	X X X X X X

PB Extra Wide Heavy-Duty Extra Outer Back Cover (Additional one week lead time)*

MODEL	FITS WIDTH	FITS HEIGHT	MSRP
HDCWW16	REFERENCE CHART ABOVE		\$200.00
SPECIFY:	WIDTH ____	<input type="checkbox"/> QTY ____	
HDCWW20	REFERENCE CHART ABOVE		\$200.00
SPECIFY:	WIDTH ____	<input type="checkbox"/> QTY ____	

* Made to order and non-refundable

MATRX® PB DEEP HEAVY-DUTY BACK: Bariatric – Positioning Back

- Easy Set Hardware fits ¾", 7/8", 1" diameter back canes
- See Page 13 for other no-charge hardware options
- 7" Contour Depth
- Weight capacity up to 500 lbs
- Backs include reinforcement plate for headrest

- Standard Meshtex Fabric upholstery



PB Deep Heavy-Duty Back



QUANTITY	MODEL	FITS WIDTH	FITS HEIGHT	MSRP
<input type="checkbox"/> QTY ____	DHD2016	20",21"	16"	\$935.00
<input type="checkbox"/> QTY ____	DHD2020	20"-21"	20"	\$935.00
<input type="checkbox"/> QTY ____	DHD2216	22", 23"	16"	\$935.00
<input type="checkbox"/> QTY ____	DHD2220	22", 23"	20"	\$935.00
<input type="checkbox"/> QTY ____	DHD2416	24", 25"	16"	\$935.00
<input type="checkbox"/> QTY ____	DHD2420	24", 25"	20"	\$935.00

Deep Heavy-Duty Extra Outer Back Cover

MODEL	FITS WIDTH	FITS HEIGHT	MSRP
DDC	REFER TO BACK SIZES ABOVE		\$110.00
SPECIFY:	WIDTH ____	HEIGHT ____	<input type="checkbox"/> QTY ____

PB Extra Wide Deep Heavy-Duty Back (Additional one week lead time)*



MODEL	FITS WIDTH	FITS HEIGHT	MSRP
DHDWW16	REFERENCE CHART BELOW		\$1535.00
SPECIFY:	WIDTH ____	<input type="checkbox"/> QTY ____	
DHDWW20	REFERENCE CHART BELOW		\$1535.00
SPECIFY:	WIDTH ____	<input type="checkbox"/> QTY ____	

W I D T H (WW)	
HEIGHT	
16"	25" X 26" X 27" X 28" X 29" X 30"
20"	X X X X X X X

Extra Wide Deep Heavy-Duty Extra Outer Back Cover (Additional one week lead time)*

MODEL	FITS WIDTH	FITS HEIGHT	MSRP
DDCWW16	REFER TO BACK SIZES ABOVE		\$210.00
SPECIFY:	WIDTH ____	HEIGHT ____	<input type="checkbox"/> QTY ____
DDCWW20	REFER TO BACK SIZES ABOVE		\$210.00
SPECIFY:	WIDTH ____	HEIGHT ____	<input type="checkbox"/> QTY ____

* Made to order and non-refundable

MATRX® ELITE BACK: Standard – Positioning Back

- Easy Set Hardware fits ¾", 7/8", 1" diameter back canes
- See Page 12 for other no-charge hardware options
- 3" Contour Depth
- Weight capacity up to 300 lbs
- All Elite Backs 15" and wider x 14" and higher are pre-drilled for headrest attachment, and include headrest adapter plate
- Extra outer cover includes foam

Standard Meshtex Fabric upholstery



Elite Back



QUANTITY	MODEL	FITS WIDTH	FITS HEIGHT	MSRP
<input type="checkbox"/> QTY ___	PBE1210	12"-13"	10"	\$625.00
<input type="checkbox"/> QTY ___	PBE1212	12"-13"	12"	\$625.00
<input type="checkbox"/> QTY ___	PBE1214	12"-13"	14"	\$625.00
<input type="checkbox"/> QTY ___	PBE1216	12"-13"	16"	\$625.00
<input type="checkbox"/> QTY ___	PBE1218	12"-13"	18"	\$625.00
<input type="checkbox"/> QTY ___	PBE1310	13"-14"	10"	\$625.00
<input type="checkbox"/> QTY ___	PBE1312	13"-14"	12"	\$625.00
<input type="checkbox"/> QTY ___	PBE1314	13"-14"	14"	\$625.00
<input type="checkbox"/> QTY ___	PBE1316	13"-14"	16"	\$625.00
<input type="checkbox"/> QTY ___	PBE1318	13"-14"	18"	\$625.00
<input type="checkbox"/> QTY ___	PBE1410	14"-15"	10"	\$625.00
<input type="checkbox"/> QTY ___	PBE1412	14"-15"	12"	\$625.00
<input type="checkbox"/> QTY ___	PBE1414	14"-15"	14"	\$625.00
<input type="checkbox"/> QTY ___	PBE1416	14"-15"	16"	\$625.00
<input type="checkbox"/> QTY ___	PBE1418	14"-15"	18"	\$625.00
<input type="checkbox"/> QTY ___	PBE1420	14"-15"	20"	\$625.00
<input type="checkbox"/> QTY ___	PBE1510	15"-16"	10"	\$625.00
<input type="checkbox"/> QTY ___	PBE1512	15"-16"	12"	\$625.00
<input type="checkbox"/> QTY ___	PBE1514	15"-16"	14"	\$625.00
<input type="checkbox"/> QTY ___	PBE1516	15"-16"	16"	\$625.00
<input type="checkbox"/> QTY ___	PBE1518	15"-16"	18"	\$625.00
<input type="checkbox"/> QTY ___	PBE1520	15"-16"	20"	\$625.00
<input type="checkbox"/> QTY ___	PBE1610	16"-17"	10"	\$625.00
<input type="checkbox"/> QTY ___	PBE1612	16"-17"	12"	\$625.00
<input type="checkbox"/> QTY ___	PBE1614	16"-17"	14"	\$625.00
<input type="checkbox"/> QTY ___	PBE1616	16"-17"	16"	\$625.00
<input type="checkbox"/> QTY ___	PBE1618	16"-17"	18"	\$625.00
<input type="checkbox"/> QTY ___	PBE1620	16"-17"	20"	\$625.00
<input type="checkbox"/> QTY ___	PBE1710	17"-18"	10"	\$625.00
<input type="checkbox"/> QTY ___	PBE1712	17"-18"	12"	\$625.00
<input type="checkbox"/> QTY ___	PBE1714	17"-18"	14"	\$625.00
<input type="checkbox"/> QTY ___	PBE1716	17"-18"	16"	\$625.00
<input type="checkbox"/> QTY ___	PBE1718	17"-18"	18"	\$625.00
<input type="checkbox"/> QTY ___	PBE1720	17"-18"	20"	\$625.00
<input type="checkbox"/> QTY ___	PBE1810	18"-19"	10"	\$625.00
<input type="checkbox"/> QTY ___	PBE1812	18"-19"	12"	\$625.00
<input type="checkbox"/> QTY ___	PBE1814	18"-19"	14"	\$625.00
<input type="checkbox"/> QTY ___	PBE1816	18"-19"	16"	\$625.00
<input type="checkbox"/> QTY ___	PBE1818	18"-19"	18"	\$625.00
<input type="checkbox"/> QTY ___	PBE1820	18"-19"	20"	\$625.00

More sizes are available on the next page

INVACARE® MATRX® BACKS – INTERNATIONAL PRICE LIST/ORDER FORM

MATRX® ELITE BACK: Standard – Positioning Back

- Easy Set Hardware fits ¾", 7/8", 1" diameter back canes
- See Page 12 for other no-charge hardware options
- 3" Contour Depth
- Weight capacity up to 300 lbs
- All Elite Backs 15" and wider x 14" and higher are pre-drilled for headrest attachment, and include headrest adapter plate
- Extra outer cover includes foam

Standard Meshtex Fabric upholstery



Elite Back Continues



QUANTITY	MODEL	FITS WIDTH	FITS HEIGHT	MSRP
<input type="checkbox"/> QTY ____	PBE1910	17"-18"	16"	\$625.00
<input type="checkbox"/> QTY ____	PBE1912	17"-18"	18"	\$625.00
<input type="checkbox"/> QTY ____	PBE1916	17"-18"	20"	\$625.00
<input type="checkbox"/> QTY ____	PBE1918	18"-19"	10"	\$625.00
<input type="checkbox"/> QTY ____	PBE1920	18"-19"	12"	\$625.00
<input type="checkbox"/> QTY ____	PBE2010	18"-19"	12"	\$625.00
<input type="checkbox"/> QTY ____	PBE2012	18"-19"	16"	\$625.00
<input type="checkbox"/> QTY ____	PBE2016	18"-19"	10"	\$625.00
<input type="checkbox"/> QTY ____	PBE2018	18"-19"	12"	\$625.00
<input type="checkbox"/> QTY ____	PBE2020	18"-19"	12"	\$625.00

Elite Extra Back Cover (Refer to back sizes for proper fit)

MODEL	FITS WIDTH	FITS HEIGHT	MSRP
PBEC	REFER TO BACK SIZES ABOVE		\$135.00
SPECIFY:	WIDTH ____	HEIGHT ____	<input type="checkbox"/> QTY ____

MATRX® ELITE FULL BACK: Standard – Positioning Back

- Easy Set Hardware fits ¾", 7/8", 1" diameter back canes
- See Page 12 for other no-charge hardware options
- Back Shell made without scapular cutout and includes SuperSoft foam
- Weight capacity up to 300 lbs
- Pre-drilled for headrest attachment, and includes headrest adapter plate
- Extra outer cover includes foam

Standard Meshtex Fabric upholstery



Elite Full Back



QUANTITY	MODEL	FITS WIDTH	FITS HEIGHT	MSRP
<input type="checkbox"/> QTY ____	PBE1620-FB	16"-17"	20"	\$625.00
<input type="checkbox"/> QTY ____	PBE1720-FB	17"-18"	20"	\$625.00
<input type="checkbox"/> QTY ____	PBE1820-FB	18"-19"	20"	\$625.00
<input type="checkbox"/> QTY ____	PBE1920-FB	19"-20"	20"	\$625.00
<input type="checkbox"/> QTY ____	PBE2020-FB	20"-21"	20"	\$625.00

Elite Full Extra Back Cover (Refer to back sizes above for proper fit)

MODEL	FITS WIDTH	FITS HEIGHT	MSRP
PBEC-FB	REFER TO BACK SIZES ABOVE		\$135.00
SPECIFY:	WIDTH ____	HEIGHT ____	<input type="checkbox"/> QTY ____

MATRX® ELITE TR BACK: Standard – Positioning Back

- Easy Set Hardware fits 3/4", 7/8", 1" diameter back canes
- See Page 12 for other no-charge hardware options
- Back Shell made without scapular cutout and includes SuperSoft foam. Recommended for Tilt/Recline applications
- 3" Contour Depth
- Weight capacity up to 300 lbs
- Pre-drilled for headrest attachment, includes headrest adapter plate
- Extra outer cover includes foam

- Standard Meshtex Fabric upholstery



Elite TR Back



QUANTITY	MODEL	FITS WIDTH	FITS HEIGHT	MSRP
<input type="checkbox"/> QTY ____	PBE1616-TR	16"-17"	16"	\$650.00
<input type="checkbox"/> QTY ____	PBE1618-TR	16"-17"	18"	\$650.00
<input type="checkbox"/> QTY ____	PBE1620-TR	16"-17"	20"	\$650.00
<input type="checkbox"/> QTY ____	PBE1716-TR	17"-18"	16"	\$650.00
<input type="checkbox"/> QTY ____	PBE1718-TR	17"-18"	18"	\$650.00
<input type="checkbox"/> QTY ____	PBE1720-TR	17"-18"	20"	\$650.00
<input type="checkbox"/> QTY ____	PBE1816-TR	18"-19"	16"	\$650.00
<input type="checkbox"/> QTY ____	PBE1818-TR	18"-19"	18"	\$650.00
<input type="checkbox"/> QTY ____	PBE1820-TR	18"-19"	20"	\$650.00
<input type="checkbox"/> QTY ____	PBE1916-TR	19"-20"	16"	\$650.00
<input type="checkbox"/> QTY ____	PBE1918-TR	19"-20"	18"	\$650.00
<input type="checkbox"/> QTY ____	PBE1920-TR	19"-20"	20"	\$650.00
<input type="checkbox"/> QTY ____	PBE2016-TR	20"-21"	16"	\$650.00
<input type="checkbox"/> QTY ____	PBE2018-TR	20"-21"	18"	\$650.00
<input type="checkbox"/> QTY ____	PBE2020-TR	20"-21"	20"	\$650.00

Elite TR Extra Back Cover (Refer to back sizes above for proper fit)

MODEL	FITS WIDTH	FITS HEIGHT	MSRP
PBEC-TR	REFER TO BACK SIZES ABOVE		\$135.00
SPECIFY:	WIDTH ____	HEIGHT ____	<input type="checkbox"/> QTY ____

INVACARE® MATRX® BACKS – INTERNATIONAL PRICE LIST/ORDER FORM

MATRX® ELITE DEEP BACK: Standard – Positioning Back

- Easy Set Hardware fits ¾", 7/8", 1" diameter back canes
- See Page 12 for other no-charge hardware options
- 6" Contour Depth
- Weight capacity up to 300 lbs
- All Elite Backs 15" and wider x 16" and higher are pre-drilled for headrest attachment, includes headrest adapter plate
- Extra outer cover includes foam

- Standard Meshtex Fabric upholstery



Elite Deep Back



QUANTITY	MODEL	FITS WIDTH	FITS HEIGHT	MSRP
<input type="checkbox"/> QTY ___	EDB1412	14"-15"	12"	\$760.00
<input type="checkbox"/> QTY ___	EDB1416	14"-15"	16"	\$760.00
<input type="checkbox"/> QTY ___	EDB1418	14"-15"	18"	\$760.00
<input type="checkbox"/> QTY ___	EDB1420	14"-15"	20"	\$760.00
<input type="checkbox"/> QTY ___	EDB1510	15"-16"	10"	\$760.00
<input type="checkbox"/> QTY ___	EDB1512	15"-16"	12"	\$760.00
<input type="checkbox"/> QTY ___	EDB1516	15"-16"	16"	\$760.00
<input type="checkbox"/> QTY ___	EDB1518	15"-16"	18"	\$760.00
<input type="checkbox"/> QTY ___	EDB1520	15"-16"	20"	\$760.00
<input type="checkbox"/> QTY ___	EDB1610	16"-17"	10"	\$760.00
<input type="checkbox"/> QTY ___	EDB1612	16"-17"	12"	\$760.00
<input type="checkbox"/> QTY ___	EDB1616	16"-17"	16"	\$760.00
<input type="checkbox"/> QTY ___	EDB1618	16"-17"	18"	\$760.00
<input type="checkbox"/> QTY ___	EDB1620	16"-17"	20"	\$760.00
<input type="checkbox"/> QTY ___	EDB1710	17"-18"	10"	\$760.00
<input type="checkbox"/> QTY ___	EDB1712	17"-18"	12"	\$760.00
<input type="checkbox"/> QTY ___	EDB1716	17"-18"	16"	\$760.00
<input type="checkbox"/> QTY ___	EDB1718	17"-18"	18"	\$760.00
<input type="checkbox"/> QTY ___	EDB1720	17"-18"	20"	\$760.00
<input type="checkbox"/> QTY ___	EDB1810	18"-19"	10"	\$760.00
<input type="checkbox"/> QTY ___	EDB1812	18"-19"	12"	\$760.00
<input type="checkbox"/> QTY ___	EDB1816	18"-19"	16"	\$760.00
<input type="checkbox"/> QTY ___	EDB1818	18"-19"	18"	\$760.00
<input type="checkbox"/> QTY ___	EDB1820	18"-19"	20"	\$760.00
<input type="checkbox"/> QTY ___	EDB1910	19"-20"	10"	\$760.00
<input type="checkbox"/> QTY ___	EDB1912	19"-20"	12"	\$760.00
<input type="checkbox"/> QTY ___	EDB1916	19"-20"	16"	\$760.00
<input type="checkbox"/> QTY ___	EDB1918	19"-20"	18"	\$760.00
<input type="checkbox"/> QTY ___	EDB1920	19"-20"	20"	\$760.00
<input type="checkbox"/> QTY ___	EDB2010	20"-21"	10"	\$760.00
<input type="checkbox"/> QTY ___	EDB2012	20"-21"	12"	\$760.00
<input type="checkbox"/> QTY ___	EDB2016	20"-21"	16"	\$760.00
<input type="checkbox"/> QTY ___	EDB2018	20"-21"	18"	\$760.00
<input type="checkbox"/> QTY ___	EDB2020	20"-21"	20"	\$760.00

Elite Deep Extra Back Cover (Refer to back sizes above for proper fit)

MODEL	FITS WIDTH	FITS HEIGHT	MSRP
EDBC	REFER TO BACK SIZES ABOVE		\$135.00
SPECIFY:	WIDTH ___	HEIGHT ___	<input type="checkbox"/> QTY ___

MATRX® ELITE HEAVY-DUTY BACK: Bariatric – Positioning Back

- Easy Set Hardware fits 3/4", 7/8", 1" diameter back canes
- See Page 13 for other no-charge hardware options
- 3" Contour Depth
- Weight capacity up to 500 lbs
- Pre-drilled for headrest attachment, includes headrest adapter plate
- Extra outer cover includes foam

- Standard Meshtex Fabric upholstery



Elite Heavy-Duty Back



QUANTITY	MODEL	FITS WIDTH	FITS HEIGHT	MSRP
<input type="checkbox"/> QTY ____	PBE2116-HD	20" - 21"	16"	\$795.00
<input type="checkbox"/> QTY ____	PBE2120-HD	20" - 21"	20"	\$795.00
<input type="checkbox"/> QTY ____	PBE2216-HD	22" - 23"	16"	\$795.00
<input type="checkbox"/> QTY ____	PBE2220-HD	22" - 23"	20"	\$795.00
<input type="checkbox"/> QTY ____	PBE2316-HD	23" - 24"	16"	\$795.00
<input type="checkbox"/> QTY ____	PBE2320-HD	23" - 24"	20"	\$795.00
<input type="checkbox"/> QTY ____	PBE2416-HD	24" - 25"	16"	\$795.00
<input type="checkbox"/> QTY ____	PBE2420-HD	24" - 25"	20"	\$795.00

Elite Heavy-Duty Extra Outer Back Cover (Refer to back sizes above for proper fit)

MODEL	FITS WIDTH	FITS HEIGHT	MSRP
PBEC-HD	REFER TO BACK SIZES ABOVE		\$135.00
SPECIFY:	WIDTH ____	HEIGHT ____	<input type="checkbox"/> QTY ____

Elite Extra Wide Heavy-Duty Back (Additional one week lead time)*



MODEL	FITS WIDTH	FITS HEIGHT	MSRP
PBEWW16-HD	REFERENCE CHART BELOW		\$1395.00
SPECIFY:	WIDTH ____	<input type="checkbox"/> QTY ____	
PBEWW20-HD	REFERENCE CHART BELOW		\$1395.00
SPECIFY:	WIDTH ____	<input type="checkbox"/> QTY ____	

W I D T H (WW)	
HEIGHT	
16"	X X X X X X X
20"	X X X X X X X

Duty Extra Outer Back Cover (Additional one week lead time)*

MODEL	FITS WIDTH	FITS HEIGHT	MSRP
PBECWW16	REFERENCE CHART ABOVE		ADD\$235.00
SPECIFY:	WIDTH ____	<input type="checkbox"/> QTY ____	
PBECWW20	REFERENCE CHART ABOVE		ADD\$235.00
SPECIFY:	WIDTH ____	<input type="checkbox"/> QTY ____	

* Made to order and none refundable

MATRIX® ELITE LATERALS: Available for Matrix Elite, Elite Deep & Elite Heavy-Duty Backs

Fixed - Available for Matrix Elite, Elite Deep & Elite Heavy-Duty Backs



QUANTITY	MODEL	DESCRIPTION	MSRP
<input type="checkbox"/> QTY ___	PBLS-L	Standard Fixed Hardware-Left or Right (each)	\$110.00
<input type="checkbox"/> QTY ___	PBLS-R	<i>Pad 5" High x 7.5" Long</i>	
<input type="checkbox"/> QTY ___	PBSS-L	Small Fixed Hardware-Left or Right (each)	\$110.00
<input type="checkbox"/> QTY ___	PBSS-R	<i>Pad 3" High x 5.5" Long</i>	
<input type="checkbox"/> QTY ___	PBOL-L	Offset Lateral Pad w/Fixed Hardware (each)	\$110.00
<input type="checkbox"/> QTY ___	PBOL-R	<i>Pad 5" High x 7.5" Long</i>	

Swing-Away Hardware* - Available only for Matrix Elite Backs



QUANTITY	MODEL	DESCRIPTION	MSRP
<input type="checkbox"/> QTY ___	SAHEB	Swing-Away Lateral Hardware Only (each)	\$225.00

Offset, Swing-Away Hardware* Available only for Matrix Elite Backs

Allows 3" variance (higher or lower) pad position relative to attachment slot.



QUANTITY	MODEL	DESCRIPTION	MSRP
<input type="checkbox"/> QTY ___	OSAH-L	Offset, Swing-Away Lateral Hardware Only (each)	\$225.00
<input type="checkbox"/> QTY ___	OSAH-R	<i>Left or Right available</i>	

Pads for Offset, Swing-Away and Swing-Away Hardware



QUANTITY	MODEL	DESCRIPTION	MSRP
<input type="checkbox"/> QTY ___	EPAD-XS	Extra Small Pad: 3.25" High x 4.25" Long (each)	\$110.00
<input type="checkbox"/> QTY ___	EPAD-S	Small Pad: 3.5" High x 5.25" Long (each)	\$110.00
<input type="checkbox"/> QTY ___	EPAD-M	Medium Pad: 3.75" High x 5.75" Long (each)	\$110.00
<input type="checkbox"/> QTY ___	EPAD-L	Large Pad: 4.25" High x 6.5" Long (each)	\$110.00

MATRIX® MINI LATERALS: Available for Matrix Mini Backs

Swing-Away Hardware* - Available only for Matrix Mini Backs



QUANTITY	MODEL	DESCRIPTION	MSRP
<input type="checkbox"/> QTY ___	SAHMINI	Swing-Away Lateral Hardware Only (each)	\$225.00

Pads for Swing-Away Hardware









QUANTITY	MODEL	DESCRIPTION	MSRP
<input type="checkbox"/> QTY ___	EPAD-XS	Extra Small Black Pad: 3.25" High x 4.25" Long (each)	\$110.00
<input type="checkbox"/> QTY ___	EPAD-XSP	Extra Small Pink Pad: 3.25" High x 4.25" Long (each)	\$110.00
<input type="checkbox"/> QTY ___	EPAD-XSB	Extra Small Blue Pad: 3.25" High x 4.25" Long (each)	\$110.00

HARDWARE UPGRADES: Available for Standard Matrix PB and Elite Backs

Hardware Options

No-charge substitutions for Easy Set hardware are available & must be specified at time of ordering back

	NO-CHARGE SUBSTITUTION MODEL #	HARDWARE ONLY MODEL #	MSRP
	Easy Set Hardware (pair) - Fits 3/4", 7/8", 1" back canes (Standard on PB & Elite Backs)	<input type="checkbox"/> PBMH\$225.00
	Easy Set Hardware (pair) - Fits 1 1/8" back canes	<input type="checkbox"/> PBMH-QUN/C
		<input type="checkbox"/> PBMH-Q\$225.00
	Easy Set Hardware (pair) - Fits 1 1/4" back canes	<input type="checkbox"/> PBMX-XUN/C
		<input type="checkbox"/> PBMX-X\$225.00
	Heavy-Duty Easy Set Hardware (pair) - Fits 3/4", 7/8", 1" back canes (Upcharge on Standard PB & Elite Backs)	<input type="checkbox"/> PBMH-QU\$125.00
	Extended Range Mounting Hardware (pair) - Provides additional 2" Depth adjustment (Available for Standard PB & Elite Backs) NOT compatible with SMH/SMH-U	<input type="checkbox"/> ERMHUN/C
		<input type="checkbox"/> ERMH\$285.00
	Split Mounting Hardware (pair) - (Available for Standard PB & Elite Backs) For avoidance of obstacles on Frame	<input type="checkbox"/> SMHUN/C
		<input type="checkbox"/> SMH\$250.00
	Fixed Space Saver Mounting Hardware (pair) Permanent mount recommended for Power Wheelchair application or when portability is not required	
	For Matrix PB & Elite Standard or Heavy-Duty Backs	<input type="checkbox"/> FSSHUN/C
		<input type="checkbox"/> FSSH\$200.00
	For Formula Tilt	<input type="checkbox"/> FSSH-FTUN/C
	<input type="checkbox"/> FSSH-FT\$200.00	
	For Formula Tilt/Recline	<input type="checkbox"/> FSSH-FTRUN/C
	<input type="checkbox"/> FSSH-FTR\$200.00	
	ESR Mounting Hardware Allows PB & Elite Backs to be mounted to Motion Concepts TRx Power Positioning with ESR	<input type="checkbox"/> ESRMH\$260.00

HARDWARE UPGRADES: Available for Bariatric, Matrx Heavy-Duty PB and Elite Backs

Hardware Options

No-charge substitutions for Heavy-Duty Easy Set hardware are available & must be specified at time of ordering back



NO-CHARGE SUBSTITUTION MODEL #	HARDWARE ONLY MODEL #	MSRP
Heavy-Duty Easy Set Hardware (pair) - Fits 3/4", 7/8", 1" back canes (Standard on PB & Elite Heavy-Duty Backs)		
	<input type="checkbox"/> HDMH	\$325.00
Easy Set Hardware (pair) - Fits 1 1/8" back canes		
<input type="checkbox"/> PBMH-QU	N/C	<input type="checkbox"/> PBMH-Q
		\$225.00
Easy Set Hardware (pair) - Fits 1 1/4" back canes		
<input type="checkbox"/> PBMX-XU	N/C	<input type="checkbox"/> PBMX-X
		\$225.00



Fixed Space Saver Mounting Hardware (pair) Permanent mount recommended for Power Wheelchair application or when portability is not required For Matrx PB & Elite Standard or Heavy-Duty Backs		
<input type="checkbox"/> FSSHU	N/C	<input type="checkbox"/> FSSH
		\$200.00
For Formula Tilt		
<input type="checkbox"/> FSSH-FTU	N/C	<input type="checkbox"/> FSSH-FT
		\$200.00
For Formula Tilt/Recline		
<input type="checkbox"/> FSSH-FTRU	N/C	<input type="checkbox"/> FSSH-FTR
		\$200.00



ESR Mounting Hardware Allows Heavy-Duty PB & Elite Backs to be mounted to Motion Concepts TRx Power Positioning with ESR		
	<input type="checkbox"/> ESRMH-HD	\$310.00

HARDWARE UPGRADES: Available for Pediatric, Matrx Mini Backs


Options:



NO-CHARGE SUBSTITUTION MODEL #	HARDWARE ONLY MODEL #	MSRP
Mini Set Mounting Hardware (pair)- Fits 3/4", 7/8", 1" back canes (Standard on Mini Backs)		
	<input type="checkbox"/> MSMH	\$225.00

BACK OPTIONS & MODIFICATIONS

Modifications

	QUANTITY	MODEL	COMPONENT	MSRP
	<input type="checkbox"/> QTY ____	MAC	Miscellaneous Modification <i>Special Request - Approval required</i>	\$100.00
			SPECIFY: _____ _____ _____	
	<input type="checkbox"/> QTY ____	WSC	Wipeable Slipcover <i>Moisture resistant, wipeable cover with a zipper Available for Matrix PB or Elite Back</i>	\$150.00
			SPECIFY: <input type="checkbox"/> PB <input type="checkbox"/> PB Deep <input type="checkbox"/> PB Heavy-Duty <input type="checkbox"/> PB Deep Heavy-Duty <input type="checkbox"/> Elite <input type="checkbox"/> Elite Full <input type="checkbox"/> Elite Tilt/Recline <input type="checkbox"/> Elite Deep <input type="checkbox"/> Elite Heavy-Duty	
			SPECIFY: WIDTH ____ HEIGHT ____ <input type="checkbox"/> QTY ____	
	<input type="checkbox"/> QTY ____	WSCU	Wipeable Slipcover - ordered with back <i>Moisture resistant, wipeable cover with a zipper Available for Matrix PB or Elite Back</i>	\$100.00
	<input type="checkbox"/> QTY ____	EBSS	Elite Back/Cover - SuperSoft Foam <i>SuperSoft Foam - 1" soft over 1" firm HR foam</i>	\$235.00
			SPECIFY: <input type="checkbox"/> Elite <input type="checkbox"/> Elite Deep <input type="checkbox"/> Elite Heavy-Duty	
			SPECIFY: WIDTH ____ HEIGHT ____ <input type="checkbox"/> QTY ____	
	<input type="checkbox"/> QTY ____	PBVF	PB/Elite Back/Cover Modification <i>Visco Foam over HR Foam</i>	\$235.00
			SPECIFY: <input type="checkbox"/> PB <input type="checkbox"/> PB Deep <input type="checkbox"/> PB Heavy-Duty <input type="checkbox"/> PB Deep Heavy-Duty <input type="checkbox"/> Elite <input type="checkbox"/> Elite Full <input type="checkbox"/> Elite Tilt/Recline <input type="checkbox"/> Elite Deep <input type="checkbox"/> Elite Heavy-Duty	
			SPECIFY: WIDTH ____ HEIGHT ____ <input type="checkbox"/> QTY ____	
	<input type="checkbox"/> QTY ____	EMP	Extended Mounting Pins <i>Allows Matrix PR or Elite Backs to be mounted on wheelchairs that are 2" wider than the back model width --only for backs 20" wide and 300 lbs weight capacity. Not compatible with FSSH.</i>	\$60.00
	<input type="checkbox"/> QTY ____	CHIEB	Chest Harness Interface Kit - for Elite Back <i>Chest Harness is not included - mounting kit only</i>	\$95.00
	<input type="checkbox"/> QTY ____	CHIPB	Chest Harness Interface Kit - for PB Back <i>Chest Harness is not included - mounting kit only</i>	\$45.00
	<input type="checkbox"/> QTY ____	EHRA-S	Small Headrest Adapter - for Elite Back <i>For Attachment of headrests to Elite 12" wide backs</i>	\$35.00
	<input type="checkbox"/> QTY ____	EHRA-M	Medium Headrest Adapter - for Elite Back <i>For Attachment of headrests to Elite Backs 13" - 17" wide backs</i>	\$35.00
	<input type="checkbox"/> QTY ____	EHRA-L	Large Headrest Adapter - for Elite Back <i>For Attachment of headrests to Elite Heavy-Duty 18" - 22" wide backs</i>	\$35.00

BACK OPTIONS & MODIFICATIONS

Modifications

QUANTITY	MODEL	COMPONENT	MSRP
<input type="checkbox"/> QTY _____	STARTEX-F	Startex - Fabric Side up <i>Available for PB and Elite Backs</i>	\$50.00
SPECIFY:	<input type="checkbox"/> PB <input type="checkbox"/> PB Deep <input type="checkbox"/> PB Heavy-Duty <input type="checkbox"/> PB Deep Heavy-Duty <input type="checkbox"/> Elite <input type="checkbox"/> Elite Full <input type="checkbox"/> Elite Tilt/Recline <input type="checkbox"/> Elite Deep <input type="checkbox"/> Elite Heavy-Duty		
	SPECIFY: WIDTH _____	HEIGHT _____	<input type="checkbox"/> QTY _____
<input type="checkbox"/> QTY _____	STARTEX-S	Startex - Smooth Side up <i>Available for PB and Elite Backs</i>	\$50.00
SPECIFY:	<input type="checkbox"/> PB <input type="checkbox"/> PB Deep <input type="checkbox"/> PB Heavy-Duty <input type="checkbox"/> PB Deep Heavy-Duty <input type="checkbox"/> Elite <input type="checkbox"/> Elite Full <input type="checkbox"/> Elite Tilt/Recline <input type="checkbox"/> Elite Deep <input type="checkbox"/> Elite Heavy-Duty		
	SPECIFY: WIDTH _____	HEIGHT _____	<input type="checkbox"/> QTY _____
<input type="checkbox"/> QTY _____	NSF	Non-Standard Fabric Selection <i>Available for all backs; excluding Elite & Mini models</i>	\$75.00
	SPECIFY: _____		

	SPECIFY: WIDTH _____	HEIGHT _____	<input type="checkbox"/> QTY _____

INVACARE® MATRX® BACKS – INTERNATIONAL PRICE LIST/ORDER FORM

MATRX® MX1 BACK: Standard – Positioning Back

NOTE: Light weight, carbon fiber back with Spacer mesh/Vinyl cover and foam. Available with 4 mounting hardware options: Mono Mount (RM), Fixed Dual Mount (CF), Quick Release Dual Mount (CR) & SuperLite (SL)

- 3" Contour Depth
- Weight capacity up to 250 lbs
- Extra outer cover includes foam

- Space mesh fabric with durable reinforced vinyl panels



Spacer Mesh/Vinyl
Outer Cover

MX1 Back with Mono Mount Hardware

Includes: adjustable, non-removable rigidizer mount hardware for ¾" cane diameter.



QUANTITY	MODEL	FITS WIDTH	FITS HEIGHT	MSRP
<input type="checkbox"/> QTY ___	IMX11409-RM	14"	09"	\$775.00
<input type="checkbox"/> QTY ___	IMX11412-RM	14"	12"	\$805.00
<input type="checkbox"/> QTY ___	IMX11416-RM	14"	16"	\$835.00
<input type="checkbox"/> QTY ___	IMX11509-RM	15"	09"	\$775.00
<input type="checkbox"/> QTY ___	IMX11512-RM	15"	12"	\$805.00
<input type="checkbox"/> QTY ___	IMX11516-RM	15"	16"	\$835.00
<input type="checkbox"/> QTY ___	IMX11609-RM	16"	09"	\$775.00
<input type="checkbox"/> QTY ___	IMX11612-RM	16"	12"	\$805.00
<input type="checkbox"/> QTY ___	IMX11616-RM	16"	16"	\$835.00
<input type="checkbox"/> QTY ___	IMX11709-RM	17"	09"	\$775.00
<input type="checkbox"/> QTY ___	IMX11712-RM	17"	12"	\$805.00
<input type="checkbox"/> QTY ___	IMX11716-RM	17"	16"	\$835.00
<input type="checkbox"/> QTY ___	IMX11809-RM	18"	09"	\$775.00
<input type="checkbox"/> QTY ___	IMX11812-RM	18"	12"	\$805.00
<input type="checkbox"/> QTY ___	IMX11816-RM	18"	16"	\$835.00

MX1 Back with Fixed Dual Mount Hardware

Includes: adjustable, non-removable hardware for 1" cane diameter.



QUANTITY	MODEL	FITS WIDTH	FITS HEIGHT	MSRP
<input type="checkbox"/> QTY ___	IMX11409-CF	14"	09"	\$825.00
<input type="checkbox"/> QTY ___	IMX11412-CF	14"	12"	\$855.00
<input type="checkbox"/> QTY ___	IMX11416-CF	14"	16"	\$885.00
<input type="checkbox"/> QTY ___	IMX11509-CF	15"	09"	\$825.00
<input type="checkbox"/> QTY ___	IMX11512-CF	15"	12"	\$855.00
<input type="checkbox"/> QTY ___	IMX11516-CF	15"	16"	\$885.00
<input type="checkbox"/> QTY ___	IMX11609-CF	16"	09"	\$825.00
<input type="checkbox"/> QTY ___	IMX11612-CF	16"	12"	\$855.00
<input type="checkbox"/> QTY ___	IMX11616-CF	16"	16"	\$885.00
<input type="checkbox"/> QTY ___	IMX11709-CF	17"	09"	\$825.00
<input type="checkbox"/> QTY ___	IMX11712-CF	17"	12"	\$855.00
<input type="checkbox"/> QTY ___	IMX11716-CF	17"	16"	\$885.00
<input type="checkbox"/> QTY ___	IMX11809-CF	18"	09"	\$825.00
<input type="checkbox"/> QTY ___	IMX11812-CF	18"	12"	\$855.00
<input type="checkbox"/> QTY ___	IMX11816-CF	18"	16"	\$885.00

MATRX® MX1 BACK: Standard – Positioning Back

MX1 Back with Quick Release Dual Mount Hardware

Includes: adjustable, removable quick release hardware for 3/4", 7/8" & 1" cane diameter.



QUANTITY	MODEL	FITS WIDTH	FITS HEIGHT	MSRP
<input type="checkbox"/> QTY ____	IMX11409-CR	14"	09"	\$850.00
<input type="checkbox"/> QTY ____	IMX11412-CR	14"	12"	\$880.00
<input type="checkbox"/> QTY ____	IMX11416-CR	14"	16"	\$910.00
<input type="checkbox"/> QTY ____	IMX11509-CR	15"	09"	\$850.00
<input type="checkbox"/> QTY ____	IMX11512-CR	15"	12"	\$880.00
<input type="checkbox"/> QTY ____	IMX11516-CR	15"	16"	\$910.00
<input type="checkbox"/> QTY ____	IMX11609-CR	16"	09"	\$850.00
<input type="checkbox"/> QTY ____	IMX11612-CR	16"	12"	\$880.00
<input type="checkbox"/> QTY ____	IMX11616-CR	16"	16"	\$910.00
<input type="checkbox"/> QTY ____	IMX11709-CR	17"	09"	\$850.00
<input type="checkbox"/> QTY ____	IMX11712-CR	17"	12"	\$880.00
<input type="checkbox"/> QTY ____	IMX11716-CR	17"	16"	\$910.00
<input type="checkbox"/> QTY ____	IMX11809-CR	18"	09"	\$850.00
<input type="checkbox"/> QTY ____	IMX11812-CR	18"	12"	\$880.00
<input type="checkbox"/> QTY ____	IMX11816-CR	18"	16"	\$910.00

MX1 Back with SuperLite Mount Hardware

Includes: non-adjustable, non-removable hardware for 1" cane diameter.







QUANTITY	MODEL	FITS WIDTH	FITS HEIGHT	MSRP
<input type="checkbox"/> QTY ____	IMX11409-SL	14"	09"	\$750.00
<input type="checkbox"/> QTY ____	IMX11412-SL	14"	12"	\$780.00
<input type="checkbox"/> QTY ____	IMX11416-SL	14"	16"	\$810.00
<input type="checkbox"/> QTY ____	IMX11509-SL	15"	09"	\$750.00
<input type="checkbox"/> QTY ____	IMX11512-SL	15"	12"	\$780.00
<input type="checkbox"/> QTY ____	IMX11516-SL	15"	16"	\$810.00
<input type="checkbox"/> QTY ____	IMX11609-SL	16"	09"	\$750.00
<input type="checkbox"/> QTY ____	IMX11612-SL	16"	12"	\$780.00
<input type="checkbox"/> QTY ____	IMX11616-SL	16"	16"	\$810.00
<input type="checkbox"/> QTY ____	IMX11709-SL	17"	09"	\$750.00
<input type="checkbox"/> QTY ____	IMX11712-SL	17"	12"	\$780.00
<input type="checkbox"/> QTY ____	IMX11716-SL	17"	16"	\$810.00
<input type="checkbox"/> QTY ____	IMX11809-SL	18"	09"	\$750.00
<input type="checkbox"/> QTY ____	IMX11812-SL	18"	12"	\$780.00
<input type="checkbox"/> QTY ____	IMX11816-SL	18"	16"	\$810.00

MX1 Extra Back Cover (Refer to back sizes above for proper fit)

MODEL	FITS WIDTH	FITS HEIGHT	MSRP
MX1C	REFER TO BACK SIZES ABOVE		\$170.00
SPECIFY:	WIDTH ____	HEIGHT ____	<input type="checkbox"/> QTY ____

MATRIX MX1 HARDWARE UPGRADES:

MX1 Hardware Options

	QUANTITY	MODEL	COMPONENT	MSRP
	<input type="checkbox"/> QTY ____	RGDMT-9	Mono Mount Mounting Hardware for 9"H Back <i>Includes hardware for ¾" cane diameter</i>	\$195.00
	<input type="checkbox"/> QTY ____	RGDMT-12	Mono Mount Mounting Hardware for 12"H Back <i>Includes hardware for ¾" cane diameter</i>	\$195.00
	<input type="checkbox"/> QTY ____	RGDMT-16	Mono Mount Mounting Hardware for 16"H Back <i>Includes hardware for ¾" cane diameter</i>	\$195.00
	<input type="checkbox"/> QTY ____	FBCMT	Dual Mount Mounting Hardware (pair) <i>Includes hardware for 1" cane diameter</i>	\$245.00
	<input type="checkbox"/> QTY ____	RBCMT	Quick Release Dual Mount Mounting Hardware (pair) <i>Includes hardware for 7/8", ¾" and 1" cane diameter</i>	\$270.00
	<input type="checkbox"/> QTY ____	BCO-78	Band Clamp Only (each) <i>Clamp for 7/8" cane diameter</i>	\$8.00
	<input type="checkbox"/> QTY ____	BCO-34	Band Clamp Only (each) <i>Clamp for ¾" cane diameter</i>	\$8.00
	<input type="checkbox"/> QTY ____	BCO-1	Band Clamp Only (each) <i>Clamp for 1" cane diameter</i>	\$8.00
	<input type="checkbox"/> QTY ____	SLMT-01	SuperLite Mounting Hardware (pair) <i>Includes hardware for 1" cane diameter</i>	\$150.00
	<input type="checkbox"/> QTY ____	SLMT-LW	SuperLite Mounting Hardware (pair) <i>For Lasher Wheelchair Includes hardware for 1¼" cane diameter</i>	\$150.00

MATRIX® GENERA BACK: Standard – Positioning Back

NOTE: Fully adjustable, quick release mounting hardware to fit ¾", 7/8" and 1" diameter back canes.


- 3" Contour Depth
- Weight capacity up to 250 lbs

- Breathable cover with vinyl backing



Outer Cover

Genera Back

	QUANTITY	MODEL	FITS WIDTH	FITS HEIGHT	MSRP
	<input type="checkbox"/> QTY ____	GB1616	16"	16"	\$395.00
	<input type="checkbox"/> QTY ____	GB1620	16"	20"	\$395.00
	<input type="checkbox"/> QTY ____	GB1816	18"	16"	\$395.00
	<input type="checkbox"/> QTY ____	GB1820	18"	20"	\$395.00
	<input type="checkbox"/> QTY ____	GB2016	20"	16"	\$395.00
	<input type="checkbox"/> QTY ____	GB2020	20"	20"	\$395.00





Genera Extra Back Cover (Refer to back sizes above for proper fit)

	MODEL	FITS WIDTH	FITS HEIGHT	MSRP
	GBC	REFER TO BACK SIZES ABOVE		\$95.00
SPECIFY:	WIDTH ____	DEPTH ____	<input type="checkbox"/> QTY ____	

MATRIX POSITIONING ACCESSORIES





Headrest Supports

Includes: Standard Hardware, ¼" shim and knob style mounting plate; exception Onyx Pad

	QUANTITY	MODEL	COMPONENT	MSRP
	<input type="checkbox"/> QTY ____	MSSH ^P	Matrix Super Soft Headrest Pad - 10" <i>Comfortable HR foam, Startex Cover</i>	\$325.00
	<input type="checkbox"/> QTY ____	OPSH ¹	Onyx Pad <i>Auto-style headrest with Onyx Leatherette cover 8" x 8" x 4"</i>	\$350.00
	<input type="checkbox"/> QTY ____	SPSH	Standard Pad <i>Self skinned polyurethane foam 6.5" x 7"</i>	\$295.00
	<input type="checkbox"/> QTY ____	STSH-N	Stealth NINO Pad - 6" <i>Self skinned polyurethane foam 6" x 3"</i>	\$350.00
	<input type="checkbox"/> QTY ____	STSH-XS	Stealth XS Pad - 8" <i>Self skinned polyurethane foam 8" x 6"</i>	\$350.00
	<input type="checkbox"/> QTY ____	STSH-S	Stealth S Pad - 10" <i>Self skinned polyurethane foam 10" x 6"</i>	\$350.00
	<input type="checkbox"/> QTY ____	STSH-L	Stealth L Pad - 14" <i>Self skinned polyurethane foam 14" x 6"</i>	\$350.00

1 Comes standard with Fixed mounting hardware & fixed mounting clamp

Headrest Hardware

	QUANTITY	MODEL	COMPONENT	MSRP
	<input type="checkbox"/> QTY ____	MAHU	Multi-Axis Hardware <i>Depth & Height Adjustable through horizontal slide bar. Dual pivot ball for pad orientation</i>	\$210.00
	<input type="checkbox"/> QTY ____	WBHU	Stealth World's Best Hardware <i>Depth & Height Adjustable through dual lonk system. Dual pivot ball for pad orientation</i>	\$210.00
	<input type="checkbox"/> QTY ____	FXMT	Fixed Mounting Clip <i>Fixed, Non-removable, Angle adjustable headrest mounting clump</i>	\$120.00
	<input type="checkbox"/> QTY ____	SHIM	¼" Shim Plate <i>Insert for headrest mounting</i>	\$8.00

MATRX ACCESSORIES

Accessories



QUANTITY	MODEL	COMPONENT	MSRP
<input type="checkbox"/> QTY ____	PFLAP-XS	X-Small Privacy Flap <i>10" W x 16" H - fits 12"-14" wide backs</i>	\$40.00
<input type="checkbox"/> QTY ____	PFLAP-S	Small Privacy Flap <i>13" W x 16" H - fits 15"-16" wide backs</i>	\$40.00
<input type="checkbox"/> QTY ____	PFLAP-M	Medium Privacy Flap <i>15" W- 16" H - fits 17"-18" wide backs</i>	\$40.00
<input type="checkbox"/> QTY ____	PFLAP-L	Large Privacy Flap <i>17"W - 16" H - fits 19"-20" wide backs</i>	\$40.00
<input type="checkbox"/> QTY ____	PFLAP-XL	X-Large Privacy Flap <i>19" W- 16" H - fits 21"-22" wide backs</i>	\$40.00
<input type="checkbox"/> QTY ____	PFLAP-XXL	XX-Large Privacy Flap <i>21" W- 16" H - fits 23"-24" wide backs</i>	\$40.00
<input type="checkbox"/> QTY ____	PBCR	Posture Back Cane Replacement (pair) <i>1" Diameter back cane replacement kit</i>	\$115.00
<input type="checkbox"/> QTY ____	PBCI ³	Posture Back Cane Clamp Inserts (pair) <i>Insert for use with 3/4" back banes</i>	N/C

