





INVACARE® MATRX®

ELITE

(10" HIGH MODELS)
BACK SPECIFICATIONS



BACK SHELL DIMENSIONS

OVERALL DIMENSIONS (CONTACT AREA)

MODEL	BACK SHELL DIMENSIONS				OVERALL DIMENSIONS (CONTACT AREA)													
	A WIDTH		B HEIGHT		C LATERAL CONTOUR HEIGHT		D HEIGHT		E INNER MAXIMUM WIDTH		F OUTER MAXIMUM WIDTH		G HEIGHT		WEIGHT WITHOUT HARDWARE		WEIGHT CAPACITY	
	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	in	cm	lb	kg	lb	kg
PBE1210	9¾	25	10¾	27.5	8	20	2	5	8½	21.5	11½	28	12	31	1.8	0.8	300	136
PBE1310	10¼	26.5	10¾	27.5	8	20	2	5	9	23	12	31	12	31	2.0	0.9	300	136
PBE1410	12½	32	10¼	26.5	6	15.5	3.5	9	11	28	14	36	12	31	2.2	1.0	300	136
PBE1510	13½	34.5	10¼	26.5	6	15.5	3.5	9	12	31	15	38	12	31	2.4	1.1	300	136
PBE1610	14½	37	10¼	26.5	6	15.5	3.5	9	13	33	16	41	12	31	2.5	1.1	300	136
PBE1710	15½	39.5	10¼	26.5	6	15.5	3.5	9	14	36	17	43	12	31	2.7	1.2	300	136
PBE1810	16½	42	10¼	26.5	6	15.5	3.5	9	15	38	18	46	12	31	2.9	1.3	300	136
PBE1910	17½	44.5	10¼	26.5	6	15.5	3.5	9	16	41	19	48	12	31	3.0	1.4	300	136
PBE2010	18½	47	10¼	26.5	6	15.5	3.5	9	17	43	20	51	12	31	3.2	1.5	300	136

Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd.
E. Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Motion Concepts
www.motionconcepts.com

USA
700 Ensminger Rd.,
Suite 112
Tonawanda, NY 14150
(888) 433-6818
Fax: (888) 433-6834

Canada
84 Citation Dr., Unit 1-7
Concord, ON L4K 3C1
(866) 748-7943
(905) 695-0134
Fax: 905-695-0138

© 2014 Invacare Corporation. All rights reserved.
Trademarks are identified by the symbols™ and®.
All trademarks are owned by or licensed to Invacare
Corporation unless otherwise noted. Specifications
are subject to change without notification.
Form 12-572 Rev. 06/14 140442