





INVACARE® MATRX® ELITE FULL BACK SPECIFICATIONS



| | BACK SHELL DIMENSIONS | | | | OVERALL DIMENSIONS (CONTACT AREA) | | | | | | | | | |
|-----------|-----------------------|------|-------------|------|-----------------------------------|----|--------------------------------|----|-------------|----|-------------------------------|-----|--------------------|-----|
| | A WIDTH | | B HEIGHT | | C INNER MAXIMUM WIDTH | | D OUTER MAXIMUM WIDTH | | E HEIGHT | | WEIGHT WITHOUT HARDWARE | | WEIGHT CAPACITY | |
| MODEL | in | cm | in | cm | in | cm | in | cm | in | cm | lb | kg | lb | kg |
| PBEFB1620 | 14½ | 37 | 19½ | 49.5 | 15 | 38 | 15 | 38 | 21 | 53 | 5.4 | 2.4 | 300 | 136 |
| PBEFB1720 | 15½ | 39.5 | 19½ | 49.5 | 16 | 41 | 16 | 41 | 21 | 53 | 5.6 | 2.5 | 300 | 136 |
| PBEFB1820 | 16½ | 42 | 19½ | 49.5 | 17 | 43 | 17 | 43 | 21 | 53 | 5.8 | 2.6 | 300 | 136 |
| PBEFB1920 | 17½ | 44.5 | 19½ | 49.5 | 18 | 46 | 18 | 46 | 21 | 53 | 6.0 | 2.7 | 300 | 136 |
| PBEFB2020 | 18½ | 47 | 19½ | 49.5 | 19 | 48 | 19 | 48 | 21 | 53 | 6.2 | 2.8 | 300 | 136 |

Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd.
E. Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Motion Concepts
www.motionconcepts.com

USA
700 Ensminger Rd.,
Suite 112
Tonawanda, NY 14150
(888) 433-6818
Fax: (888) 433-6834

Canada
84 Citation Dr., Unit 1-7
Concord, ON L4K 3C1
(866) 748-7943
(905) 695-0134
Fax: 905-695-0138

© 2014 Invacare Corporation. All rights reserved.
Trademarks are identified by the symbols™ and®.
All trademarks are owned by or licensed to Invacare
Corporation unless otherwise noted. Specifications
are subject to change without notification.
Form 12-573 Rev. 06/14 140443