





# INVACARE® MATRX® ELITE TR BACK SPECIFICATIONS



		BACK SHELL DIMENSIONS				OVERALL DIMENSIONS (CONTACT AREA)									
NOMINAL WIDTH	NOMINAL HEIGHT	A		B		C		D		E		CONTOUR DEPTH (loaded)		WEIGHT W/O HARDWARE	
		inches	cm	inches	cm	inches	cm	inches	cm	inches	cm	inches	cm	lb	kg
15	16	13.5	34.3	15.5	39.4	15	38.1	11	27.9	17	43.2	3	7.6	3.3	1.5
15	18	13.5	34.3	17.5	44.5	15	38.1	11	27.9	19	48.3	3	7.6	3.6	1.6
15	20	13.5	34.3	19.5	49.5	15	38.1	11	27.9	21	53.3	3	7.6	3.9	1.8
16	16	14.5	36.8	15.5	39.4	16	40.6	12	30.5	17	43.2	3	7.6	3.5	1.6
16	18	14.5	36.8	17.5	44.5	16	40.6	12	30.5	19	48.3	3	7.6	3.8	1.7
16	20	14.5	36.8	19.5	49.5	16	40.6	12	30.5	21	53.3	3	7.6	4.1	1.9
17	16	15.5	39.4	15.5	39.4	17	43.2	13	33.0	17	43.2	3	7.6	3.7	1.7
17	18	15.5	39.4	17.5	44.5	17	43.2	13	33.0	19	48.3	3	7.6	4.0	1.8
17	20	15.5	39.4	19.5	49.5	17	43.2	13	33.0	21	53.3	3	7.6	4.4	2.0
18	16	16.5	41.9	15.5	39.4	18	45.7	14	35.6	17	43.2	3	7.6	3.8	1.7
18	18	16.5	41.9	17.5	44.5	18	45.7	14	35.6	19	48.7	3	7.6	4.2	1.9
18	20	16.5	41.9	19.5	49.5	18	45.7	14	35.6	21	53.3	3	7.6	4.5	2.0
19	16	17.5	44.5	15.5	39.4	19	48.3	15	38.1	17	43.2	3	7.6	4.0	1.8
19	18	17.5	44.5	17.5	44.5	19	48.3	15	38.1	19	48.3	3	7.6	4.4	2.0
19	20	17.5	44.5	19.5	49.5	19	48.3	15	38.1	21	53.3	3	7.6	4.8	2.2
20	16	18.5	47.0	15.5	39.4	20	50.8	16	40.6	17	43.2	3	7.6	4.2	1.9
20	18	18.5	47.0	17.5	44.5	20	50.8	16	40.6	19	48.3	3	7.6	4.7	2.1
20	20	18.5	47.0	19.5	49.5	20	50.8	16	40.6	21	53.3	3	7.6	5.1	2.3

**WEIGHT CAPACITY:** 300 lb. / 136 kg.



**Invacare Corporation** www.invacare.com **USA** One Invacare Way, Elyria, Ohio, 44035-4190. (800) 333-6900.  
**Canada** 570 Matheson Blvd. E., Unit 8, Mississauga, Ontario, L4Z 4G4 Canada. (800) 668-5324. **Motion Concepts**  
 www.motionconcepts.com. **USA** 700 Ensinger Rd., Suite 112, Tonawanda, NY 14150. (888) 433-6818 Fax: (888) 433-  
 6834. **Canada** 84 Citation Dr., Unit 1-7, Concord, ON L4K 3C1. (866) 748-7943. (905) 695-0134 Fax: 905-695-0138. © 2016  
 Invacare Corporation. All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by  
 or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.  
 Form 16-091 160510