

# MOTION CONCEPTS

## Airmesh Calf Strap



### Package Contents:

- (1) 36" (91 cm) Airmesh Calf Strap
- (2) Velcro™ Anchoring Strips

### Mounting Location Options:

The Sports Calf Panel will fit around 16" (41 cm) - 18" (46 cm) widths

- ◆ Behind the calf
- ◆ In front of the lower leg
- ◆ Above the knee



### Installation & Fitting Instructions:

1. Determine the desired mounting location.
2. Wrap the Velcro™ Strips around each footrest hangar bracket tube where the Airmesh Calf Strap is to be located.
3. Wrap the flexible part of the Airmesh Calf Strap over each Velcro™ Strip (this prevents sliding), with the padded side facing the body.

**Cleaning:** Machine wash warm (90° F, 32° C) with mild detergent. Do not use bleach. Line dry.

**Warranty:** All Motion Concepts products carry an unconditional guarantee against defects in workmanship and materials. Items found to be defective will be repaired or replaced at the discretion of Motion Concepts.

84 Citation Drive, Units 1-7, Concord, Ontario L4K 3C1  
Toll Free Tel: 1.866.748.7943 Fax: 905.695.0138  
700 Ensminger Rd. Suite 112, Tonawanda, New York 14150  
Toll Free Tel: 1.888.433.6818 Toll Free Fax: 1.888.433.6834



# MOTION CONCEPTS

## Airmesh Calf Strap



### Package Contents:

- (1) 36" (91 cm) Airmesh Calf Strap
- (2) Velcro™ Anchoring Strips

### Mounting Location Options:

The Sports Calf Panel will fit around 16" (41 cm) - 18" (46 cm) widths

- ◆ Behind the calf
- ◆ In front of the lower leg
- ◆ Above the knee



### Installation & Fitting Instructions:

1. Determine the desired mounting location.
2. Wrap the Velcro™ Strips around each footrest hangar bracket tube where the Airmesh Calf Strap is to be located.
3. Wrap the flexible part of the Airmesh Calf Strap over each Velcro™ Strip (this prevents sliding), with the padded side facing the body.

**Cleaning:** Machine wash warm (90° F, 32° C) with mild detergent. Do not use bleach. Line dry.

**Warranty:** All Motion Concepts products carry an unconditional guarantee against defects in workmanship and materials. Items found to be defective will be repaired or replaced at the discretion of Motion Concepts.

84 Citation Drive, Units 1-7, Concord, Ontario L4K 3C1  
Toll Free Tel: 1.866.748.7943 Fax: 905.695.0  
700 Ensminger Rd. Suite 112, Tonawanda, New York 14150  
Toll Free Tel: 1.888.433.6818 Toll Free Fax: 1.888.433.6834

