

INVACARE® MATRX®
ELITE™ TR
BACK



FREEDOM DESIGNS® P.R.O. CG™ Tilt-In-Space Wheelchair with Elite™ TR Back*

The **Invacare Matrx Elite TR Back** is specially designed for use in manual or power tilt-only or tilt/recline applications. It includes an extra layer of SuperSoft HR foam for improved immersion and comfort, and incorporates full support across the shoulder area. Bariatric Elite TR Back available with 500 lb. capacity.

Call to order

1.800.333.6900
invacare.com

1.888.433.6818
motionconcepts.com



matrx®
SEATING SERIES

Features and Functional Benefits

Cover – Oversized cover featuring SuperSoft HR foam layered over firm HR foam for comfort and protection. Attractive Meshtex cover is water-resistant and breathable.

Hardware – “Easy Set” hardware installs and adjusts in minutes – all adjustments can be made with user seated in the wheelchair. Other hardware styles may be substituted for unique mounting applications; Split, Fixed Space Saver, Extended Range, as well as many versions specific to power wheelchair applications including ESR Mount.

Back Shell – Stylish, durable aluminum shell includes contoured support in shoulder area, slotted for lateral support hardware attachment.

Adjustable Lumbo-Sacral Support – Included.

Headrest Adapter – Required for easy, durable and safe mounting of headrest hardware and to provide offset for clearance of foam overlap.

Positioning Options

Lateral Supports – Full range of fixed, offset and swing-away laterals available for order separately.

Fixed – 2.5" of vertical adjustment on back pan. Mounts inside back shell without interfering with depth/angle adjustment. Available in standard 5" H x 7.5" L and pediatric 3" H x 5.5" L sizes.

Offset Fixed – Provides an additional 3" of vertical position (higher or lower) relative to mounting slots. Available 5" H x 7.5" L.

Swing-Away – 2.5" of vertical adjustment on back pan. Mounts inside back shell without interfering with depth/angle adjustment. Easy push-button release, four pad sizes (see below).

Offset Swing-Away – Provides an additional 3" of vertical position (higher or lower). Mounts inside back shell without interfering with depth/angle adjustment. Easy push-button release, four pad sizes (see below).

Lateral Pads for Swing-Away Laterals – XS: 3.25" H x 4.25" L, S: 3.5" H x 5.25" L, M: 3.75" H x 5.75" L, L: 4.25" H x 6.5" L

HCPCS Code: E2620



1. Fixed Lateral Pads
2. Offset Fixed Lateral Pads
3. Swing-Away Laterals
4. Offset, Swing-Away Laterals
5. 3" Contour Depth
6. Privacy Flap



Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd. E.
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Motion Concepts
www.motionconcepts.com
info@motionconcepts.com

USA
700 Ensinger Rd., Suite 112
Tonawanda, NY 14150
(888) 433-6818
Fax: (888) 433-6834

Canada
84 Citation Dr., Unit 1-7
Concord, ON L4K 3C1
(866) 748-7943
(905) 695-0134
Fax: 905-695-0138

Model #	Height	Width	Depth** Adjustment	Angle Adjustment	Product Weight	Weight Capacity	Limited Warranty
PBE TR	16", 18", 20"	15" - 20"	4"	40° range +/- 20°	3.5 lb. and up	300 lb.	2 Year
PBE TRHD	20"	21" - 24"	4"	40° range +/- 20°	8.5 lb. and up	500 lb.	2 Year

*FREEDOM DESIGNS® P.R.O. CG™ Tilt-In-Space Wheelchair weight capacity limit is 250 lbs.

**Limitation of depth/angle adjustment if used on chairs narrower than model width

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program and should consult an attorney or other advisor to discuss specific situations in further detail.



© 2016
Invacare Corporation.
All rights reserved.
Trademarks are identified
by the symbols ™ and
©. All trademarks are
owned by or licensed to
Invacare Corporation unless
otherwise noted. Matrix
is a registered trademark
of Motion Concepts.
Specifications are subject
to change without
notification.
Form No. 15-035
Rev. 09/16
160958