


ELITE E2

BACK SPECIFICATIONS



E2S BACK SIZE	A BACK SHELL LENGTH	B BACK CANE WIDTH	C BACK DEPTH	D BACK SHELL WIDTH	E INSIDE FOAM WIDTH	F BACK FOAM LENGTH	G APEX OF LATERAL	H FOAM THICKNESS	I BACK SHELL DEPTH
1406	6"	13"-15"	3.5"	12.3"	10.5"	7.5"	3"	1.25"	3.5"
1408	8"	13"-15"	3.5"	12.3"	10.5"	9.5"	4"	1.25"	3.5"
1410	10"	13"-15"	3.5"	12.3"	10.5"	11.5"	6"	1.25"	3.5"
1412	12"	13"-15"	3.5"	12.3"	10.5"	13.5"	7"	1.25"	3.5"
1414	14"	13"-15"	3.5"	12.3"	10.5"	15.5"	7.7"	1.25"	3.5"
1416	16"	13"-15"	3.5"	12.3"	10.5"	17.5"	8.25"	1.25"	3.5"
1418	18"	13"-16"	3.5"	13.3"	11.5"	19.5"	8.25"	1.25"	3.5"
1420	20"	13"-16"	3.5"	13.3"	11.5"	21.5"	8.25"	1.25"	3.5"
1506	6"	14"-16"	3.5"	13.3"	11.5"	7.5"	3"	1.25"	3.5"
1508	8"	14"-16"	3.5"	13.3"	11.5"	9.5"	4"	1.25"	3.5"
1510	10"	14"-16"	3.5"	13.3"	11.5"	11.5"	6"	1.25"	3.5"
1512	12"	14"-16"	3.5"	13.3"	11.5"	13.5"	7"	1.25"	3.5"
1514	14"	14"-16"	3.5"	13.3"	11.5"	15.5"	7.5"	1.25"	3.5"
1516	16"	14"-16"	3.5"	13.3"	11.5"	17.5"	8.25"	1.25"	3.5"
1518	18"	14"-17"	3.5"	13.3"	11.5"	19.5"	8.25"	1.25"	3.5"
1520	20"	14"-17"	3.5"	13.3"	11.5"	21.5"	8.25"	1.25"	3.5"
1606	6"	15"-17"	3.5"	14.3"	12.5"	7.5"	3"	1.25"	3.5"
1608	8"	15"-17"	3.5"	14.3"	12.5"	9.5"	4"	1.25"	3.5"
1610	10"	15"-17"	3.5"	14.3"	12.5"	11.5"	6"	1.25"	3.5"
1612	12"	15"-17"	3.5"	14.3"	12.5"	13.5"	7"	1.25"	3.5"
1614	14"	15"-17"	3.5"	14.3"	12.5"	15.5"	7.5"	1.25"	3.5"
1616	16"	15"-17"	3.5"	14.3"	12.5"	17.5"	8.25"	1.25"	3.5"
1618	18"	15"-18"	3.5"	14.3"	12.5"	19.5"	8.25"	1.25"	3.5"
1620	20"	15"-18"	3.5"	14.3"	12.5"	21.5"	8.25"	1.25"	3.5"

MORE SIZE SPECIFICATIONS ON THE BACK

INVACARE® MATRIX®
ELITE E2
BACK SPECIFICATIONS

E2S BACK SIZE	A BACK SHELL HEIGHT	B BACK CANE WIDTH	C BACK DEPTH	D BACK SHELL WIDTH	E INSIDE FOAM WIDTH	F BACK FOAM HEIGHT	G APEX OF LATERAL	H FOAM THICKNESS	I BACK SHELL DEPTH
1706	6"	16-18"	3.5"	15.3"	13.5"	7.5"	3"	1.25"	3.5"
1708	8"	16-18"	3.5"	15.3"	13.5"	9.5"	4"	1.25"	3.5"
1710	10"	16-18"	3.5"	15.3"	13.5"	11.5"	6"	1.25"	3.5"
1712	12"	16-18"	3.5"	15.3"	13.5"	13.5"	7"	1.25"	3.5"
1714	14"	16-18"	3.5"	15.3"	13.5"	15.5"	7.5"	1.25"	3.5"
1716	16"	16-18"	3.5"	15.3"	13.5"	17.5"	8.25"	1.25"	3.5"
1718	18"	16-19"	3.5"	15.3"	13.5"	19.5"	8.25"	1.25"	3.5"
1720	20"	16-19"	3.5"	15.3"	13.5"	21.5"	8.25"	1.25"	3.5"
1806	6"	17-19"	3.5"	16.3"	14.5"	7.5"	3"	1.25"	3.5"
1808	8"	17-19"	3.5"	16.3"	14.5"	9.5"	4"	1.25"	3.5"
1810	10"	17-19"	3.5"	16.3"	14.5"	11.5"	6"	1.25"	3.5"
1812	12"	17-19"	3.5"	16.3"	14.5"	13.5"	7"	1.25"	3.5"
1814	14"	17-19"	3.5"	16.3"	14.5"	15.5"	7.5"	1.25"	3.5"
1816	16"	17-19"	3.5"	16.3"	14.5"	17.5"	8.25"	1.25"	3.5"
1818	18"	17-20"	3.5"	16.3"	14.5"	19.5"	8.25"	1.25"	3.5"
1820	20"	17-20"	3.5"	16.3"	14.5"	21.5"	8.25"	1.25"	3.5"
1906	6"	18-20"	3.5"	17.3"	15.5"	7.5"	3"	1.25"	3.5"
1908	8"	18-20"	3.5"	17.3"	15.5"	9.5"	4"	1.25"	3.5"
1910	10"	18-20"	3.5"	17.3"	15.5"	11.5"	6"	1.25"	3.5"
1912	12"	18-20"	3.5"	17.3"	15.5"	13.5"	7"	1.25"	3.5"
1914	14"	18-20"	3.5"	17.3"	15.5"	15.5"	7.5"	1.25"	3.5"
1916	16"	18-20"	3.5"	17.3"	15.5"	17.5"	8.25"	1.25"	3.5"
1918	18"	18-21"	3.5"	17.3"	15.5"	19.5"	8.25"	1.25"	3.5"
1920	20"	18-21"	3.5"	17.3"	15.5"	21.5"	8.25"	1.25"	3.5"
2006	6"	19-21"	3.5"	18.3"	16.5"	7.5"	3"	1.25"	3.5"
2008	8"	19-21"	3.5"	18.3"	16.5"	9.5"	4"	1.25"	3.5"
2010	10"	19-21"	3.5"	18.3"	16.5"	11.5"	6"	1.25"	3.5"
2012	12"	19-21"	3.5"	18.3"	16.5"	13.5"	7"	1.25"	3.5"
2014	14"	19-21"	3.5"	18.3"	16.5"	15.5"	7.5"	1.25"	3.5"
2016	16"	19-21"	3.5"	18.3"	16.5"	17.5"	8.25"	1.25"	3.5"
2018	18"	19-22"	3.5"	18.3"	16.5"	19.5"	8.25"	1.25"	3.5"
2020	20"	19-22"	3.5"	18.3"	16.5"	21.5"	8.25"	1.25"	3.5"

Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Motion Concepts
www.motionconcepts.com

USA
700 Ensminger Rd. Suite 112
Tonawanda, NY 14150
(888) 433-6818
Fax: (888) 433-6834

Canada
84 Citation Dr.,
Unit 1-12
Concord, ON L4K 3C1
(866) 748-7943
(905) 695-0134
Fax: (905) 695-0138

©2020 Motion Concepts.
All rights reserved. Trademarks are
identified by the symbols ™ and ®. All
trademarks are owned by or licensed to
Motion Concepts unless otherwise noted.
Specifications are subject to change
without notification.
Form No. TRD0602 REV A

