




## Armpad Options



### 1. Multi-position arm pads

Soft, gel-like feel provides extra comfort and protection for arms  
Highly durable skinned polyurethane resists cracking and peeling  
Triple track extrusion allows adjustability of mounting position and attachment location for other accessories  
3" wide; desk length 10" or full length 14"

### 2. Ergonomic arm trough/hand pad kit

Molded, skinned polyurethane fits the natural arm curvature  
Unique cut-out at elbow provides users with flexibility of various arm positions  
Triple track extrusion allows a variety of mounting positions  
Complete with flat hand pad or interfaces with other manufacturer's hand pieces

### 3. Gel arm pads

Fabric covered polymer gel arm pads  
Provides extra protection for highly sensitive arms  
Available in 3 models: Std. Narrow: 2" wide x 12" long  
Std. Wide: 3.5" wide x 12" long  
Long Wide: 3.5" wide x 14" long





### 4. Multi-Axis Armrest Mount

Allows infinite adjustment of angle and rotation of arm pad  
Interfaces with multi-position arm pad

### 5. Joystick Mounting Bar for Multi-Axis Armrest Mount

Available for 3/4" or 7/8" joystick bar size

## Head Support Accessories



### 1. Otto Bock™ Head Support

Easy to clean self skinning polyurethane foam  
Anatomically shaped to provide lateral flexion control  
7" wide x 5" high or 9" wide x 6" high

### 2. Standard

Easy to clean self skinning polyurethane foam  
Gentle lateral contouring  
6.5" wide x 7" high

### 3. Stealth Comfort Plus

Soft, comfortable layered foam on coated metal shell  
Removable, washable cover  
Four pad sizes: 8", 10" or 14" wide x 6" high and 6" wide x 3.75" high

### 4. Head Support Kit: Onyx auto style kit

Onyx leatherette cover & Auto-style contouring  
8" wide x 8" high x 4" thick  
Dual pivot ball for pad orientation  
Limited adjustability – recommended for Contour or High backs  
Height adjustable  
Includes fixed mount clamp

### 5. Standard Hardware

Comes standard with Otto Bock™, Standard and Stealth Comfort Plus Pads





### 6. Multi-Axis

Upgrade option for Otto Bock, Standard and Stealth Comfort Plus pads  
Dual pivot linkage for pad orientation  
Depth adjustable through horizontal slide bar  
Height adjustable  
Includes mounting clamp for attachment to back with knob style release for headrest removal

### 7. Stealth "World's Best"

Upgrade option for Otto Bock, Standard and Stealth Comfort Plus Pads  
Dual pivot balls for pad orientation  
Depth adjustable through dual link system  
Height adjustable  
Includes mounting clamp for attachment to back with latch style release for headrest removal

## Other Seating and Positioning Accessories



### 1. Air Mesh Calf Strap

Available in 4 sizes: M, L, MD, LG For 18" – 22" wide chairs  
 For 15" – 18" wide chairs L: 3.5" high  
 M: 3.5" high LG: 9" high  
 MD: 9" high

### 2. Padded Lap Belts

2" wide belts, comfortable padding  
 Center pull tightening, push release buckle  
 2 or 4 point attachment

### 3. Cushion Rigidizer

Lightweight plastic seat insert fits between seat rails for level support and eliminates hammock effect of seat upholstery. Securely attached to bottom of cushion inside cushion cover

### 4. Drop Base

Provides stable support base for the cushion  
 Lowers seat-to-floor height which aids in foot propulsion, opens up seat to back angle, allows easy lateral transfers and improves basic functionality when in a seated position  
 Width adjustable: Three sizes fit wheelchairs from 14" to 22" wide  
 Height and angle adjustable: up to 2" of drop provided  
 Portable quick-release hardware allows easy removal of base from wheelchair

### 5. Basic Headrest

Padded to provide comfortable occipital support with horizontal and vertical adjustable hardware and outer cover easily removes for laundering

### 6. Adjustable Headrests

Contoured head and neck supports available in four styles  
 Dual-ball swivel hardware for multi-plane adjustment  
 Durable steel reinforced polyurethane pads - easily wipe clean

### 7. Full Tray - Slide-on Mount

Break-resistant clear polyester resin tray surface  
 Width adjustable mounting hardware securely fastens tray to chair arm pads  
 Molded-in rim helps keep items on tray

### 8. Full Tray - Hook and Loop Mount

Lightweight and easy to handle with break-resistant clear polycarbonate tray surface  
 Easily attaches to chair armrest with hook and loop straps  
 Removable trim allows thorough cleaning

### 9. Half Tray - Swingaway

Break-resistant clear polycarbonate tray surface  
 Tray securely bolts to wheelchair armrest and is durable  
 Tray easily swings away when not in use and molded-in rim helps keep items on tray

### 10. Foam-In-Place Kit

Easy-to-use Foam-in-Place kit provides custom foam mold insert in minutes  
 Provides additional support and positioning for users who need a more custom fit  
 Custom-molded insert can be easily removed for modification or remolded to meet the changing needs of the user (Can be used with Invacare® Matrx® Personal Back™10, Personal Back™ 10 Plus, PB and PB Deep Backs)

### 11. Belts and Straps

Three sizes of shoulder supports fit child, junior and adult  
 Three sizes of padded lap belts fit child, junior and adult; pinch release buckle  
 Push-button style release provides secure attachment  
 Three sizes of four-point lap belt to fit child, junior and adult; push button release buckle

### 12. Privacy Flaps

The Privacy Flaps interface with the Invacare Matrx PB, Invacare Matrx PB Deep, Invacare Matrx Elite, Invacare Matrx Elite Deep and Invacare Matrx MX1 backs, up to 20" widths. They are designed to cover the gap between the bottom of the backrest and the seat pan and attach with hook and loop fasteners. Four standard sizes are available to fit backs up to 20" wide. Privacy flaps for backs 21" and wider are available as options in the back configurator