

INVACARE® MATRX®

**PS**  
CUSHION

---



---

The **Invacare® Matrx® PS Cushion** is an exceptionally lightweight cushion designed to provide superior positioning, stability and comfort. The design features an anatomically contoured shape with a unique waffled ischial relief recess and an innovative reversible outer cover. It is a great choice for a very wide range of individual needs, including pediatric, geriatric, bariatric and highly active.

---

**Features and Functional Benefits**

**Reversible outer cover** - Innovative reversible Dartex® outer cover may be used smooth or fabric side out. Dartex® fabric is moisture proof and breathable. Hook and loop fastener strips

**Inner cover** - Moisture resistant. Easily wiped clean. Zippered cover allows access for cushion modifications

**Cushion Foam** - Durable high resilient foam prevents "bottoming out". Anatomically molded shape for optimal positioning & comfort. Rear radius for sacral support. Waterfall front edge. Made with Ultra-Fresh™ for antibacterial and odor protection

**Ischial relief area** - Waffled recess provides immersion of bony prominences. Helps maintain pelvic position and prevent sliding

**Custom Modifications** - Leg length discrepancies, width/depth reductions. Foam, inner and outer covers modified to your specifications

**Cushion Ridgizer** - Lightweight, slim profile solid seat insert. Contoured to eliminate "hammock effect" of sling upholstery

**Pelvic Obliquity Kit** - Easily installed on Invacare Matrix PS cushions. Accommodates .5" or 1" obliquities

**HCPCS Code: E2605/E2606**



- 1. Custom Modifications
- 2. Ischial Relief Area
- 3. Cushion Ridgizer
- 4. Pelvic Obliquity Kit



Model #	Width *	Depth**	Height Abductor/Adductor	Height at Leg Trough	Product Weight	Weight Capacity	Limited Warranty
PS (Pediatric)	10", 12", 14"	10", 12", 14", 16"	3.75"/3.5"	2.5"	1 lb. and up	150 lb.	2 Years
PS (Standard)	16", 18", 20"	16", 18", 20"	4"/3.5"	2.5"	2 lb. and up	300 lb.	2 Years
PS (Heavy-Duty)	21" - 30"	16" - 30"	4"/3.5"	2.5"	7 lb. and up	600 lb.	2 Years

\* Actual width is 1/2" narrower than stated width  
 \*\* Actual Depth is 1" longer than stated depth

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program and should consult an attorney or other advisor to discuss specific situations in further detail.